

Single-Element ICP-MS Kit B (9 Elements)

Kit Contains Individual 100 ml Bottles of Each Listed Element

Elements	Concentrations	Elements	Concentrations	Elements	Concentrations
Antimony (Sb)	10 mg/L	Niobium (Nb)	10 mg/L	Tin (Sn)	10 mg/L
Hafnium (Hf)	10 mg/L	Tantalum (Ta)	10 mg/L	Tungsten (W)	10 mg/L
Molybdenum (Mo)	10 mg/L	Tellurium (Te)	10 mg/L	Zirconium (Zr)	10 mg/L

Partnumber	Volume	Matrix	
LK1-04360101	100 ml	2% HNO ₃ + Tr HF	€ 132,00

Crown your labwork