



Single-Element ICP-MS Kit

Single-Element ICP-MS Kit C (7 Elements)

Kit Contains Individual 100 ml Bottles of Each Listed Element

Elements	Concentrations	Elements	Concentrations	Elements	Concentrations
Gold (Au)	10 mg/L	Palladium (Pd)	10 mg/L	Ruthenium (Ru)	10 mg/L
Iridium (Ir)	10 mg/L	Platinum (Pt)	10 mg/L		
Osmium (Os)	10 mg/L	Rhodium (Rh)	10 mg/L		

Partnumber	Volume	Matrix
LK1-04370101	100 ml	2% HCl

Crown your labwork