

Ice Pack Required? ☐

Quote Specification Sheet

Page 1.000

| | | | | |
|---|---------------------------|----------------------------|--|---|
| Catalog No: LK2- CB100815-04 | Storage: ≤ 6 °C | Solvent: Acetone | Container: 2 ml Ampule, Amber Glass | Description: Custom Semivolatile Mix, 20-3627, 100 mg/L, 2 x 2 ml |
|---|---------------------------|----------------------------|--|---|

*Estimated Days to Expiration: 180

| Compound | CAS No | Concentration |
|-----------------------|------------|----------------|
| Bromophos ethyl | 4824-78-6 | 100 +/- 5 mg/L |
| Bromophos methyl | 2104-96-3 | 100 +/- 5 mg/L |
| chlorpyrifos | 2921-88-2 | 100 +/- 5 mg/L |
| coumaphos | 56-72-4 | 100 +/- 5 mg/L |
| diazinon | 333-41-5 | 100 +/- 5 mg/L |
| dichlorvos | 62-73-7 | 100 +/- 5 mg/L |
| dimethoate | 60-51-5 | 100 +/- 5 mg/L |
| disulfoton | 298-04-4 | 100 +/- 5 mg/L |
| ethion | 563-12-2 | 100 +/- 5 mg/L |
| fenitrothion | 122-14-5 | 100 +/- 5 mg/L |
| fenthion | 55-38-9 | 100 +/- 5 mg/L |
| Guthion® | 86-50-0 | 100 +/- 5 mg/L |
| Guthion Ethyl | 2642-71-9 | 100 +/- 5 mg/L |
| malathion | 121-75-5 | 100 +/- 5 mg/L |
| methidathion | 950-37-8 | 100 +/- 5 mg/L |
| methyl parathion | 298-00-0 | 100 +/- 5 mg/L |
| parathion | 56-38-2 | 100 +/- 5 mg/L |
| phosdrin™ (mevinphos) | 7786-34-7 | 100 +/- 5 mg/L |
| pyrazophos | 13457-18-6 | 100 +/- 5 mg/L |

*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval: _____

Date: _____

Quote Specification Sheet

Page 2.000

Catalog No: LK2-CB100815

| Compound | CAS No | Concentration |
|------------|------------|----------------|
| triazophos | 24017-47-8 | 100 +/- 5 mg/L |

*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval: _____ Date: _____

Utrechtseweg 5, 1213TK Hilversum, The Netherlands,
T. +31 84 875 63 44 : F. +31 84 875 63 44 : VAT: NL854948429B01
www.labkings.com