

Ice Pack Required? ☐

Quote Specification Sheet

Page 1.000

Catalog No: LK2- CB100837 -5PCK	Storage: ≤ -10 °C	Solvent: Acetone	Container: 1 mL Ampule, Amber Glass	Description: Custom Semivolatile Mix, 22-3461, 100 mg/L, 5 x 1 ml
---	-----------------------------	----------------------------	--	---

*Estimated Days to Expiration: 385

Compound	CAS No	Concentration
aldrin	309-00-2	100 +/- 5 mg/l
a-BHC	319-84-6	100 +/- 5 mg/l
b-BHC	319-85-7	100 +/- 5 mg/l
d-BHC	319-86-8	100 +/- 5 mg/l
α-BHC	58-89-9	100 +/- 5 mg/l
cis-chlordane	5103-71-9	100 +/- 5 mg/l
o,p'-DDD	53-19-0	100 +/- 5 mg/L
p,p'-DDD	72-54-8	100 +/- 5 mg/l
o,p'-DDE	3424-82-6	100 +/- 5 mg/L
p,p'-DDE	72-55-9	100 +/- 5 mg/l
o,p'-DDT	789-02-6	100 +/- 5 mg/L
p,p'-DDT	50-29-3	100 +/- 5 mg/l
dieldrin	60-57-1	100 +/- 5 mg/l
endosulfan I	959-98-8	100 +/- 5 mg/l
endosulfan sulfate	1031-07-8	100 +/- 5 mg/l
endrin	72-20-8	100 +/- 5 mg/l
heptachlor	76-44-8	100 +/- 5 mg/l
heptachlor epoxide isomer A	28044-83-9	100 +/- 5 mg/l
heptachlor epoxide isomer B	1024-57-3	100 +/- 5 mg/l
isobenzan	297-78-9	100 +/- 5 mg/l
isodrin	465-73-6	100 +/- 5 mg/l
trans-chlordane	5103-74-2	100 +/- 5 mg/l

*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval:

Date: