

Ice Pack Required? ☐

Quote Specification Sheet

Page 1.000

Catalog No:	Storage:	Solvent:	Container:	Description:
LK2- CB100856 -5PCK	≤ -10 °C	P/T Methanol	1 ml Ampule, Amber Glass	Custom Volatile Mix, 40-3376, 8,000 mg/L, 5 x 1 ml

*Estimated Days to Expiration: TBD

Compound	CAS No	Concentration
1.1-Dimethylcyclohexane	590-66-9	8000 +/- 400 mg/L
1.2-diethylbenzene	135-01-3	8000 +/- 400 mg/L
1.3-diethylbenzene	141-93-5	8000 +/- 400 mg/L
acenaphthene	83-32-9	8000 +/- 400 mg/L
allylbenzene	300-57-2	8000 +/- 400 mg/L
benzene	71-43-2	8000 +/- 400 mg/L
sec-butylbenzene	135-98-8	8000 +/- 400 mg/L
cumene	98-82-8	8000 +/- 400 mg/L
cyclohexane	110-82-7	8000 +/- 400 mg/L
cyclopentane	287-92-3	8000 +/- 400 mg/L
n-decane (C10)	124-18-5	8000 +/- 400 mg/L
2.3-dimethylbutane	79-29-8	8000 +/- 400 mg/L
2.2-dimethylpentane	590-35-2	8000 +/- 400 mg/L
dodecane (C12)	112-40-3	8000 +/- 400 mg/L
ethylbenzene	100-41-4	8000 +/- 400 mg/L
2-ethyltoluene	611-14-3	8000 +/- 400 mg/L
3-ethyltoluene	620-14-4	8000 +/- 400 mg/L
4-ethyltoluene	622-96-8	8000 +/- 400 mg/L
heptane (C7)	142-82-5	8000 +/- 400 mg/L
n-hexane (C6)	110-54-3	8000 +/- 400 mg/L
2-methylbutane	78-78-4	8000 +/- 400 mg/L
methyl cyclohexane	108-87-2	8000 +/- 400 mg/L
methylcyclopentane	96-37-7	8000 +/- 400 mg/L

*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval:

Date:

Quote Specification Sheet

Page 2.000

Catalog No: LK2-CB100856

Compound	CAS No	Concentration			
2-methylpentane	107-83-5	8000	+/-	400	mg/L
3-methylpentane	96-14-0	8000	+/-	400	mg/L
a-methylstyrene	98-83-9	8000	+/-	400	mg/L
naphthalene	91-20-3	8000	+/-	400	mg/L
n-hexadecane (C16)	544-76-3	8000	+/-	400	mg/L
nonane (C9)	111-84-2	8000	+/-	400	mg/L
octane (C8)	111-65-9	8000	+/-	400	mg/L
n-pentane (C5)	109-66-0	8000	+/-	400	mg/L
1,2,3-trimethylbenzene	526-73-8	8000	+/-	400	mg/L
styrene	100-42-5	8000	+/-	400	mg/L
toluene	108-88-3	8000	+/-	400	mg/L
trans-1,3-pentadiene (stabilized with TBC)	2004-70-8	8000	+/-	400	mg/L
1,2,4-trimethylbenzene	95-63-6	8000	+/-	400	mg/L
1,3,5-trimethylbenzene	108-67-8	8000	+/-	400	mg/L
m-xylene	108-38-3	8000	+/-	400	mg/L
o-xylene	95-47-6	8000	+/-	400	mg/L
p-xylene	106-42-3	8000	+/-	400	mg/L

*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval:

Date:

Utrechtseweg 5, 1213TK Hilversum, The Netherlands,
T. +31 84 875 63 44 : F. +31 84 875 63 44 : VAT: NL854948429B01
www.labking.com