

Ice Pack Required? ☐

Quote Specification Sheet

Page 1.000

Catalog No: LK2- CB100857	Storage: 2°C - 8°C	Solvent: Methanol	Container: 1 ml Ampule, Amber Glass	Description: Custom Organotin Mix, 4-1, 1000 mg/L, 5 x 1 ml
--	------------------------------	-----------------------------	--	---

*Estimated Days to Expiration: 730

Compound	CAS No	Concentration
n-butyltin trichloride as n-butyltin	1118-46-3	1000 +/- 50 mg/L
di-n-butyltin dichloride as di-n-butyltin	683-18-1	1000 +/- 50 mg/L
tri-n-butyltin chloride as tri-n-butyltin	1461-22-9	1000 +/- 50 mg/L
triophenyltin chloride as triphenyltin	639-58-7	1000 +/- 50 mg/L

*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval:

Date: