

Ice Pack Required? ☒

Quote Specification Sheet

Page 1.000

Catalog No:	Storage:	Solvent:	Container:	Description:
LK2- CB232803-SS -10PCK	≤ -10 °C	P/T Methanol	1 mL Ampule, Amber Glass	Custom Volatile Mix, Second Source, 45-1282, 2,000 mg/L, 10 x 1 ml

Days to Expiration: 730

Compound	CAS No	Concentration
bromotrichloromethane	75-62-7	2000 +/- 100 mg/L
2,3-dichloro-1-propene	78-88-6	2000 +/- 100 mg/L
3-chlorotoluene	108-41-8	2000 +/- 100 mg/L
allyl chloride	107-05-1	2000 +/- 100 mg/L
tert-amyl methyl ether (TAME)	994-05-8	2000 +/- 100 mg/L
bromobenzene	108-86-1	2000 +/- 100 mg/L
bromochloromethane	74-97-5	2000 +/- 100 mg/L
bromodichloromethane	75-27-4	2000 +/- 100 mg/L
bromomethane	74-83-9	2000 +/- 100 mg/L
n-butylbenzene	104-51-8	2000 +/- 100 mg/L
sec-butylbenzene	135-98-8	2000 +/- 100 mg/L
tert-butylbenzene	98-06-6	2000 +/- 100 mg/L
chloroethane	75-00-3	2000 +/- 100 mg/L
chloromethane	74-87-3	2000 +/- 100 mg/L
chloroprene	126-99-8	2000 +/- 100 mg/L
2-chlorotoluene	95-49-8	2000 +/- 100 mg/L
4-chlorotoluene	106-43-4	2000 +/- 100 mg/L
cyclohexane	110-82-7	2000 +/- 100 mg/L
n-decane (C10)	124-18-5	2000 +/- 100 mg/L

We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval: _____

Date: _____

Quote Specification Sheet

Page 2.000

Catalog No: LK2-

Compound	CAS No	Concentration
dibromochloromethane	124-48-1	2000 +/- 100 mg/L
1,2-dibromo-3-chloropropane	96-12-8	2000 +/- 100 mg/L
1,2-dibromoethane	106-93-4	2000 +/- 100 mg/L
dibromomethane	74-95-3	2000 +/- 100 mg/L
dichlorodifluoromethane	75-71-8	2000 +/- 100 mg/L
2,2-dichloropropane	594-20-7	2000 +/- 100 mg/L
1,1-dichloropropylene	563-58-6	2000 +/- 100 mg/L
cis-1,3-dichloropropylene	10061-01-5	2000 +/- 100 mg/L
trans-1,3-dichloropropylene	10061-02-6	2000 +/- 100 mg/L
dodecane (C12)	112-40-3	2000 +/- 100 mg/L
fluorobenzene	462-06-6	2000 +/- 100 mg/L
heptane (C7)	142-82-5	2000 +/- 100 mg/L
hexachlorobutadiene	87-68-3	2000 +/- 100 mg/L
n-hexane (C6)	110-54-3	2000 +/- 100 mg/L
isooctane	540-84-1	2000 +/- 100 mg/L
isopropyl ether	108-20-3	2000 +/- 100 mg/L
4-isopropyltoluene	99-87-6	2000 +/- 100 mg/L
a-methylstyrene	98-83-9	2000 +/- 100 mg/L
nonane (C9)	111-84-2	2000 +/- 100 mg/L
octane (C8)	111-65-9	2000 +/- 100 mg/L
pentachloroethane	76-01-7	2000 +/- 100 mg/L
n-pentane (C5)	109-66-0	2000 +/- 100 mg/L
1,1,1,2-tetrachloroethane	630-20-6	2000 +/- 100 mg/L
1,1,2,2-tetrachloroethane	79-34-5	2000 +/- 100 mg/L
trichlorofluoromethane	75-69-4	2000 +/- 100 mg/L
1,2,3-trichloropropane	96-18-4	2000 +/- 100 mg/L
n-undecane (C11)	1120-21-4	2000 +/- 100 mg/L

We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval: _____

Date: _____

Tappersweg 59, 2031 ET Haarlem, The Netherlands,
T. +31 84 875 63 44 : F. +31 84 875 63 44 : VAT: NL854948429B01
www.labkings.com