

Ice Pack Required? ☐

Quote Specification Sheet

Page 1.000

Catalog No:	Storage:	Solvent:	Container:	Description:
LK2- CB232806 -5PCK	≤ -10 °C	Acetone	1 mL Ampule, Amber Glass	Custom OCB Pesticide Mix, 38-2438, 100 mg/L, 5 x 1 ml

Days to Expiration:550

Compound	CAS No	Concentration
heptachlor epoxide isomer A	28044-83-9	100 +/- 5 mg/L
1,2,3,5-tetrachlorobenzene	634-90-2	100 +/- 5 mg/L
aldrin	309-00-2	100 +/- 5 mg/L
a-BHC	319-84-6	100 +/- 5 mg/L
b-BHC	319-85-7	100 +/- 5 mg/L
d-BHC	319-86-8	100 +/- 5 mg/L
g-BHC	58-89-9	100 +/- 5 mg/L
cis-chlordane	5103-71-9	100 +/- 5 mg/L
trans-chlordane	5103-74-2	100 +/- 5 mg/L
1-chloronaphthalene	90-13-1	100 +/- 5 mg/L
2-chloronaphthalene	91-58-7	100 +/- 5 mg/L
o,p'-DDD	53-19-0	100 +/- 5 mg/L
p,p'-DDD	72-54-8	100 +/- 5 mg/L
o,p'-DDE	3424-82-6	100 +/- 5 mg/L
p,p'-DDE	72-55-9	100 +/- 5 mg/L
o,p'-DDT	789-02-6	100 +/- 5 mg/L
p,p'-DDT	50-29-3	100 +/- 5 mg/L
dieldrin	60-57-1	100 +/- 5 mg/L
endosulfan I	959-98-8	100 +/- 5 mg/L

We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval: _____

Date: _____

Quote Specification Sheet

Page 2.000

Catalog No: LK2-CB232806

Compound	CAS No	Concentration
endosulfan II	33213-65-9	100 +/- 5 mg/L
endosulfan sulfate	1031-07-8	100 +/- 5 mg/L
endrin	72-20-8	100 +/- 5 mg/L
heptachlor	76-44-8	100 +/- 5 mg/L
heptachlor epoxide isomer B	1024-57-3	100 +/- 5 mg/L
hexachlorobenzene	118-74-1	100 +/- 5 mg/L
hexachlorobutadiene	87-68-3	100 +/- 5 mg/L
isobenzan	297-78-9	100 +/- 5 mg/L
isodrin	465-73-6	100 +/- 5 mg/L
methoxychlor	72-43-5	100 +/- 5 mg/L
pentachlorobenzene	608-93-5	100 +/- 5 mg/L
pentachloronitrobenzene	82-68-8	100 +/- 5 mg/L
propachlor	1918-16-7	100 +/- 5 mg/L
tecnazene	117-18-0	100 +/- 5 mg/L
1,2,4,5-tetrachlorobenzene	95-94-3	100 +/- 5 mg/L
1,2,3,4-tetrachlorobenzene	634-66-2	100 +/- 5 mg/L
1,2,3-trichlorobenzene	87-61-6	100 +/- 5 mg/L
1,2,4-trichlorobenzene	120-82-1	100 +/- 5 mg/L
1,3,5-trichlorobenzene	108-70-3	100 +/- 5 mg/L

We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval: _____

Date: _____

Tappersweg 59, 2031 ET Haarlem, The Netherlands,
T. +31 84 875 63 44 : F. +31 84 875 63 44 : VAT: NL854948429B01
www.labkings.com