

Ice Pack Required? ☐

## Quote Specification Sheet

Page 1.000

<b>Catalog No:</b> LK2- CB257406	<b>Storage:</b> ≤ -10 °C	<b>Solvent:</b> P/T Methanol	<b>Container:</b> 2 mL Ampule, Amber Glass	<b>Description:</b> Custom Volatile Mixture, 12-3692, 1000 mg/L, 2 mL
--	-----------------------------	---------------------------------	--	--

\*Estimated Days to Expiration: TBD

Compound	CAS No	Concentration			
benzene	71-43-2	1000	+/-	50	mg/L
bromodichloromethane	75-27-4	1000	+/-	50	mg/L
bromoform	75-25-2	1000	+/-	50	mg/L
carbon tetrachloride	56-23-5	1000	+/-	50	mg/L
chloroform	67-66-3	1000	+/-	50	mg/L
cis-1.2-dichloroethylene	156-59-2	1000	+/-	50	mg/L
dibromochloromethane	124-48-1	1000	+/-	50	mg/L
1.2-dichloroethane	107-06-2	1000	+/-	50	mg/L
methylene chloride	75-09-2	1000	+/-	50	mg/L
tetrachloroethylene	127-18-4	1000	+/-	50	mg/L
1.1.1-trichloroethane	71-55-6	1000	+/-	50	mg/L
trichloroethylene	79-01-6	1000	+/-	50	mg/L

\*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval:

Date: