

Ice Pack Required? ☐

Quote Specification Sheet

Page 1.000

Catalog No: LK3- CB100805 -5PCK	Storage: ≤ 6 °C	Solvent: CS2:Methanol 9:1	Container: 1 ml Ampule, Amber Glass	Description: Custom VOC Mix, 68-1780, 4000:8000 mg/L, 5 x 1 mL
---	---------------------------	-------------------------------------	--	--

Days to Expiration:180

Compound	CAS No	Concentration
allyl chloride	107-05-1	8000 +/- 400 mg/L
benzene	71-43-2	8000 +/- 400 mg/L
bromochloromethane	74-97-5	8000 +/- 400 mg/L
bromodichloromethane	75-27-4	8000 +/- 400 mg/L
bromoform	75-25-2	8000 +/- 400 mg/L
tert-butyl ethyl ether (ETBE)	637-92-3	8000 +/- 400 mg/L
carbon tetrachloride	56-23-5	8000 +/- 400 mg/L
chlorobenzene	108-90-7	8000 +/- 400 mg/L
chloroethane	75-00-3	8000 +/- 400 mg/L
chloroform	67-66-3	8000 +/- 400 mg/L
chloromethane	74-87-3	8000 +/- 400 mg/L
2-chloropropane	75-29-6	8000 +/- 400 mg/L
2-chlorotoluene	95-49-8	8000 +/- 400 mg/L
3-chlorotoluene	108-41-8	8000 +/- 400 mg/L
4-chlorotoluene	106-43-4	8000 +/- 400 mg/L
cis-1,2-dichloroethylene	156-59-2	8000 +/- 400 mg/L
cumene	98-82-8	8000 +/- 400 mg/L
cyclohexane	110-82-7	8000 +/- 400 mg/L
cyclopentane	287-92-3	8000 +/- 400 mg/L

We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval: _____

Date: _____

Quote Specification Sheet

Page 2.000

Catalog No: LK3-CB100805

Compound	CAS No	Concentration
n-decane (C10)	124-18-5	8000 +/- 400 mg/L
dibromochloromethane	124-48-1	8000 +/- 400 mg/L
1,2-dibromoethane	106-93-4	8000 +/- 400 mg/L
dibromomethane	74-95-3	8000 +/- 400 mg/L
1,2-dichlorobenzene	95-50-1	8000 +/- 400 mg/L
1,3-dichlorobenzene	541-73-1	8000 +/- 400 mg/L
1,4-dichlorobenzene	106-46-7	8000 +/- 400 mg/L
1,1-dichloroethane	75-34-3	8000 +/- 400 mg/L
1,2-dichloroethane	107-06-2	8000 +/- 400 mg/L
1,1-dichloroethylene	75-35-4	8000 +/- 400 mg/L
trans-1,2-dichloroethylene	156-60-5	8000 +/- 400 mg/L
1,2-dichloropropane	78-87-5	8000 +/- 400 mg/L
1,3-dichloropropane	142-28-9	8000 +/- 400 mg/L
2,3-dichloro-1-propene	78-88-6	8000 +/- 400 mg/L
cis-1,3-dichloropropylene	10061-01-5	8000 +/- 400 mg/L
trans-1,3-dichloropropylene	10061-02-6	8000 +/- 400 mg/L
dodecane (C12)	112-40-3	8000 +/- 400 mg/L
ethylbenzene	100-41-4	8000 +/- 400 mg/L
2-ethyltoluene	611-14-3	8000 +/- 400 mg/L
3-ethyltoluene	620-14-4	8000 +/- 400 mg/L
4-ethyltoluene	622-96-8	8000 +/- 400 mg/L
heptane (C7)	142-82-5	8000 +/- 400 mg/L
hexachloroethane	67-72-1	8000 +/- 400 mg/L
n-hexane (C6)	110-54-3	8000 +/- 400 mg/L
isooctane	540-84-1	8000 +/- 400 mg/L
4-isopropyltoluene	99-87-6	8000 +/- 400 mg/L
methyl cyclohexane	108-87-2	8000 +/- 400 mg/L
methylcyclopentane	96-37-7	8000 +/- 400 mg/L
methylene chloride	75-09-2	8000 +/- 400 mg/L
2-methylpentane	107-83-5	8000 +/- 400 mg/L

We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval: _____

Date: _____

Utrechtseweg 5, 1213TK Hilversum, The Netherlands,
T. +31 84 875 63 44 : F. +31 84 875 63 44 : VAT: NL854948429B01
www.labkings.com

Quote Specification Sheet

Page 3.000

Catalog No: LK3-CB100805

Compound	CAS No	Concentration
3-methylpentane	96-14-0	8000 +/- 400 mg/L
a-methylstyrene	98-83-9	8000 +/- 400 mg/L
methyl t-butyl ether	1634-04-4	8000 +/- 400 mg/L
naphthalene	91-20-3	8000 +/- 400 mg/L
nonane (C9)	111-84-2	8000 +/- 400 mg/L
octane (C8)	111-65-9	8000 +/- 400 mg/L
pentachloroethane	76-01-7	8000 +/- 400 mg/L
n-pentane (C5)	109-66-0	8000 +/- 400 mg/L
n-propylbenzene	103-65-1	8000 +/- 400 mg/L
1,2,3-trimethylbenzene	526-73-8	8000 +/- 400 mg/L
styrene	100-42-5	8000 +/- 400 mg/L
1,1,1,2-tetrachloroethane	630-20-6	8000 +/- 400 mg/L
1,1,2,2-tetrachloroethane	79-34-5	8000 +/- 400 mg/L
tetrachloroethylene	127-18-4	8000 +/- 400 mg/L
toluene	108-88-3	8000 +/- 400 mg/L
1,2,3-trichlorobenzene	87-61-6	8000 +/- 400 mg/L
1,2,4-trichlorobenzene	120-82-1	8000 +/- 400 mg/L
1,3,5-trichlorobenzene	108-70-3	8000 +/- 400 mg/L
1,1,1-trichloroethane	71-55-6	8000 +/- 400 mg/L
1,1,2-trichloroethane	79-00-5	8000 +/- 400 mg/L
trichloroethylene	79-01-6	8000 +/- 400 mg/L
1,2,3-trichloropropane	96-18-4	8000 +/- 400 mg/L
1,2,4-trimethylbenzene	95-63-6	8000 +/- 400 mg/L
1,3,5-trimethylbenzene	108-67-8	8000 +/- 400 mg/L
n-undecane (C11)	1120-21-4	8000 +/- 400 mg/L
vinyl chloride	75-01-4	8000 +/- 400 mg/L
m-xylene	108-38-3	8000 +/- 400 mg/L
o-xylene	95-47-6	8000 +/- 400 mg/L
p-xylene	106-42-3	8000 +/- 400 mg/L

We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval: _____

Date: _____

Utrechtseweg 5, 1213TK Hilversum, The Netherlands,
T. +31 84 875 63 44 : F. +31 84 875 63 44 : VAT: NL854948429B01
www.labkings.com