

Ice Pack Required?

## Quote Specification Sheet

Page 1.000

<b>Catalog No:</b> LK2- CB101835 -100	<b>Storage:</b> ≤ -10 °C	<b>Solvent:</b> Acetone	<b>Container:</b> 1 mL Ampule, Amber Glass	<b>Description:</b> Custom Pesticide Mixture 16-8883 100 mg/L 1 ml
------------------------------------------------	-----------------------------	----------------------------	--------------------------------------------------	-----------------------------------------------------------------------

\*Estimated Days to Expiration: 385

Compound	CAS No	Concentration
azoxvstrobin	131860-33-8	100 +/- 5 mg/L
bupirimate	41483-43-6	100 +/- 5 mg/L
2-chloro-4.6-diamino-1.3.5-triazine	3397-62-4	100 +/- 5 mg/L
chlorovrifos-methyl	5598-13-0	100 +/- 5 mg/L
Command (clomazone)	81777-89-1	100 +/- 5 mg/L
cvcloxdim	101205-02-1	100 +/- 5 mg/L
2.4-dichlorophenoxyacetic acid	94-75-7	100 +/- 5 mg/L
fluopicolide	239110-15-7	100 +/- 5 mg/L
imazamox	114311-32-9	100 +/- 5 mg/L
linuron	330-55-2	100 +/- 5 mg/L
MCPA acid	94-74-6	100 +/- 5 mg/L
penconazol	66246-88-6	100 +/- 5 mg/L
Propiconazol (mixture of isomers)	60207-90-1	100 +/- 5 mg/L
sulcotrione	99105-77-8	100 +/- 5 mg/L
Terbumeton	33693-04-8	100 +/- 5 mg/L
triclopr	55335-06-3	100 +/- 5 mg/L

\*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval:

Date: