

Ice Pack Required? ☐

Quote Specification Sheet

Page 1.000

Catalog No:	Storage:	Solvent:	Container:	Description:
LK3- CB101801 -2PCK	≤ -10 °C	Acetone	1 mL Ampule, Amber Glass	Custom Pesticide Mix, 45-8165, 1000 mg/L, 2 x 1 ml

*Estimated Days to Expiration: 730

Compound	CAS No	Concentration
alachlor	15972-60-8	1000 +/- 50 mg/l
a-BHC	319-84-6	1000 +/- 50 mg/l
b-BHC	319-85-7	1000 +/- 50 mg/l
d-BHC	319-86-8	1000 +/- 50 mg/l
α-BHC	58-89-9	1000 +/- 50 mg/l
cis-chlordane	5103-71-9	1000 +/- 50 mg/l
chlorfenvinphos (mixture of z and e isomers)	470-90-6	1000 +/- 50 mg/l
chlorpyrifos	2921-88-2	1000 +/- 50 mg/l
o,p'-DDD	53-19-0	1000 +/- 50 mg/l
p,p'-DDD	72-54-8	1000 +/- 50 mg/l
o,p'-DDE	3424-82-6	1000 +/- 50 mg/l
p,p'-DDE	72-55-9	1000 +/- 50 mg/l
o,p'-DDT	789-02-6	1000 +/- 50 mg/l
p,p'-DDT	50-29-3	1000 +/- 50 mg/l
desethyl atrazine	6190-65-4	1000 +/- 50 mg/l
diazinon	333-41-5	1000 +/- 50 mg/l
2,6-dichlorobenzamide	2008-58-4	1000 +/- 50 mg/l
dieldrin	60-57-1	1000 +/- 50 mg/l
dimepiperate	61432-55-1	1000 +/- 50 mg/l
dimethenamid	87674-68-8	1000 +/- 50 mg/l
endosulfan I	959-98-8	1000 +/- 50 mg/l
endosulfan II	33213-65-9	1000 +/- 50 mg/l
endrin	72-20-8	1000 +/- 50 mg/l

Compounds will precipitate at slightly less than room temperature. Warm solution and agitate gently in ampule, checking for complete solubility prior to use.

*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval:

Date:

Quote Specification Sheet

Page 2.000

Catalog No: LK3-CB101801

Compound	CAS No	Concentration			
heptachlor	76-44-8	1000	+/-	50	mg/l
heptachlor epoxide isomer B	1024-57-3	1000	+/-	50	mg/l
hexazinon	51235-04-2	1000	+/-	50	mg/l
isodrin	465-73-6	1000	+/-	50	mg/l
kepone	143-50-0	1000	+/-	50	mg/l
metalaxyl	57837-19-1	1000	+/-	50	mg/l
metolachlor	51218-45-2	1000	+/-	50	mg/l
mirex	2385-85-5	1000	+/-	50	mg/l
molinate	2212-67-1	1000	+/-	50	mg/l
oxadiazon	19666-30-9	1000	+/-	50	mg/l
oxadixyl	77732-09-3	1000	+/-	50	mg/l
pendimethalin	40487-42-1	1000	+/-	50	mg/l
Pirimicarb	23103-98-2	1000	+/-	50	mg/l
pretilachlor	51218-49-6	1000	+/-	50	mg/l
Propanil	709-98-8	1000	+/-	50	mg/l
terbuthylazine-desethyl	30125-63-4	1000	+/-	50	mg/l
thiobencarb	28249-77-6	1000	+/-	50	mg/l
Thiocarbazyl	36756-79-3	1000	+/-	50	mg/l
trans-chlordane	5103-74-2	1000	+/-	50	mg/l
trifluralin	1582-09-8	1000	+/-	50	mg/l
vinclozolin	50471-44-8	1000	+/-	50	mg/l
aldrin	309-00-2	1000	+/-	50	mg/l

Compounds will precipitate at slightly less than room temperature. Warm solution and agitate gently in ampule, checking for complete solubility prior to use.

*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval:

Date:

Utrechtseweg 5, 1213TK Hilversum, The Netherlands,
T. +31 84 875 63 44 : F. +31 84 875 63 44 : VAT: NL854948429B01
www.labking.com