



LabKings B.V.
Wateringweg 17
2031 EK Haarlem
The Netherlands

T. +31 84 875 63 44
F. +31 84 875 63 44
info@labkings.com
www.labkings.com

Partnumber	LK1-CB724515-500	Total	Bottling	Price per
Volume	500ml	500 mL	1 x 500ml	€ 96,00
Name Mix	ICP-MS Setup/Stability/Masscal Solution			
Shelf life	18 months			
Matrix	0.5% HCl			

Element	Conc.µg/L
1 Ba	10
2 Al	1
3 Cd	1
4 Ce	1
5 Cr	1
6 Cu	1
7 In	1
8 Mg	1
9 Mn	1
10 Pb	1
11 Rh	1
12 Th	1
13	
14	
15	
16	
17	
18	
19	
20	
21	
22	
23	
24	
25	
26	
27	
28	
29	
30	
31	
32	
33	
34	
35	
36	
37	
38	
39	
40	

BIC: ABNANL2A

IBAN: NL69ABNAO474142884

VAT: NL854948429B01

Crown your labwork