

Partnumber	LK1-CB724517-500	Total	Bottling	Price per
Volume	500ml	500 mL	1 x 500ml	€ 138,00
Name Mix	ICP-MS Setup Solution			
Shelf life	18 months			
Matrix	1% HNO ₃ , 500ml			
Element	Conc.µg/L			
1 Be	1			
2 Ce	1			
3 Fe	1			
4 In	1			
5 Li	1			
6 Mg	1			
7 Pb	1			
8 U	1			
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				
31				
32				
33				
34				
35				
36				
37				
38				
39				
40				

BIC: ABNANL2A

IBAN: NL69ABNAO474142884

VAT: NL854948429B01

Crown your labwork