



LabKings B.V.
Wateringweg 17
2031 EK Haarlem
The Netherlands

T. +31 84 875 63 44
F. +31 84 875 63 44
info@labkings.com
www.labkings.com

Partnumber	LK1-CB724519-250	Total	Bottling	Price per
Volume	250ml	250 mL	1 x 250ml	€ 178,00
Name Mix	ICP-MS Stability Solution			
Shelf life	18 months			
Matrix	2% HNO3			

Element	Conc.µg/L
1 Al	200
2 Ba	200
3 Ce	200
4 Co	200
5 Cu	200
6 In	200
7 Li	200
8 Mg	200
9 Mn	200
10 Ni	200
11 Pb	200
12 Tb	200
13 U	200
14 Zn	200
15	
16	
17	
18	
19	
20	
21	
22	
23	
24	
25	
26	
27	
28	
29	
30	
31	
32	
33	
34	
35	
36	
37	
38	
39	
40	

BIC: ABNANL2A

IBAN: NL69ABNAO474142884

VAT: NL854948429B01

Crown your labwork