



LabKings B.V.
Wateringweg 17
2031 EK Haarlem
The Netherlands

T. +31 84 875 63 44
F. +31 84 875 63 44
info@labkings.com
www.labkings.com

Partnumber	LK1-CB724520-500	Total	Bottling	Price per
Volume	500ml	500 mL	1 x 500ml	€ 166,00
Name Mix	ICP-MS Stability Solution			
Shelf life	18 months			
Matrix	1% HNO3			
Element		Conc.µg/L		
1	Co	10		
2	Cu	10		
3	In	10		
4	Se	10		
5	Cd	1		
6	Cr	1		
7	Fe	1		
8	Mg	1		
9	Pb	1		
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				
31				
32				
33				
34				
35				
36				
37				
38				
39				
40				

BIC: ABNANL2A

IBAN: NL69ABNAO474142884

VAT: NL854948429B01

Crown your labwork