

Item	Model	Description	Cat. No.	"Unit Price (Euro) "																																				
Pipette																																								
1	Autoclavable Digital Micro Pipette Nichipet Premium LT	New evolution of "Fatigue Free" Lighter model of Nichipet Premium LT! Enhanced UV resistance.																																						
		<table><tr><th>Volume Range (ul)</th><th>Increments (ul)</th><th></th><th></th></tr><tr><td>0.1-2</td><td>0.002</td><td>LK9-00NLT-2</td><td>€ 199,75</td></tr><tr><td>0.5-10</td><td>0.01</td><td>LK9-00NLT-10</td><td>€ 176,00</td></tr><tr><td>2-20</td><td>0.02</td><td>LK9-00NLT-20</td><td>€ 176,00</td></tr><tr><td>10-100</td><td>0.1</td><td>LK9-00NLT-100</td><td>€ 176,00</td></tr><tr><td>20-200</td><td>0.2</td><td>LK9-00NLT-200</td><td>€ 176,00</td></tr><tr><td>100-1000</td><td>1.0</td><td>LK9-00NLT-1000</td><td>€ 176,00</td></tr></table>	Volume Range (ul)	Increments (ul)			0.1-2	0.002	LK9-00NLT-2	€ 199,75	0.5-10	0.01	LK9-00NLT-10	€ 176,00	2-20	0.02	LK9-00NLT-20	€ 176,00	10-100	0.1	LK9-00NLT-100	€ 176,00	20-200	0.2	LK9-00NLT-200	€ 176,00	100-1000	1.0	LK9-00NLT-1000	€ 176,00										
Volume Range (ul)	Increments (ul)																																							
0.1-2	0.002	LK9-00NLT-2	€ 199,75																																					
0.5-10	0.01	LK9-00NLT-10	€ 176,00																																					
2-20	0.02	LK9-00NLT-20	€ 176,00																																					
10-100	0.1	LK9-00NLT-100	€ 176,00																																					
20-200	0.2	LK9-00NLT-200	€ 176,00																																					
100-1000	1.0	LK9-00NLT-1000	€ 176,00																																					
		*If white body color selected, please put (white) at the end of Cat. No.																																						
		One instrument with one piece of appropriate tip and grease for maintenance																																						
Pipette Starter Package		Most usable 3 Nichipet Premium LT models package.																																						
S-PAC-A	CE	NLT-10 , NLT-100, NLT-1000 (1 set each)	LK9-00SP-AL	€ 485,00																																				
		BMT2-SSR, BMT2-SGR, BMT2-LGR (1 rack each) NTO (Tube Opener) 3pcs																																						
S-PAC-B	CE	NLT-20, NLT-200, NLT-1000 (1 set each)	LK9-00SP-BL	€ 485,00																																				
		BMT2-SGRx2, BMT2-LGR x1 rack NTO (Tube Opener) 3pcs																																						
S-PAC-H	CE	NLT-2, NLT-10, NLT-100 (1 set each)	LK9-00SP-HL	€ 510,00																																				
		BMT2-UTWRx2, BMT2-SSR, BMT2-SGR(1 rack each) NTO (Tube Opener) 3pcs																																						
		*If white body color selected, please put (white) at the end of Cat. No.																																						
2	Autoclavable Digital Micro Pipette Nichipet Premium	Fatigue-free operation with enlarged finger rest. Enhanced UV resistance.																																						
		<table><tr><th>Volume Range (ul)</th><th>Increments (ul)</th><th></th><th></th></tr><tr><td>0.1-2</td><td>0.002</td><td>LK9-00NPP-2</td><td>€ 195,00</td></tr><tr><td>0.5-10</td><td>0.01</td><td>LK9-00NPP-10</td><td>€ 165,00</td></tr><tr><td>2-20</td><td>0.02</td><td>LK9-00NPP-20</td><td>€ 165,00</td></tr><tr><td>10-100</td><td>0.1</td><td>LK9-00NPP-100</td><td>€ 165,00</td></tr><tr><td>20-200</td><td>0.2</td><td>LK9-00NPP-200</td><td>€ 165,00</td></tr><tr><td>100-1000</td><td>1.0</td><td>LK9-00NPP-1000</td><td>€ 165,00</td></tr><tr><td>1000-5000</td><td>10.0</td><td>LK9-00NPP-5000</td><td>€ 195,00</td></tr><tr><td>1000-10000</td><td>10.0</td><td>LK9-00NPP-10000</td><td>€ 195,00</td></tr></table>	Volume Range (ul)	Increments (ul)			0.1-2	0.002	LK9-00NPP-2	€ 195,00	0.5-10	0.01	LK9-00NPP-10	€ 165,00	2-20	0.02	LK9-00NPP-20	€ 165,00	10-100	0.1	LK9-00NPP-100	€ 165,00	20-200	0.2	LK9-00NPP-200	€ 165,00	100-1000	1.0	LK9-00NPP-1000	€ 165,00	1000-5000	10.0	LK9-00NPP-5000	€ 195,00	1000-10000	10.0	LK9-00NPP-10000	€ 195,00		
Volume Range (ul)	Increments (ul)																																							
0.1-2	0.002	LK9-00NPP-2	€ 195,00																																					
0.5-10	0.01	LK9-00NPP-10	€ 165,00																																					
2-20	0.02	LK9-00NPP-20	€ 165,00																																					
10-100	0.1	LK9-00NPP-100	€ 165,00																																					
20-200	0.2	LK9-00NPP-200	€ 165,00																																					
100-1000	1.0	LK9-00NPP-1000	€ 165,00																																					
1000-5000	10.0	LK9-00NPP-5000	€ 195,00																																					
1000-10000	10.0	LK9-00NPP-10000	€ 195,00																																					
		One instrument with one piece of appropriate tip.																																						
		Tube Opener (3pcs./Bag)	LK9-00NTO-001	€ 15,00																																				
Pipette Starter Package		Most usable 3 Nichipet Premium models package.																																						
S-PAC-A	CE	NPP-10, NPP-100, NPP-1000 (1 set each)	LK9-00SP-AP	€ 455,00																																				
		BMT2-SSR, BMT2-SGR, BMT2-LGR (1 rack each) NTO (Tube Opener)																																						
S-PAC-B	CE	NPP-20, NPP-200, NPP-1000 (1 set each)	LK9-00SP-BP	€ 460,00																																				
		BMT2-SGRx2, BMT2-LGR (1 rack each) NTO (Tube Opener)																																						
S-PAC-H	CE	NPP-2, NPP-10, NPP-100 (1 set each)	LK9-00SP-HP	€ 482,00																																				
		BMT2-UTRx2, BMT2-SSR, BMT2-SGR(1 rack each) NTO (Tube Opener)																																						

3 Nichipet EXII

CE

Enhanced UV resistance.



Volume Range (ul)	Increments (ul)		
0.1-2	0.002	LK9-00NPX2-2	€ 143,00
0.5-10	0.01	LK9-00NPX2-10	€ 122,00
2-20	0.02	LK9-00NPX2-20	€ 122,00
10-100	0.1	LK9-00NPX2-100	€ 122,00
20-200	0.2	LK9-00NPX2-200	€ 122,00
100-1000	1.0	LK9-00NPX2-1000	€ 122,00
1000-5000	10.0	LK9-00NPX2-5000	€ 142,00
1000-10000	10.0	LK9-00NPX2-10000	€ 142,00

One instrument with one piece of appropriate tip.

Pipette Starter Package

S-PAC-A

CE

Most usable 3 Nichipet EXII models package.
NPX2-10, NPX2-100, NPX2-1000 (1 set each)
BMT2-SSR, BMT2-SGR, BMT2-LGR (1 rack each)
NTO (Tube Opener)

LK9-00SP-A2

€ 360,00

S-PAC-B

CE

NPX2-20, NPX2-200, NPX2-1000 (1 set each)
BMT2-SGR×2, BMT2-LGR
NTO (Tube Opener)

LK9-00SPBN2

€ 360,00

S-PAC-H

CE

NPX2-2, NPX2-10, NPX2-100 (1 set each)
BMT2-UTR×2, BMT2-SSR, BMT2-SGR
NTO (Tube Opener)

LK9-00SPH2

€ 386,00

4 A Value Added Nichipet EXII

Nichipet EX PlusII

CE

New mechanism added Nichipet EXII with popular features.
Advanced durability against organic solvent dispensing.
Installing Perfluoro rubber O-ring as standard option.



Volume Range (ul)		
0.1-2	LK9-00NPLO2-2	€ 166,00
0.5-10	LK9-00NPLO2-10	€ 158,00
2-20	LK9-00NPLO2-20	€ 158,00
10-100	LK9-00NPLO2-100	€ 158,00
20-200	LK9-00NPLO2-200	€ 158,00
100-1000	LK9-00NPLO2-1000	€ 158,00
1000-5000	LK9-00NPLO2-5000	€ 166,00
1000-10000	LK9-00NPLO2-10000	€ 166,00

One instrument with one piece of appropriate tip.

Pipette Stand

Rotary Stand: 6 pipettes can be stored.
6 pipettes can be stored

LK9-00MLT-STD2

€ 34,00

LK9-00BM-STD

€ 98,88

5 Light & Ergonomic Pipette

Nichiryo Le

CE

Suitable light weight as a tool for the daily work.
Less probability against RSI (repetitive strain injury) eventually.



Volume Range (ul)		
0.5-10	LK9-00NLE-10	€ 79,00
2-20	LK9-00NLE-20	€ 79,00
10-100	LK9-00NLE-100	€ 79,00
20-200	LK9-00NLE-200	€ 79,00
100-1000	LK9-00NLE-1000	€ 79,00

Pipette Holder

Stand to store one pipette

LK9-00NLE-STD

€ 7,50

6 Autoclavable Multi Channel Pipette

Nichipet EXII MULTI

CE

Four sizes cover a wide volume range (0.5-300ul).
Head revolves 360° to allow for comfortable dispensing.
Lighter model of Nichipet EXII MULTI



8-Channel

Volume Range (ul)	Increments (ul)		
0.5-10	0.01	LK9-00NPM-8V	€ 439,00
5-100	0.1	LK9-00NPM-8S	€ 389,00
20-200	0.2	LK9-00NPM-8L	€ 389,00
30-300	1.0	LK9-00NPM-8K	€ 439,00

12-Channel

Volume Range (ul)	Increments (ul)		
0.5-10	0.01	LK9-00NPM-12V	€ 550,00
5-100	0.1	LK9-00NPM-12S	€ 526,00
20-200	0.2	LK9-00NPM-12L	€ 526,00
30-300	1.0	LK9-00NPM-12K	€ 550,00

Stand 2 pipettes can be stored

LK9-007000S

€ 111,00

Tray Polypropylene Tray to pool reagent (50pcs./box)

LK9-0088-TY

€ 94,00

7 Triple Volume Micro Pipette
Nichipet V II

Built-in tip ejector being available in three different models.
Simply turn the thumb knob to change the volume of the

Volume Setting (ul)

10/20/50	LK9-00NPV2-S	€ 142,00
50/100/200	LK9-00NPV2-M	€ 142,00
200/500/1000	LK9-00NPV2-L	€ 142,00

One instrument with one piece of appropriate tip.

8 Fixed Volume
Nichipet F II

Smooth movement of the plunger and light weight.
9 different standard volumes ranging from 10ul to 1000ul.

LK9-00NPF2-10	€ 113,00
LK9-00NPF2-20	€ 113,00
LK9-00NPF2-25	€ 113,00
LK9-00NPF2-30	€ 113,00
LK9-00NPF2-50	€ 113,00
LK9-00NPF2-100	€ 113,00
LK9-00NPF2-200	€ 113,00
LK9-00NPF2-500	€ 113,00
LK9-00NPF2-1000	€ 113,00

9 Disposable Tips
Premium Tip

Universal fit.
5ML and 10ML tips designed to fit Nichiryo pipettes exclusively.

Bulk Tips

Volume Range (ul)	Color	Q'ty/Bag		
0.1-10	clear	1000	LK9-00BMT2-UT	€ 46,00
0.5-10	clear	1000	LK9-00BMT2-SS	€ 60,00
2-200	clear	1000	LK9-00BMT2-SE	€ 46,00
2-200	clear	1000	LK9-00BMT2-SG	€ 46,00
50-300	clear	1000	LK9-00BMT2-K	€ 46,00
100-1000	clear	1000	LK9-00BMT2-LE	€ 46,00
100-1000	clear	1000	LK9-00BMT2-LG	€ 46,00
1000-5000	clear	200	LK9-00BMT2-X	€ 30,00
1000-10000	clear	200	LK9-00BMT2-Z	€ 36,00

Racked Tips

Volume Range (ul)	Color	Q'ty/Bag		
0.1-10	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-UTWR	€ 76,00
0.5-10	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-SSR	€ 99,63
2-200	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-SGR	€ 76,00
2-200	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-SER	€ 76,00
50-300	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-KR	€ 76,00
100-1000	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-LGR	€ 76,00
100-1000	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-LER	€ 76,00
1000-5000	clear	100 (50pcs.x2 ^c / _s)	LK9-00BMT2-XR	€ 44,00
1000-10000	clear	80 (40pcs.x2 ^c / _s)	LK9-00BMT2-ZR	€ 46,00

Racked Tips (sterilized)

Volume Range (ul)	Color	Q'ty/Bag		
0.1-10	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-UTWRB	€ 94,63
0.5-10	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-SSRB	€ 99,87
2-200	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-SGRB	€ 90,00
2-200	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-SERB	€ 90,00
50-300	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-KRB	€ 90,00
100-1000	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-LGRB	€ 90,00
100-1000	clear	960 (96pcs.x10 ^c / _s)	LK9-00BMT2-LERB	€ 90,00
1000-5000	clear	100 (50pcs.x2 ^c / _s)	LK9-00BMT2-XRB	€ 50,00
1000-10000	clear	80 (40pcs.x2 ^c / _s)	LK9-00BMT2-ZRB	€ 57,00

Rack only

for BMT2-UT	5pcs	LK9-00ETS-UTRK	€ 21,00
for BMT2-SE	5pcs	LK9-00ETS-SERK	€ 21,00
for BMT2-SG	5pcs	LK9-00ETS-SGRK	€ 21,00
for BMT2-LE	5pcs	LK9-00ETS-LERK	€ 21,00
for BMT2-LG	5pcs	LK9-00ETS-LGRK	€ 21,00
for BMT2-X	2pcs	LK9-00RK-X	€ 21,00
for BMT2-Z	2pcs	LK9-00RK-Z	€ 21,00

9-1 Filtered Disposable Tips
Premium Filter Tip

Universal fit. Ideal for molecular and microbiology applications.
Prevents transferring volatile or hazardous samples.

Bulk Tips

Volume Range (ul)	Color	Q'ty/Bag		
0.1-10	clear	1000	LK9-00BMT2-FUT	€ 97,00
0.5-10	clear	1000	LK9-00BMT2-FSS	€ 97,00
2-20	clear	1000	LK9-00BMT2-FSG	€ 97,00
10-100	clear	1000	LK9-00BMT2-FSE	€ 97,00
20-200	clear	1000	LK9-00BMT2-FAG	€ 97,00
100-1000	clear	1000	LK9-00BMT2-FLG	€ 126,00
100-1000	clear	1000	LK9-00BMT2-FLE	€ 126,00

Racked Tips (sterilized)

Volume Range (ul)	Color	Q'ty/Bag		
0.1-10	clear	960 (96pcs.x10 ⁶ /s)	LK9-00BMT2-FUTRB	€ 138,00
0.5-10	clear	960 (96pcs.x10 ⁶ /s)	LK9-00BMT2-FSSRB	€ 138,00
2-20	clear	960 (96pcs.x10 ⁶ /s)	LK9-00BMT2-FSGRB	€ 138,00
10-100	clear	960 (96pcs.x10 ⁶ /s)	LK9-00BMT2-FSERB	€ 138,00
20-200	clear	960 (96pcs.x10 ⁶ /s)	LK9-00BMT2-FAGRB	€ 166,00
100-1000	clear	1000 (100pcs.x10 ⁶ /s)	LK9-00BMT2-FLGRB	€ 207,00
100-1000	clear	1000 (100pcs.x10 ⁶ /s)	LK9-00BMT2-FLERB	€ 207,00

10 **Economy-Tip**

Bulk Tips

0.1-10	clear	1000	LK9-00ETS-UT	€ 14,00
2-200	clear	1000	LK9-00ETS-SG	€ 14,00
2-200	clear	1000	LK9-00ETS-SE	€ 14,00
100-1000	clear	1000	LK9-00ETS-LG	€ 17,00
100-1000	clear	1000	LK9-00ETS-LE	€ 17,00

10-1 **Economy-Filter Tip**

DNA Rnase/Dnase, 120°C Autoclavable, Non-Pyrogenic Endotoxin,
Low Cost Trace Metals, Medical Grade Resin and Packing

Racked Tips (sterilized)

Volume Range (ul)	Color	Q'ty/Bag		
0.1-10	clear	960 (96pcs.x10 ⁶ /s)	LK9-00ETS-FUTRB	€ 82,00
0.5-10	clear	960 (96pcs.x10 ⁶ /s)	LK9-00ETS-FSSRB	€ 90,00
2-20	clear	960 (96pcs.x10 ⁶ /s)	LK9-00ETS-FSGRB	€ 90,00
10-100	clear	960 (96pcs.x10 ⁶ /s)	LK9-00ETS-FSERB	€ 90,00
20-200	clear	960 (96pcs.x10 ⁶ /s)	LK9-00ETS-FAGRB	€ 90,00
100-1000	clear	576 (96pcs.x6 ⁷ /s)	LK9-00ETS-FLGRB	€ 60,00

11 Replaceable Rack Tips

Smart Pack

Easy & quick tip setting. Eco-friendly & low-waste.

Volume Range (ul)	Color	Q'ty/Bag		
0.1-10	clear	960pcs(96pcs×10decks)	LK9-00ETS-JUTW	€ 40,00
0.1-10	clear	9600pcs (96pcs×10decks×10boxes)	LK9-00ETS-JUTW10	€ 390,00
2-200	clear	960pcs(96pcs×5decks×2boxes)	LK9-00ETS-JSG	€ 30,00
2-200	clear	9600pcs (96pcs×5decks×2×10boxes	LK9-00ETS-JSG10	€ 290,00
100-1000	clear	960pcs (96pcs×5decks×2)	LK9-00ETS-JLG	€ 36,00
100-1000	clear	9600pcs (96pcs×5decks×2×10boxes	LK9-00ETS-JLG10	€ 339,00

12 Glass Tip Pipette
Nichipet Eco

Combining the accuracy and precision of the micropipette
Glass tips: the most reliable and stable material for this research.

Volume Range (ul)	Increments (ul)		
20-200	0.2	LK9-00NPC-200	€ 228,00
100-1000	1.0	LK9-00NPC-1000	€ 228,00
1000-5000	10.0	LK9-00NPC-5000	€ 228,00
1000-10000	50.0	LK9-00NPC-10ML	€ 228,00

Comes with a PTFE nozzle and ceramic piston.

Glass Tips

Volume Range (ul)	Q'ty/Box		
20-200	10	LK9-00GLT-S	€ 84,00
100-1000	10	LK9-00GLT-L	€ 84,00
1000-5000	10	LK9-00GLT-X	€ 84,00
1000-10000	10	LK9-00GLT-Z	€ 84,00

13 Variable Volume Macro Pipette
MODEL 3100

A unique macro pipette with excellent accuracy and precision.
Easy operation: any volume up to 10ml is achieved with one stroke.

Volume Range (ml)	Increments (ml)		
0.1-1.0	0.1	LK9-0031N-1000	€ 173,00
0.4-2.0	0.1	LK9-0031N-2000	€ 173,00
1.0-5.0	0.5	LK9-0031N-5000	€ 173,00
2.0-10.0	0.5	LK9-0031N-10000	€ 173,00

One instrument with one appropriate tip, one piston installation tool, and silicone grease

Disposable Tips

Volume Range (ml)	Color	Q'ty/Bag		
0.1-1.0	red	1000	LK9-00MPT-1	€ 52,00
0.4-2.0	yellow	200	LK9-00MPT-2	€ 33,00
1.0-5.0	green	200	LK9-00BMT2-X	€ 33,00
2.0-10.0	clear	200	LK9-00BMT2-Z	€ 40,00

14 Ultra Macro Pipette
MODEL 5100

Two models, each with three standard volume settings.
Volume change done by placing a volume space in the pipette.

Volume Setting (ml)		
15.0-25.0	LK9-005100S	€ 250,00
30.0-50.0	LK9-005100L	€ 304,00

One instrument with 3 disposable tips

Disposable Tips

Capacity (ml)	Q'ty/Bag		
25.0	5	LK9-0051TS	€ 16,00
50.0	5	LK9-0051TL	€ 16,00

Dispenser

15 Ergonomically Designed Repetitive Dispenser
NICHIMATE STEPPER

Up to 40 different volumes from 1ul to 5000ul available. LK9-00NSTP

Syringes	Capacity (ml)	Q'ty/Box		
	0.05	100	LK9-00NS-V	€ 172,00
	0.5	100	LK9-00NS-S	€ 84,00
	1.25	100	LK9-00NS-M	€ 72,00
	2.5	100	LK9-00NS-N	€ 72,00
	5.0	100	LK9-00NS-L	€ 72,00
	12.5	100	LK9-00NS-X	€ 72,00
	25.0	50	LK9-00NS-Y	€ 63,00
	50.0	25	LK9-00NS-Z	€ 60,00

*One box of Y and Z syringes come with one Adaptor: NS-ADP

Sterilized Syringes	Capacity (ml)	Q'ty/Box		
	0.5	100	LK9-00NS-SB	€ 128,00
	1.25	100	LK9-00NS-MB	€ 128,00
	2.5	100	LK9-00NS-NB	€ 128,00
	5.0	100	LK9-00NS-LB	€ 128,00
	12.5l	100	LK9-00NS-XB	€ 128,00

Adaptor	Attachment for 25.0ml Syringes(00-NS-Y)	LK9-00NS-ADPY	€ 14,00
	Attachment for 50.0ml Syringes(00-NS-Z)	LK9-00NS-ADPZ	€ 14,00

Stand	Stands to hold two dispensers	LK9-00NS-STD	€ 102,00
-------	-------------------------------	--------------	----------

16 Repetitive Syringe Dispenser
MODEL 8100

35 ways of dispensing in the volume range between 1 and 5000l LK9-008100 € 205,00
Fast change of syringe and volume adjustment dial

Syringe	Volume (ml)	Q'ty/Bag		
	0.06	100	LK9-00SG-V2	€ 103,00
	0.6	100	LK9-00SG-S ^{*1}	€ 80,00
	1.5	100	LK9-00SG-M ^{*2}	€ 80,00
	3.0	100	LK9-00SG-N3	€ 80,00
	6.0	100	LK9-00SG-L	€ 80,00
	15.0	50	LK9-00SG-Y	€ 52,00
	60.0	50	LK9-00SG-Z ^{*3}	€ 99,76

^{*1,*2} SG-S and SG-M should be equipped with a Nichiryo pipette tip (type TS) on a nozzle.

^{*3} SG-Z with an adaptor to mount the syringe to the handle

Adapter	Adaptor for SG-Z: 60ml syringe	LK9-00SG-ZA	€ 5,00
---------	--------------------------------	-------------	--------

Stand	Stand to hold two dispensers with syringes	LK9-008100S	€ 105,00
-------	--	-------------	----------

17	8-Channel Repetitive Syringe Dispenser MODEL 8800	Syringe type 8-channel dispenser	LK9-008800	€ 430,00
		Syringe 520ul syringe with a silicone greased O-ring, autoclavable (100pcs)	LK9-0088-SAR	€ 80,00
		1300ul syringe with a silicone greased O-ring, autoclavable (100pcs)	LK9-0088-SBR	€ 80,00
	Accessories	Cartridge A for 520ul syringe, autoclavable (10pcs./box)	LK9-0088-CA	€ 29,00
		Cartridge B for 1300ul syringe, autoclavable (10pcs./box)	LK9-0088-CB	€ 29,00
		Tray to pool reagent (50 pcs./box)	LK9-0088-TY	€ 93,00
	Stand	Stand to hold two dispensers	LK9-008800S	€ 80,00
18	Digital Hand-held Dispenser MODEL 2100DG	A digitally volume reading handy dispenser Dispensing from most liquid container is available.		
		<u>Volume Range (ml)</u> <u>Increments (ml)</u>		
		0.2-1.0 0.01	LK9-0021DG-1	€ 214,00
		0.4-2.0 0.01	LK9-0021DG-2	€ 214,00
		1.0-5.0 0.02	LK9-0021DG-5	€ 214,00
		2.0-10.0 0.1	LK9-0021DG-10	€ 214,00
19	Pipette Controller Pipette Mate NEO CE	The new model following Pipette Mate Pro Charger available in Euro, UK, US, and Japanese versions	LK9-00PMNEO	€ 290,00
	Parts and disposables	Head Fitting	LK9-00PM-N005	€ 52,00
		Tube Holder(Silicone Pipette Holder)	LK9-00PM-N013	€ 31,00
		Filter set (0.20um) <Clear> (5pcs.)	LK9-00PM-N100	€ 56,00
		*Please contact us about the price for parts and disposables.		
20	Volume Bottle Dispenser ACCUPENSER JUNIOR	The most accurate dispenser Five models cover a range from 0.2ml to 30.0ml.		
		<u>Volume Range (ml)</u> <u>Increments (ml)</u>		
		0.2-1.0 0.05	LK9-00AJ-1	€ 183,00
		0.4-2.0 0.1	LK9-00AJ-2	€ 183,00
		1.0-5.0 0.1	LK9-00AJ-5	€ 183,00
		2.0-10.0 0.25	LK9-00AJ-10	€ 183,00
		One instrument with one amber glass bottle, one funnel, one silicone rubber, and extension dispensing tube with glass tip		
21	Variable Volume Bottle-top Dispenser DISPET	High quality dispenser		
		<u>Volume Range (ml)</u> <u>Increments (ml)</u>		
		0.2-1.0 0.05	LK9-00DP-1	€ 123,00
		0.4-2.0 0.1	LK9-00DP-2	€ 123,00
		1.0-5.0 0.1	LK9-00DP-5	€ 123,00
		2.0-10.0 0.25	LK9-00DP-10	€ 123,00
		*If you need reagent bottles, please indicate "B" after the cat. No. ex) DP-1B		
	Adapter	<u>Size</u>		
		18mm	LK9-00DP-1118	€ 2,75
		19mm	LK9-00DP-1119	€ 2,75
		22mm	LK9-00DP-1122	€ 2,75
		24mm	LK9-00DP-1124	€ 2,75
		25mm	LK9-00DP-1125	€ 2,75
		28mm	LK9-00DP-1128	€ 2,75
		33mm	LK9-00DP-1133	€ 2,75
		38mm	LK9-00DP-1138	€ 2,75
		45mm	LK9-00DP-1145	€ 2,75
22	Bottle-top Dispenser DISPET EX II	Autoclavable bottle-top dispenser Liquid return function in order to minimize reagent loss.		
		<u>Volume Range (ml)</u> <u>Increments (ml)</u>		
		0.2-2.0 0.05	LK9-00DPX2-20	€ 285,00
		0.5-5.0 0.10	LK9-00DPX2-50	€ 285,00
		1.0-10.0 0.20	LK9-00DPX2-100	€ 285,00
		2.5-25.0 0.50	LK9-00DPX2-250	€ 375,00
		5.0-50.0 1.00	LK9-00DPX2-500	€ 405,00
		10.0-100.0 2.00	LK9-00DPX2-1000	€ 520,00
	Adapter	<u>Size</u>		
		28mm	LK9-00DPX-003A	€ 20,00
		32mm	LK9-00DPX-003B	€ 20,00
		38mm	LK9-00DPX-003F ^{*1}	€ 20,00
		38mm	LK9-00DPX-003C ^{*2}	€ 20,00
		40mm	LK9-00DPX-003D	€ 20,00
		45mm	LK9-00DPX-003E	€ 20,00
		^{*1} For DPX2-25, 50, 100		
		^{*2} For DPX2-250, 500, 1000		

23 Bottle-top Dispenser
DISPET Mini

Bottle-top dispenser
Liquid return function in order to minimize reagent loss.

Volume Range (ul)

100-250	LK9-00DPM-250	€ 194,00
500-1000	LK9-00DPM-1000	€ 194,00
1000-2000	LK9-00DPM-2000	€ 194,00

Centrifuge

24 **Nichimate Flash** Mini Centrifuge LK9-00NCF € 334,00

25 **Centrifuge Tube** Centrifugal strength: Up to 20000G (High strength model)
Volume (ml) Q'ty/Bag
15 400pcs(50pics×8bags) LK9-00ETS-CT-15 € 90,00
50 200pcs(20pics×10bags) LK9-00ETS-CT-50 € 56,00

26 **Micro Tube** Reducing contamination during one hand operation
Volume (ml) Q'ty/Bag
1.5 1,000pcs LK900-ETS-MC-15 € 19,00
2.0 1,000pcs LK900-ETS-MC-20 € 25,00

BIC: ABNANL2A

IBAN: NL69ABNA0474142884

VAT: NL854948429B01

Crown your labwork