

Partnumber	LK1-CB724526-100	Total	Bottling	Price per
Volume	100ml	100 mL	1 x 100ml	€ 140,00
Name Mix	ICP-OES ICAL Solution			
Shelf life	18 months			
Matrix	2% HNO3 / 2% HCl			

Element	Conc.µg/mL
1 S	50
2 Ce	10
3 Cu	10
4 Eu	10
5 Fe	10
6 In	10
7 K	10
8 Ni	10
9 P	10
10 Si	10
11 Ti	10
12 V	10
13 Y	10
14 Zr	10
15 Mn	5
16 Mo	5
17 Na	5
18 Sc	5
19 Be	2
20 Li	2
21 Sr	2
22 Ca	1
23	
24	
25	
26	
27	
28	
29	
30	
31	
32	
33	
34	
35	
36	
37	
38	
39	
40	

BIC: ABNANL2A

IBAN: NL69ABNAO474142884

VAT: NL854948429B01

Crown your labwork