

Date Received: _____

Certificate of Analysis Rev 0 Page 1

Catalog No.: LK3-CB528508 **Lot No.:** 276017 **Storage:** ≤ -10 °C **Solvent:** Acetone **Exp. Date:** 5-Apr-2018 **Description:** Custom Pesticide Mix, 51-2738, 1000 mg/L, 2 x 1 ml -2PCK

Compound	CAS No.	Purity (%)	Neat Material Lot No.	Concentration
heptachlor	76-44-8	99.8	261.3.4P	1008 +/- 5.43 mg/L
endosulfan I	959-98-8	99.5	266.4.2.3P	985.1 +/- 5.39 mg/L
mirex	2385-85-5	99	352.4.2.2P	1005 +/- 7.3 mg/L
o,p'-DDE	3424-82-6	99.5	396.421.1P	995 +/- 7.27 mg/L
o,p'-DDT	789-02-6	99	395.3.15P	999.8 +/- 7.28 mg/L
thiobencarb	28249-77-6	99.9	641.8.3P	994.1 +/- 11.27 mg/L
vinclozolin	50471-44-8	99.5	1065.3.1P	985.1 +/- 5.84 mg/L
2,6-dichlorobenzamide	2008-58-4	98.4	2743.7.1P	984 +/- 11.15 mg/L
cis-chlordane	5103-71-9	99.5	264.4.3.2P	990.1 +/- 5.4 mg/L
trans-chlordane	5103-74-2	99.5	265.4.2.4P	1005 +/- 5.41 mg/L
diazinon	333-41-5	99	578.3.10P	999.8 +/- 7.28 mg/L
metalaxyl	57837-19-1	98.6	1055.286.1P	990.9 +/- 11.22 mg/L
oxadixyl	77732-09-3	99.2	1651.3.5.1P	992 +/- 7.25 mg/L
pendimethalin	40487-42-1	99	1254.3.1.1P	985.1 +/- 7.22 mg/L
Pirimicarb	23103-98-2	98.5	1636.8.1.1P	990.9 +/- 16.22 mg/L
Propanil	709-98-8	99.6	645.8.1.1P	991.1 +/- 11.24 mg/L
pretilachlor	51218-49-6	97	2495.7.1P	994 +/- 5.29 mg/L
dimepiperate	61432-55-1	99.6	2494.8.2P	986.1 +/- 11.19 mg/L
alachlor	15972-60-8	99	552.8.3P	988 +/- 12.73 mg/L
aldrin	309-00-2	99.2	258.3.5P	987.1 +/- 7.23 mg/L
ametryne	834-12-8	98.9	913.3.3P	986.9 +/- 16.49 mg/L
atrazine	1912-24-9	99	337.7.2P	990 +/- 5.37 mg/L
a-BHC	319-84-6	99.3	237.3.3P	997.9 +/- 5.39 mg/L

Certified By: _____

Madeline Luby
Madeline Luby

All weights are traceable through N. I. S. T. Test No. 822/264157-00.
Concentration (correct for purity) and uncertainty (95% confidence) values listed are determined gravimetrically.

This standard was manufactured by a laboratory accredited to ISO/IEC 17025:2005 (certificate number 3031.01) by the American Association of Laboratory Accreditation (A2LA). The manufacturer's quality system is audited and registered by NSF-ISR to ISO 9001:2008 (certificate number IZ391-IS4).

Certificate of Analysis

Page 2


Catalog No.: LK3-CB528508-2PCK

Lot No.: 276017

Expiration Date: 5-Apr-2018

Compound	CAS No.	Purity (%)	Neat Material Lot No.	Concentration	
b-BHC	319-85-7	98.5	259.8.2.3P	994.8 +/- 5.35	mg/L
g-BHC	58-89-9	99	238.1.2.1P	999.8 +/- 11.31	mg/L
chlorfenvinphos (mixture of z and E isomers)	470-90-6	95	925.286.1.1P	992.3 +/- 11.15	mg/L
p,p'-DDD	72-54-8	97	234.1.1P	989.2 +/- 15.72	mg/L
p,p'-DDE	72-55-9	99.5	236.7.1.1P	995 +/- 5.4	mg/L
p,p'-DDT	50-29-3	98.3	233.1.1P	997.6 +/- 11.27	mg/L
desethyl atrazine	6190-65-4	99.5	1017.8.1.1.1P	990.1 +/- 11.23	mg/L
dieldrin	60-57-1	91.8	240.1.1P	999.8 +/- 11.16	mg/L
endosulfan II	33213-65-9	99.6	253.4.2.6P	1006 +/- 5.42	mg/L
endrin	72-20-8	97.5	255.3.11P	1004 +/- 7.24	mg/L
heptachlor epoxide isomer B	1024-57-3	99.5	262.3.5P	1005 +/- 5.41	mg/L
hexazinon	51235-04-2	99.5	1176.3.3P	999.9 +/- 11.32	mg/L
isodrin	465-73-6	100	242.1.2.1P	1025 +/- 11.56	mg/L
metolachlor	51218-45-2	97.1	909.8.4.1P	976.2 +/- 15.96	mg/L
molinate	2212-67-1	99.3	1005.3.11.3P	997.9 +/- 11.3	mg/L
oxadiazon	19666-30-9	99.5	561.286.1P	990.1 +/- 11.23	mg/L
chlorpyrifos	2921-88-2	99.75	543.3.3P	997.5 +/- 5.41	mg/L
simazine	122-34-9	98.7	916.3.3P	987 +/- 7.21	mg/L
terbuthylazine	5915-41-3	98.8	684.8.2.2P	986 +/- 16.14	mg/L
terbuthylazine-desethyl	30125-63-4	99.6	1642.201.1.1.1P	991.1 +/- 7.26	mg/L
Thiocarbaryl	36756-79-3	98	1674.421.1P	989.7 +/- 7.2	mg/L
kepone	143-50-0	98.5	285.282.1.1P	994.8 +/- 5.35	mg/L
d-BHC	319-86-8	98.6	260.4.5.1P	995.8 +/- 5.36	mg/L
trifluralin	1582-09-8	99.1	335.8.1.1P	986.1 +/- 11.18	mg/L
dimethenamid	87674-68-8	99.6	2037.8.2.1P	991.1 +/- 11.24	mg/L
prometryn	7287-19-6	99.5	329.3.3P	999.9 +/- 11.32	mg/L
propazine	139-40-2	98.8	334.3.4P	983.1 +/- 5.36	mg/L
o,p'-DDD	53-19-0	99.8	235.115.2.1P	993.1 +/- 5.41	mg/L

Certified By:



Madeline Luby

All weights are traceable through N. I. S. T. Test No. 822/264157-00.
Concentration (correct for purity) and uncertainty (95% confidence) values listed are determined gravimetrically.

This standard was manufactured by a laboratory accredited to ISO/IEC 17025:2005 (certificate number 3031.01) by the American Association of Laboratory Accreditation (A2LA). The manufacturer's quality system is audited and registered by NSF-ISR to ISO 9001:2008 (certificate number IZ391-IS4).