

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-031402	PCB Congener Mix 2-02	decachlorobiphenyl	2051-24-3	100	Isooctane	1 ml, Amber Glass	€ 103,50
		2,4,6-trichlorobiphenyl (BZ# 30)	35693-92-6	100			
LK3-001104	Aroclor 1016 + 1260 MDL Solution, 100 ug/L				Methanol	1 ml, Amber Glass	€ 55,50
LK3-001107	OEXT PCB Spike, 10 mg/L , 2-1,				Hexane	25 ml, BR, Amber Glass	€ 55,50
LK3-001108	Aroclor 1016 + 1260 Solution, 1,000 mg/L				Isooctane	1 ml, Amber Glass	€ 67,50
LK3-001111	Aroclor 1016 + 1260 Solution, 10,000 mg/L				P/T Methanol	1 ml, Amber Glass	€ 55,50
		Aroclor 1016	12674-11-2				
		Aroclor 1260	11096-82-5				
LK3-011201	PCB Congener Mix (BZ #'s 103, 2,2',3,3',4,5,5',6-octachlorobiphenyl (BZ# 198)		68194-17-2	80	Isooctane	1 ml, Amber Glass	€ 63,00
		2,2',4,5',6-pentachlorobiphenyl (BZ# 103)	60145-21-3	80			
		2,2',3,4,4',5,6,6'-octachlorobiphenyl (BZ# 204)	74472-52-9	80			
LK3-030301	PCB Congener Mix 3-1	2-chlorobiphenyl (BZ# 1)	2051-60-7	100	Isooctane	1 ml, Amber Glass	€ 55,50
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	100			
		2,2',3,3',4,5,5',6'-octachlorobiphenyl (BZ# 210)	52663-75-9	100			
LK3-225001	PCB Intermediate Standard, 3	Aroclor 1016	12674-11-2	10	Hexane	25 ml, BR, Amber Glass	€ 55,50
		Aroclor 1260	11096-82-5	10			
		decachlorobiphenyl (BZ 209)	2051-24-3	1			
LK3-011703-5pck	Custom PCB Mix, 4	2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	100	Isooctane	1 ml, Amber Glass	€ 355,50
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
		2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	100			
		2,3,3',5,5'-pentachlorobiphenyl (BZ# 111)	39635-32-0	100			
LK3-103302	Custom PCB Conginer Mix, 4-1	2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	1000	Isooctane	1 ml, Amber Glass	€ 40,50
		2,2',3,3',4,4',5,6'-octachlorobiphenyl (BZ# 196)	42740-50-1	1000			
		2,2',3,4,4',5,5',6-octachlorobiphenyl (BZ# 203)	52663-76-0	1000			
		2,2',3,3',4,5,5',6,6'-nonachlorobiphenyl (BZ# 208)	52663-77-1	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-102302	Custom PCB Mix, 10 mg/L ,	4-13			Isooctane	1 ml, Amber Glass	€ 55,50
LK3-102301	Custom PCB Mix, 100 mg/L ,	4-13				1 ml, Amber Glass	€ 72,00
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3				
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7				
		2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4				
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9				
LK3-031805	PCB Mix, 4-15				Isooctane	1 ml, Amber Glass	€ 48,00
LK3-031806	PCB Mix, 4-15				Nonane	1 ml, Amber Glass	€ 48,00
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	10			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	10			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	10			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	10			
LK3-247807	Custom PCB Mix, 5-34				Isooctane	1.2 ml, Amber Glass	€ 48,00
LK3-247807-5pck	Custom PCB Mix, 5-34					1.2 ml, Amber Glass	€ 198,00
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	2			
		3,3',4,5'-tetrachlorobiphenyl (BZ# 79)	41464-48-6	2			
		2,2',3,4',6-pentachlorobiphenyl (BZ# 91)	68194-05-8	2			
		2,2',3,3',4,5,6'-heptachlorobiphenyl (BZ# 174)	38411-25-5	2			
		2,2',3,3',4,4',5,6'-octachlorobiphenyl (BZ# 196)	42740-50-1	2			
LK3-050101	PCB Congener Mix, 10 mg/L ,	6-1			Isooctane	1 ml, Amber Glass	€ 72,00
LK3-050102	PCB Congener Mix, 100 mg/L ,	6-1			Acetone	1 ml, Amber Glass	€ 118,50
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8				
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8				
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4				
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6				
		2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1				
		2,2',3,4,4',6,6'-heptachlorobiphenyl (BZ# 184)	74472-48-3				

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-107501	PCB Kit, Aroclor, 200 mg/L				Methanol	1 ml each, Amber	€ 135,00
LK3-107505	PCB Kit, Aroclor, 1,000 mg/L				Isooctane	1 ml each, Amber	€ 189,00
		Aroclor 1016	12674-11-2				
		Aroclor 1232	11141-16-5				
		Aroclor 1248	12672-29-6				
		Aroclor 1260	11096-82-5				
		Aroclor 1221	11104-28-2				
		Aroclor 1242	53469-21-9				
		Aroclor 1254	11097-69-1				
LK3-016301	PCB Congener Mix-7, 100 mg/L	2,3,3',4,4',6-hexachlorobiphenyl (BZ# 158)	74472-42-7	100	Isooctane	1 ml, Amber Glass	€ 118,50
		2,3,4,4',5,6-hexachlorobiphenyl (BZ# 166)	41411-63-6	100			
		2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1	100			
		2,2',3,4,4',6,6'-heptachlorobiphenyl (BZ# 184)	74472-48-3	100			
		2,3,4,4'-tetrachlorobiphenyl (BZ# 60)	33025-41-1	100			
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	100			
		2,2',3,4',5-pentachlorobiphenyl (BZ# 90)	68194-07-0	100			
LK3-243001	Method 525, 525.1 PCB Mix, 500 mg/L	1,2-chlorobiphenyl (BZ# 1)	2051-60-7	500	Hexane	1 ml, Amber Glass	€ 237,00
		2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	500			
		2,2',3,3',4,4',6-heptachlorobiphenyl (BZ# 171)	52663-71-5	500			
		2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	500			
		2,2',3,3',4,5',6'-octachlorobiphenyl (BZ# 200)	40186-71-8	500			
		2,2',3',4,6-pentachlorobiphenyl (BZ# 98)	60233-25-2	500			
		2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	500			
LK3-031605	PCB Congener Mix, 10 mg/L , 7-14				Isooctane	1 ml, Amber Glass	€ 31,50
LK3-031603	PCB Congener Mix, 500 mg/L , 7-14					1 ml, Amber Glass	€ 1.221,00
LK3-031601	PCB Congener Mix, 1,000 mg/L , 7-14					1 ml, Amber Glass	€ 394,50
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5				
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3				
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2				

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6				
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2				
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1				
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3				
LK3-031302	PCB Congener Mix 7-8	2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100	Isooctane	1 ml, Amber Glass	€ 142,50
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	100			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	100			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	100			
LK3-031705	PCB Mix, 8-14	2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	10	Isooctane	1 ml, Amber Glass	€ 55,50
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	10			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	10			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	10			
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	10			
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	10			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	10			
LK3-304005	EPA 525, 525.1 PCB Mix, 8-65				Hexane	1 ml, Amber Glass	€ 135,00
LK3-304006	EPA 525, 525.1 PCB Mix, 8-65				Acetone	1 ml, Amber Glass	€ 135,00
		2-chlorobiphenyl (BZ# 1)	2051-60-7	500			
		2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	500			
		2,2',3,3',4,4',6-heptachlorobiphenyl (BZ# 171)	52663-71-5	500			
		2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	500			
		2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 200)	40186-71-8	500			
		2,2',3',4,6-pentachlorobiphenyl (BZ# 98)	60233-25-2	500			
		2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	500			
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	500			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-012101	Method 680 PCB Congener Calibration Mix (Various concentrations)				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 198,00
LK3-012102	Method 680 PCB Congener Calibration Mix, 500 mg/L					1 ml, Amber Glass	€ 394,50
		2-chlorobiphenyl (BZ# 1)	2051-60-7	100			
		2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	100			
		2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	100			
		2,2',4,6-tetrachlorobiphenyl (BZ# 50)	62796-65-0	200			
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	200			
		2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	200			
		2,2',3,4',5,6,6'-heptachlorobiphenyl (BZ# 188)	74487-85-7	300			
		2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 200)	40186-71-8	300			
		decachlorobiphenyl	2051-24-3	500			
LK3-031502	PCB Congener Mix 9-6	2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100	Isooctane	1 ml, Amber Glass	€ 181,50
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	100			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	100			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	100			
		decachlorobiphenyl (BZ 209)	2051-24-3	100			
		2,4,6-trichlorobiphenyl (BZ# 30)	35693-92-6	100			
LK3-243502	Custom PCB Mix, 100 mg/L,	9-24			Tol.:Isooct. 1:1	1 ml, Amber Glass	€ 118,50
LK3-243501	Custom PCB Mix, 500 mg/L,	9-24				1 ml, Amber Glass	€ 355,50
		2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	100			
		2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	100			
		2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	100			
		2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	100			
		2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	100			
		2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	100			
		2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1	100			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	100			
LK3-031201	PCB Congener Mix 10-1	3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	100	Isooctane	1 ml, Amber Glass	€ 198,00
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	100			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	100			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	100			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	100			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	100			
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	100			
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	100			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	31508-00-6	100			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	100			
LK3-102203	Custom PCB Mix, 10-5	2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	5	Hexane	1 ml, Amber Glass	€ 55,50
		decachlorobiphenyl (BZ 209)	2051-24-3	5			
		2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	5			
		2,2',3,4,5,5',6-heptachlorobiphenyl (BZ# 185)	52712-05-7	5			
		2,2',3,3',6,6'-hexachlorobiphenyl (BZ# 136)	38411-22-2	10			
		2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	10			
		2,3',4,5',6-pentachlorobiphenyl (BZ# 121)	56558-18-0	10			
		2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	10			
		2-chlorobiphenyl (BZ# 1)	2051-60-7	100			
		3,3'-dichlorobiphenyl (BZ# 11)	2050-67-1	100			
LK3-244001	Custom PCB Mix, 11-1	3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	100	Hexane	1 ml, Amber Glass	€ 213,00
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	100			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	100			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	100			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	100			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	100			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	100			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	100			
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	100			
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	100			
LK3-030603	PCB Congener Mix 1, 12-1	2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10	Heptane	1 ml, Amber Glass	€ 63,00
		2,2',3,4,4',5-hexachlorobiphenyl (BZ# 137)	35694-06-5	10			
		2,2',3,4,5,6'-hexachlorobiphenyl (BZ# 143)	68194-15-0	10			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	10			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	10			
		2,2',3,5-tetrachlorobiphenyl (BZ# 43)	70362-46-8	10			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	10			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10			
		2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	10			
LK3-237006	PCB Congener Mix, 2 mg/L , 12-1				Nonane	1 ml, Amber Glass	€ 72,00
LK3-237001	PCB Congener Mix, 2 mg/L , 12-1				Isooctane	1 ml, Amber Glass	€ 55,50
LK3-237002	PCB Congener Mix, 10 mg/L , 12-1					1 ml, Amber Glass	€ 55,50
LK3-237003	PCB Congener Mix, 100 mg/L , 12-1					1 ml, Amber Glass	€ 237,00
LK3-237004	PCB Congener Mix, 200 mg/L , 12-1					1 ml, Amber Glass	€ 355,50
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3				
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4				
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4				
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0				
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6				
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3				
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8				
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4				
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7				

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6				
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6				
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9				
LK3-030901	PCB Congener Solution, 10 mg/L , 13-3				Nonane	1 ml, Amber Glass	€ 55,50
LK3-030903	PCB Congener Solution, 10 mg/L , 13-3				Isooctane	1 ml, Amber Glass	€ 55,50
LK3-03090201	PCB Congener Solution, 100 mg/L , 13-4					1 ml, Amber Glass	€ 72,00
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5				
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3				
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3				
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4				
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2				
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6				
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8				
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3				
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2				
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1				
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4				
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6				
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3				
LK3-062101	PCB Congener Solution 14-1	3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	10	Isooctane	1 ml, Amber Glass	€ 63,00
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	10			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	10			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	10			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	10			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	10			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	10			
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	10			
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	10			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	10			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	10			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	10			
LK3-237103	Custom PCB Calibration Mix, 14-22	2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	10	Isooctane	1 ml, Amber Glass	€ 87,00
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10			
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	10			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	10			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			
		2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	10			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	10			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
		2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	10			
		decachlorobiphenyl (BZ 209)	2051-24-3	10			
LK3-247306	Custom PCB Mix, 2 mg/L, 14-22				Nonane	1 ml, Amber Glass	€ 72,00
LK3-247303	Custom PCB Mix, 100 mg/L, 14-22				Isooctane	1 ml, Amber Glass	€ 237,00
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	100			
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	100			
		2',3,4-trichlorobiphenyl (BZ# 33)	38444-86-9	100			
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	100			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	100			
		2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	100			
		2,4,4',5-tetrachlorobiphenyl (BZ# 74)	32690-93-0	100			
		3,3',4,5'-tetrachlorobiphenyl (BZ# 79)	41464-48-6	100			
		2,2',3,4',6-pentachlorobiphenyl (BZ# 91)	68194-05-8	100			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)	52712-04-6	100			
		2,2',3,3',4,5,6'-heptachlorobiphenyl (BZ# 174)	38411-25-5	100			
		2,2',3,3',4,4',5,6'-octachlorobiphenyl (BZ# 196)	42740-50-1	100			
		2,2',3,4,4',5,5',6-octachlorobiphenyl (BZ# 203)	52663-76-0	100			
LK3-247001	PCB Congener Mix 20, 15-1				Isooctane	1 ml, Amber Glass	€ 87,00
LK3-247002	PCB Congener Mix 20, 15-1				Cyclohexane	1 ml, Amber Glass	€ 55,50
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	10			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	10			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	10			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	10			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	10			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	10			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	10			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	10			
LK3-714005	PCB Mix, 16-37	2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10	Isooctane	5 ml, Amber Glass	€ 394,50
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	10			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			
		2,2',3,4',5,6-hexachlorobiphenyl (BZ# 149)	38380-04-0	10			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	10			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	10			
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	10			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
		2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	10			
		decachlorobiphenyl (BZ 209)	2051-24-3	10			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10			
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	10			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	10			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
LK3-020301	PCB Congener Solution 17-1	2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10	Isooctane	1 ml, Amber Glass	€ 63,00
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	10			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	10			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	10			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	10			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	10			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	10			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	10			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	10			
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	10			
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	10			
		2,2',3,4,4',5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	10			
LK3-237505	PCB Congener Mix, 100 mg/L , 18-1				Methylene Chloric	1 ml, Amber Glass	€ 300,00
LK3-237502	PCB Congener Mix, 10 mg/L , 18-1				Isooctane	1 ml, Amber Glass	€ 87,00
LK3-237503	PCB Congener Mix, 100 mg/L , 18-1					1 ml, Amber Glass	€ 300,00
LK3-237504	PCB Congener Mix, 200 mg/L , 18-1					1 ml, Amber Glass	€ 394,50
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3				
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4				

			Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
			2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4				
			2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0				
			2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6				
			2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3				
			3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8				
			2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4				
			2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7				
			2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6				
			3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6				
			2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9				
			2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5				
			2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3				
			2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2				
			2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2				
			2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1				
			2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3				
LK3-030001	PCB Congener Mix	18-1	2,4-dichlorobiphenyl (BZ# 7)	33284-50-3	100	Isooctane	1 ml, Amber Glass	€ 315,00
			2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	100			
			2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-1	100			
			2,2',3,5-tetrachlorobiphenyl (BZ# 43)	70362-46-8	100			
			2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
			2,3,4,4'-tetrachlorobiphenyl (BZ# 60)	33025-41-1	100			
			2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
			2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	100			
			2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	100			
			2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			
			2,2',3,4,4',5-hexachlorobiphenyl (BZ# 137)	35694-06-5	100			
			2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
			2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	100			
			2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
			2,2',3,4,5,5',6-heptachlorobiphenyl (BZ# 185)	52712-05-7	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	100			
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	100			
		decachlorobiphenyl	2051-24-3	100			
LK3-031001	PCB Congener Mix, 100 mg/L ,	18-2			Isooctane	1 ml, Amber Glass	€ 315,00
LK3-031101	PCB Congener Mix, 100 mg/L ,	18-3	34883-43-7	100		1 ml, Amber Glass	€ 315,00
		2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7	100			
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	100			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-1	100			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	100			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	100			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	100			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	100			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			
		2,2',3,4,4',5-hexachlorobiphenyl (BZ# 137)	35694-06-5	100			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	100			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	100			
		2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	100			
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	100			
		decachlorobiphenyl	2051-24-3	100			
LK3-060304	PCB Congener Mix, 1 mg/L ,	18-4			Isooctane	1 ml, Amber Glass	€ 198,00
LK3-060303	PCB Congener Mix, 10 mg/L ,	18-4				1 ml, Amber Glass	€ 268,50
LK3-060302	PCB Congener Mix, 100 mg/L ,	18-4				1 ml, Amber Glass	€ 394,50
LK3-060301	PCB Congener Mix, 250 mg/L ,	18-4			Tol.:Hex. 1:1	1 ml, Amber Glass	€ 709,50
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	1			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	1			

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	1			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	1			
		2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	1			
		2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	1			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	1			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	1			
		2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	1			
		2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	1			
		2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	1			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	1			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	1			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	1			
		2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1	1			
		2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	1			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	1			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	1			
LK3-050603	Custom PCB Mix, 19-02	2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	100	Hexane	1 ml, Amber Glass	€ 324,00
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	100			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	100			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	100			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	100			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	100			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	100			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	100			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	100			
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	100			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	100			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	100			
LK3-050001	PCB Congener Mix, 10 mg/L ,	20-1			Isooctane	1 ml, Amber Glass	€ 118,50
LK3-050002	PCB Congener Mix, 100 mg/L ,	20-1			Acetone	1 ml, Amber Glass	€ 315,00
		2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7				
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2				
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5				
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5				
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3				
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0				
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3				
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2				
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4				
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6				
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8				
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3				
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2				
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1				
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6				
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3				
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0				
		2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2				
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9				
		decachlorobiphenyl (BZ 209)	2051-24-3				
LK3-050405-5pck	Custom PCB Mix, 2 mg/L ,	20-02			Isooctane	1 ml, Amber Glass	€ 472,50
LK3-050402	Custom PCB Mix, 200 mg/L ,	20-02				1 ml, Amber Glass	€ 433,50
LK3-050406	Custom PCB Mix, 100 mg/L ,	20-02			Hexane	1 ml, Amber Glass	€ 355,50
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5				

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3				
3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3				
3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4				
2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2				
2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4				
2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0				
2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6				
2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3				
3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8				
2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3				
2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2				
2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1				
2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4				
2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7				
2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6				
3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6				
2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6				
2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3				
2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9				

LK3-050503	Custom PCB Mix, 21-03			Isooctane	1 ml, Amber Glass	€ 135,00
LK3-050503-5pck	Custom PCB Mix, 21-03				5 x 1 ml, Amber	€ 535,50

2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10
2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10
3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	10
3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	10
2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10
2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	10
2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	10
2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10
2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	10

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	10			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	10			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	10			
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	10			
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	10			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	10			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	10			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	10			
		2,4,6-trichlorobiphenyl (BZ# 30)	35693-92-6	10			
LK3-219703	Custom 8082 PCB Mix, 21-45	2-chlorobiphenyl (BZ# 1)	2051-60-7	100	Isooctane	1 ml, Amber Glass	€ 355,50
		2,4,5,6-tetrachloro-m-xylene	877-09-8	100			
		2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	100			
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	100			
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	100			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	100			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	100			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	100			
		2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	100			
		2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	100			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
		2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)	52712-04-6	100			
		2,2',3,4,4',5-hexachlorobiphenyl (BZ# 137)	35694-06-5	100			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	100			
		2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1	100			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	100			
		decachlorobiphenyl (BZ 209)		100			
LK3-057003-5pck	PCB Congener Mix 22-39	2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7	100	Hexane	1 ml, Amber Glass	€ 828,00
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	100			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	100			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	100			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	100			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	100			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	100			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	100			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	100			
		2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	100			
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	100			
		decachlorobiphenyl (BZ 209)	2051-24-3	100			
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	100			
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	100			
		2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1	100			
		2,2',3,4,4',6,6'-heptachlorobiphenyl (BZ# 184)	74472-48-3	100			
LK3-030202	PCB Congener Mix 25-1, 8 ug/	2-chlorobiphenyl (BZ# 1)	2051-60-7	8	Isooctane	1 ml, Amber Glass	€ 72,00
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	8			
		2,2',3,3',4,5,5',6'-octachlorobiphenyl (BZ# 210)	52663-75-9	8			
		2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7	8			
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	8			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-1	8			

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	8			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	8			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	8			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	8			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	8			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	8			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	8			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	8			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	8			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	8			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	8			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	8			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	8			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	8			
		2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	8			
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	8			
		decachlorobiphenyl	2051-24-3	8			
		2,2',3,4',5,6,6'-heptachlorobiphenyl (BZ# 188)	74487-85-7	8			
		2,2',4,6,6'-pentachlorobiphenyl (BZ# 104)	56558-16-8	8			
LK3-030402	PCB Congener Mix, 10 mg/L ,	25-1			Isooctane	1 ml, Amber Glass	€ 198,00
LK3-030203	PCB Congener Mix, 200 mg/L ,	25-1				1 ml, Amber Glass	€ 72,00
LK3-030201	PCB Congener Mix, 2,000 mg/L ,	25-1				1 ml, Amber Glass	€ 87,00
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	10			
		2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7	10			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	10			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
		2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1	10			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	10			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	10			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	10			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	10			
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	10			
		2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	10			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	10			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	10			
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	10			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	10			
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	10			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	10			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	10			
		2,2',3,4,4',6,6'-heptachlorobiphenyl (BZ# 184)	74472-48-3	10			
LK3-219203	PCB Congener Mix, 20 ug/L (ppb) , 26-5				Nonane	1 ml, Amber Glass	€ 118,50
LK3-219204	PCB Congener Mix, 200 ug/L (ppb) , 26-5					1 ml, Amber Glass	€ 157,50
LK3-219205	PCB Congener Mix, 4,000 ug/L (ppb) , 26-5					1 ml, Amber Glass	€ 198,00
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2				
		2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7				
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5				
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6				
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3				
		2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1				
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0				
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3				
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2				
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1				
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4				
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6				

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9				
		2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2				
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2				
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4				
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6				
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8				
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8				
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5				
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8				
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3				
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0				
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3				
		2,2',3,4,4',6,6'-heptachlorobiphenyl (BZ# 184)	74472-48-3				
		decachlorobiphenyl (BZ 209)	2051-24-3				
LK3-011301	PCB Congener Solution 27-1	2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	2	Nonane	1 ml, Amber Glass	€ 198,00
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	2			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	2			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	2			
		2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	2			
		2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	2			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	2			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	2			
		2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	2			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	2			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	2			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	2			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	2			
		2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	2			
		2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	2			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	2			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	2			

			Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
			2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	2			
			2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	2			
			3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	2			
			2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	2			
			2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	2			
			2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	2			
			2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1	2			
			2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	2			
			2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	2			
			2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	2			
LK3-011001	PCB Congener Solution	28-1	2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	100	Isooctane	1 ml, Amber Glass	€ 441,00
			2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
			3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	100			
			3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	100			
			2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	100			
			2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	100			
			2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
			2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	100			
			2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	100			
			2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	100			
			2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	100			
			2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	100			
			3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	100			
			2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
			2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	100			
			2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	100			
			2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	100			
			2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
			2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	100			
			2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	100			
			2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	100			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	100			
		2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	100			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
		2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1	100			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	100			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	100			
LK3-011601	PCB Congener Solution, 28-1,	1,2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	100	Isooctane	1 ml, Amber Glass	€ 394,50
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	100			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	100			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	100			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	100			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	100			
		2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7	100			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	100			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	100			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	100			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	100			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	100			
		2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	100			
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	100			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	100			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	100			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	100			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	100			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	100			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	100			
		decachlorobiphenyl (BZ 209)	2051-24-3	100			
LK3-011101	PCB Congener Solution, 10 mg/L , 29-1					1 ml, Amber Glass	€ 103,50
LK3-011103	PCB Congener Solution, 100 mg/L , 29-1					1 ml, Amber Glass	€ 472,50
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5		Isooctane		
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3				
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3				
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4				
		2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6				
		2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7				
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2				
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4				
		2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9				
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0				
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6				
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3				
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8				
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2				
		2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8				
		2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0				
		2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5				
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1				
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4				
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7				
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6				
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6				
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6				
		2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4				
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3				

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1				
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0				
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9				
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3				
LK3-102003	Custom PCB Conginer Mix, 29-1, 2 mg/L				Isooctane	1 ml, Amber Glass	€ 118,50
LK3-102002	Custom PCB Conginer Mix, 29-1, 20 mg/L					1 ml, Amber Glass	€ 205,50
		4-chlorobiphenyl (BZ# 3)	2051-62-9	2			
		4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	2			
		2,3,3'-trichlorobiphenyl (BZ# 20)	38444-84-7	2			
		3,4,4'-trichlorobiphenyl (BZ# 37)	38444-90-5	2			
		2,2',4,5-tetrachlorobiphenyl (BZ# 48)	70362-47-9	2			
		2,3,4,4'-tetrachlorobiphenyl (BZ# 60)	33025-41-1	2			
		2',3,4,5-tetrachlorobiphenyl (BZ# 76)	70362-48-0	2			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	2			
		2,2',3,3',6-pentachlorobiphenyl (BZ# 84)	52663-60-2	2			
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	2			
		2,2',3,4,6'-pentachlorobiphenyl (BZ# 89)	73575-57-2	2			
		2,2',4,6,6'-pentachlorobiphenyl (BZ# 104)	56558-16-8	2			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	2			
		2,3,4,4',6-pentachlorobiphenyl (BZ# 115)	74472-38-1	2			
		2,3',4,4',6-pentachlorobiphenyl (BZ# 119)	56558-17-9	2			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	2			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	2			
		2,2',3,3',4,6'-hexachlorobiphenyl (BZ# 132)	38380-05-1	2			
		2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)	52712-04-6	2			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	2			
		2,3,3',4',5',6-hexachlorobiphenyl (BZ# 164)	74472-45-0	2			
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	2			
		2,3',4,4',5',6-hexachlorobiphenyl (BZ# 168)	59291-65-5	2			
		2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	2			
		2,2',3,4,4',6,6'-heptachlorobiphenyl (BZ# 184)	74472-48-3	2			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	2			
		2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 200)	40186-71-8	2			
		2,2',3,3',5,5',6,6'-octachlorobiphenyl (BZ# 202)	2136-99-4	2			
		2,2',3,3',4,4',5,6,6'-nonachlorobiphenyl (BZ# 207)	52663-79-3	2			
LK3-036103-5pck	PCB Congener Check Solution	29-1 (2 mg/L) PCB IS Only			Isooctane	1 ml, Amber Glass	€ 198,00
LK3-036198-5pck	PCB Congener Calibration Set	29-1 (0.2, 1, 2, 5, 15 mg/L) PCB IS Only				5 x 1 ml, Amber	€ 591,00
		4-chlorobiphenyl (BZ# 3)	2051-62-9	2			
		4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	2			
		2,3,3'-trichlorobiphenyl (BZ# 20)	38444-84-7	2			
		3,4,4'-trichlorobiphenyl (BZ# 37)	38444-90-5	2			
		2,2',4,5-tetrachlorobiphenyl (BZ# 48)	70362-47-9	2			
		2,3,4,4'-tetrachlorobiphenyl (BZ# 60)	33025-41-1	2			
		2',3,4,5-tetrachlorobiphenyl (BZ# 76)	70362-48-0	2			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)		2			
		2,2',3,3',6-pentachlorobiphenyl (BZ# 84)	52663-60-2	2			
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	2			
		2,2',3,4,6'-pentachlorobiphenyl (BZ# 89)	73575-57-2	2			
		2,2',4,6,6'-pentachlorobiphenyl (BZ# 104)	56558-16-8	2			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	2			
		2,3,4,4',6-pentachlorobiphenyl (BZ# 115)	74472-38-1	2			
		2,3',4,4',6-pentachlorobiphenyl (BZ# 119)	56558-17-9	2			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	2			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	2			
		2,2',3,3',4,6'-hexachlorobiphenyl (BZ# 132)	38380-05-1	2			
		2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)	52712-04-6	2			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	2			
		2,3,3',4',5',6-hexachlorobiphenyl (BZ# 164)	74472-45-0	2			
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	2			
		2,3',4,4',5',6-hexachlorobiphenyl (BZ# 168)	59291-65-5	2			
		2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	2			
		2,2',3,4,4',6,6'-heptachlorobiphenyl (BZ# 184)	74472-48-3	2			

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	2			
		2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 200)	40186-71-8	2			
		2,2',3,3',5,5',6,6'-octachlorobiphenyl (BZ# 202)	2136-99-4	2			
		2,2',3,3',4,4',5,6,6'-nonachlorobiphenyl (BZ# 207)	52663-79-3	2			
		13C12-4-chlorobiphenyl (BZ# 3)	2051-62-9	0,075			
		13C12-4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	0,075			
		13C12-2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	0,075			
		13C12-2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	0,075			
		13C12-2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	0,075			
		13C12-2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	0,075			
		13C12-2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	0,075			
		13C12-2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	0,075			
		13C12-2,2',3,3',4,5,5',6,6'-nonachlorobiphenyl (BZ# 2)	52663-77-1	0,075			
		13C12-decachlorobiphenyl	2051-24-3	0,075			
LK3-030801	PCB Congener Solution 32-1	2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7	100	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 567,00
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	100			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	100			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	100			
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	100			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	100			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	100			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	100			
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	100			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	100			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	100			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	100			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	100			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	100			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	100			
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	100			
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	100			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	100			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	100			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
		2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1	100			
		2,2',3,4,4',6,6'-heptachlorobiphenyl (BZ# 184)	74472-48-3	100			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	100			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	100			
		2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	100			
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	100			
		decachlorobiphenyl	2051-24-3	100			
LK3-030701	PCB Congener Solution, 10 mg/L , 32-2				Isooctane	1 ml, Amber Glass	€ 229,50
LK3-030703	PCB Congener Solution, 100 mg/L , 32-2					1 ml, Amber Glass	€ 433,50
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5				
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3				
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3				
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4				
		2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6				
		2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7				
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2				
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4				
		2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9				
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0				
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6				
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3				
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8				
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2				

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8				
		2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0				
		2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5				
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1				
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4				
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7				
		2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6				
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6				
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6				
		2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4				
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3				
		2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1				
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0				
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9				
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3				
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2				
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5				
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3				
LK3-241703	Custom PCB Congener Mix, 40-	2,2',3,3',4-pentachlorobiphenyl (BZ# 82)	52663-62-4	10	Isooctane	1 ml, Amber Glass	€ 355,50
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	10			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	10			
		2,2',3,3',5,6,6'-heptachlorobiphenyl (BZ# 179)	52663-64-6	10			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	10			
		2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1	10			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	10			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	10			
		2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	10			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	10			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
		2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	10			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	10			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	10			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	10			
		2,3,3',4,4',6-hexachlorobiphenyl (BZ# 158)	74472-42-7	10			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	10			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	10			
		2,3,4,4',5,6-hexachlorobiphenyl (BZ# 166)	41411-63-6	10			
		2,3,4,4'-tetrachlorobiphenyl (BZ# 60)	33025-41-1	10			
		2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7	10			
		2,4,4',5-tetrachlorobiphenyl (BZ# 74)	32690-93-0	10			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	10			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	10			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	10			
		3,4,4'-trichlorobiphenyl (BZ# 37)	38444-90-5	10			
		2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	10			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	10			
		2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	10			
		2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	10			
		2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	10			
		2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	10			
		2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	10			
		2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	10			
LK3-247203-5pck	PCB Congener Solution	42-3 2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	100	Isooctane	1 ml, Amber Glass	€ 2.512,50
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	100			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	100			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	100			
2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	100			
2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	100			
2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	100			
2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	100			
2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	100			
2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	100			
3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	100			
2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	100			
2,2',3,4',5,6-hexachlorobiphenyl (BZ# 149)	38380-04-0	100			
2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	100			
2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	100			
2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	100			
2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	100			
3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	100			
2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	100			
2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	100			
2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1	100			
2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	100			
2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	100			
2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	100			
2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	100			
2',3,4-trichlorobiphenyl (BZ# 33)	38444-86-9	100			
2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	100			
2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	100			
2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	100			
2,4,4',5-tetrachlorobiphenyl (BZ# 74)	32690-93-0	100			
3,3',4,5'-tetrachlorobiphenyl (BZ# 79)	41464-48-6	100			

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,4',6-pentachlorobiphenyl (BZ# 91)	68194-05-8	100			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			
		2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)	52712-04-6	100			
		2,2',3,3',4,5,6'-heptachlorobiphenyl (BZ# 174)	38411-25-5	100			
		2,2',3,3',4,4',5,6'-octachlorobiphenyl (BZ# 196)	42740-50-1	100			
		2,2',3,4,4',5,5',6-octachlorobiphenyl (BZ# 203)	52663-76-0	100			
LK3-241501	Custom PCB Congener Mix, 42-	2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	100	Isooctane	1 ml, Amber Glass	€ 591,00
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	100			
		2,4,6-trichlorobiphenyl (BZ# 30)	35693-92-6	100			
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	100			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	100			
		2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	100			
		2,4,4',5-tetrachlorobiphenyl (BZ# 74)	32690-93-0	100			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	100			
		3,3',4,5'-tetrachlorobiphenyl (BZ# 79)	41464-48-6	100			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	100			
		2,2',3,4',6-pentachlorobiphenyl (BZ# 91)	68194-05-8	100			
		2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	100			
		2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	100			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	100			
		2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	100			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	100			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	100			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	100			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	100			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	100			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
		2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)	52712-04-6	100			
		2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	100			

	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
	2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	100			
	2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	100			
	2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
	2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	100			
	2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	100			
	2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	100			
	3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	100			
	2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	100			
	2,2',3,3',4,5,6-heptachlorobiphenyl (BZ# 174)	38411-25-5	100			
	2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	100			
	2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
	2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1	100			
	2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	100			
	2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	100			
	2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	100			
	2,2',3,4,4',5,5',6-octachlorobiphenyl (BZ# 203)	52663-76-0	100			
	2,3,3',4',5,6-hexachlorobiphenyl (BZ# 163)	74472-44-9	100			
LK3-224504-5pck Custom PCB Mix, 47-3	3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	2	Isooctane	1 ml, Amber Glass	€ 1.095,00
	2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	2			
	2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	2			
	2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	2			
	2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	2			
	2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	2			
	2-chlorobiphenyl (BZ# 1)	2051-60-7	2			
	2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	2			
	2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	2			
	2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	2			
	2,4,6-trichlorobiphenyl (BZ# 30)	35693-92-6	2			
	2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	2			
	2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	2			
	2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	2			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	2			
2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	2			
2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	2			
2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	2			
2,2',3,4,4',5-hexachlorobiphenyl (BZ# 137)	35694-06-5	2			
2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	2			
2,2',3,4',5,6-hexachlorobiphenyl (BZ# 149)	38380-04-0	2			
2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	2			
2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	2			
2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1	2			
2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	2			
2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	2			
2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	2			
2,2',3,3',4,4',5,6'-octachlorobiphenyl (BZ# 196)	42740-50-1	2			
2,2',3,4,4',5,5',6-octachlorobiphenyl (BZ# 203)	52663-76-0	2			
decachlorobiphenyl (BZ 209)	38380-07-3	2			
2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	2			
2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	2			
2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	2			
3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	2			
2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	2			
2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	2			
2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	2			
3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	2			
2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	2			
2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	2			
2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	2			
2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	2			
3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	2			
2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	2			
2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	2			
2,3,4',5,6-pentachlorobiphenyl (BZ# 117)	68194-11-6	2			

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,3,3',4,5,5',6-heptachlorobiphenyl (BZ# 192)	74472-51-8	2			
LK3-224001	Custom PCB Congener Mix, 48-	2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	10	Nonane	5 ml, Amber Glass	€ 709,50
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10			
		2,4,6-trichlorobiphenyl (BZ# 30)	35693-92-6	10			
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	10			
		3,3',5-trichlorobiphenyl (BZ# 36)	38444-87-0	10			
		2,2',3,3'-tetrachlorobiphenyl (BZ# 40)	38444-93-8	10			
		2,2',3,5-tetrachlorobiphenyl (BZ# 43)	70362-46-8	10			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	10			
		2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	10			
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	10			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	10			
		2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	10			
		2,4,4',5-tetrachlorobiphenyl (BZ# 74)	32690-93-0	10			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	10			
		3,3',4,5'-tetrachlorobiphenyl (BZ# 79)	41464-48-6	10			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	52663-72-6	10			
		2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	10			
		2,2',3,4',6-pentachlorobiphenyl (BZ# 91)	68194-05-8	10			
		2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	10			
		2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	10			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10			
		2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	10			
		2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	10			
		2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	10			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	10			
		3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	10			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	10			
		2,2',3,4,4',5-hexachlorobiphenyl (BZ# 137)	35694-06-5	10			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			
		2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)	52712-04-6	10			
		2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	10			
		2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	10			
		2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	10			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	10			
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	10			
		3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	10			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	10			
		2,2',3,3',4,5,6'-heptachlorobiphenyl (BZ# 174)	38411-25-5	10			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
		2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1	10			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	10			
		2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	10			
		2,2',3,3',4,4',5,6'-octachlorobiphenyl (BZ# 196)	42740-50-1	10			
		2,2',3,3',4,5,6,6'-octachlorobiphenyl (BZ# 199)	52663-73-7	10			
		2,2',3,4,4',5,5',6-octachlorobiphenyl (BZ# 203)	52663-76-0	10			
LK3-241001	Custom PCB Congener Mix, 48-	2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	10	Nonane	5 ml, Amber Glass	€ 709,50
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10			
		2,4,6-trichlorobiphenyl (BZ# 30)	35693-92-6	10			
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	10			
		2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	10			
		2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	10			
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	10			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	10			
		2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	10			
		2,4,4',5-tetrachlorobiphenyl (BZ# 74)	32690-93-0	10			
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	10			
		3,3',4,5'-tetrachlorobiphenyl (BZ# 79)	41464-48-6	10			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	10			
2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	10			
2,2',3,4',6-pentachlorobiphenyl (BZ# 91)	68194-05-8	10			
2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	10			
2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	10			
2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10			
2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	10			
2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	10			
2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	10			
2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	10			
3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	10			
2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	10			
2,2',3,4,4',5-hexachlorobiphenyl (BZ# 137)	35694-06-5	10			
2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			
2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)	52712-04-6	10			
2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	10			
2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	10			
2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	10			
2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	10			
2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	10			
2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	10			
3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	10			
2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	10			
2,2',3,3',4,5,6'-heptachlorobiphenyl (BZ# 174)	38411-25-5	10			
2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	10			
2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1	10			
2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	10			
2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	10			
2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	10			

	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
	2,2',3,3',4,4',5,6'-octachlorobiphenyl (BZ# 196)	42740-50-1	10			
	2,2',3,3',4,5,5',6'-octachlorobiphenyl (BZ# 201)	52663-75-9	10			
	2,2',3,4,4',5,5',6-octachlorobiphenyl (BZ# 203)	52663-76-0	10			
LK3-036098-5pck	PCB Congener Calibration Set 69-2 (0.2, 1, 2, 5, 15 mg/L, PCB IS Only)			Isooctane	1 ml, Amber Glass	€ 630,00
	2-chlorobiphenyl (BZ# 1)	2051-60-7	15			
	3-chlorobiphenyl (BZ# 2)	2051-61-8	15			
	2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	15			
	2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7	15			
	2,5-dichlorobiphenyl (BZ# 9)	34883-39-1	15			
	3,4-dichlorobiphenyl (BZ# 12)	2974-92-7	15			
	2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	15			
	2,3',5-trichlorobiphenyl (BZ# 26)	38444-81-4	15			
	2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	15			
	2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	15			
	2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	15			
	2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	15			
	2,2',3,6-tetrachlorobiphenyl (BZ# 45)	70362-45-7	15			
	2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	15			
	2,2',4,6-tetrachlorobiphenyl (BZ# 50)	62796-65-0	15			
	2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	15			
	2,3,3',4'-tetrachlorobiphenyl (BZ# 56)	41464-43-1	15			
	2,3,4,5-tetrachlorobiphenyl (BZ# 61)	33284-53-6	15			
	2,3,4',5-tetrachlorobiphenyl (BZ# 63)	74472-34-7	15			
	2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	15			
	2,3',4,6-tetrachlorobiphenyl (BZ# 69)	60233-24-1	15			
	2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	15			
	2,4,4',5-tetrachlorobiphenyl (BZ# 74)	32690-93-0	15			
	3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	15			
	2,2',3,3',4-pentachlorobiphenyl (BZ# 82)	52663-62-4	15			
	2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	15			
	2,2',3,4,6-pentachlorobiphenyl (BZ# 88)	55215-17-3	15			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
2,2',3,4',5-pentachlorobiphenyl (BZ# 90)	68194-07-0	15			
2,2',3,5,5'-pentachlorobiphenyl (BZ# 92)	52663-61-3	15			
2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	15			
2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	15			
2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	15			
2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	15			
2,3,3',4,5-pentachlorobiphenyl (BZ# 106)	70424-69-0	15			
2,3,3',4,5-pentachlorobiphenyl (BZ# 107)	70424-68-9	15			
2,3,3',4,5'-pentachlorobiphenyl (BZ# 108)	70362-41-3	15			
2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	15			
2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	15			
3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	15			
2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	15			
2,2',3,3',4,5'-hexachlorobiphenyl (BZ# 130)	52663-66-8	15			
2,2',3,3',4,6'-hexachlorobiphenyl (BZ# 132)	38380-05-1	15			
2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	15			
2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	15			
2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	15			
2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	15			
2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	15			
2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	15			
2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	15			
2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	15			
2,3,3',4,4',6-hexachlorobiphenyl (BZ# 158)	74472-42-7	15			
2,3,3',4,5,5'-hexachlorobiphenyl (BZ# 159)	39635-35-3	15			
2,3,3',4',5,6-hexachlorobiphenyl (BZ# 163)	74472-44-9	15			
2,3,3',5,5',6-hexachlorobiphenyl (BZ# 165)	74472-46-1	15			
3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	15			
2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	15			
2,2',3,3',4,5,5'-heptachlorobiphenyl (BZ# 172)	52663-74-8	15			
2,2',3,3',4,5,6'-heptachlorobiphenyl (BZ# 174)	38411-25-5	15			
2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	15			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1	15			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	15			
		2,2',3,4',5,6,6'-heptachlorobiphenyl (BZ# 188)	74487-85-7	15			
		2,3,3',4,4',5,6-heptachlorobiphenyl (BZ# 190)	41411-64-7	15			
		2,3,3',4',5,5',6-heptachlorobiphenyl (BZ# 193)	69782-91-8	15			
		2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	15			
		2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	15			
		2,2',3,3',4,5,5',6'-octachlorobiphenyl (BZ# 201)	52663-75-9	15			
		2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	15			
		decachlorobiphenyl	2051-24-3	0,075			
		13C12-4-chlorobiphenyl (BZ# 3)	2051-62-9	0,075			
		13C12-4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	0,075			
		13C12-2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	0,075			
		13C12-2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	0,075			
		13C12-2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	0,075			
		13C12-2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	0,075			
		13C12-2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	0,075			
		13C12-2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	0,075			
		13C12-2,2',3,3',4,5,5',6,6'-nonachlorobiphenyl (BZ# 2)	52663-77-1	0,075			
		13C12-decachlorobiphenyl	2051-24-3	0,075			
LK3-036601	PCB Congener Custom Mix 94	2-chlorobiphenyl (BZ# 1)	2051-60-7	20	Isooctane	5 ml, BR, Amber Glass	€ 2.757,00
		3-chlorobiphenyl (BZ# 2)	2051-61-8	20			
		2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	20			
		2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7	20			
		2,5-dichlorobiphenyl (BZ# 9)	34883-39-1	20			
		3,4-dichlorobiphenyl (BZ# 12)	2974-92-7	20			
		2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	20			
		2,3',5-trichlorobiphenyl (BZ# 26)	38444-81-4	20			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	20			
		2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	20			
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	20			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	20			
2,2',3,6-tetrachlorobiphenyl (BZ# 45)	70362-45-7	20			
2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	20			
2,2',4,6-tetrachlorobiphenyl (BZ# 50)	62796-65-0	20			
2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	20			
2,3,3',4'-tetrachlorobiphenyl (BZ# 56)	41464-43-1	20			
2,3,4,5-tetrachlorobiphenyl (BZ# 61)	33284-53-6	20			
2,3,4',5-tetrachlorobiphenyl (BZ# 63)	74472-34-7	20			
2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	20			
2,3',4,6-tetrachlorobiphenyl (BZ# 69)	60233-24-1	20			
2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	20			
2,4,4',5-tetrachlorobiphenyl (BZ# 74)	32690-93-0	20			
3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	20			
2,2',3,3',4-pentachlorobiphenyl (BZ# 82)	52663-62-4	20			
2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	20			
2,2',3,4,6-pentachlorobiphenyl (BZ# 88)	55215-17-3	20			
2,2',3,4',5-pentachlorobiphenyl (BZ# 90)	68194-07-0	20			
2,2',3,5,5'-pentachlorobiphenyl (BZ# 92)	52663-61-3	20			
2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	20			
2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	20			
2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	20			
2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	20			
2,3,3',4,5-pentachlorobiphenyl (BZ# 106)	70424-69-0	20			
2,3,3',4',5-pentachlorobiphenyl (BZ# 107)	70424-68-9	20			
2,3,3',4,5'-pentachlorobiphenyl (BZ# 108)	70362-41-3	20			
2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	20			
2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	20			
3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	20			
2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	20			
2,2',3,3',4,5'-hexachlorobiphenyl (BZ# 130)	52663-66-8	20			
2,2',3,3',4,6'-hexachlorobiphenyl (BZ# 132)	38380-05-1	20			
2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	20			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	20			
2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	20			
2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	20			
2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	20			
2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	20			
2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	20			
2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	20			
2,3,3',4,4',6-hexachlorobiphenyl (BZ# 158)	74472-42-7	20			
2,3,3',4,5,5'-hexachlorobiphenyl (BZ# 159)	39635-35-3	20			
2,3,3',4',5,6-hexachlorobiphenyl (BZ# 163)	74472-44-9	20			
2,3,3',5,5',6-hexachlorobiphenyl (BZ# 165)	74472-46-1	20			
3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	20			
2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	20			
2,2',3,3',4,5,5'-heptachlorobiphenyl (BZ# 172)	52663-74-8	20			
2,2',3,3',4,5,6'-heptachlorobiphenyl (BZ# 174)	38411-25-5	20			
2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	20			
2,2',3,4,4',5,6-heptachlorobiphenyl (BZ# 183)	52663-69-1	20			
2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	20			
2,2',3,4',5,6,6'-heptachlorobiphenyl (BZ# 188)	74487-85-7	20			
2,3,3',4,4',5,6-heptachlorobiphenyl (BZ# 190)	41411-64-7	20			
2,3,3',4',5,5',6-heptachlorobiphenyl (BZ# 193)	69782-91-8	20			
2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	20			
2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	20			
2,2',3,3',4,5,5',6'-octachlorobiphenyl (BZ# 201)	52663-75-9	20			
2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	20			
decachlorobiphenyl	2051-24-3	20			
4-chlorobiphenyl (BZ# 3)	2051-62-9	20			
4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	20			
2,3,3'-trichlorobiphenyl (BZ# 20)	38444-84-7	20			
3,4,4'-trichlorobiphenyl (BZ# 37)	38444-90-5	20			
2,2',4,5-tetrachlorobiphenyl (BZ# 48)	70362-47-9	20			
2,3,4,4'-tetrachlorobiphenyl (BZ# 60)	33025-41-1	20			

Compound		CAS#	Concentr. mg/L	Matrix	Container	List Price
2',3,4,5-tetrachlorobiphenyl (BZ# 76)		70362-48-0	20			
3,4,4',5-tetrachlorobiphenyl (BZ# 81)			20			
2,2',3,3',6-pentachlorobiphenyl (BZ# 84)		52663-60-2	20			
2,2',3,4,6'-pentachlorobiphenyl (BZ# 89)		73575-57-2	20			
2,2',4,6,6'-pentachlorobiphenyl (BZ# 104)		56558-16-8	20			
2,3,4,4',5-pentachlorobiphenyl (BZ# 114)		74472-37-0	20			
2,3,4,4',6-pentachlorobiphenyl (BZ# 115)		74472-38-1	20			
2,3',4,4',6-pentachlorobiphenyl (BZ# 119)		56558-17-9	20			
2',3,4,4',5-pentachlorobiphenyl (BZ# 123)		65510-44-3	20			
2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)		52712-04-6	20			
2,3,3',4',5,6-hexachlorobiphenyl (BZ# 164)		74472-45-0	20			
2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)		52663-72-6	20			
2,3',4,4',5,6-hexachlorobiphenyl (BZ# 168)		59291-65-5	20			
2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)		52663-70-4	20			
2,2',3,4,4',6,6'-heptachlorobiphenyl (BZ# 184)		74472-48-3	20			
2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)		39635-31-9	20			
2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 200)		40186-71-8	20			
2,2',3,3',5,5',6,6'-octachlorobiphenyl (BZ# 202)		2136-99-4	20			
2,2',3,3',4,4',5,6,6'-nonachlorobiphenyl (BZ# 207)		52663-79-3	20			
LK3-102102	Custom PCB Conginer Mix, 96-:					
	2-chlorobiphenyl (BZ# 1)	2051-60-7	1000	Isooctane	1 ml, Amber Glass	€ 828,00
	4-chlorobiphenyl (BZ# 3)	2051-62-9	1000			
	2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	1000			
	3,4-dichlorobiphenyl (BZ# 12)	2974-92-7	1000			
	4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	1000			
	2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	1000			
	2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	1000			
	3,4,4'-trichlorobiphenyl (BZ# 37)	38444-90-5	1000			
	2,2',3,5'-tetrachlorobiphenyl (BZ# 44)	41464-39-5	1000			
	2,2',3,6-tetrachlorobiphenyl (BZ# 45)	70362-45-7	1000			
	2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	1000			
	2,2',4,6-tetrachlorobiphenyl (BZ# 50)	62796-65-0	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	1000			
2,3,3',4'-tetrachlorobiphenyl (BZ# 56)	41464-43-1	1000			
2,3,4,5-tetrachlorobiphenyl (BZ# 61)	33284-53-6	1000			
2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	1000			
2',3,4,5-tetrachlorobiphenyl (BZ# 76)	70362-48-0	1000			
3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	1000			
2,2',3,3',4-pentachlorobiphenyl (BZ# 82)	52663-62-4	1000			
2,2',3,4,5'-pentachlorobiphenyl (BZ# 87)	38380-02-8	1000			
2,2',3,4,6-pentachlorobiphenyl (BZ# 88)	55215-17-3	1000			
2,2',3,5,5'-pentachlorobiphenyl (BZ# 92)	52663-61-3	1000			
2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	1000			
2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	1000			
2,2',4,6,6'-pentachlorobiphenyl (BZ# 104)	56558-16-8	1000			
2,3,3',4,4'-pentachlorobiphenyl (BZ# 105)	32598-14-4	1000			
2,3,3',4',5-pentachlorobiphenyl (BZ# 107)	70424-68-9	1000			
2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	1000			
2,3,4,4',5-pentachlorobiphenyl (BZ# 114)	74472-37-0	1000			
2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	1000			
3,3',4,4',5-pentachlorobiphenyl (BZ# 126)	57465-28-8	1000			
2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	1000			
2,2',3,3',4,5'-hexachlorobiphenyl (BZ# 130)	52663-66-8	1000			
2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	1000			
2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	1000			
2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	1000			
2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	1000			
2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	1000			
2,3,3',4,5,5'-hexachlorobiphenyl (BZ# 159)	39635-35-3	1000			
2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	1000			
2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1	1000			
2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	1000			
2,2',3,4',5,6,6'-heptachlorobiphenyl (BZ# 188)	74487-85-7	1000			
2,4'-dichlorobiphenyl (BZ# 8)	34883-43-7	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
2,5-dichlorobiphenyl (BZ# 9)	34883-39-1	1000			
2,3,3'-trichlorobiphenyl (BZ# 20)	38444-84-7	1000			
2,3',5-trichlorobiphenyl (BZ# 26)	38444-81-4	1000			
2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	1000			
2,3,4,4'-tetrachlorobiphenyl (BZ# 60)	33025-41-1	1000			
2,3,4',5-tetrachlorobiphenyl (BZ# 63)	74472-34-7	1000			
2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	1000			
2,2',3,4,6'-pentachlorobiphenyl (BZ# 89)	73575-57-2	1000			
2,2',4,5',6-pentachlorobiphenyl (BZ# 103)	60145-21-3	1000			
2,3,3',4,5-pentachlorobiphenyl (BZ# 106)	70424-69-0	1000			
2,3,4,4',6-pentachlorobiphenyl (BZ# 115)	74472-38-1	1000			
2,3',4,4',6-pentachlorobiphenyl (BZ# 119)	56558-17-9	1000			
2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	1000			
2,3,3',4,4',6-hexachlorobiphenyl (BZ# 158)	74472-42-7	1000			
2,3,3',4',5',6-hexachlorobiphenyl (BZ# 164)	74472-45-0	1000			
2,3,3',5,5',6-hexachlorobiphenyl (BZ# 165)	74472-46-1	1000			
2,3',4,4',5,5'-hexachlorobiphenyl (BZ# 167)	52663-72-6	1000			
2,3',4,4',5',6-hexachlorobiphenyl (BZ# 168)	59291-65-5	1000			
3,3',4,4',5,5'-hexachlorobiphenyl (BZ# 169)	32774-16-6	1000			
2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	1000			
2,2',3,3',4,5,5'-heptachlorobiphenyl (BZ# 172)	52663-74-8	1000			
2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	1000			
2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	1000			
2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	1000			
2,2',3,3',4,4',5,6-octachlorobiphenyl (BZ# 195)	52663-78-2	1000			
2,2',3,3',4,5,5',6-octachlorobiphenyl (BZ# 198)	68194-17-2	1000			
2,2',3,3',4,5,5',6'-octachlorobiphenyl (BZ# 201)	52663-75-9	1000			
2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 200)	40186-71-8	1000			
2,2',3,3',4,4',5,5',6-nonachlorobiphenyl (BZ# 206)	40186-72-9	1000			
2,2',3,3',4,4',5,6,6'-nonachlorobiphenyl (BZ# 207)	52663-79-3	1000			
3-chlorobiphenyl (BZ# 2)	2051-61-8	1000			
2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	1000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',4,5-tetrachlorobiphenyl (BZ# 48)	70362-47-9	1000			
		2,3',4,6-tetrachlorobiphenyl (BZ# 69)	60233-24-1	1000			
		3,4,4',5-tetrachlorobiphenyl (BZ# 81)	70362-50-4	1000			
		2,2',3,3',6-pentachlorobiphenyl (BZ# 84)	52663-60-2	1000			
		2,2',3,4',5-pentachlorobiphenyl (BZ# 90)	68194-07-0	1000			
		2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	1000			
		2,3,3',4,5'-pentachlorobiphenyl (BZ# 108)	70362-41-3	1000			
		2',3,4,4',5-pentachlorobiphenyl (BZ# 123)	65510-44-3	1000			
		2,2',3,3',4,6'-hexachlorobiphenyl (BZ# 132)	38380-05-1	1000			
		2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)	52712-04-6	1000			
		2,4,4',5-tetrachlorobiphenyl (BZ# 74)	32690-93-0	1000			
		2,3,3',4,4',5'-hexachlorobiphenyl (BZ# 157)	69782-90-7	1000			
		2,3,3',4',5,6-hexachlorobiphenyl (BZ# 163)	74472-44-9	1000			
		2,2',3,3',4,5,6'-heptachlorobiphenyl (BZ# 174)	38411-25-5	1000			
		2,2',3,4,4',6,6'-heptachlorobiphenyl (BZ# 184)	74472-48-3	1000			
		2,3,3',4,4',5,6-heptachlorobiphenyl (BZ# 190)	41411-64-7	1000			
		2,3,3',4',5,5',6-heptachlorobiphenyl (BZ# 193)	69782-91-8	1000			
		2,2',3,3',5,5',6,6'-octachlorobiphenyl (BZ# 202)	2136-99-4	1000			
		decachlorobiphenyl (BZ 209)	2051-24-3	1000			
		2,3,3',4,4',5-hexachlorobiphenyl (BZ# 156)	38380-08-4	1000			
LK3-200501	Custom Pesticide Mix 1, 3-1	1,2,4,5-tetrachlorobenzene	95-94-3	100	Acetone	1 ml, Amber Glass	€ 40,50
		pentachlorobenzene	608-93-5	100			
		propachlor	1918-16-7	100			
LK3-038003-4pck	Custom Pesticide Mix 2-20	Captafol	2425-06-1	100	Acetonitrile	1 ml, Amber Glass	€ 189,00
		dinoseb	88-85-7	100			
LK3-201001	Custom Pesticide Mix 2, 4-1	pentachloroanisole	1825-21-4	100	Acetone	1 ml, Amber Glass	€ 48,00
		hexachlorobenzene	118-74-1	100			
		mirex	2385-85-5	100			
		trifluralin	1582-09-8	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-101701	Custom Pesticide Mix #2, 7-7, Various Concentrations				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 55,50
LK3-101701-5pck	Custom Pesticide Mix #2, 7-7, Various Concentrations					5 x 1 ml, Amber	€ 220,50
		endrin	72-20-8	80			
		g-BHC	58-89-9	20			
		heptachlor	76-44-8	30			
		heptachlor epoxide isomer B	1024-57-3	40			
		methoxychlor	72-43-5	200			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	40			
		decachlorobiphenyl	2051-24-3	80			
LK3-204101	Custom Pesticide Solution NPII	phorate	298-02-2	1000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 55,50
		fenchlorphos	299-84-3	1000			
		tetrachlorvinphos (Rabon)	961-11-5	5000			
LK3-204001	Custom Pesticide Solution NPII	phorate	298-02-2	2	Hexane	1 ml, Amber Glass	€ 40,50
		fenchlorphos	299-84-3	2			
		tetrachlorvinphos (Rabon)	961-11-5	10			
		1,3-dimethyl-2-nitrobenzene	81-20-9	25			
LK3-251002	Custom Pesticide Mix, 3-2	trichlorfon	52-68-6	1000	Methanol	1 ml, Amber Glass	€ 40,50
		dimethoate	60-51-5	1000			
		malathion	121-75-5	1000			
LK3-201501	Custom Pesticide Mix 3, 5-1	oxychlordane	27304-13-8	100	Acetone	1 ml, Amber Glass	€ 276,00
		cis-nonachlor	5103-73-1	100			
		trans-nonachlor	39765-80-5	100			
		a-chlordane	5103-71-9	100			
		g-chlordane	5103-74-2	100			
LK3-101401	Custom Pesticide Mix #3, 6-6, Various Concentrations				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 40,50
LK3-101401-5pck	Custom Pesticide Mix #3, 6-6, Various Concentrations					1 ml, Amber Glass	€ 126,00

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		p,p'-DDE	72-55-9	1			
		p,p'-DDT	50-29-3	2			
		p,p'-DDD	72-54-8	2			
		endrin	72-20-8	2			
		endrin aldehyde	7421-93-4	2			
		endrin ketone	53494-70-5	2			
LK3-005101	Custom Pesticide Mix, 3-34	p,p'-DDT	50-29-3	50	Hexane	4.5 ml, Vial w/ TFE Cap	€ 72,00
		p,p'-DDD	72-54-8	50			
		p,p'-DDE	72-55-9	50			
LK3-074201	Custom Pesticides Mix, 3-65	mirex	2385-85-5	2000	Cyclohexane	1 ml, Amber Glass	€ 55,50
		heptachlor	76-44-8	2000			
		aldrin	309-00-2	2000			
LK3-016201	Custom Pesticide Solution	EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	5	Hexane	1 ml, Amber Glass	€ 118,50
		metolachlor	51218-45-2	5			
		atrazine	1912-24-9	5			
		tributyl phosphate	126-73-8	10			
		triphenyl phosphate	115-86-6	10			
LK3-109501	Custom Pesticide Solution, 4-1	alloxydim-sodium	66003-55-2	100	Acetonitrile	1 ml, Amber Glass	€ 142,50
		clethodim	99129-21-2	100			
		sethoxydim	74051-80-2	100			
		tralkoxydim	87820-88-0	100			
LK3-239202	Custom Pesticide Mix, 4-3	chlorpyrifos	2921-88-2	200	Toluene	1 ml, Amber Glass	€ 72,00
		desethyl atrazine	6190-65-4	200			
		molinate	2212-67-1	200			
		metolachlor	51218-45-2	200			
LK3-201201	Custom Pesticide Mix, 4-6	trifluralin	1582-09-8	100	Acetone	1 ml, Amber Glass	€ 55,50

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		alachlor	15972-60-8	100			
		simazine	122-34-9	100			
		atrazine	1912-24-9	100			
LK3-102804	Custom Pesticide Mix, 4-12	hexachlorobenzene	118-74-1	10	Methylene Chloric	5 ml, Amber Glass	€ 118,50
		a-BHC	319-84-6	10			
		g-BHC	58-89-9	10			
		heptachlor	76-44-8	10			
LK3-205303	Pesticide Mix, 4-25	bentazon	25057-89-0	100	Acetone	1 ml, Amber Glass	€ 48,00
		MCPA acid	94-74-6	100			
		MCPA acid	7085-19-0	100			
		2,4-dichlorophenoxyacetic acid	94-75-7	100			
LK3-203203-5pck	Custom Pesticide Mix, 5-03	chlorpyrifos	2921-88-2	100	Acetone	1 ml, Amber Glass	€ 198,00
		dimethenamid	87674-68-8	100			
		metolachlor	51218-45-2	100			
		Oxadiazon	19666-30-9	100			
		Terbuthylazine	5915-41-3	100			
LK3-203101	Custom Pesticide Solution NPII	vernolate	1929-77-7	1000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 87,00
		atrazine	1912-24-9	1000			
		chlorpyrifos	2921-88-2	1000			
		dimethoate	60-51-5	1000			
		Alachlor	15972-60-8	5000			
LK3-045603	Custom Pesticide Mix, 5-37	chlormephos	24934-91-6	100	Acetone: Isooctan	1 ml, Amber Glass	€ 79,50
		Fonophos	944-22-9	100			
		Tolclofos-Methyl	57018-04-9	100			
		methyl parathion	298-00-0	100			
		parathion	56-38-2	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-203001	Custom Pesticide Solution NPII	vernalate	1929-77-7	2	Hexane	1 ml, Amber Glass	€ 48,00
		atrazine	1912-24-9	2			
		chlorpyrifos	2921-88-2	2			
		dimethoate	60-51-5	2			
		Alachlor	15972-60-8	10			
		1,3-dimethyl-2-nitrobenzene	81-20-9	25			
LK3-223103	Custom Pesticide Solution with Surrogates, 6-1				Isooctane	1 ml, Amber Glass	€ 72,00
		p,p'-DDT	50-29-3	100			
		p,p'-DDD	72-54-8	80			
		o,p'-DDD	53-19-0	80			
		p,p'-DDE	72-55-9	40			
		o,p'-DDT	789-02-6	80			
		o,p'-DDE	3424-82-6	40			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	40			
		1,3-dicyanobenzene	626-17-5	80			
LK3-223203	Custom Pesticide Solution with Surrogates, 6-1				Isooctane	1 ml, Amber Glass	€ 72,00
		p,p'-DDT	50-29-3	100			
		p,p'-DDD	72-54-8	80			
		o,p'-DDD	53-19-0	80			
		p,p'-DDE	72-55-9	40			
		o,p'-DDT	789-02-6	80			
		o,p'-DDE	3424-82-6	40			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	40			
		decachlorobiphenyl (BZ 209)	626-17-5	80			
LK3-223001	Custom Pesticide Solution, 10 mg/L, 6-1				Cyclohexane	1 ml, Amber Glass	€ 31,50
LK3-223002	Custom Pesticide Solution, 500 mg/L, 6-2				Methanol	1 ml, Amber Glass	€ 40,50
		p,p'-DDT	50-29-3	10			
		p,p'-DDD	72-54-8	10			
		o,p'-DDD	53-19-0	10			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		p,p'-DDE	72-55-9	10			
		o,p'-DDT	789-02-6	10			
		o,p'-DDE	3424-82-6	10			
LK3-239001	Custom Pesticide Solution, 6-2	atrazine	1912-24-9	10	Acetone	1 ml, Amber Glass	€ 79,50
		desethyl atrazine	6190-65-4	10			
		Alachlor	15972-60-8	10			
		Simazine	122-34-9	10			
		Terbuthylazine	5915-41-3	10			
		metazachlor	67129-08-2	10			
LK3-044501-5pck	Custom Pesticide Solution, 6-2	fenoxon	80-38-6	100	Isooctane	1 ml, Amber Glass	€ 591,00
		captan	133-06-2	100			
		Folpet	133-07-3	100			
		Bifenthrin	82657-04-3	100			
		acetamiprid	135410-20-7	100			
		famoxadone	131807-57-3	100			
LK3-203403	Custom Pesticide Mix, 6-22				Acetone	1 ml, Amber Glass	€ 55,50
LK3-203403-5pck	Custom Pesticide Mix, 6-22					5 x 1 ml, Amber	€ 220,50
		dimethoate	60-51-5	100			
		atrazine	1912-24-9	100			
		propazine	139-40-2	100			
		Guthion®	86-50-0	100			
		pirimiphos methyl	29232-93-7	100			
		fenthion	55-38-9	100			
LK3-233203-5pck	*Custom Pesticide Mix, 6-24	flufenacet	142459-58-3	100	Acetone	1 ml, Amber Glass	€ 496,50
		isoxaflutole	141112-29-0	100			
		propachlor	1918-16-7	100			
		Propanil	709-98-8	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Oxadiazon	19666-30-9	100			
		Isofenphos	25311-71-1	100			
LK3-201601	Custom Pesticide Mix, 6-26	o,p'-DDD	53-19-0	1000	Methanol	1 ml, Amber Glass	€ 552,00
		o,p'-DDE	3424-82-6	1000			
		o,p'-DDT	789-02-6	1000			
		g-chlordane	5103-74-2	1000			
		a-chlordane	5103-71-9	1000			
		oxychlordane	27304-13-8	1000			
LK3-008401	*Custom Pesticide Mix, 6-27	Dicamba	1918-00-9	1000	Acetone	1 ml, Amber Glass	€ 87,00
		picloram	1918-02-1	1000			
		2,4-dichlorophenoxyacetic acid	94-75-7	1000			
		2,4,5-trichlorophenoxyacetic acid	93-76-5	1000			
		dinoseb	88-85-7	1000			
		diclofop-methyl	51338-27-3	1000			
LK3-203603	Pesticide Mix, 100 mg/L, 6-35				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 55,50
LK3-203603-5pck	Pesticide Mix, 100 mg/L, 6-35					5 x 1 ml, Amber	€ 220,50
LK3-203601	Pesticide Mix, 1,000 mg/L, 6-35					1 ml, Amber Glass	€ 103,50
LK3-203601-5pck	Pesticide Mix, 1,000 mg/L, 6-35					5 x 1 ml, Amber	€ 409,50
		p,p'-DDT	50-29-3	100			
		hexachlorobenzene	118-74-1	100			
		aldrin	309-00-2	100			
		pirimiphos methyl	29232-93-7	100			
		methyl parathion	298-00-0	100			
		chloropyrifos-methyl	5598-13-0	100			
LK3-324003-4pck	Pesticide Mix, 6-78	Baythroid (mixture of isomers)	68359-37-5	100	Acetonitrile	1 ml, Amber Glass	€ 276,00
		lambda cyhalothrin	91465-08-6	100			
		cypermethrin	52315-07-8	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		deltamethrin	52918-63-5	100			
		esfenvalerate	66230-04-4	100			
		fenvalerate	51630-58-1	100			
LK3-324501-4pck	Custom Pesticide Mix, 6-79	2,4,5-T	93-76-5	1000	Acetonitrile	1 ml, Amber Glass	€ 465,00
		2,4-dichlorophenoxyacetic acid	94-75-7	1000			
		dichlorprop (2,4-DP)	120-36-5	1000			
		MCPA acid	94-74-6	1000			
		MCPB	94-81-5	1000			
		MCPB (Mecoprop)	93-65-2	1000			
LK3-003204-5pck	Pesticide Mix 6-99	chlorpyrifos	2921-88-2	200	Acetone	1 ml, Amber Glass	€ 552,00
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	200			
		fonofos	944-22-9	200			
		phorate	298-02-2	200			
		terbufos	13071-79-9	200			
		triallate	2303-17-5	200			
LK3-202101	Custom Pesticide Solution NPI, methyl parathion	methyl parathion	298-00-0	500	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 118,50
		dichlorvos	62-73-7	2500			
		diazinon	333-41-5	500			
		disulfoton	298-04-4	500			
		malathion	121-75-5	500			
		Carbophenothion	786-19-6	500			
		Guthion®	86-50-0	500			
LK3-303003	Pesticide Mix, 7-106	methyl parathion	298-00-0	100	Methanol	1 ml, Amber Glass	€ 55,50
		parathion	56-38-2	100			
		malathion	121-75-5	100			
		dimethoate	60-51-5	100			
		dichlorvos	62-73-7	100			
		trichlorfon	52-68-6	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Demeton O&S (24.4% O & 67.2% S)	8065-48-3	100			
LK3-305505	Pesticide Mix, 7-12	mirex	2385-85-5	250	Methanol	2 ml, Amber Glass	€ 355,50
		oxychlordane	27304-13-8	250			
		chlorpyrifos	2921-88-2	250			
		cis-nonachlor	5103-73-1	250			
		trans-nonachlor	39765-80-5	250			
		hexachlorobutadiene	87-68-3	250			
		hexachloroethane	67-72-1	250			
LK3-046001-5pck	Custom Pesticide Solution, 7-4	fipronil	120068-37-3	100	Isooctane:Toluene	1 ml, Amber Glass	€ 552,00
		Systhane (Myclobutanil)	88671-89-0	100			
		endosulfan sulfate	1031-07-8	100			
		esfenvalerate	66230-04-4	100			
		tetradifon	116-29-0	100			
		pyridaben	96489-71-3	100			
		indoxacarb	173584-44-6	100			
LK3-046501-5pck	Custom Pesticide Solution, 7-6	fenhexamid	126833-17-8	100	Isooctane:Toluene	1 ml, Amber Glass	€ 615,00
		Permethrin (mixture of isomers)	52645-53-1	100			
		allethrin	584-79-2	100			
		tetramethrin	7696-12-0	100			
		Cypermethrin	52315-07-8	100			
		fenvalerate	51630-58-1	100			
		Deltamethrin	52918-63-5	100			
LK3-239903	Custom Pesticide Solution, 100 mg/L, 7-26				Methanol	1 ml, Amber Glass	€ 72,00
LK3-239902	Custom Pesticide Solution, 1,000 mg/L, 7-26					1 ml, Amber Glass	€ 103,50
LK3-239902-5pck	Custom Pesticide Solution, 1,000 mg/L, 7-26					5 x 1 ml, Amber	€ 496,50
		alachlor	15972-60-8	1000			
		desethyl atrazine	6190-65-4	1000			
		terbuthylazine	5915-41-3	1000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		terbuthylazine-desethyl	30125-63-4	1000			
		Cyanazine	21725-46-2	1000			
		metolachlor	51218-45-2	1000			
		pendimethalin	40487-42-1	1000			
LK3-046703	*Custom Pesticide Solution, 7-26				Acetone	1 ml, Amber Glass	€ 87,00
LK3-046703-5pck	*Custom Pesticide Solution, 7-26					5 x 1 ml, Amber	€ 346,50
		Bifenthrin	82657-04-3	100			
		Baythroid (mixture of isomers)	68359-37-5	100			
		esfenvalerate	66230-04-4	100			
		fenvalerate	51630-58-1	100			
		lambda cyhalothrin	91465-08-6	100			
		Permethrin	52645-53-1	100			
		Cypermethrin	52315-07-8	500			
LK3-239701-5pck	*Custom Pesticide Mix, 7-10	pebulate	1114-71-2	1000	Acetone	1 ml, Amber Glass	€ 748,50
		phorate	298-02-2	1000			
		propyzamide (pronamide)	23950-58-5	1000			
		disulfoton	298-04-4	1000			
		triallate	2303-17-5	1000			
		methyl parathion	298-00-0	1000			
		parathion	56-38-2	1000			
LK3-213501	Custom Pesticide Mix, 7-27	parathion	56-38-2	1000	Methanol	1 ml, Amber Glass	€ 103,50
		Abate®	3383-96-8	1000			
		2,4-dichlorophenol	120-83-2	1000			
		methyl parathion	298-00-0	1000			
		malathion	121-75-5	1000			
		pentachlorophenol	87-86-5	1000			
		2,3,4,6-tetrachlorophenol	58-90-2	1000			
LK3-305705	Custom Pesticide Mix, 7-37	endrin ketone	53494-70-5	10	Isooctane	1 ml, Amber Glass	€ 72,00

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		a-chlordane	5103-71-9	10			
		g-chlordane	5103-74-2	10			
		oxychlordane	27304-13-8	10			
		endosulfan sulfate	1031-07-8	10			
		o,p'-DDD	53-19-0	10			
		o,p'-DDE	3424-82-6	10			
LK3-051806	Custom Pesticide Mix, 7-56				Methanol	1 ml, Amber Glass	€ 87,00
LK3-051806-5pck	Custom Pesticide Mix, 7-56					5 x 1 ml, Amber	€ 355,50
		bis(2-ethylhexyl)phthalate	117-81-7	1000			
		pentachlorophenol	87-86-5	1000			
		2,6-Dichlorobenzonitrile (Dichlobenil)	1194-65-6	1000			
		2,6-dichlorobenzamide	2008-58-4	1000			
		epichlorohydrin	106-89-8	1000			
		malathion	121-75-5	1000			
		chlorpyrifos	2921-88-2	1000			
LK3-320103-4pck	Custom Pesticide mix, 7-70	Isfenphos	25311-71-1	100	Acetone	1 ml, Amber Glass	€ 355,50
		ethion	563-12-2	100			
		parathion	56-38-2	100			
		fensulfothion	115-90-2	100			
		bendiocarb	22781-23-3	100			
		phosphamidon	13171-21-6	100			
		Phenthoate	2597-03-7	100			
LK3-003301	Pesticide Mix 7-99				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 198,00
LK3-003301-5pck	Pesticide Mix 7-99					5 x 1 ml, Amber	€ 780,00
		2-phenylphenol	90-43-7	1000			
		terbuthylazine-desethyl	30125-63-4	1000			
		terbuthylazine	5915-41-3	1000			
		heptachlor	76-44-8	1000			
		fenthion	55-38-9	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endosulfan I	959-98-8	1000			
		o,p'-DDT	789-02-6	1000			
LK3-208001	Custom Pesticide Mix, 8-1	methamidophos	10265-92-6	100	Hexane	1 ml, Amber Glass	€ 63,00
		phosmet	732-11-6	100			
		methidathion	950-37-8	50			
		atrazine	1912-24-9	250			
		bladex	21725-46-2	250			
		molinate	2212-67-1	250			
		Simazine	122-34-9	250			
		acetochlor	34256-82-1	250			
LK3-202001	Custom Pesticide Solution NPI,	methyl parathion	298-00-0	2	Hexane	1 ml, Amber Glass	€ 55,50
		dichlorvos	62-73-7	10			
		diazinon	333-41-5	2			
		disulfoton	298-04-4	2			
		malathion	121-75-5	2			
		Carbophenothion	786-19-6	2			
		Guthion®	86-50-0	2			
		1,3-dimethyl-2-nitrobenzene	81-20-9	25			
LK3-208101	Custom Pesticide Mix, 8-2, Various Concentrations				Hexane	1 ml, Amber Glass	€ 63,00
LK3-594501-4pck	Custom Pesticide Mix, 8-128				Methanol	1 ml, Amber Glass	€ 378,00
		methamidophos	10265-92-6	100			
		phosmet	732-11-6	50			
		methidathion	950-37-8	100			
		atrazine	1912-24-9	250			
		bladex	21725-46-2	250			
		molinate	2212-67-1	250			
		simazine	122-34-9	250			
		thiobencarb	28249-77-6	250			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-048101	Custom Pesticide Solution, 8-7				Acetone	1 ml, Amber Glass	€ 198,00
LK3-048101-5pck	Custom Pesticide Solution, 8-7					5 x 1 ml, Amber	€ 709,50
		Methiocarb	2032-65-7	100			
		formothion	2540-82-1	100			
		Carbaryl	63-25-2	100			
		flusilazole	85509-19-9	100			
		Benalaxyl	71626-11-4	100			
		etoxazole	153233-91-1	100			
		Phosalone	2310-17-0	100			
		Pyrazophos	13457-18-6	100			
LK3-033904-4pck	Custom Pesticide Mix, 8-28	dimethomorph (mixture of isomers)	110488-70-5	200	Acetonitrile	1 ml, Amber Glass	€ 220,50
		pyrimethanil	53112-28-0	200			
		Metalaxyl	57837-19-1	200			
		Vinclozolin	50471-44-8	200			
		tebuconazol	107534-96-3	200			
		cyprodinil	121552-61-2	200			
		penconazol	66246-88-6	200			
		chlorpyrifos	2921-88-2	200			
LK3-205903-4pck	*Custom Pesticide Solution, 100 mg/L, 8-40				Methanol	4 x 1 ml, Amber	€ 220,50
LK3-205901-4pck	*Custom Pesticide Solution, 1,000 mg/L, 8-40					4 x 1 ml, Amber	€ 346,50
		2,4,5-T	93-76-5				
		2,4,5-TP (Silvex)	93-72-1				
		2,4-dichlorophenoxyacetic acid	94-75-7				
		brominal	1689-84-5				
		Dicamba	1918-00-9				
		diclofop-methyl	51338-27-3				
		MCPA acid	94-74-6				
		picloram	1918-02-1				
LK3-051405	Custom Pesticide Solution, 8-4 chlorpyrifos		2921-88-2	500	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 150,00

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		dieldrin	60-57-1	500			
		endosulfan (alpha + beta)	115-29-7	500			
		heptachlor	76-44-8	500			
		heptachlor epoxide isomer A	28044-83-9	500			
		heptachlor epoxide isomer B	1024-57-3	500			
		hexachlorobenzene	118-74-1	500			
		methoxychlor	72-43-5	500			
LK3-114202	Custom Pesticide Mix, 8-41	ametryne	834-12-8	2000	Acetone	1 ml, Amber Glass	€ 213,00
		disulfoton	298-04-4	2000			
		ethoprophos (prophos)	13194-48-4	2000			
		Fenarimol	60168-88-9	2000			
		fluridone	59756-60-4	2000			
		mevinphos (Phosdrin (TM))	7786-34-7	2000			
		napropamide	15299-99-7	2000			
		norflurazon	27314-13-2	2000			
LK3-114701	Custom Pesticide Mix, 8-44				M-t-BE	1 ml, Amber Glass	€ 135,00
LK3-114701-5pck	Custom Pesticide Mix, 8-44					5 x 1 ml, Amber	€ 261,00
		Benfluralin (Benefin)	1861-40-1	1000			
		Isopropalin	33820-53-0	1000			
		Oxadiazon	19666-30-9	1000			
		Oxyfluorfen	42874-03-3	1000			
		pendimethalin	40487-42-1	1000			
		profluralin	26399-36-0	1000			
		propachlor	1918-16-7	1000			
		trifluralin	1582-09-8	1000			
LK3-045501-5pck	Custom Pesticide Solution, 9-5	triazophos	24017-47-8	100	Isooctane:Toluene	1 ml, Amber Glass	€ 937,50
		pyrimethanil	53112-28-0	100			
		Metalaxyl	57837-19-1	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		fenpropidin	67306-00-7	100			
		fludioxonil	131341-86-1	100			
		oxadixyl	77732-09-3	100			
		tebufenpyrad	119168-77-3	100			
		fenbuconazole	114369-43-6	100			
		cyproconazole (mixture of isomers)	113096-99-4	100			
		bitertanole (mixture of isomers)	55179-31-2	100			
LK3-047001-5pck	Custom Pesticide Solution, 9-4	bupirimate	41483-43-6	100	Isooctane:Toluene	1 ml, Amber Glass	€ 937,50
		quinoxifen	124495-18-7	100			
		nuarimol	63284-71-9	100			
		danitol	39515-41-8	100			
		lambda cyhalothrin	91465-08-6	100			
		acrinathrin [ISO]	101007-06-1	100			
		Baythroid (mixture of isomers)	68359-37-5	100			
		flucythrinate (mixture of isomers)	70124-77-5	100			
		tau-fluvalinate	102851-06-9	100			
LK3-048001-5pck	Custom Pesticide Solution, 9-4	Methiocarb	2032-65-7	100	Isooctane:Toluene	1 ml, Amber Glass	€ 828,00
		formothion	2540-82-1	100			
		Carbaryl	63-25-2	100			
		fenchlorphos (Ronnell)	299-84-3	100			
		flusilazole	85509-19-9	100			
		Benalaxyl	71626-11-4	100			
		etoxazole	153233-91-1	100			
		phosalone	2310-17-0	100			
		Pyrazophos	13457-18-6	100			
LK3-046601	Pesticide Solution, 9-7				Acetone	1 ml, Amber Glass	€ 229,50
LK3-046601-5pck	Custom Pesticide Solution, 9-7					5 x 1 ml, Amber	€ 913,50
		fenhexamid	126833-17-8	100			
		allethrin	584-79-2	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		tetramethrin	7696-12-0	100			
		Permethrin	52645-53-1	100			
		Cypermethrin	52315-07-8	100			
		fenvalerate	51630-58-1	100			
		Deltamethrin	52918-63-5	100			
		Baythroid (mixture of isomers)	68359-37-5	100			
		trifluralin	1582-09-8	100			
LK3-239801	Pesticide Mix, 9-9	p-phenylenediamine	106-50-3	2000	Methylene Chloric	1 ml, Amber Glass	€ 213,00
		3,3'-dimethylbenzidine	119-93-7	2000			
		1,4-naphthoquinone	130-15-4	2000			
		di-allate	2303-16-4	2000			
		propyzamide (pronamide)	23950-58-5	2000			
		isodrin	465-73-6	2000			
		chlorobenzilate	510-15-6	2000			
		kepone	143-50-0	2000			
		1,2-dinitrobenzene	528-29-0	2000			
LK3-233701	*Custom Pesticide Mix, 9-12	simazine	122-34-9	500	Methylene Chloric	1 ml, Amber Glass	€ 72,00
		Cyanazine	21725-46-2	1000			
		atrazine	1912-24-9	1000			
		prometon	1610-18-0	1000			
		prometryn	7287-19-6	1000			
		pendimethalin	40487-42-1	1000			
		chlorpyrifos	2921-88-2	1000			
		diazinon	333-41-5	1000			
		Isofenphos	25311-71-1	1000			
LK3-016401	Custom Pesticide Mix, 9-14	pendimethalin	40487-42-1	200	Methanol	1 ml, Amber Glass	€ 103,50
		sulprofos (Bolstar®)	35400-43-2	200			
		trifluralin	1582-09-8	200			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	200			
		diphenamid	957-51-7	200			
		tribufos	78-48-8	200			
		tributyl phosphate	126-73-8	200			
		triphenyl phosphate	115-86-6	200			
		ethion	563-12-2	200			
LK3-208304	Custom Pesticide Mix, 9-15	benthiocarb	28249-77-6	250	Acetone	1 ml, Amber Glass	€ 79,50
		carbophenothion	786-19-6	250			
		dioxathion (mixture of isomers)	78-34-2	250			
		diphenamid	957-51-7	250			
		phosmet	732-11-6	250			
		phosalone	2310-17-0	250			
		prometon	1610-18-0	250			
		simazine	122-34-9	250			
		methidathion	950-37-8	250			
LK3-043103	Custom Pesticide Mix, 9-23				Acetone	1 ml, Amber Glass	€ 55,50
LK3-043103-5pck	Custom Pesticide Mix, 9-23					1 ml, Amber Glass	€ 276,00
		profenofos	41198-08-7	100			
		triazophos	24017-47-8	100			
		Pirimiphos-ethyl	23505-41-1	100			
		pirimiphos methyl	29232-93-7	100			
		quinalphos	13593-03-8	100			
		Pyridaphenthion	119-12-0	100			
		Omethoate	1113-02-6	100			
		sulfotep	3689-24-5	100			
		tetrachlorvinphos (Rabon)	961-11-5	100			
LK3-253305-4pck	Custom Pesticide Mix, 9-31	Carbaryl	63-25-2	10	Acetone	1 ml, Amber Glass	€ 567,00
		chlorpyrifos	2921-88-2	10			
		dimethoate	60-51-5	10			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endosulfan sulfate	1031-07-8	10			
		fenitrothion	122-14-5	10			
		fenthion	55-38-9	10			
		fludioxonil	131341-86-1	10			
		methidathion	950-37-8	10			
		Procymidone	32809-16-8	10			
LK3-004201	P1, EPA Pesticide Standard, 5, 10, 20 & 50 mg/L, 10-1				Tol.:Hex.:Ethyl Ac	1 ml, Amber Glass	€ 79,50
LK3-004204	Custom Pesticide Standard, 100, 200, 400 & 1,600 mg/L, 10-1				Acetone	1.2 ml, Amber Glass	€ 55,50
LK3-004205	Custom Pesticide Standard, 250, 500, 1,000 & 4,000 mg/L, 10-1					1.2 ml, Amber Glass	€ 72,00
LK3-004206	Custom Pesticide Standard, 500, 1,000, 2,000 & 8,000 mg/L, 10-1					1.2 ml, Amber Glass	€ 87,00
		g-BHC	58-89-9	5			
		aldrin	309-00-2	10			
		dieldrin	60-57-1	10			
		endosulfan I	959-98-8	10			
		heptachlor	76-44-8	10			
		heptachlor epoxide	1024-57-3	10			
		p,p'-DDT	50-29-3	20			
		endosulfan II	33213-65-9	20			
		endrin aldehyde	7421-93-4	25			
		methoxychlor	72-43-5	50			
LK3-043501-5pck	Custom Pesticide Solution, 10 ethoprophos (prophos)		13194-48-4	100	Isooctane:Toluene	1 ml, Amber Glass	€ 622,50
		Ethoxyquin	91-53-2	100			
		monocrotophos	6923-22-4	100			
		Fonophos	944-22-9	100			
		phosphamidon	13171-21-6	100			
		malathion	121-75-5	100			
		Bromophos methyl	2104-96-3	100			
		Bromophos ethyl	4824-78-6	100			
		phosmet	732-11-6	100			
		iprovalicarb	140923-17-7	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-047501-5pck	Custom Pesticide Solution, 10	triadimefon	43121-43-3	100	Isooctane	1 ml, Amber Glass	€ 1.063,50
		tetraconazole	112281-77-3	100			
		penconazol	66246-88-6	100			
		clazolinate	84332-86-5	100			
		Hexaconazole	79983-71-4	100			
		diclobutrazol (mixture of stereo isomers)	75736-33-3	100			
		alpha-cypermethrin	67375-30-8	100			
		Propiconazol (mixture of isomers)	60207-90-1	100			
		Difenoconazole (mixture of isomers)	119446-68-3	100			
		dimethomorph (mixture of isomers)	110488-70-5	100			
LK3-245504	Custom Pesticide Standard Mix B, 100, 200 & 400 mg/L, 10-17				Acetone	1.2 ml, Amber Glass	€ 55,50
LK3-245505	Custom Pesticide Standard Mix B, 200, 500 & 1,000 mg/L, 10-17					1.2 ml, Amber Glass	€ 72,00
LK3-245506	Custom Pesticide Standard Mix B, 500, 1,000 & 2,000 mg/L, 10-17					1.2 ml, Amber Glass	€ 87,00
		a-BHC	319-84-6	100			
		d-BHC	319-86-8	200			
		a-chlordane	5103-71-9	200			
		g-chlordane	5103-74-2	200			
		p,p'-DDT	50-29-3	200			
		endrin	72-20-8	400			
		p,p'-DDD	72-54-8	400			
		endosulfan sulfate	1031-07-8	400			
		endrin ketone	53494-70-5	400			
		b-BHC	319-85-7	400			
LK3-073403	Custom Pesticide Solution, 10	chlorpyrifos-methyl	5598-13-0	100	Ethyl Acetate	1 ml, Amber Glass	€ 166,50
		parathion	56-38-2	100			
		malathion	121-75-5	100			
		methyl parathion	298-00-0	100			
		flucythrinate (mixture of isomers)	70124-77-5	100			
		g-BHC	58-89-9	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endrin	72-20-8	100			
		o,p'-DDD	53-19-0	100			
		p,p'-DDT	50-29-3	100			
		o,p'-DDT	789-02-6	100			
LK3-028601	Custom Pesticide Solution, 10-4				Acetone	1 ml, Amber Glass	€ 150,00
LK3-028601-5pck	Custom Pesticide Solution, 10-4					5 x 1 ml, Amber	€ 748,50
		butylate	2008-41-5	1000			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000			
		molinate	2212-67-1	1000			
		vernolate	1929-77-7	1000			
		thiodicarb	59669-26-0	1000			
		pebulate	1114-71-2	1000			
		carbofuran phenol	1563-38-8	2000			
		promecarb	2631-37-0	1000			
		carbosulfan	55285-14-8	1000			
		triallate	2303-17-5	1000			
LK3-045203-5pck	Custom Pesticide Mix, 10-21	fenvalerate	51630-58-1	100	Acetone	1 ml, Amber Glass	€ 315,00
		nuarimol	63284-71-9	100			
		Quinalphos	13593-03-8	100			
		o,p'-DDD	53-19-0	100			
		endosulfan II	33213-65-9	100			
		endosulfan I	959-98-8	100			
		captan	133-06-2	100			
		Propiconazol (mixture of isomers)	60207-90-1	100			
		triadimefon	43121-43-3	100			
		pendimethalin	40487-42-1	100			
LK3-103003	Custom Pesticide Mix, 10-24				Hexane	1 ml, Amber Glass	€ 55,50
LK3-103003-4pck	Custom Pesticide Mix, 10-24					4 x 1 ml, Amber	€ 220,50

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		g-BHC	58-89-9	100			
		endosulfan II	33213-65-9	100			
		endosulfan I	959-98-8	100			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	100			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	100			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	100			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	100			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	100			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	100			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	100			
LK3-073103	Custom Pesticide Mix, 10-25	chloropyrifos-methyl	5598-13-0	100	Ethyl Acetate	1 ml, Amber Glass	€ 166,50
		parathion	56-38-2	100			
		malathion	121-75-5	100			
		methyl parathion	298-00-0	100			
		g-BHC	58-89-9	100			
		endrin	72-20-8	100			
		p,p'-DDD	72-54-8	100			
		p,p'-DDT	50-29-3	100			
		o,p'-DDT	789-02-6	100			
		flucythrinate (mixture of isomers)	70124-77-5	100			
LK3-053204	*Custom Pesticide Mix, 100 mg/L, 10-5				Methanol	1 ml, Amber Glass	€ 63,00
LK3-013501	*Custom Pesticide Mix, 100 mg/L, 10-5				Acetonitrile	1 ml, Amber Glass	€ 87,00
		aldicarb	116-06-3	100			
		aldicarb sulfone	1646-88-4	100			
		aldicarb sulfoxide	1646-87-3	100			
		Carbaryl	63-25-2	100			
		carbofuran	1563-66-2	100			
		3-hydroxycarbofuran	16655-82-6	100			
		Methiocarb	2032-65-7	100			
		methomyl	16752-77-5	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		oxamyl	23135-22-0	100			
		propoxur	114-26-1	100			
LK3-051701	*Custom Pesticide Mix, 10-56 1ml				Methanol	1 ml, Amber Glass	€ 103,50
LK3-051701-4pck	*Custom Pesticide Mix, 10-56 4 X 1ml					4 x 1 ml, Amber	€ 283,50
		bis(2-ethylhexyl)phthalate	117-81-7	1000			
		pentachlorophenol	87-86-5	1000			
		2,6-Dichlorobenzonitrile (Dichlobenil)	1194-65-6	1000			
		2,6-dichlorobenzamide	2008-58-4	1000			
		epichlorohydrin	106-89-8	1000			
		malathion	121-75-5	1000			
		chlorpyrifos	2921-88-2	1000			
		epoxiconazole (ISO)	133855-98-8	1000			
		dimethoate	60-51-5	1000			
		diazinon	333-41-5	1000			
LK3-250003	Custom Pesticide Mix, 10 mg/L, 11-1				Acetonitrile	1 ml, Amber Glass	€ 118,50
LK3-250001	Custom Pesticide Mix, 1,000 mg/L, 11-1					1 ml, Amber Glass	€ 370,50
		chloridazon	1698-60-8				
		bentazon	25057-89-0				
		2,4-D	94-75-7				
		MCPA acid	94-74-6				
		diuron	330-54-1				
		isoproturon	34123-59-6				
		Triclopyr	55335-06-3				
		MCPP acid	7085-19-0				
		linuron	330-55-2				
		chlorfenvinphos (mixture of z and E isomers)	470-90-6				
		trifluralin	1582-09-8				
LK3-105302	Custom Pesticide Solution, 50,100 & 500 mg/L, 11-2				Isooctane	1 ml, Amber Glass	€ 79,50
LK3-105301	Custom Pesticide Solution, 200,400 & 2,000 mg/L, 11-2					1 ml, Amber Glass	€ 79,50

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		a-BHC	319-84-6	50			
		decachlorobiphenyl	2051-24-3	100			
		dieldrin	60-57-1	100			
		endosulfan I	959-98-8	50			
		endrin	72-20-8	100			
		g-BHC	58-89-9	50			
		heptachlor	76-44-8	50			
		methoxychlor	72-43-5	500			
		p,p'-DDD	72-54-8	100			
		p,p'-DDT	50-29-3	100			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	50			
LK3-044001-5pck	Custom Pesticide Solution, 11	Dichloran	99-30-9	100	Isooctane	1 ml, Amber Glass	€ 756,00
		Vinclozolin	50471-44-8	100			
		Dichlofluanid	1085-98-9	100			
		tolyfluanid	731-27-1	100			
		procymidone	32809-16-8	100			
		Imazalil	35554-44-0	100			
		kresoxim methyl	143390-89-0	100			
		trifloxystrobin	141517-21-7	100			
		iprodione	36734-19-7	100			
		Fenarimol	60168-88-9	100			
		azoxystrobin	131860-33-8	100			
LK3-033503	Custom Pesticide Mix, 11-14				Acetone	1 ml, Amber Glass	€ 87,00
LK3-033503-5pck	Custom Pesticide Mix, 11-14					5 x 1 ml, Amber	€ 346,50
		hexachlorobenzene	118-74-1	100			
		o,p'-DDE	3424-82-6	100			
		g-BHC	58-89-9	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDD	72-54-8	100			
		Kelthane ® (Dicofol)	115-32-2	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		penconazol	66246-88-6	100			
		Procymidone	32809-16-8	100			
		p,p'-DDT	50-29-3	100			
		acrinathrin [ISO]	101007-06-1	100			
		Baythroid (mixture of isomers)	68359-37-5	100			
LK3-052203-4pck	*Custom Pesticide Mix, 100 mg/L , 11-19				Methanol	4 x 1 ml, Amber	€ 283,50
LK3-052201-4pck	*Custom Pesticide Mix, 1,000 mg/L , 11-19					4 x 1 ml, Amber	€ 409,50
		diazinon	333-41-5				
		dimethoate	60-51-5				
		fenitrothion	122-14-5				
		chlorpyrifos	2921-88-2				
		malathion	121-75-5				
		flamprop-methyl (metavane)	52756-25-9				
		methyl parathion	298-00-0				
		parathion	56-38-2				
		Propanil	709-98-8				
		triallate	2303-17-5				
		trifluralin	1582-09-8				
LK3-253001	*Custom Pesticide Mix, 11-20				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 229,50
LK3-253002	*Custom Pesticide Mix, 11-21				Acetone	1 ml, Amber Glass	€ 229,50
		Acephate	30560-19-1	100			
		Guthion Ethyl	2642-71-9	100			
		diazinon	333-41-5	100			
		dimethoate	60-51-5	100			
		fenitrothion	122-14-5	100			
		fenthion	55-38-9	100			
		malathion	121-75-5	100			
		monocrotophos	6923-22-4	100			
		phosalone	2310-17-0	100			
		profenofos	41198-08-7	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		quinalphos	13593-03-8	100			
LK3-114102	Custom Pesticide Mix, 11-23	ametryne	834-12-8	2000	Acetone	1 ml, Amber Glass	€ 276,00
		disulfoton	298-04-4	2000			
		ethoprophos (prophos)	13194-48-4	2000			
		Fenarimol	60168-88-9	2000			
		fluridone	59756-60-4	2000			
		hexazinon	51235-04-2	2000			
		mevinphos (Phosdrin (TM))	7786-34-7	2000			
		napropamide	15299-99-7	2000			
		norflurazon	27314-13-2	2000			
		prometon	1610-18-0	2000			
		prometryn	7287-19-6	2000			
LK3-043703	Custom Pesticide Mix, 11-24	g-BHC	58-89-9	100	Hexane	1 ml, Amber Glass	€ 166,50
		endrin	72-20-8	100			
		o,p'-DDD	53-19-0	100			
		p,p'-DDT	50-29-3	100			
		o,p'-DDT	789-02-6	100			
		methyl parathion	298-00-0	100			
		malathion	121-75-5	100			
		coumaphos	56-72-4	100			
		hexachlorobenzene	118-74-1	100			
		isodrin	465-73-6	100			
		chloropyrifos-methyl	5598-13-0	100			
LK3-048501-5pck	Custom Pesticide Mix, 11-4	propoxur	114-26-1	100	Isooctane:Toluene	1 ml, Amber Glass	€ 937,50
		ethiofencarb	29973-13-5	100			
		methyl parathion	298-00-0	100			
		fenitrothion	122-14-5	100			
		mepanipyrim	110235-47-7	100			
		ethion	563-12-2	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		tebuconazol	107534-96-3	100			
		Pyridaphenthion	119-12-0	100			
		benfuracarb	82560-54-1	100			
		furathiocarb	65907-30-4	100			
		etaconazole	60207-93-4	100			
LK3-246201	Custom Pesticide Mix, 11-6				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 103,50
LK3-246201-5pck	Custom Pesticide Mix, 11-6					5 x 1 ml, Amber	€ 409,50
		dioxathion	78-34-2	1000			
		trichlorfon	52-68-6	1000			
		thionazine (zinophos)	297-97-2	1000			
		Carbophenothion	786-19-6	1000			
		ethion	563-12-2	1000			
		chloropyrifos-methyl ester	5598-13-0	1000			
		famphur	52-85-7	1000			
		Terbufos	13071-79-9	1000			
		fenitrothion	122-14-5	1000			
		dichlofenthion	97-17-6	1000			
		Fonophos	944-22-9	1000			
LK3-043001-5pck	Custom Pesticide Mix, 12-2	methamidophos	10265-92-6	100	Isooctane	1 ml, Amber Glass	€ 622,50
		acephate	30560-19-1	100			
		Omethoate	1113-02-6	100			
		disulfoton	298-04-4	100			
		chloropyrifos-methyl	5598-13-0	100			
		pirimiphos methyl	29232-93-7	100			
		chlorpyrifos	2921-88-2	100			
		Pirimiphos-ethyl	23505-41-1	100			
		buprofezin	69327-76-0	100			
		fenoxycarb	72490-01-8	100			
		Guthion®	86-50-0	100			
		Guthion Ethyl	2642-71-9	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-071104	Custom Pesticide Mix, 12-6	EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	10	Methylene Chloric	1 ml, Amber Glass	€ 94,50
		butylate	2008-41-5	10			
		vernolate	1929-77-7	10			
		pebulate	1114-71-2	10			
		molinate	2212-67-1	10			
		Cycloate	1134-23-2	10			
		atrazine	1912-24-9	10			
		Fonophos	944-22-9	10			
		fenitrothion	122-14-5	10			
		parathion	56-38-2	10			
		carbophenothion methyl	953-17-3	10			
		Carbophenothion	786-19-6	10			
LK3-073703	Custom Pesticide Mix, 12-11				Acetone	1 ml, Amber Glass	€ 94,50
LK3-073703-5pck	Custom Pesticide Mix, 12-11					5 x 1 ml, Amber	€ 378,00
		aldrin	309-00-2	100			
		isodrin	465-73-6	100			
		endrin	72-20-8	100			
		dieldrin	60-57-1	100			
		b-BHC	319-85-7	100			
		alachlor	15972-60-8	100			
		chlordane (mixture of isomers)	57-74-9	100			
		flucythrinate (mixture of isomers)	70124-77-5	100			
		Fenarimol	60168-88-9	100			
		chlorothalonil	1897-45-6	100			
		Bifenthrin	82657-04-3	100			
		Systhane (TM)	88671-89-0	100			
LK3-203803	Custom Pesticide Mix, 12-22, 100 mg/L				Hexane	1 ml, Amber Glass	€ 189,00
		a-BHC	319-84-6	100			
		b-BHC	319-85-7	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		o,p'-DDE	3424-82-6	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDD	72-54-8	100			
		chlorpyrifos	2921-88-2	100			
		chlordan (Mix of Isomers)	57-74-9	100			
		pirimiphos methyl	29232-93-7	100			
		dieldrin	60-57-1	100			
		aldrin	309-00-2	100			
		disulfoton	298-04-4	100			
		sulfotep	3689-24-5	100			
LK3-253903-5pck	Custom Pesticide Mix, 100 mg/L, 12-31				Acetone	5 x 1 ml, Amber	€ 630,00
LK3-253901	Custom Pesticide Mix, 1,000 mg/L, 12-31					1 ml, Amber Glass	€ 229,50
		Guthion Ethyl	2642-71-9				
		Guthion®	86-50-0				
		fenthion	55-38-9				
		malathion	121-75-5				
		fenitrothion	122-14-5				
		Pyrazophos	13457-18-6				
		pirimiphos methyl	29232-93-7				
		Pirimiphos-ethyl	23505-41-1				
		methyl parathion	298-00-0				
		parathion	56-38-2				
		chlorfenvinphos (mixture of z and E isomers)	470-90-6				
		phosalone	2310-17-0				
LK3-046403-4pck	Custom Pesticide Mix, 12-32	Permethrin	52645-53-1	100	Methanol	1 ml, Amber Glass	€ 465,00
		allethrin	584-79-2	100			
		Bifenthrin	82657-04-3	100			
		Cypermethrin	52315-07-8	100			
		danitol	39515-41-8	100			
		Baythroid (mixture of isomers)	68359-37-5	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		esfenvalerate	66230-04-4	100			
		fenvalerate	51630-58-1	100			
		tau-fluvalinate	102851-06-9	100			
		Deltamethrin	52918-63-5	100			
		lambda cyhalothrin	91465-08-6	100			
		prallethrin	23031-36-9	100			
LK3-105202	Custom Pesticide Solution, 50 & 100 mg/L, 13-2				Isooctane	1 ml, Amber Glass	€ 87,00
LK3-105201	Custom Pesticide Solution, 200 & 400 mg/L, 13-2					1 ml, Amber Glass	€ 87,00
		a-chlordane	5103-71-9	50			
		aldrin	309-00-2	50			
		b-BHC	319-85-7	50			
		d-BHC	319-86-8	50			
		decachlorobiphenyl	2051-24-3	100			
		endosulfan II	33213-65-9	100			
		endosulfan sulfate	1031-07-8	100			
		endrin aldehyde	7421-93-4	100			
		endrin ketone	53494-70-5	100			
		g-chlordane	5103-74-2	50			
		endosulfan I	959-98-8	50			
		p,p'-DDE	72-55-9	100			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	50			
LK3-239502	Custom Pesticide Solution, 13 alachlor		15972-60-8	500	Methylene Chloric	1 ml, Amber Glass	€ 111,00
		butachlor	23184-66-9	500			
		metolachlor	51218-45-2	500			
		metribuzin	21087-64-9	500			
		simazine	122-34-9	500			
		bladex	21725-46-2	500			
		propachlor	1918-16-7	500			
		propazine	139-40-2	500			
		trifluralin	1582-09-8	500			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		pendimethalin	40487-42-1	500			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	500			
		Fonophos	944-22-9	500			
		bromacil	314-40-9	500			
LK3-079601	Custom Pesticide Mix, 13-7				Acetone	1 ml, Amber Glass	€ 103,50
LK3-07960105	Custom Pesticide Mix, 13-7					5 ml, Amber Glass	€ 457,50
		g-BHC	58-89-9	1000			
		heptachlor	76-44-8	1000			
		aldrin	309-00-2	1000			
		heptachlor epoxide isomer B	1024-57-3	1000			
		endosulfan I	959-98-8	1000			
		dieldrin	60-57-1	1000			
		p,p'-DDE	72-55-9	1000			
		o,p'-DDD	53-19-0	1000			
		endrin	72-20-8	1000			
		endosulfan II	33213-65-9	1000			
		p,p'-DDD	72-54-8	1000			
		o,p'-DDT	789-02-6	1000			
		p,p'-DDT	50-29-3	1000			
LK3-106004	Pesticide Mix, 13-8, Various Concentrations, 5 & 10 mg/L				Hexane	1 ml, Amber Glass	€ 55,50
		b-BHC	319-85-7	5			
		d-BHC	319-86-8	5			
		aldrin	309-00-2	5			
		a-chlordane	5103-71-9	5			
		g-chlordane	5103-74-2	5			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	5			
		heptachlor epoxide isomer B	1024-57-3	5			
		p,p'-DDE	72-55-9	10			
		endosulfan sulfate	1031-07-8	10			
		endrin aldehyde	7421-93-4	10			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endosulfan II	33213-65-9	10			
		endrin ketone	53494-70-5	10			
		decachlorobiphenyl (BZ 209)	2051-24-3	10			
LK3-047703	Custom Pesticide Mix, 13-10				Acetone	1 ml, Amber Glass	€ 103,50
LK3-047703-5pck	Custom Pesticide Mix, 13-10					5 x 1 ml, Amber	€ 409,50
		trifluralin	1582-09-8	100			
		a-BHC	319-84-6	100			
		propyzamide (pronamide)	23950-58-5	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		Vinclozolin	50471-44-8	100			
		heptachlor	76-44-8	100			
		chlorthal-dimethyl (dacthal)	1861-32-1	100			
		heptachlor epoxide isomer A	28044-83-9	100			
		endosulfan sulfate	1031-07-8	100			
		Folpet	133-07-3	100			
		alpha-cypermethrin	67375-30-8	100			
		Hexaconazole	79983-71-4	100			
		phosalone	2310-17-0	100			
LK3-052605	Custom Pesticide Solution, 13	dimethoate	60-51-5	200	Methylene Chloric	1 ml, Amber Glass	€ 213,00
		pendimethalin	40487-42-1	200			
		trifluralin	1582-09-8	200			
		metolachlor	51218-45-2	200			
		atrazine	1912-24-9	200			
		propazine	139-40-2	200			
		propachlor	1918-16-7	200			
		triallate	2303-17-5	200			
		acetochlor	34256-82-1	200			
		metribuzin	21087-64-9	200			
		dimethenamid	87674-68-8	200			
		Command	81777-89-1	200			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		alachlor	15972-60-8	200			
LK3-045303	Custom Pesticide Mix, 13-17				Ethyl Acetate	1 ml, Amber Glass	€ 213,00
LK3-045303-4pck	Custom Pesticide Mix, 13-17					4 x 1 ml, Amber	€ 765,00
		coumaphos	56-72-4	100			
		Tolclofos-Methyl	57018-04-9	100			
		disulfoton	298-04-4	100			
		sulfotep	3689-24-5	100			
		Phenthoate	2597-03-7	100			
		Fenarimol	60168-88-9	100			
		nuarimol	63284-71-9	100			
		Vinclozolin	50471-44-8	100			
		fenvalerate	51630-58-1	100			
		tau-fluvalinate	102851-06-9	100			
		ditalimfos	5131-24-8	100			
		fenamiphos	22224-92-6	100			
		Isofenphos	25311-71-1	100			
LK3-206106-5pck	Custom Pesticide Mix, 13-42	a-BHC	319-84-6	0,4	Isooctane:Toluene	5 ml, Amber Glass	€ 787,50
		2,4,5,6-tetrachloro-m-xylene	877-09-8	0,8			
		b-BHC	319-85-7	0,8			
		d-BHC	319-86-8	0,8			
		aldrin	309-00-2	0,8			
		a-chlordane	5103-71-9	0,8			
		g-chlordane	5103-74-2	0,8			
		p,p'-DDE	72-55-9	0,8			
		endrin	72-20-8	1,6			
		p,p'-DDD	72-54-8	1,6			
		decachlorobiphenyl (BZ 209)	2051-24-3	1,6			
		endosulfan sulfate	1031-07-8	2			
		endrin ketone	53494-70-5	2			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-359504-4pck	Custom Pesticide Mix, 13-67	dichlorvos	62-73-7	10	Acetone	1 ml, Amber Glass	€ 472,50
		ethoprophos (prophos)	13194-48-4	10			
		monocrotophos	6923-22-4	10			
		dimethoate	60-51-5	10			
		chloropyrifos-methyl	5598-13-0	10			
		Tolclofos-Methyl	57018-04-9	10			
		methyl parathion	298-00-0	10			
		malathion	121-75-5	10			
		parathion	56-38-2	10			
		methidathion	950-37-8	10			
		triazophos	24017-47-8	10			
		Pyridaphenthion	119-12-0	10			
		phosalone	2310-17-0	10			
LK3-102401	Pesticide Mix, 14 Compounds	Alachlor	15972-60-8	200	Acetone	1 ml, Amber Glass	€ 213,00
		Benfluralin (Benefin)	1861-40-1	200			
		Captan	133-06-2	200			
		Carbophenothion	786-19-6	200			
		hexachlorobenzene	118-74-1	200			
		Kelthane ® (Dichlobenil)	115-32-2	200			
		Nitrofen	1836-75-5	200			
		o,p'-DDD	53-19-0	200			
		o,p'-DDE	3424-82-6	200			
		o,p'-DDT	789-02-6	200			
		pentachloronitrobenzene	82-68-8	200			
		Schradan (DMPA)	152-16-9	200			
		Chlorothalonil	1897-45-6	200			
		Propachlor	1918-16-7	500			
LK3-079806	Custom Pesticide Mix, 14-11	hexachlorobutadiene	87-68-3	10	Isooctane	1 ml, Amber Glass	€ 55,50
		2,6-Dichlorobenzonitrile (Dichlobenil)	1194-65-6	10			
		trifluralin	1582-09-8	10			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		g-BHC	58-89-9	10			
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10			
		endosulfan I	959-98-8	10			
		endosulfan II	33213-65-9	10			
		2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	10			
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		endosulfan sulfate	1031-07-8	10			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
LK3-119703	Custom Pesticide Mix, 14-14	Guthion Ethyl	2642-71-9	100	Acetone	1 ml, Amber Glass	€ 237,00
		Guthion®	86-50-0	100			
		diazinon	333-41-5	100			
		dimethoate	60-51-5	100			
		fenthion	55-38-9	100			
		phosalone	2310-17-0	100			
		malathion	121-75-5	100			
		fenitrothion	122-14-5	100			
		Pyrazophos	13457-18-6	100			
		pirimiphos methyl	29232-93-7	100			
		Pirimiphos-ethyl	23505-41-1	100			
		methyl parathion	298-00-0	100			
		parathion	56-38-2	100			
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	100			
LK3-047303	Custom Pesticide Mix, 14-23	a-BHC	319-84-6	100	Ethyl Acetate	1 ml, Amber Glass	€ 229,50
		b-BHC	319-85-7	100			
		o,p'-DDE	3424-82-6	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDD	72-54-8	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		heptachlor	76-44-8	100			
		aldrin	309-00-2	100			
		dieldrin	60-57-1	100			
		heptachlor epoxide isomer A	28044-83-9	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		chlorpyrifos	2921-88-2	100			
		chlordane (Mix of Isomers)	57-74-9	100			
		fenitrothion	122-14-5	100			
		pirimiphos methyl	29232-93-7	100			
LK3-215501	Custom Pesticide Mix, 14-26	Alachlor	15972-60-8	1000	Acetone	1 ml, Amber Glass	€ 111,00
		brominal	1689-84-5	1000			
		diazinon	333-41-5	1000			
		phorate	298-02-2	1000			
		trifluralin	1582-09-8	1000			
		atrazine	1912-24-9	1000			
		chlorpyrifos	2921-88-2	1000			
		dimethoate	60-51-5	1000			
		Terbufos	13071-79-9	1000			
		prometryn	7287-19-6	1000			
		desethyl atrazine	6190-65-4	1000			
		bladex	21725-46-2	1000			
		metolachlor	51218-45-2	1000			
		2,4,6-trichlorophenol	88-06-2	1000			
LK3-215002	Custom Pesticide Solution, 2 & 10 mg/L, 15-1				Methanol	1 ml, Amber Glass	€ 79,50
LK3-215001	Custom Pesticide Solution, 20 & 100 mg/L, 15-1					1 ml, Amber Glass	€ 118,50
		hexachlorocyclopentadiene	77-47-4	10			
		trifluralin	1582-09-8	2			
		propachlor	1918-16-7	10			
		hexachlorobenzene	118-74-1	2			
		Chlorothalonil	1897-45-6	2			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		a-chlordane	5103-71-9	2			
		heptachlor	76-44-8	2			
		kepone	143-50-0	10			
		mirex	2385-85-5	2			
		alachlor	15972-60-8	10			
		trans-nonachlor	39765-80-5	2			
		cis-nonachlor	5103-73-1	2			
		o,p'-DDD	53-19-0	2			
		o,p'-DDE	3424-82-6	2			
LK3-108501	Custom Pesticide Solution, 15	Tebuthiuron	34014-18-1	100	Acetonitrile	1 ml, Amber Glass	€ 198,00
		Monuron	150-68-5	100			
		diuron	330-54-1	100			
		monolinuron	1746-81-2	100			
		linuron	330-55-2	100			
		chlorbromuron	13360-45-7	100			
		Siduron	1982-49-6	100			
		fluometuron	2164-17-2	100			
		metobromuron	3060-89-7	100			
		neburon	555-37-3	100			
		metoxuron	19937-59-8	100			
		isoproturon	34123-59-6	100			
		chlortoluron	15545-48-9	100			
		chloroxuron	1982-47-4	100			
		fenuron	101-42-8	100			
LK3-215305	Pesticide Mix, 1 mg/L, 15-13				Nonane	1.2 ml, Amber Glass	€ 118,50
LK3-215304	Pesticide Mix, 10 mg/L, 15-13					1.2 ml, Amber Glass	€ 150,00
		hexachlorobenzene	118-74-1	10			
		a-BHC	319-84-6	1			
		b-BHC	319-85-7	1			
		g-BHC	58-89-9	1			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		aldrin	309-00-2	1			
		dieldrin	60-57-1	1			
		endrin	72-20-8	1			
		g-chlordane	5103-74-2	1			
		oxychlordane	27304-13-8	1			
		trans-nonachlor	39765-80-5	1			
		heptachlor	76-44-8	1			
		heptachlor epoxide isomer B	1024-57-3	1			
		p,p'-DDT	50-29-3	1			
		p,p'-DDE	72-55-9	1			
		p,p'-DDD	72-54-8	1			
LK3-215101	Custom Pesticide Solution, 15	hexachlorocyclopentadiene	77-47-4	100	Methanol	1 ml, Amber Glass	€ 87,00
		trifluralin	1582-09-8	20			
		propachlor	1918-16-7	100			
		hexachlorobenzene	118-74-1	20			
		Chlorothalonil	1897-45-6	20			
		a-chlordane	5103-71-9	20			
		g-chlordane	5103-74-2	20			
		kepone	143-50-0	100			
		mirex	2385-85-5	20			
		alachlor	15972-60-8	100			
		trans-nonachlor	39765-80-5	20			
		cis-nonachlor	5103-73-1	20			
		o,p'-DDD	53-19-0	20			
		o,p'-DDE	3424-82-6	20			
		o,p'-DDT	789-02-6	20			
LK3-101803	Pesticide Mix, 15-19	2-hydroxyethyl methacrylate	868-77-9	100	Methylene Chloric	1 ml, Amber Glass	€ 118,50
		carbofuran	1563-66-2	100			
		diazinon	333-41-5	100			
		Carbaryl	63-25-2	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		phorate	298-02-2	100			
		dimethoate	60-51-5	100			
		Abate®	3383-96-8	100			
		Terbufos	13071-79-9	100			
		malathion	121-75-5	100			
		triallate	2303-17-5	100			
		methyl parathion	298-00-0	100			
		chlorpyrifos	2921-88-2	100			
		parathion	56-38-2	100			
		dichlorvos	62-73-7	100			
		Guthion®	86-50-0	100			
LK3-246502	*Custom Pesticide Mix, 16-1	chlorpyrifos	2921-88-2	50	Hexane	1 ml, Amber Glass	€ 135,00
		dimethoate	60-51-5	50			
		methyl parathion	298-00-0	50			
		phorate	298-02-2	50			
		ethoprophos (prophos)	13194-48-4	50			
		prothiofos	34643-46-4	50			
		pendimethalin	40487-42-1	50			
		sulprofos (Bolstar®)	35400-43-2	50			
		trifluralin	1582-09-8	50			
		diazinon	333-41-5	50			
		tribufos	78-48-8	100			
		Guthion®	86-50-0	100			
		mevinphos (Phosdrin (TM))	7786-34-7	100			
		tributyl phosphate	126-73-8	100			
		triphenyl phosphate	115-86-6	100			
		dibrom	300-76-5	200			
LK3-218001	Custom Pesticide Solution, 10 mg/L, 16-1				Methanol	1 ml, Amber Glass	€ 103,50
LK3-218002	Custom Pesticide Solution, 100 mg/L, 16-1					1 ml, Amber Glass	€ 103,50
LK3-218005	Custom Pesticide Solution, 100 mg/L, 16-1				Acetone	1 ml, Amber Glass	€ 103,50

alachlor	15972-60-8
aldrin	309-00-2
a-BHC	319-84-6
b-BHC	319-85-7
g-BHC	58-89-9
a-chlordane	5103-71-9
g-chlordane	5103-74-2
o,p'-DDE	3424-82-6
p,p'-DDE	72-55-9
o,p'-DDD	53-19-0
p,p'-DDD	72-54-8
o,p'-DDT	789-02-6
p,p'-DDT	50-29-3
dieldrin	60-57-1
endrin	72-20-8
atrazine	1912-24-9

LK3-236503	Custom Pesticide Mix, 16-3			Methanol	1 ml, Amber Glass	€ 198,00
LK3-236504	Custom Pesticide Mix, 16-3			Hexane	1 ml, Amber Glass	€ 198,00
	a-BHC	319-84-6	100			
	b-BHC	319-85-7	100			
	g-BHC	58-89-9	100			
	hexachlorobenzene	118-74-1	100			
	heptachlor	76-44-8	100			
	aldrin	309-00-2	100			
	heptachlor epoxide isomer B	1024-57-3	100			
	o,p'-DDE	3424-82-6	100			
	o,p'-DDD	53-19-0	100			
	o,p'-DDT	789-02-6	100			
	p,p'-DDE	72-55-9	100			
	p,p'-DDD	72-54-8	100			
	p,p'-DDT	50-29-3	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		dieldrin	60-57-1	100			
		endrin	72-20-8	100			
		methoxychlor	72-43-5	100			
LK3-214502	Pesticide Mix, 16-8				Acetone	1 ml, Amber Glass	€ 87,00
LK3-214502-5pck	Pesticide Mix, 16-8					5 x 1 ml, Amber	€ 346,50
		alachlor	15972-60-8	2000			
		atrazine	1912-24-9	2000			
		bladex	21725-46-2	2000			
		bromacil	314-40-9	2000			
		endrin	72-20-8	2000			
		endrin aldehyde	7421-93-4	2000			
		endrin ketone	53494-70-5	2000			
		g-BHC	58-89-9	2000			
		heptachlor	76-44-8	2000			
		heptachlor epoxide isomer B	1024-57-3	2000			
		hexachlorobenzene	118-74-1	2000			
		hexachlorocyclopentadiene	77-47-4	2000			
		methoxychlor	72-43-5	2000			
		metolachlor	51218-45-2	2000			
		metribuzin	21087-64-9	2000			
		trifluralin	1582-09-8	2000			
LK3-246302	Custom Pesticide Solution, 16	aldicarb	116-06-3	2000	Acetone	1 ml, Amber Glass	€ 355,50
		Carbaryl	63-25-2	2000			
		carbofuran	1563-66-2	2000			
		chlorpyrifos	2921-88-2	2000			
		diazinon	333-41-5	2000			
		dichlorvos	62-73-7	2000			
		dimethoate	60-51-5	2000			
		malathion	121-75-5	2000			
		Guthion®	86-50-0	2000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		methyl parathion	298-00-0	2000			
		mevinphos (Phosdrin (TM))	7786-34-7	2000			
		parathion	56-38-2	2000			
		phorate	298-02-2	2000			
		Terbufos	13071-79-9	2000			
		triallate	2303-17-5	2000			
		Abate®	3383-96-8	2000			
LK3-003104-5pck	Pesticide Mix 16-50	acetochlor	34256-82-1	200	Acetone	1 ml, Amber Glass	€ 922,50
		alachlor	15972-60-8	200			
		atrazine	1912-24-9	200			
		bladex	21725-46-2	200			
		desethyl atrazine	6190-65-4	200			
		atrazine desisopropyl	1007-28-9	200			
		dimethenamid	87674-68-8	200			
		Ethalfuralin	55283-68-6	200			
		metolachlor	51218-45-2	200			
		metribuzin	21087-64-9	200			
		pendimethalin	40487-42-1	200			
		prometon	1610-18-0	200			
		propachlor	1918-16-7	200			
		propazine	139-40-2	200			
		simazine	122-34-9	200			
		trifluralin	1582-09-8	200			
LK3-236701	Custom Pesticide Mix, 17-3	p,p'-DDD	72-54-8	100	Hexane	1 ml, Amber Glass	€ 72,00
		p,p'-DDE	72-55-9	100			
		p,p'-DDT	50-29-3	100			
		dieldrin	60-57-1	100			
		endrin	72-20-8	100			
		o,p'-DDD	53-19-0	100			
		a-BHC	319-84-6	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		b-BHC	319-85-7	100			
		g-BHC	58-89-9	100			
		d-BHC	319-86-8	100			
		aldrin	309-00-2	100			
		o,p'-DDT	789-02-6	100			
		isodrin	465-73-6	100			
		endosulfan I	959-98-8	100			
		hexachlorobenzene	118-74-1	100			
		pentachlorobenzene	608-93-5	100			
		o,p'-DDE	3424-82-6	100			
LK3-236203	Custom Pesticide Solution, 17	aldrin	309-00-2	100	Methanol: Toluene	1 ml, Amber Glass	€ 276,00
		a-BHC	319-84-6	100			
		a-chlordane	5103-71-9	100			
		dieldrin	60-57-1	100			
		endosulfan II	33213-65-9	100			
		endosulfan sulfate	1031-07-8	100			
		endrin	72-20-8	100			
		g-BHC	58-89-9	100			
		g-chlordane	5103-74-2	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		hexachlorobenzene	118-74-1	100			
		methoxychlor	72-43-5	100			
		mirex	2385-85-5	100			
		p,p'-DDD	72-54-8	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDT	50-29-3	100			
LK3-021301	Custom Pesticides Mix, 17-35	heptachlor	76-44-8	2000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 276,00
		aldrin	309-00-2	2000			
		dieldrin	60-57-1	2000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		p,p'-DDT	50-29-3	2000			
		p,p'-DDD	72-54-8	2000			
		p,p'-DDE	72-55-9	2000			
		g-BHC	58-89-9	2000			
		d-BHC	319-86-8	2000			
		b-BHC	319-85-7	2000			
		a-BHC	319-84-6	2000			
		a-chlordane	5103-71-9	2000			
		g-chlordane	5103-74-2	2000			
		endosulfan I	959-98-8	2000			
		endosulfan II	33213-65-9	2000			
		endrin	72-20-8	2000			
		heptachlor epoxide isomer B	1024-57-3	2000			
		methoxychlor	72-43-5	2000			
LK3-037903-4pck	Pesticide Mix 17-47	Guthion®	86-50-0	100	Acetonitrile	1 ml, Amber Glass	€ 828,00
		Guthion Ethyl	2642-71-9	100			
		Bromophos ethyl	4824-78-6	100			
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	100			
		coumaphos	56-72-4	100			
		tribufos	78-48-8	100			
		diazinon	333-41-5	100			
		dimethoate	60-51-5	100			
		malathion	121-75-5	100			
		methamidophos	10265-92-6	100			
		monocrotophos	6923-22-4	100			
		parathion	56-38-2	100			
		methyl parathion	298-00-0	100			
		mevinphos (Phosdrin (TM))	7786-34-7	100			
		propetamphos	31218-83-4	100			
		profenofos	41198-08-7	100			
		quinalphos	13593-03-8	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-246001	Custom Pesticide Mix, 18-1	chlorpyrifos	2921-88-2	10	Isooctane	1 ml, Amber Glass	€ 135,00
		chloropyrifos-methyl ester	5598-13-0	10			
		Guthion®	86-50-0	10			
		phorate	298-02-2	10			
		fenitrothion	122-14-5	10			
		dimethoate	60-51-5	10			
		parathion	56-38-2	10			
		methyl parathion	298-00-0	10			
		trichlorfon	52-68-6	10			
		diazinon	333-41-5	10			
		fenthion	55-38-9	10			
		dichlorvos	62-73-7	10			
		malathion	121-75-5	10			
		Fonophos	944-22-9	10			
		tetraethyldithio pyrophosphate	3689-24-5	10			
		Terbufos	13071-79-9	10			
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	10			
		monocrotophos	6923-22-4	10			
LK3-214002	Custom Pesticide Solution, 100 mg/L , 18-1				Methanol	1 ml, Amber Glass	€ 157,50
LK3-214001	Custom Pesticide Solution, 200 mg/L , 18-1					1 ml, Amber Glass	€ 198,00
		alachlor	15972-60-8	100			
		atrazine	1912-24-9	100			
		bromacil	314-40-9	100			
		butachlor	23184-66-9	100			
		diazinon	333-41-5	100			
		dimethoate	60-51-5	100			
		hexachlorobenzene	118-74-1	100			
		hexachlorocyclopentadiene	77-47-4	100			
		metolachlor	51218-45-2	100			
		metribuzin	21087-64-9	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		molinate	2212-67-1	100			
		pentachlorophenol	87-86-5	100			
		prometon	1610-18-0	100			
		prometryn	7287-19-6	100			
		propachlor	1918-16-7	100			
		simazine	122-34-9	100			
		thiobencarb	28249-77-6	100			
		acetochlor	34256-82-1	100			
LK3-246603	*Custom Pesticide Mix, 18-2	EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	500	Acetone	1 ml, Amber Glass	€ 198,00
		butylate	2008-41-5	500			
		propachlor	1918-16-7	500			
		trifluralin	1582-09-8	500			
		prometon	1610-18-0	500			
		simazine	122-34-9	500			
		atrazine	1912-24-9	500			
		propazine	139-40-2	500			
		Terbufos	13071-79-9	500			
		acetochlor	34256-82-1	500			
		metribuzin	21087-64-9	500			
		alachlor	15972-60-8	500			
		Ametryne	834-12-8	500			
		malathion	121-75-5	500			
		metolachlor	51218-45-2	500			
		bladex	21725-46-2	500			
		pendimethalin	40487-42-1	500			
		butachlor	23184-66-9	500			
LK3-239302	*Custom Pesticide Mix, 18-6				Toluene	1 ml, Amber Glass	€ 118,50
LK3-239302-5pck	*Custom Pesticide Mix, 18-6					1 ml, Amber Glass	€ 591,00
		alachlor	15972-60-8	200			
		atrazine	1912-24-9	200			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		desethyl atrazine	6190-65-4	200			
		chlorpyrifos	2921-88-2	200			
		g-BHC	58-89-9	200			
		linuron	330-55-2	200			
		Metalaxyl	57837-19-1	200			
		metolachlor	51218-45-2	200			
		molinate	2212-67-1	200			
		Oxadiazon	19666-30-9	200			
		Oxadixyl	77732-09-3	200			
		pendimethalin	40487-42-1	200			
		prometon	1610-18-0	200			
		propyzamide (pronamide)	23950-58-5	200			
		simazine	122-34-9	200			
		Terbumeton	33693-04-8	200			
		terbuthylazine	5915-41-3	200			
		Terbuthylazine-desethyl	30125-63-4	200			
LK3-045803	Custom Pesticide Mix, 18-8				Acetone	1 ml, Amber Glass	€ 103,50
LK3-045803-5pck	Custom Pesticide Mix, 18-8					5 x 1 ml, Amber	€ 513,00
		fenthion	55-38-9	100			
		Fonophos	944-22-9	100			
		formothion	2540-82-1	100			
		Heptenophos	23560-59-0	100			
		malathion	121-75-5	100			
		methidathion	950-37-8	100			
		paraoxon ethyl	311-45-5	100			
		methyl parathion	298-00-0	100			
		phorate	298-02-2	100			
		phosalone	2310-17-0	100			
		phosmet	732-11-6	100			
		fenchlorphos	299-84-3	100			
		fenitrothion	122-14-5	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		methamidophos	10265-92-6	100			
		mevinphos (Phosdrin (TM))	7786-34-7	100			
		parathion	56-38-2	100			
		Isofenphos	25311-71-1	100			
		phosphamidon	13171-21-6	100			
LK3-052303	Custom Pesticide Mix, 18-9				Acetone	1 ml, Amber Glass	€ 135,00
LK3-052303-5pck	Custom Pesticide Mix, 18-9					5 x 1 ml, Amber Glass	€ 513,00
		chlorpyrifos-methyl	5598-13-0	100			
		Guthion Ethyl	2642-71-9	100			
		Bromophos methyl	2104-96-3	100			
		acephate	30560-19-1	100			
		buprofezin	69327-76-0	100			
		chlormephos	24934-91-6	100			
		carbophenothion	786-19-6	100			
		dichlorvos	62-73-7	100			
		Guthion®	86-50-0	100			
		coumaphos	56-72-4	100			
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	100			
		fenamiphos	22224-92-6	100			
		dimethoate	60-51-5	100			
		ethion	563-12-2	100			
		Bromophos ethyl	4824-78-6	100			
		chlorpyrifos	2921-88-2	100			
		diazinon	333-41-5	100			
		disulfoton	298-04-4	100			
LK3-076105-5pck	Custom Pesticide Solution, 18	chinomethionate	2439-01-2	100	Methanol	1 ml, Amber Glass	€ 780,00
		esfenvalerate	66230-04-4	100			
		flucythrinate (mixture of isomers)	70124-77-5	100			
		tau-fluvalinate	102851-06-9	100			
		oxadiazon	19666-30-9	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Oxyfluorfen	42874-03-3	100			
		penconazol	66246-88-6	100			
		phosalone	2310-17-0	100			
		Pyrazophos	13457-18-6	100			
		pyridaben	96489-71-3	100			
		tralomethrin	66841-25-6	100			
		buprofezin	69327-76-0	100			
		malathion	121-75-5	100			
		phosmet	732-11-6	100			
		Pirimicarb	23103-98-2	100			
		benomyl	17804-35-2	100			
		thiophanate methyl	23564-05-8	100			
		Methiocarb	2032-65-7	100			
LK3-084703-5pck	Custom Pesticide Mix, 18-47	methyl parathion	298-00-0	100	Methanol	1 ml, Amber Glass	€ 780,00
		hexachlorocyclopentadiene	77-47-4	100			
		2,4,5-T	93-76-5	100			
		2,4,5-TP (Silvex)	93-72-1	100			
		2,4-DB	94-82-6	100			
		aldicarb	116-06-3	100			
		atrazine	1912-24-9	100			
		bendiocarb	22781-23-3	100			
		bladex	21725-46-2	100			
		isoproturon	34123-59-6	100			
		chlortoluron	15545-48-9	100			
		chlorpyrifos	2921-88-2	100			
		dalapon	75-99-0	100			
		dinoseb	88-85-7	100			
		methoxychlor	72-43-5	100			
		MCPA acid	94-74-6	100			
		oxamyl	23135-22-0	100			
		picloram	1918-02-1	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-249505	Pesticide Mix, 19-4	dichlorvos	62-73-7	200	Acetone:Hexane 1 2	ml, Amber Glass	€ 276,00
		mevinphos (Phosdrin (TM))	7786-34-7	200			
		Demeton O&S (29.0% O & 62.6% S)	8065-48-3	200			
		ethoprophos (prophos)	13194-48-4	200			
		sulfotep	3689-24-5	200			
		diazinon	333-41-5	200			
		fenchlorphos (Ronnell)	299-84-3	200			
		fenthion	55-38-9	200			
		methyl parathion	298-00-0	200			
		malathion	121-75-5	200			
		prothiofos	34643-46-4	200			
		parathion	56-38-2	200			
		tetrachlorvinphos (Rabon)	961-11-5	200			
		sulprofos (Bolstar®)	35400-43-2	200			
		fensulfothion	115-90-2	200			
		EPN	2104-64-5	200			
		Guthion®	86-50-0	200			
		coumaphos	56-72-4	200			
		tributylphosphoro-trithioite (Merphos)	150-50-5	200			
LK3-249706	Pesticide Mix, 20-14	Guthion®	86-50-0	200	Methanol	1 ml, Amber Glass	€ 261,00
		sulprofos (Bolstar®)	35400-43-2	200			
		chlorpyrifos	2921-88-2	200			
		coumaphos	56-72-4	200			
		Demeton O&S (28.2% O & 65.2% S)	8065-48-3	200			
		diazinon	333-41-5	200			
		dichlorvos	62-73-7	200			
		disulfoton	298-04-4	200			
		ethoprophos (prophos)	13194-48-4	200			
		fensulfothion	115-90-2	200			
		fenthion	55-38-9	200			

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		tributylphosphoro-trithioite (Merphos)	150-50-5	200			
		dibrom	300-76-5	200			
		phorate	298-02-2	200			
		fenchlorphos (Ronne)	299-84-3	200			
		tetrachlorvinphos (Stirophos)	961-11-5	200			
		prothiofos	34643-46-4	200			
		trichloronate	327-98-0	200			
		mevinphos (Phosdrin (TM))	7786-34-7	200			
		methyl parathion	298-00-0	200			
LK3-105901	Pesticide Mix, 20-3	alachlor	15972-60-8	100	Acetone	1 ml, Amber Glass	€ 79,50
		aldrin	309-00-2	100			
		atrazine	1912-24-9	100			
		benzo[a]pyrene	50-32-8	100			
		dioctyl phthalate	117-81-7	100			
		bis(2-ethylhexyl) adipate	103-23-1	100			
		p,p'-DDD	72-54-8	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDT	50-29-3	100			
		a-chlordane	5103-71-9	100			
		g-chlordane	5103-74-2	100			
		endrin	72-20-8	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		hexachlorobenzene	118-74-1	100			
		hexachlorocyclopentadiene	77-47-4	100			
		g-BHC	58-89-9	100			
		methoxychlor	72-43-5	100			
		simazine	122-34-9	100			
		trans-nonachlor	39765-80-5	100			
LK3-214305	Pesticide Mix, 21-9	alachlor	15972-60-8	500	Acetone	1 ml, Amber Glass	€ 237,00

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		aldrin	309-00-2	500			
		atrazine	1912-24-9	500			
		g-BHC	58-89-9	500			
		bromacil	314-40-9	500			
		butachlor	23184-66-9	500			
		diazinon	333-41-5	500			
		hexachlorobenzene	118-74-1	500			
		hexachlorocyclopentadiene	77-47-4	500			
		methoxychlor	72-43-5	500			
		metolachlor	51218-45-2	500			
		molinate	2212-67-1	500			
		prometryn	7287-19-6	500			
		dieldrin	60-57-1	500			
		endrin	72-20-8	500			
		heptachlor	76-44-8	500			
		heptachlor epoxide isomer B	1024-57-3	500			
		propachlor	1918-16-7	500			
		simazine	122-34-9	500			
		benthiocarb	28249-77-6	500			
		trifluralin	1582-09-8	500			
LK3-017001	Custom Pesticide Mix, 21-11	Guthion®	86-50-0	200	Acetone	1 ml, Amber Glass	€ 174,00
		chlorpyrifos	2921-88-2	200			
		diazinon	333-41-5	200			
		dimethoate	60-51-5	200			
		disulfoton	298-04-4	200			
		malathion	121-75-5	200			
		methidathion	950-37-8	200			
		molinate	2212-67-1	200			
		methyl parathion	298-00-0	200			
		phorate	298-02-2	200			
		benthiocarb	28249-77-6	200			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		atrazine	1912-24-9	200			
		tributyl phosphate	126-73-8	200			
		triphenyl phosphate	115-86-6	200			
		dichlorvos	62-73-7	200			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	200			
		parathion	56-38-2	200			
		pendimethalin	40487-42-1	200			
		trifluralin	1582-09-8	200			
		Demeton O&S (29.0% O & 62.6% S)	8065-48-3	200			
		simazine	122-34-9	200			
LK3-316905-4pck	Custom Pesticide Mix, 21-25	2,4,5-T	93-76-5	500	Acetone	1 ml, Amber Glass	€ 937,50
		2,4-D	94-75-7	500			
		2,4-DB	94-82-6	500			
		chlorthiamid	1918-13-4	500			
		diphenylamine	122-39-4	500			
		EPN	2104-64-5	500			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	500			
		MCPA acid	94-74-6	500			
		MCPB	94-81-5	500			
		MCPA-methyl ester	2436-73-9	500			
		metolcarb	1129-41-5	500			
		o,p'-DDD	53-19-0	500			
		o,p'-DDE	3424-82-6	500			
		o,p'-DDT	789-02-6	500			
		p,p'-DDD	72-54-8	500			
		p,p'-DDE	72-55-9	500			
		p,p'-DDT	50-29-3	500			
		a-BHC	319-84-6	500			
		b-BHC	319-85-7	500			
		g-BHC	58-89-9	500			
		d-BHC	319-86-8	500			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-036701	*Custom Pesticide Solution, 22-1				Isooctane	1 ml, Amber Glass	€ 355,50
		methoxychlor	72-43-5	100			
		a-chlordane	5103-71-9	100			
		g-chlordane	5103-74-2	100			
		oxychlordane	27304-13-8	100			
		trans-nonachlor	39765-80-5	100			
		cis-nonachlor	5103-73-1	100			
		o,p'-DDE	3424-82-6	100			
		o,p'-DDD	53-19-0	100			
		o,p'-DDT	789-02-6	100			
		hexachlorobenzene	118-74-1	100			
		heptachlor epoxide	1024-57-3	100			
		mirex	2385-85-5	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDD	72-54-8	100			
		p,p'-DDT	50-29-3	100			
		a-BHC	319-84-6	100			
		b-BHC	319-85-7	100			
		g-BHC	58-89-9	100			
		aldrin	309-00-2	100			
		eldrin	153-18-4	100			
		heptachlor	76-44-8	100			
		dieldrin	60-57-1	100			
LK3-239401	Pesticide Mix, 22-3	acetochlor	34256-82-1	1000	Acetone	1 ml, Amber Glass	€ 276,00
		alachlor	15972-60-8	1000			
		atrazine	1912-24-9	1000			
		desethyl atrazine	6190-65-4	1000			
		Atrazine desisopropyl	1007-28-9	1000			
		dimethenamid	87674-68-8	1000			
		chlorpyrifos	2921-88-2	1000			

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
Cyanazine	21725-46-2	1000			
EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000			
Ethalfuralin	55283-68-6	1000			
Fonophos	944-22-9	1000			
metolachlor	51218-45-2	1000			
metribuzin	21087-64-9	1000			
pendimethalin	40487-42-1	1000			
phorate	298-02-2	1000			
propachlor	1918-16-7	1000			
prometon	1610-18-0	1000			
propazine	139-40-2	1000			
simazine	122-34-9	1000			
Terbufos	13071-79-9	1000			
triallate	2303-17-5	1000			
trifluralin	1582-09-8	1000			

LK3-006503	Pesticide Mix, 22-6	aldrin	309-00-2	100	Isooctane	1 ml, Amber Glass	€ 103,50
		o,p'-DDD	53-19-0	100			
		p,p'-DDD	72-54-8	100			
		o,p'-DDE	3424-82-6	100			
		p,p'-DDE	72-55-9	100			
		o,p'-DDT	789-02-6	100			
		p,p'-DDT	50-29-3	100			
		dieldrin	60-57-1	100			
		endosulfan I	959-98-8	100			
		endosulfan II	33213-65-9	100			
		endosulfan sulfate	1031-07-8	100			
		endrin	72-20-8	100			
		endrin aldehyde	7421-93-4	100			
		a-BHC	319-84-6	100			
		g-BHC	58-89-9	100			
		d-BHC	319-86-8	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		heptachlor epoxide isomer A	28044-83-9	100			
		hexachlorobenzene	118-74-1	100			
		isodrin	465-73-6	100			
		pentachlorobenzene	608-93-5	100			
LK3-037803-4pck	Pesticide Mix 22-37	aldrin	309-00-2	100	Acetonitrile	1 ml, Amber Glass	€ 922,50
		o,p'-DDD	53-19-0	100			
		p,p'-DDD	72-54-8	100			
		o,p'-DDE	3424-82-6	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDT	50-29-3	100			
		o,p'-DDT	789-02-6	100			
		dieldrin	60-57-1	100			
		endosulfan I	959-98-8	100			
		endosulfan II	33213-65-9	100			
		endrin	72-20-8	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		hexachlorobenzene	118-74-1	100			
		a-BHC	319-84-6	100			
		b-BHC	319-85-7	100			
		g-BHC	58-89-9	100			
		d-BHC	319-86-8	100			
		methoxychlor	72-43-5	100			
		mirex	2385-85-5	100			
		trifluralin	1582-09-8	100			
		isodrin	465-73-6	100			
LK3-239601	Custom Pesticide Mix, 23-3	EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	500	Acetone	1 ml, Amber Glass	€ 229,50
		butylate	2008-41-5	500			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		propachlor	1918-16-7	500			
		trifluralin	1582-09-8	500			
		prometon	1610-18-0	500			
		simazine	122-34-9	500			
		atrazine	1912-24-9	500			
		propazine	139-40-2	500			
		Terbufos	13071-79-9	500			
		acetochlor	34256-82-1	500			
		metribuzin	21087-64-9	500			
		alachlor	15972-60-8	500			
		Ametryne	834-12-8	500			
		malathion	121-75-5	500			
		metolachlor	51218-45-2	500			
		bladex	21725-46-2	500			
		pendimethalin	40487-42-1	500			
		butachlor	23184-66-9	500			
		Ethalfuralin	55283-68-6	500			
		phorate	298-02-2	500			
		triallate	2303-17-5	500			
		Fonophos	944-22-9	500			
		bromacil	314-40-9	500			
LK3-043201-5pck	Custom Pesticide Mix, 23-4	endosulfan (alpha + beta)	115-29-7	500	Acetone	1 ml, Amber Glass	€ 1.410,00
		fenobucarb	3766-81-2	500			
		diazinon	333-41-5	500			
		carbofuran	1563-66-2	500			
		Propanil	709-98-8	500			
		alachlor	15972-60-8	500			
		bentazon	25057-89-0	500			
		pendimethalin	40487-42-1	500			
		butachlor	23184-66-9	500			
		metolachlor	51218-45-2	500			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		molinate	2212-67-1	500			
		chlorothalonil	1897-45-6	500			
		iprobenfos	26087-47-8	500			
		Metalaxyl	57837-19-1	500			
		procymidone	32809-16-8	500			
		buprofezin	69327-76-0	500			
		methyl parathion	298-00-0	500			
		isoprothiolane	50512-35-1	500			
		beam™	41814-78-2	500			
		Hexaconazole	79983-71-4	500			
		cadusafos	95465-99-9	500			
		ethoprophos (prophos)	13194-48-4	500			
		Tolclofos-Methyl	57018-04-9	500			
LK3-006403	Custom Pesticide Mix, 23-6	a-BHC	319-84-6	100	Isooctane	1 ml, Amber Glass	€ 118,50
		b-BHC	319-85-7	100			
		aldrin	309-00-2	100			
		d-BHC	319-86-8	100			
		dieldrin	60-57-1	100			
		p,p'-DDT	50-29-3	100			
		g-BHC	58-89-9	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDD	72-54-8	100			
		endrin	72-20-8	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		endosulfan I	959-98-8	100			
		endosulfan II	33213-65-9	100			
		endosulfan sulfate	1031-07-8	100			
		endrin aldehyde	7421-93-4	100			
		o,p'-DDD	53-19-0	100			
		o,p'-DDE	3424-82-6	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		heptachlor epoxide isomer A	28044-83-9	100			
		hexachlorobenzene	118-74-1	100			
		isodrin	465-73-6	100			
		pentachlorobenzene	608-93-5	100			
		o,p'-DDT	789-02-6	100			
LK3-421006-10	Custom Pesticide Mix, 23-16	2,4,6-trichlorophenol	88-06-2	200	Acetone	5 ml, Amber Glass	€ 1.095,00
		pentachlorophenol	87-86-5	200			
		bentazon	25057-89-0	200			
		2,4-dichlorophenoxyacetic acid	94-75-7	200			
		Deltamethrin	52918-63-5	200			
		dimethoate	60-51-5	200			
		parathion	56-38-2	200			
		malathion	121-75-5	200			
		methyl parathion	298-00-0	200			
		dichlorvos	62-73-7	200			
		p,p'-DDD	72-54-8	200			
		p,p'-DDE	72-55-9	200			
		p,p'-DDT	50-29-3	200			
		o,p'-DDT	789-02-6	200			
		a-BHC	319-84-6	200			
		b-BHC	319-85-7	200			
		g-BHC	58-89-9	200			
		d-BHC	319-86-8	200			
		heptachlor	76-44-8	200			
		chlorothalonil	1897-45-6	200			
		chlorpyrifos	2921-88-2	200			
		atrazine	1912-24-9	200			
		bis(2-ethylhexyl)phthalate	117-81-7	200			
LK3-016501	Custom Pesticide Mix, 23-17	Guthion®	86-50-0	200	Acetone	1 ml, Amber Glass	€ 198,00
		chlorpyrifos	2921-88-2	200			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		diazinon	333-41-5	200			
		dimethoate	60-51-5	200			
		disulfoton	298-04-4	200			
		malathion	121-75-5	200			
		methidathion	950-37-8	200			
		molinate	2212-67-1	200			
		methyl parathion	298-00-0	200			
		phorate	298-02-2	200			
		phosmet	732-11-6	200			
		benthiocarb	28249-77-6	200			
		atrazine	1912-24-9	200			
		bladex	21725-46-2	200			
		tributyl phosphate	126-73-8	200			
		triphenyl phosphate	115-86-6	200			
		Demeton O&S (29.0% O & 62.6% S)	8065-48-3	200			
		dichlorvos	62-73-7	200			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	200			
		parathion	56-38-2	200			
		pendimethalin	40487-42-1	200			
		trifluralin	1582-09-8	200			
		simazine	122-34-9	200			
LK3-214704	Custom Pesticide Mix, 23-22				Acetone	1 ml, Amber Glass	€ 150,00
LK3-214704-4pck	Custom Pesticide Mix, 23-22					4 x 1 ml, Amber	€ 622,50
		aldrin	309-00-2	500			
		atrazine	1912-24-9	500			
		dieldrin	60-57-1	500			
		endrin	72-20-8	500			
		heptachlor	76-44-8	500			
		heptachlor epoxide isomer B	1024-57-3	500			
		g-BHC	58-89-9	500			
		methoxychlor	72-43-5	500			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		metolachlor	51218-45-2	500			
		metribuzin	21087-64-9	500			
		simazine	122-34-9	500			
		trifluralin	1582-09-8	500			
		chlorpyrifos	2921-88-2	500			
		diazinon	333-41-5	500			
		dimethoate	60-51-5	500			
		malathion	121-75-5	500			
		parathion	56-38-2	500			
		phorate	298-02-2	500			
		Terbufos	13071-79-9	500			
		carbofuran	1563-66-2	500			
		o,p'-DDT	789-02-6	500			
		p,p'-DDT	50-29-3	500			
		a-BHC	319-84-6	500			
LK3-210303	Pesticide Mix, 24-7	a-BHC	319-84-6	100	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 198,00
		b-BHC	319-85-7	100			
		d-BHC	319-86-8	100			
		g-BHC	58-89-9	100			
		heptachlor	76-44-8	100			
		aldrin	309-00-2	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		dieldrin	60-57-1	100			
		endosulfan I	959-98-8	100			
		endosulfan II	33213-65-9	100			
		endrin	72-20-8	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDD	72-54-8	100			
		endosulfan sulfate	1031-07-8	100			
		p,p'-DDT	50-29-3	100			
		o,p'-DDT	789-02-6	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		hexachlorobenzene	118-74-1	100			
		a-chlordane	5103-71-9	100			
		methoxychlor	72-43-5	100			
		mirex	2385-85-5	100			
		g-chlordane	5103-74-2	100			
		endrin aldehyde	7421-93-4	100			
		o,p'-DDE	3424-82-6	100			
		endrin ketone	53494-70-5	100			
LK3-037704-4pck	Pesticide Mix 24-26	Carbaryl	63-25-2	100	Acetonitrile	1 ml, Amber Glass	€ 937,50
		p,p'-DDE	72-55-9	100			
		o,p'-DDT	789-02-6	100			
		hexachlorobenzene	118-74-1	100			
		mirex	2385-85-5	100			
		trifluralin	1582-09-8	100			
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	100			
		coumaphos	56-72-4	100			
		Baythroid (mixture of isomers)	68359-37-5	100			
		cypermethrin	52315-07-8	100			
		deltamethrin	52918-63-5	100			
		diazinon	333-41-5	100			
		dicrotophos	141-66-2	100			
		dimethoate	60-51-5	100			
		esfenvalerate	66230-04-4	100			
		malathion	121-75-5	100			
		methamidophos	10265-92-6	100			
		monocrotophos	6923-22-4	100			
		parathion	56-38-2	100			
		methyl parathion	298-00-0	100			
		mevinphos (Phosdrin (TM))	7786-34-7	100			
		profenofos	41198-08-7	100			
		quinalphos	13593-03-8	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		isodrin	465-73-6	100			
LK3-248701	Custom Pesticide Mix, 1,000 mg/L , 25-7				Acetone	1 ml, Amber Glass	€ 450,00
LK3-248705	Custom Pesticide Mix, 500 mg/L , 25-7				Toluene	1 ml, Amber Glass	€ 394,50
		isoproturon	34123-59-6				
		diuron	330-54-1				
		pentachlorobenzene	608-93-5				
		trifluralin	1582-09-8				
		hexachlorobenzene	118-74-1				
		a-BHC	319-84-6				
		b-BHC	319-85-7				
		g-BHC	58-89-9				
		d-BHC	319-86-8				
		chlorpyrifos-methyl	5598-13-0				
		chlorpyrifos	2921-88-2				
		alachlor	15972-60-8				
		aldrin	309-00-2				
		isodrin	465-73-6				
		chlorfenvinphos (mixture of z and E isomers)	470-90-6				
		o,p'-DDE	3424-82-6				
		endosulfan I	959-98-8				
		endosulfan II	33213-65-9				
		endosulfan sulfate	1031-07-8				
		dieldrin	60-57-1				
		p,p'-DDT	50-29-3				
		o,p'-DDD	53-19-0				
		endrin	72-20-8				
		simazine	122-34-9				
		atrazine	1912-24-9				
LK3-085506-2pck	Custom Pesticide Mix, 25-14	o,p'-DDD	53-19-0	100	Methanol	1 ml, Amber Glass	€ 237,00
		o,p'-DDE	3424-82-6	100			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		o,p'-DDT	789-02-6	100			
		p,p'-DDD	72-54-8	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDT	50-29-3	100			
		a-BHC	319-84-6	100			
		a-chlordane	5103-71-9	100			
		aldrin	309-00-2	100			
		b-BHC	319-85-7	100			
		chlordene	3734-48-3	100			
		d-BHC	319-86-8	100			
		dieldrin	60-57-1	100			
		endosulfan I	959-98-8	100			
		endosulfan II	33213-65-9	100			
		endosulfan sulfate	1031-07-8	100			
		endrin	72-20-8	100			
		g-chlordane	5103-74-2	100			
		g-BHC	58-89-9	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer A	28044-83-9	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		hexachlorobenzene	118-74-1	100			
		methoxychlor	72-43-5	100			
		Perthane (TM)	72-56-0	100			
LK3-178603-5pck	Custom Pesticide Mix, 26-6	dichlorvos	62-73-7	100	Acetone:Hexane 1 1	ml, Amber Glass	€ 780,00
		propham	122-42-9	100			
		Heptenophos	23560-59-0	100			
		ethoprophos (prophos)	13194-48-4	100			
		phorate	298-02-2	100			
		Fonophos	944-22-9	100			
		propyzamide (pronamide)	23950-58-5	100			
		chlorothalonil	1897-45-6	100			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		formothion	2540-82-1	100			
		chloropyrifos-methyl	5598-13-0	100			
		Metalaxyl	57837-19-1	100			
		aldrin	309-00-2	100			
		chlorpyrifos	2921-88-2	100			
		cyprodinil	121552-61-2	100			
		triadimenol	55219-65-3	100			
		methidathion	950-37-8	100			
		tetrachlorvinphos (Rabon)	961-11-5	100			
		fenamiphos	22224-92-6	100			
		buprofezin	69327-76-0	100			
		oxadixyl	77732-09-3	100			
		quinoxifen	124495-18-7	100			
		trifloxystrobin	141517-21-7	100			
		propargite	2312-35-8	100			
		phosmet	732-11-6	100			
		tebufenpyrad	119168-77-3	100			
		Fenarimol	60168-88-9	100			
LK3-070904	Pesticide Mix, 27-2	mevinphos (Phosdrin (TM))	7786-34-7	500	Acetone	1 ml, Amber Glass	€ 292,50
		ethoprophos (prophos)	13194-48-4	500			
		thionazine (zinophos)	297-97-2	500			
		sulfotep	3689-24-5	500			
		phorate	298-02-2	500			
		atrazine	1912-24-9	500			
		diazinon	333-41-5	500			
		simazine	122-34-9	500			
		disulfoton	298-04-4	500			
		dichlofenthion	97-17-6	500			
		dimethoate	60-51-5	500			
		fenchlorphos (Ronnell)	299-84-3	500			
		methyl parathion	298-00-0	500			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		trichloronate	327-98-0	500			
		chlorpyrifos	2921-88-2	500			
		malathion	121-75-5	500			
		parathion	56-38-2	500			
		fenthion	55-38-9	500			
		prothiofos	34643-46-4	500			
		tetrachlorvinphos (Stirophos)	961-11-5	500			
		ethion	563-12-2	500			
		sulprofos (Bolstar®)	35400-43-2	500			
		fensulfothion	115-90-2	500			
		famphur	52-85-7	500			
		EPN	2104-64-5	500			
		Guthion®	86-50-0	500			
		coumaphos	56-72-4	500			
LK3-663203	Custom Pesticide Mix, 27-5	fenitrothion	122-14-5	100	Ethyl Acetate: Ace 1 ml, Amber Glass		€ 237,00
		malathion	121-75-5	100			
		metolachlor	51218-45-2	100			
		molinate	2212-67-1	100			
		parathion	56-38-2	100			
		tributyl phosphate	126-73-8	100			
		trifluralin	1582-09-8	100			
		alachlor	15972-60-8	100			
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	100			
		chlorpyrifos	2921-88-2	100			
		diazinon	333-41-5	100			
		dimethoate	60-51-5	100			
		methyl parathion	298-00-0	100			
		desethyl atrazine	6190-65-4	100			
		dichlofenthion	97-17-6	100			
		fenchlorphos (Ronnel)	299-84-3	100			
		simazine	122-34-9	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		terbuthylazine	5915-41-3	100			
		terbutryne	886-50-0	100			
		Bromophos ethyl	4824-78-6	100			
		Bromophos methyl	2104-96-3	100			
		ethion	563-12-2	100			
		Omethoate	1113-02-6	100			
		Propanil	709-98-8	100			
		atrazine	1912-24-9	100			
		propazine	139-40-2	100			
		Guthion Ethyl	2642-71-9	100			
LK3-105603	Custom Pesticide Solution, 30-1				Acetone	1 ml, Amber Glass	€ 157,50
LK3-105603-5pck	Custom Pesticide Solution, 30-1					5 x 1 ml, Amber	€ 630,00
		p,p'-DDT	50-29-3	100			
		p,p'-DDD	72-54-8	100			
		p,p'-DDE	72-55-9	100			
		a-BHC	319-84-6	100			
		g-BHC	58-89-9	100			
		dieldrin	60-57-1	100			
		methoxychlor	72-43-5	100			
		endosulfan II	33213-65-9	100			
		endosulfan sulfate	1031-07-8	100			
		endrin	72-20-8	100			
		endrin aldehyde	7421-93-4	100			
		aldrin	309-00-2	100			
		b-BHC	319-85-7	100			
		d-BHC	319-86-8	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		a-chlordane	5103-71-9	100			
		g-chlordane	5103-74-2	100			
		endosulfan I	959-98-8	100			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		trans-nonachlor	39765-80-5	100			
		alachlor	15972-60-8	100			
		atrazine	1912-24-9	100			
		chlorobenzilate	510-15-6	100			
		chloroneb	2675-77-6	100			
		chlorothalonil	1897-45-6	100			
		chlorthal-dimethyl (dacthal)	1861-32-1	100			
		etridiazole	2593-15-9	100			
		simazine	122-34-9	100			
		cis-permethrin	54774-45-7	100			
		trans-permethrin	51877-74-8	100			
LK3-241604	Custom Pesticide Mix, 30-8	2,2',5-trichlorobiphenyl (BZ# 18)	37680-65-2	10	Isooctane	1 ml, Amber Glass	€ 276,00
		2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-5	10			
		2,4',5-trichlorobiphenyl (BZ# 31)	16606-02-3	10			
		2',3,4-trichlorobiphenyl (BZ# 33)	38444-86-9	10			
		2,2',4,5'-tetrachlorobiphenyl (BZ# 49)	41464-40-8	10			
		2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	10			
		2,3',4,4'-tetrachlorobiphenyl (BZ# 66)	32598-10-0	10			
		2,3',4',5-tetrachlorobiphenyl (BZ# 70)	32598-11-1	10			
		2,4,4',5-tetrachlorobiphenyl (BZ# 74)	32690-93-0	10			
		2,2',3,4',6-pentachlorobiphenyl (BZ# 91)	68194-05-8	10			
		2,2',3,5',6-pentachlorobiphenyl (BZ# 95)	38379-99-6	10			
		2,2',4,4',5-pentachlorobiphenyl (BZ# 99)	38380-01-7	10			
		2,2',4,5,5'-pentachlorobiphenyl (BZ# 101)	37680-73-2	10			
		2,3,3',4',6-pentachlorobiphenyl (BZ# 110)	38380-03-9	10			
		2,2',3,3',4,4'-hexachlorobiphenyl (BZ# 128)	38380-07-3	10			
		2,2',3,4,4',5'-hexachlorobiphenyl (BZ# 138)	35065-28-2	10			
		2,2',3,4,5,5'-hexachlorobiphenyl (BZ# 141)	52712-04-6	10			
		2,2',3,4',5,5'-hexachlorobiphenyl (BZ# 146)	51908-16-8	10			
		2,2',3,4',5',6-hexachlorobiphenyl (BZ# 149)	38380-04-0	10			
		2,2',3,5,5',6-hexachlorobiphenyl (BZ# 151)	52663-63-5	10			

Catalog Number Description		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	10			
		2,2',3,3',4,4',5-heptachlorobiphenyl (BZ# 170)	35065-30-6	10			
		2,2',3,3',4,5,6'-heptachlorobiphenyl (BZ# 174)	38411-25-5	10			
		2,2',3,3',4',5,6-heptachlorobiphenyl (BZ# 177)	52663-70-4	10			
		2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	10			
		2,2',3,4,4',5',6-heptachlorobiphenyl (BZ# 183)	52663-69-1	10			
		2,2',3,4',5,5',6-heptachlorobiphenyl (BZ# 187)	52663-68-0	10			
		2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	10			
		2,2',3,3',4,4',5,6'-octachlorobiphenyl (BZ# 196)	42740-50-1	10			
		2,2',3,4,4',5,5',6-octachlorobiphenyl (BZ# 203)	52663-76-0	10			
LK3-119903	Custom Pesticide Mix, 31-9	Guthion®	86-50-0	100	Methanol	1 ml, Amber Glass	€ 213,00
		captan	133-06-2	100			
		carbophenothion	786-19-6	100			
		tribufos	78-48-8	100			
		Demeton O&S (28.2% O & 65.2% S)	8065-48-3	100			
		diazinon	333-41-5	100			
		dimethoate	60-51-5	100			
		disulfoton	298-04-4	100			
		diuron	330-54-1	100			
		esfenvalerate	66230-04-4	100			
		ethion	563-12-2	100			
		malathion	121-75-5	100			
		methidathion	950-37-8	100			
		methyl parathion	298-00-0	100			
		molinate	2212-67-1	100			
		dibrom	300-76-5	100			
		napropamide	15299-99-7	100			
		parathion	56-38-2	100			
		pendimethalin	40487-42-1	100			
		Permethrin (mixture of isomers)	52645-53-1	100			
		phorate	298-02-2	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		phosalone	2310-17-0	100			
		phosmet	732-11-6	100			
		profenofos	41198-08-7	100			
		prometryn	7287-19-6	100			
		propachlor	1918-16-7	100			
		propargite	2312-35-8	100			
		propetamphos	31218-83-4	100			
		benthiocarb	28249-77-6	100			
		mevinphos (Phosdrin (TM))	7786-34-7	100			
		chlorothalonil	1897-45-6	100			
LK3-071201	Custom Pesticide Mix, 32-6	molinate	2212-67-1	1000	Acetone	1 ml, Amber Glass	€ 1.378,50
		desethyl atrazine	6190-65-4	1000			
		terbuthylazine-desethyl	30125-63-4	1000			
		simazine	122-34-9	1000			
		atrazine	1912-24-9	1000			
		propazine	139-40-2	1000			
		g-BHC	58-89-9	1000			
		terbuthylazine	5915-41-3	1000			
		dimethenamid	87674-68-8	1000			
		alachlor	15972-60-8	1000			
		ametryne	834-12-8	1000			
		prometryn	7287-19-6	1000			
		aldrin	309-00-2	1000			
		metolachlor	51218-45-2	1000			
		chlorpyrifos	2921-88-2	1000			
		Thiocarbazyl	36756-79-3	1000			
		isodrin	465-73-6	1000			
		heptachlor epoxide isomer B	1024-57-3	1000			
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	1000			
		dimepiperate	61432-55-1	1000			
		pretilachlor	51218-49-6	1000			

Compound		CAS#	Concentr. mg/L	Matrix	Container	List Price
p,p'-DDE		72-55-9	1000			
Oxadiazon		19666-30-9	1000			
endrin		72-20-8	1000			
endosulfan II		33213-65-9	1000			
p,p'-DDD		72-54-8	1000			
hexazinon		51235-04-2	1000			
chlordane (mixture of isomers)		57-74-9	1000			
a-BHC		319-84-6	1000			
b-BHC		319-85-7	1000			
dieldrin		60-57-1	1000			
p,p'-DDT		50-29-3	1000			
LK3-178403-5pck	Custom Pesticide Mix, 33-9	mevinphos (Phosdrin (TM))	7786-34-7	100	Acetone:Hexane 1 1 ml, Amber Glass	€ 937,50
		chloroneb	2675-77-6	100		
		diphenylamine	122-39-4	100		
		chlorpropham	101-21-3	100		
		dimethoate	60-51-5	100		
		pyrimethanil	53112-28-0	100		
		tefluthrin	79538-32-2	100		
		metribuzin	21087-64-9	100		
		Tolclofos-Methyl	57018-04-9	100		
		fenitrothion	122-14-5	100		
		Dichlofluanid	1085-98-9	100		
		fenthion	55-38-9	100		
		Bromophos methyl	2104-96-3	100		
		Pirimiphos-ethyl	23505-41-1	100		
		captan	133-06-2	100		
		procymidone	32809-16-8	100		
		mepanipirim	110235-47-7	100		
		Hexaconazole	79983-71-4	100		
		fludioxonil	131341-86-1	100		
		Systhane (Myclobutanil)	88671-89-0	100		

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		cyproconazol	94361-06-5	100			
		ethion	563-12-2	100			
		Propiconazol (mixture of isomers)	60207-90-1	100			
		tebuconazol	107534-96-3	100			
		iprodione	36734-19-7	100			
		Bifenthrin	82657-04-3	100			
		fenazaquin	120928-09-8	100			
		tetradifon	116-29-0	100			
		phosalone	2310-17-0	100			
		Pyrazophos	13457-18-6	100			
		ethofenprox	80844-07-1	100			
		Deltamethrin	52918-63-5	100			
		azoxystrobin	131860-33-8	100			
LK3-119203	Custom Pesticide Mix, 34-5	a-BHC	319-84-6	100	Methanol	2 ml, Amber Glass	€ 387,00
		b-BHC	319-85-7	100			
		d-BHC	319-86-8	100			
		g-BHC	58-89-9	100			
		aldrin	309-00-2	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		p,p'-DDD	72-54-8	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDT	50-29-3	100			
		o,p'-DDE	3424-82-6	100			
		alachlor	15972-60-8	100			
		atrazine	1912-24-9	100			
		Benfluralin (Benefin)	1861-40-1	100			
		bromacil	314-40-9	100			
		butachlor	23184-66-9	100			
		chlordan (Mix of Isomers)	57-74-9	100			
		chlorpropham	101-21-3	100			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		chlorpyrifos	2921-88-2	100			
		bladex	21725-46-2	100			
		chlorthal-dimethyl (dacthal)	1861-32-1	100			
		Dichloran	99-30-9	100			
		Kelthane ® (Dicofol)	115-32-2	100			
		dieldrin	60-57-1	100			
		endosulfan I	959-98-8	100			
		endosulfan II	33213-65-9	100			
		endrin	72-20-8	100			
		methoxychlor	72-43-5	100			
		metolachlor	51218-45-2	100			
		norflurazon	27314-13-2	100			
		Oxyfluorfen	42874-03-3	100			
		pentachloronitrobenzene	82-68-8	100			
		simazine	122-34-9	100			
		trifluralin	1582-09-8	100			
LK3-006803	Custom Pesticide Mix, 49-4				Acetone	1 ml, Amber Glass	€ 213,00
LK3-006803-5pck	Custom Pesticide Mix, 49-4					5 x 1 ml, Amber	€ 937,50
		hexachlorobenzene	118-74-1	100			
		o,p'-DDD	53-19-0	100			
		isodrin	465-73-6	100			
		endosulfan sulfate	1031-07-8	100			
		d-BHC	319-86-8	100			
		heptachlor	76-44-8	100			
		endosulfan I	959-98-8	100			
		p,p'-DDT	50-29-3	100			
		p,p'-DDD	72-54-8	100			
		p,p'-DDE	72-55-9	100			
		a-BHC	319-84-6	100			
		g-BHC	58-89-9	100			
		dieldrin	60-57-1	100			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
endosulfan II	33213-65-9	100			
endrin	72-20-8	100			
aldrin	309-00-2	100			
b-BHC	319-85-7	100			
heptachlor epoxide isomer B	1024-57-3	100			
propazine	139-40-2	100			
atrazine	1912-24-9	100			
alachlor	15972-60-8	100			
dimethoate	60-51-5	100			
parathion	56-38-2	100			
trifluralin	1582-09-8	100			
o,p'-DDT	789-02-6	100			
o,p'-DDE	3424-82-6	100			
malathion	121-75-5	100			
chlorpyrifos	2921-88-2	100			
Benfluralin (Benefin)	1861-40-1	100			
Oxadiazon	19666-30-9	100			
terbuthylazine	5915-41-3	100			
metolachlor	51218-45-2	100			
terbutryne	886-50-0	100			
simazine	122-34-9	100			
desethyl atrazine	6190-65-4	100			
Metalaxyl	57837-19-1	100			
Vinclozolin	50471-44-8	100			
chloropyrifos-methyl	5598-13-0	100			
hexazinon	51235-04-2	100			
pendimethalin	40487-42-1	100			
ethion	563-12-2	100			
fenthion	55-38-9	100			
Procymidone	32809-16-8	100			
terbuthylazine-desethyl	30125-63-4	100			
Oxadixyl	77732-09-3	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		fludioxonil	131341-86-1	100			
		iprodione	36734-19-7	100			
		heptachlor epoxide isomer A	28044-83-9	100			
		chlordane (Mix of Isomers)	57-74-9	100			
LK3-005303	Custom Organochlorine Pesticide Calibration Mix A, 25, 50, 100 & 400 mg/L , 12-16				Isooctane	1 ml, Amber Glass	€ 24,00
LK3-005302	Custom Organochlorine Pesticide Calibration Mix A, 50, 100, 200 & 800 mg/L , 12-16					1 ml, Amber Glass	€ 135,00
LK3-005301	Custom Organochlorine Pesticide Calibration Mix A, 100, 200, 400 & 1,600 mg/L , 12-16					1 ml, Amber Glass	€ 103,50
		g-BHC	58-89-9	25			
		aldrin	309-00-2	50			
		dieldrin	60-57-1	50			
		endosulfan I	959-98-8	50			
		heptachlor	76-44-8	50			
		heptachlor epoxide isomer B	1024-57-3	50			
		p,p'-DDT	50-29-3	100			
		endosulfan II	33213-65-9	100			
		endrin aldehyde	7421-93-4	100			
		methoxychlor	72-43-5	400			
		decachlorobiphenyl (BZ 209)	2051-24-3	100			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	100			
LK3-005403	Custom Organochlorine Pesticide Calibration Mix B, 25, 50 & 100 mg/L , 12-17				Isooctane	1 ml, Amber Glass	€ 24,00
LK3-005402	Custom Organochlorine Pesticide Calibration Mix B, 50, 100 & 200 mg/L , 12-17					1 ml, Amber Glass	€ 135,00
LK3-005401	Custom Organochlorine Pesticide Calibration Mix B, 100, 200 & 400 mg/L , 12-17					1 ml, Amber Glass	€ 135,00
		a-BHC	319-84-6	25			
		d-BHC	319-86-8	50			
		a-chlordane	5103-71-9	50			
		g-chlordane	5103-74-2	50			
		p,p'-DDE	72-55-9	50			
		decachlorobiphenyl (BZ 209)	2051-24-3	100			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	100			
		b-BHC	319-85-7	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endrin	72-20-8	100			
		p,p'-DDD	72-54-8	100			
		endosulfan sulfate	1031-07-8	100			
		endrin ketone	53494-70-5	100			
LK3-042101	Organochlorine Pesticide Solut	aldrin	309-00-2	100	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 94,50
		a-BHC	319-84-6	100			
		b-BHC	319-85-7	100			
		dieldrin	60-57-1	100			
		endrin	72-20-8	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		o,p'-DDE	3424-82-6	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDT	50-29-3	100			
		hexachlorobenzene	118-74-1	100			
		g-BHC	58-89-9	100			
		methoxychlor	72-43-5	100			
		mirex	2385-85-5	100			
		p,p'-DDD	72-54-8	100			
		o,p'-DDT	789-02-6	100			
LK3-001525	Organochlorine Pesticide Solut	p,p'-DDT	50-29-3	2000	Acetone	1 ml, Amber Glass	€ 103,50
		p,p'-DDD	72-54-8	2000			
		p,p'-DDE	72-55-9	2000			
		a-BHC	319-84-6	2000			
		g-BHC	58-89-9	2000			
		dieldrin	60-57-1	2000			
		methoxychlor	72-43-5	2000			
		endosulfan II	33213-65-9	2000			
		endosulfan sulfate	1031-07-8	2000			
		endrin	72-20-8	2000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endrin aldehyde	7421-93-4	2000			
		aldrin	309-00-2	2000			
		b-BHC	319-85-7	2000			
		d-BHC	319-86-8	2000			
		heptachlor	76-44-8	2000			
		heptachlor epoxide isomer B	1024-57-3	2000			
		endosulfan I	959-98-8	2000			
LK3-105803	Organochlorine Pesticide Mix,	19-2			Hexane	1 ml, Amber Glass	€ 72,00
LK3-10580305	Organochlorine Pesticide Mix,	19-2				5 ml, Amber Glass	€ 229,50
		a-chlordane	5103-71-9	1			
		g-chlordane	5103-74-2	1			
		chlorpyrifos	2921-88-2	1			
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	1			
		etridiazole	2593-15-9	1			
		hexachlorobenzene	118-74-1	1			
		trifluralin	1582-09-8	1			
		o,p'-DDE	3424-82-6	1			
		o,p'-DDD	53-19-0	1			
		o,p'-DDT	789-02-6	1			
		chlordene	3734-48-3	1			
		cis-nonachlor	5103-73-1	1			
		oxychlordane	27304-13-8	1			
		chlorothalonil	1897-45-6	1,5			
		Permethrin (mixture of isomers)	52645-53-1	2			
		chloroneb	2675-77-6	3			
		chlorobenzilate	510-15-6	3			
		propachlor	1918-16-7	3			
		4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	20			
LK3-106705	Custom Organochlorine Pesticide Mix,	a-chlordane	5103-71-9	100	Hexane	1 ml, Amber Glass	€ 237,00
		g-chlordane	5103-74-2	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		chlorpyrifos	2921-88-2	100			
		etridiazole	2593-15-9	100			
		hexachlorobenzene	118-74-1	100			
		trifluralin	1582-09-8	300			
		o,p'-DDD	53-19-0	100			
		o,p'-DDT	789-02-6	100			
		chlordene	3734-48-3	100			
		cis-nonachlor	5103-73-1	100			
		oxychlordane	27304-13-8	100			
		chlorothalonil	1897-45-6	100			
		Permethrin (mixture of isomers)	52645-53-1	200			
		chloroneb	2675-77-6	300			
		chlorobenzilate	510-15-6	300			
		propachlor	1918-16-7	300			
		4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	2000			
		endosulfan sulfate	1031-07-8	100			
		chlorthal (dacthal)	1861-32-1	100			
LK3-106301	Custom Organochlorine Pesticide Mix, 19-2				Hexane	1 ml, Amber Glass	€ 237,00
LK3-10630105	Custom Organochlorine Pesticide Mix, 19-2					5 ml, Amber Glass	€ 433,50
		a-chlordane	5103-71-9	100			
		g-chlordane	5103-74-2	100			
		chlorpyrifos	2921-88-2	100			
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	100			
		etridiazole	2593-15-9	100			
		hexachlorobenzene	118-74-1	100			
		trifluralin	1582-09-8	300			
		o,p'-DDD	53-19-0	100			
		o,p'-DDT	789-02-6	100			
		chlordene	3734-48-3	100			
		cis-nonachlor	5103-73-1	100			
		oxychlordane	27304-13-8	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		chlorothalonil	1897-45-6	100			
		Permethrin (mixture of isomers)	52645-53-1	200			
		chloroneb	2675-77-6	300			
		chlorobenzilate	510-15-6	300			
		propachlor	1918-16-7	300			
		4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	2000			
		endosulfan sulfate	1031-07-8	100			
LK3-010202	Organochlorine Pesticide Solution	hexachlorobenzene	118-74-1	8	Tol.:Hex. 1:1	1 ml Ampule, Amber Glass	€ 24,00
		o,p'-DDD	53-19-0	8			
		o,p'-DDE	3424-82-6	8			
		o,p'-DDT	789-02-6	8			
		mirex	2385-85-5	8			
		heptachlor	76-44-8	8			
		heptachlor epoxide	1024-57-3	8			
		aldrin	309-00-2	8			
		a-chlordane	5103-71-9	8			
		trans-nonachlor	39765-80-5	8			
		dieldrin	60-57-1	8			
		p,p'-DDE	72-55-9	8			
		p,p'-DDD	72-54-8	8			
		p,p'-DDT	50-29-3	8			
		endosulfan I	959-98-8	8			
		endosulfan II	33213-65-9	8			
		endosulfan sulfate	1031-07-8	8			
		endosulfan ether	3369-52-6	8			
		endosulfan lactone	3868-61-9	8			
		chlorpyrifos	2921-88-2	8			
LK3-011506	Organochlorine Pesticide Calibration Solution, 10 µg/L, 20				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 55,50
LK3-011503	Organochlorine Pesticide Calibration Solution, 100 µg/L, 20					1 ml, Amber Glass	€ 103,50
LK3-011501	Organochlorine Pesticide Calibration Solution, 200 µg/L, 20					1 ml, Amber Glass	€ 55,50

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-011550	Organochlorine Pesticide Calibration Solution 22, 10 mg/L	a-BHC	319-84-6			1 ml, BR, Amber	€ 55,50
		b-BHC	319-85-7				
		d-BHC	319-86-8				
		g-BHC	58-89-9				
		aldrin	309-00-2				
		heptachlor	76-44-8				
		heptachlor epoxide	1024-57-3				
		a-chlordane	5103-71-9				
		g-chlordane	5103-74-2				
		p,p'-DDD	72-54-8				
		p,p'-DDE	72-55-9				
		p,p'-DDT	50-29-3				
		dieldrin	60-57-1				
		endrin	72-20-8				
		endrin aldehyde	7421-93-4				
		endrin ketone	53494-70-5				
		endosulfan I	959-98-8				
		endosulfan II	33213-65-9				
		endosulfan sulfate	1031-07-8				
		methoxychlor	72-43-5				
		decachlorobiphenyl	2051-24-3				
		2,4,5,6-tetrachloro-m-xylene	877-09-8				
LK3-004901	Organochlorine Pesticide Soil Spike Solution, 100 mg/L , 20 + Surrogates				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 103,50
LK3-004902	Organochlorine Pesticide Soil Spike Solution, 1,000 mg/L , 20 + Surrogates					1 ml, Amber Glass	€ 118,50
		a-BHC	319-84-6				
		b-BHC	319-85-7				
		d-BHC	319-86-8				
		g-BHC	58-89-9				
		aldrin	309-00-2				
		heptachlor	76-44-8				

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		heptachlor epoxide isomer B	1024-57-3				
		a-chlordane	5103-71-9				
		g-chlordane	5103-74-2				
		p,p'-DDD	72-54-8				
		p,p'-DDE	72-55-9				
		p,p'-DDT	50-29-3				
		dieldrin	60-57-1				
		endrin	72-20-8				
		endrin aldehyde	7421-93-4				
		endrin ketone	53494-70-5				
		endosulfan I	959-98-8				
		endosulfan II	33213-65-9				
		endosulfan sulfate	1031-07-8				
		2,4,5,6-tetrachloro-m-xylene	877-09-8				
		decachlorobiphenyl	2051-24-3				
		dibutyl chlorendate	1770-80-5				
		methoxychlor	72-43-5				
LK3-010203	Organochlorine Pesticide Solution, 100 mg/L				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 103,50
LK3-010204	Organochlorine Pesticide Solution, 500 µg/L					1 ml, Amber Glass	€ 55,50
		hexachlorobenzene	118-74-1				
		o,p'-DDD	53-19-0				
		mirex	2385-85-5				
		o,p'-DDE	3424-82-6				
		o,p'-DDT	789-02-6				
		heptachlor	76-44-8				
		heptachlor epoxide	1024-57-3				
		aldrin	309-00-2				
		a-chlordane	5103-71-9				
		trans-nonachlor	39765-80-5				
		dieldrin	60-57-1				
		p,p'-DDE	72-55-9				

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
p,p'-DDD	72-54-8				
p,p'-DDT	50-29-3				
endosulfan I	959-98-8				
endosulfan II	33213-65-9				
endosulfan sulfate	1031-07-8				
endosulfan ether	3369-52-6				
endosulfan lacton	3868-61-9				
chlorpyrifos	2921-88-2				
g-BHC	58-89-9				

LK3-056001	Organochlorine Pesticide Solut	a-BHC	319-84-6	8	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 48,00
		b-BHC	319-85-7	8			
		d-BHC	319-86-8	8			
		g-BHC	58-89-9	8			
		aldrin	309-00-2	8			
		heptachlor	76-44-8	8			
		heptachlor epoxide isomer B	1024-57-3	8			
		a-chlordane	5103-71-9	8			
		g-chlordane	5103-74-2	8			
		p,p'-DDD	72-54-8	16			
		p,p'-DDE	72-55-9	16			
		p,p'-DDT	50-29-3	16			
		dieldrin	60-57-1	16			
		endrin	72-20-8	16			
		endrin aldehyde	7421-93-4	16			
		endrin ketone	53494-70-5	16			
		endosulfan I	959-98-8	8			
		endosulfan II	33213-65-9	16			
		endosulfan sulfate	1031-07-8	16			
		methoxychlor	72-43-5	80			
		tris(2,3-dibromopropyl)phosphate	126-72-7	16			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-006005	Custom Organochlorine Pesticide Mix	a-BHC	319-84-6	10	Isooctane	5 ml, Amber Glass	€ 450,00
		aldrin	309-00-2	10			
		endosulfan II	33213-65-9	10			
		b-BHC	319-85-7	10			
		a-chlordane	5103-71-9	10			
		dieldrin	60-57-1	10			
		endosulfan sulfate	1031-07-8	10			
		endrin	72-20-8	10			
		heptachlor epoxide isomer B	1024-57-3	10			
		heptachlor	76-44-8	10			
		g-BHC	58-89-9	10			
		o,p'-DDD	53-19-0	10			
		o,p'-DDE	3424-82-6	10			
		o,p'-DDT	789-02-6	10			
		oxychlordane	27304-13-8	10			
		p,p'-DDD	72-54-8	10			
		p,p'-DDT	50-29-3	10			
		trans-nonachlor	39765-80-5	10			
		g-chlordane	5103-74-2	10			
		p,p'-DDE	72-55-9	10			
		hexachlorobenzene	118-74-1	10			
LK3-105703	Organochlorine Pesticide Mix, 22-2				Hexane	1 ml, Amber Glass	€ 72,00
LK3-10570305	Organochlorine Pesticide Mix, 22-2					5 ml, Amber Glass	€ 229,50
		aldrin	309-00-2	1			
		a-BHC	319-84-6	1			
		b-BHC	319-85-7	1			
		d-BHC	319-86-8	1			
		g-BHC	58-89-9	1			
		p,p'-DDE	72-55-9	1			
		p,p'-DDD	72-54-8	1			
		p,p'-DDT	50-29-3	1			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		dieldrin	60-57-1	1			
		endosulfan I	959-98-8	1			
		endosulfan II	33213-65-9	1			
		endosulfan sulfate	1031-07-8	1			
		endrin	72-20-8	1			
		endrin aldehyde	7421-93-4	1			
		heptachlor	76-44-8	1			
		heptachlor epoxide isomer B	1024-57-3	1			
		methoxychlor	72-43-5	1			
		decachlorobiphenyl (BZ 209)	2051-24-3	1			
		trans-nonachlor	39765-80-5	1			
		endrin ketone	53494-70-5	1			
		mirex	2385-85-5	1			
		hexachlorocyclopentadiene	77-47-4	4			
LK3-106601	Custom Organochlorine Pesticide Mix, 22-3				Hexane	1 ml, Amber Glass	€ 118,50
LK3-106601-05	Custom Organochlorine Pesticide Mix, 22-3					5 ml, Amber Glass	€ 394,50
		aldrin	309-00-2	100			
		a-BHC	319-84-6	100			
		b-BHC	319-85-7	100			
		d-BHC	319-86-8	100			
		g-BHC	58-89-9	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDD	72-54-8	100			
		p,p'-DDT	50-29-3	100			
		dieldrin	60-57-1	100			
		endosulfan I	959-98-8	100			
		endosulfan II	33213-65-9	100			
		endrin	72-20-8	100			
		endrin aldehyde	7421-93-4	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		methoxychlor	72-43-5	100			
		decachlorobiphenyl (BZ 209)	2051-24-3	100			
		trans-nonachlor	39765-80-5	100			
		endrin ketone	53494-70-5	100			
		mirex	2385-85-5	100			
		hexachlorocyclopentadiene	77-47-4	500			
		o,p'-DDE	3424-82-6	100			
LK3-713505	Organochlorine Pesticide Mix,	hexachlorobenzene	118-74-1	10	Isooctane	5 ml, Amber Glass	€ 670,50
		a-BHC	319-84-6	10			
		endosulfan II	33213-65-9	10			
		endosulfan sulfate	1031-07-8	10			
		b-BHC	319-85-7	10			
		trans-nonachlor	39765-80-5	10			
		oxychlordane	27304-13-8	10			
		o,p'-DDE	3424-82-6	10			
		p,p'-DDT	50-29-3	10			
		p,p'-DDD	72-54-8	10			
		p,p'-DDE	72-55-9	10			
		o,p'-DDD	53-19-0	10			
		o,p'-DDT	789-02-6	10			
		endrin	72-20-8	10			
		aldrin	309-00-2	10			
		dieldrin	60-57-1	10			
		heptachlor	76-44-8	10			
		heptachlor epoxide isomer B	1024-57-3	10			
		g-BHC	58-89-9	10			
		a-chlordane	5103-71-9	10			
		g-chlordane	5103-74-2	10			
		hexachlorobutadiene	87-68-3	10			
		isodrin	465-73-6	10			
		endosulfan (alpha + beta)	115-29-7	10			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		HCH (BHC)	608-73-1	10			
		pentachlorobenzene	608-93-5	10			
		endosulfan I	959-98-8	10			
LK3-001202	Organochlorine Pesticide Mixt	p,p'-DDT	50-29-3	100	Acetone	1 ml, Amber Glass	€ 162,00
		p,p'-DDD	72-54-8	100			
		p,p'-DDE	72-55-9	100			
		a-BHC	319-84-6	100			
		g-BHC	58-89-9	100			
		dieldrin	60-57-1	100			
		methoxychlor	72-43-5	100			
		endosulfan II	33213-65-9	100			
		endosulfan sulfate	1031-07-8	100			
		endrin	72-20-8	100			
		endrin aldehyde	7421-93-4	100			
		aldrin	309-00-2	100			
		b-BHC	319-85-7	100			
		d-BHC	319-86-8	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		a-chlordane	5103-71-9	100			
		g-chlordane	5103-74-2	100			
		endosulfan I	959-98-8	100			
		trans-nonachlor	39765-80-5	100			
		etridiazole	2593-15-9	100			
		simazine	122-34-9	100			
		Permethrin	52645-53-1	100			
		chlorthal-dimethyl (dacthal)	1861-32-1	100			
		alachlor	15972-60-8	100			
		atrazine	1912-24-9	100			
		chlorobenzilate	510-15-6	100			
		chlorothalonil	1897-45-6	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		chloroneb	2675-77-6	100			
LK3-001701	Organochlorine Pesticide Solution in Acetone, 30-6				Acetone	1 ml, Amber Glass	€ 55,50
LK3-001701-2pck	Organochlorine Pesticide Solution in Acetone, 30-6					2 x 25 ml, BR, Amber	€ 237,00
		a-BHC	319-84-6	200			
		b-BHC	319-85-7	200			
		d-BHC	319-86-8	200			
		g-BHC	58-89-9	200			
		aldrin	309-00-2	200			
		heptachlor	76-44-8	200			
		heptachlor epoxide isomer B	1024-57-3	200			
		a-chlordane	5103-71-9	200			
		g-chlordane	5103-74-2	200			
		p,p'-DDD	72-54-8	200			
		p,p'-DDE	72-55-9	200			
		p,p'-DDT	50-29-3	200			
		dieldrin	60-57-1	200			
		endrin	72-20-8	200			
		endrin aldehyde	7421-93-4	200			
		endrin ketone	53494-70-5	200			
		endosulfan I	959-98-8	200			
		endosulfan II	33213-65-9	200			
		endosulfan sulfate	1031-07-8	200			
		methoxychlor	72-43-5	200			
		etridiazole	2593-15-9	200			
		propachlor	1918-16-7	200			
		trifluralin	1582-09-8	200			
		chlorthal-dimethyl (dacthal)	1861-32-1	200			
		chlorothalonil	1897-45-6	200			
		chlorpyrifos	2921-88-2	200			
		trans-nonachlor	39765-80-5	200			
		mirex	2385-85-5	200			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Permethrin (mixture of isomers)	52645-53-1	400			
		hexachlorocyclopentadiene	77-47-4	500			
LK3-038105	Organochlorine Pesticide Mix,	hexachlorobenzene	118-74-1	10	Isooctane	5 ml, Amber Glass	€ 661,50
		a-BHC	319-84-6	10			
		endosulfan II	33213-65-9	10			
		endosulfan sulfate	1031-07-8	10			
		b-BHC	319-85-7	10			
		trans-nonachlor	39765-80-5	10			
		oxychlordane	27304-13-8	10			
		o,p'-DDE	3424-82-6	10			
		p,p'-DDT	50-29-3	10			
		p,p'-DDD	72-54-8	10			
		p,p'-DDE	72-55-9	10			
		o,p'-DDD	53-19-0	10			
		o,p'-DDT	789-02-6	10			
		endrin	72-20-8	10			
		aldrin	309-00-2	10			
		dieldrin	60-57-1	10			
		heptachlor	76-44-8	10			
		heptachlor epoxide isomer B	1024-57-3	10			
		g-BHC	58-89-9	10			
		a-chlordane	5103-71-9	10			
		g-chlordane	5103-74-2	10			
		hexachlorobutadiene	87-68-3	10			
		isodrin	465-73-6	10			
		endosulfan (alpha + beta)	115-29-7	10			
		pentachlorobenzene	608-93-5	10			
		endosulfan I	959-98-8	10			
		dichlorvos	62-73-7	10			
		pentachlorophenol	87-86-5	10			
		pentachloroanisole	1825-21-4	10			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	10			
		chlorpyrifos	2921-88-2	10			
		cypermethrin	52315-07-8	10			
		d-BHC	319-86-8	10			
LK3-106206	Custom Organochlorine Pesticide Mix, 41-6				Hexane	1 ml, Amber Glass	€ 135,00
LK3-10620605	Custom Organochlorine Pesticide Mix, 41-6					5 ml, Amber Glass	€ 418,50
		aldrin	309-00-2	1			
		a-BHC	319-84-6	1			
		b-BHC	319-85-7	1			
		d-BHC	319-86-8	1			
		g-BHC	58-89-9	1			
		p,p'-DDE	72-55-9	1			
		p,p'-DDD	72-54-8	1			
		p,p'-DDT	50-29-3	1			
		dieldrin	60-57-1	1			
		endosulfan I	959-98-8	1			
		endosulfan II	33213-65-9	1			
		endosulfan sulfate	1031-07-8	1			
		endrin	72-20-8	1			
		endrin aldehyde	7421-93-4	1			
		heptachlor	76-44-8	1			
		heptachlor epoxide isomer B	1024-57-3	1			
		methoxychlor	72-43-5	1			
		decachlorobiphenyl (BZ 209)	2051-24-3	1			
		trans-nonachlor	39765-80-5	1			
		endrin ketone	53494-70-5	1			
		mirex	2385-85-5	1			
		hexachlorocyclopentadiene	77-47-4	4			
		a-chlordane	5103-71-9	1			
		g-chlordane	5103-74-2	1			
		chlorpyrifos	2921-88-2	1			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	1			
		etridiazole	2593-15-9	1			
		hexachlorobenzene	118-74-1	1			
		trifluralin	1582-09-8	1			
		o,p'-DDE	3424-82-6	1			
		o,p'-DDD	53-19-0	1			
		o,p'-DDT	789-02-6	1			
		chlordene	3734-48-3	1			
		cis-nonachlor	5103-73-1	1			
		oxychlordane	27304-13-8	1			
		chlorothalonil	1897-45-6	1,5			
		Permethrin (mixture of isomers)	52645-53-1	2			
		chloroneb	2675-77-6	3			
		chlorobenzilate	510-15-6	3			
		propachlor	1918-16-7	3			
		4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	20			
LK3-207301	Custom Herbicide Soution, 2-	dalapon	75-99-0	4000	M-t-BE	0.5 oz, BR, Amber	€ 118,50
		2,4-DB	94-82-6	4000			
LK3-207803	8151 Herbicide Mix, 2,000 mg/l , 2-1				Hexane	1 ml, Amber Glass	€ 55,50
LK3-207802	8151 Herbicide Mix, 20,000 mg/l , 2-1				Acetone:Hexane 1	1 ml, Amber Glass	€ 118,50
LK3-207806	8151 Herbicide Mix, 200,000 mg/l , 2-1				Methanol	1 ml, Amber Glass	€ 276,00
		MCPA acid	94-74-6				
		MCPP acid	7085-19-0				
LK3-021101	Chlorinate Herbicide Mix (3 Co	dalapon	79-99-0	1000	Methanol	1 ml, Amber Glass	€ 48,00
		MCPP acid	7085-19-0	1000			
		MCPA acid	94-74-6	1000			
LK3-231001	Custom Herbicide Mix, 3-1	Fenarimol	60168-88-9	500	Hexane	1 ml, Amber Glass	€ 87,00
		Systhane (TM)	88671-89-0	500			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Propiconazol	60207-90-1	500			
LK3-255003	Chlorinated Herbicides Mix, 3-1, 100 mg/L				Methanol	1 ml, Amber Glass	€ 25,50
LK3-255001	Chlorinated Herbicides Mix, 3-1, 1,000 mg/L					1 ml, Amber Glass	€ 40,50
		2,4-dichlorophenoxyacetic acid	94-75-7				
		2,4,5-TP (Silvex)	93-72-1				
		2,4,5-trichlorophenoxyacetic acid	93-76-5				
LK3-255103	Methylated Chlorinated Herbic 2,4-D-methyl ester		1928-38-7	100	Methanol	1 ml, Amber Glass	€ 31,50
		silvex methyl ester	4841-20-7	100			
		2,4,5-T methyl ester	1928-37-6	100			
LK3-232001	Custom Herbicide Mix, 3-2	2,4-dichlorophenoxyacetic acid	94-75-7	500	Hexane	1 ml, Amber Glass	€ 87,00
		Dicamba	1918-00-9	500			
		MCPP acid	7085-19-0	500			
LK3-233001	Custom Herbicide Solution, 5- molinate		2212-67-1	1000	M-t-BE	1 ml, Amber Glass	€ 55,50
		simazine	122-34-9	1000			
		thiobencarb	28249-77-6	1000			
		atrazine	1912-24-9	1000			
		alachlor	15972-60-8	1000			
LK3-235001	Custom Herbicide Solution, 6- bentazon		25057-89-0	200	Methanol	1 ml, Amber Glass	€ 72,00
		dinoseb	88-85-7	200			
		picloram	1918-02-1	200			
		2,4,5-TP (Silvex)	93-72-1	200			
		2,4-dichlorophenoxyacetic acid	94-75-7	200			
		dalapon	75-99-0	2000			
LK3-234001	Custom Herbicide Derivatives 5 bentazon methyl ester		Not Available	200	Hexane	1 ml, Amber Glass	€ 72,00
		dinoseb-methyl ether	6099-79-2	200			
		picloram methyl ester	14143-55-6	200			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		silvex methyl ester	4841-20-7	200			
		2,4-D-methyl ester	1928-38-7	200			
		dalapon-methyl ester	17640-02-7	2000			
LK3-234102	Custom Herbicide Derivatives S	2,4-D-methyl ester	1928-38-7	200	Hexane	1 ml, Amber Glass	€ 63,00
		2,4-DB methyl ester	18625-12-2	200			
		2,4,5-T methyl ester	1928-37-6	200			
		silvex methyl ester	4841-20-7	200			
		dicamba methyl ester	6597-78-0	200			
		dichlorprop-methyl ester		200			
		dinoseb-methyl ether	6099-79-2	200			
LK3-008501	Herbicide Mix, 7-3	2,4-D	94-75-7	200	Methanol	1 ml, Amber Glass	€ 55,50
		2,4-DB	94-82-6	200			
		2,4,5-trichlorophenoxyacetic acid	93-76-5	200			
		2,4,5-TP (Silvex)	93-72-1	200			
		Dicamba	1918-00-9	200			
		dichlorprop (2,4-DP)	120-36-5	200			
		dinoseb	88-85-7	200			
LK3-008601	Custom Herbicide Mix, 7-9				Acetone	1 ml, Amber Glass	€ 63,00
LK3-008601-4pck	Custom Herbicide Mix, 7-9					1 ml, Amber Glass	€ 213,00
		2,4-D	94-75-7	500			
		dalapon	75-99-0	500			
		Dicamba	1918-00-9	500			
		dinoseb	88-85-7	500			
		pentachlorophenol	87-86-5	500			
		picloram	1918-02-1	500			
		2,4,5-TP (Silvex)	93-72-1	500			
LK3-008307	Herbicide Acid Free Mix, 7-15	pentachlorophenol	87-86-5	100	Methanol	1 ml, Amber Glass	€ 87,00
		picloram	1918-02-1	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Dicamba	1918-00-9	100			
		2,4,5-TP (Silvex)	93-72-1	100			
		2,4,5-T	93-76-5	100			
		2,4-D	94-75-7	100			
		dalapon	75-99-0	500			
LK3-207101	Custom Herbicide Mix, Underivatized Acids, EPA Method 515, 8-1				M-t-BE	1 ml, Amber Glass	€ 118,50
		3,5-dichlorobenzoic acid	51-36-5	1000			
		Dicamba	1918-00-9	1000			
		2-methyl-4,6-dinitrophenol	534-52-1	2000			
		pentachlorophenol	87-86-5	1000			
		2,4,5-TP (Silvex)	93-72-1	1000			
		2,4,5-T	93-76-5	1000			
		picloram	1918-02-1	1000			
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	1000			
LK3-234703	Herbicide Mix, Derivatized, EPA Method 515, 8-1				Hexane:Toluene: 1:1	1 ml, Amber Glass	€ 55,50
		2,4,5-T methyl ester	1928-37-6	100			
		dicamba methyl ester	6597-78-0	100			
		2,4-D-methyl ester	1928-38-7	100			
		silvex methyl ester	4841-20-7	100			
		dinoseb-methyl ether	6099-79-2	100			
		pentachloroanisole	1825-21-4	100			
		picloram-methyl ester	14143-55-6	100			
		dalapon methyl ester	17640-02-7	100			
LK3-008201	Herbicide Solution (Mix 1), 8-1	2,4-D	94-75-7	200	Methanol	1 ml, Amber Glass	€ 72,00
		dichlorprop (2,4-DP)	120-36-5	200			
		2,4-DB	94-82-6	200			
		2,4,5-trichlorophenoxyacetic acid	93-76-5	200			
		2,4,5-TP (Silvex)	93-72-1	200			
		Dicamba	1918-00-9	200			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		dinoseb	88-85-7	200			
		pentachlorophenol	87-86-5	200			
LK3-207201	Custom Herbicide Mix, Underivatized Acids, EPA Method 515, 8-2				M-t-BE	1 ml, Amber Glass	€ 118,50
		4-nitrophenol	100-02-7	2000			
		2, 4 Dichlorophenylacetic acid	19719-28-9	4000			
		dichlorprop (2,4-DP)	120-36-5	4000			
		2,4-dichlorophenoxyacetic acid	94-75-7	4000			
		chloramben	133-90-4	4000			
		dinoseb	88-85-7	2000			
		bentazon	25057-89-0	4000			
		acifluorfen	50594-66-6	2000			
LK3-205001	Custom Acid Herbicide Mix 1, Method 515.1, 8-2				Acetone	1 ml, Amber Glass	€ 55,50
		Dicamba	1918-00-9	250			
		2,4-dichlorophenoxyacetic acid	94-75-7	250			
		pentachlorophenol	87-86-5	125			
		2,4,5-TP (Silvex)	93-72-1	250			
		2,4,5-T	93-76-5	250			
		dinoseb	88-85-7	250			
		bentazon	25057-89-0	250			
		Chlorthal (Dacthal)	1861-32-1	250			
LK3-206001	Custom Acid Herbicide Mix 2, Method 515.1, 8-3				Acetone	1 ml, Amber Glass	€ 55,50
		dichlorprop (2,4-DP)	120-36-5				
		chloramben	133-90-4				
		2,4-DB	94-82-6				
		picloram	1918-02-1				
		acifluorfen	50594-66-6				
		Triclopyr	55335-06-3				
		brominal	1689-84-5				
		dalapon	75-99-0				

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-205103	Custom Acid Herbicide Mix, Method 515.1, 8-12				M-t-BE:Acetone 9	1 ml, Amber Glass	€ 72,00
LK3-205105	Custom Acid Herbicide Mix, Method 515.1, 8-13				Methanol	1 ml, Amber Glass	€ 72,00
		Dicamba	1918-00-9	100			
		2,4-dichlorophenoxyacetic acid	94-75-7	100			
		pentachlorophenol	87-86-5	100			
		2,4,5-TP (Silvex)	93-72-1	100			
		2,4,5-T	93-76-5	100			
		dinoseb	88-85-7	100			
		picloram	1918-02-1	100			
		dalapon	75-99-0	100			
LK3-206504	Method 515 Chlorinated Free Acid Herbicide, 8-22				Methanol	1 ml, Amber Glass	€ 48,00
		acifluorfen	50594-66-6	200			
		bentazon	25057-89-0	200			
		chlormamben	133-90-4	200			
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	200			
		3,5-dichlorobenzoic acid	51-36-5	200			
		4-nitrophenol	100-02-7	200			
		pentachlorophenol	87-86-5	200			
		picloram	1918-02-1	200			
LK3-020604	Chlorinated Herbicide Methyl Derivatives, 10 & 1,000 mg/L				Hexane	1 ml, Amber Glass	€ 63,00
LK3-020601	Chlorinated Herbicide Methyl Derivatives, 100 & 10,000 mg/L				Methanol	1 ml, Amber Glass	€ 72,00
LK3-020602	Chlorinated Herbicide Methyl Derivatives, Various Concentrations				Methanol	1 ml, Amber Glass	€ 72,00
		2,4-D methyl ester	1928-38-7	10			
		2,4-DB methyl ester	18625-12-2	10			
		2,4,5-T methyl ester	1928-37-6	10			
		Silvex methyl ester	4841-20-7	10			
		Dalapon methyl ester	17640-02-7	10			
		Dicamba methyl ester	6597-78-0	10			
		Dichlorprop methyl ester	57153-17-0	10			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Dinoseb methyl ester	6099-79-2	10			
		MCPA methyl ester	2436-73-9	1000			
		MCPP methyl ester	23844-56-6	1000			
LK3-020703	EPA 8150 Herbicide Methyl Derivatives Mix, 200mg/L				Methanol	1 ml, Amber Glass	€ 60,00
LK3-020704	EPA 8150 Herbicide Methyl Derivatives Mix, Various Concentrations					1 ml, Amber Glass	€ 87,00
		Dicamba	1918-00-9	100			
		dalapon	75-99-0	100			
		dinoseb	88-85-7	100			
		2,4,5-TP (Silvex)	93-72-1	100			
		2,4-D	94-75-7	100			
		2,4-DB	94-82-6	100			
		2,4,5-T	93-76-5	100			
		dichlorprop (2,4-DP)	120-36-5	100			
		MCPA acid	94-74-6	10000			
		MCPP acid	7085-19-0	10000			
LK3-008701	Custom Herbicide Mix, 10-6	2,4-D	94-75-7	1000	Methanol	1 ml, Amber Glass	€ 79,50
		2,4-DB	94-82-6	1000			
		2,4,5-trichlorophenoxyacetic acid	93-76-5	1000			
		2,4,5-TP (Silvex)	93-72-1	1000			
		Dicamba	1918-00-9	1000			
		dichlorprop (2,4-DP)	120-36-5	1000			
		dinoseb	88-85-7	1000			
		dalapon	75-99-0	1000			
		pentachlorophenol	87-86-5	1000			
		4-nitrophenol	100-02-7	1000			
LK3-233403	Triazinic Herbicides Mix, 100 mg/L				Acetone	1 ml, Amber Glass	€ 111,00
LK3-233401	Triazinic Herbicides Mix, 1000 mg/L					1 ml, Amber Glass	€ 166,50
		atrazine	1912-24-9				
		desethyl atrazine	6190-65-4				

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Atrazine desisopropyl	1007-28-9				
		terbuthylazine	5915-41-3				
		terbuthylazine-desethyl	30125-63-4				
		propazine	139-40-2				
		hexazinon	51235-04-2				
		simazine	122-34-9				
		alachlor	15972-60-8				
		metribuzin	21087-64-9				
		Cyanazine	21725-46-2				
LK3-207702	8151 Herbicide Mix, 11-2, Various Concentrations (All in Acid Form)				Hexane	1 ml, Amber Glass	€ 87,00
		pentachlorophenol	87-86-5	100			
		2,4,5-T	93-76-5	200			
		2,4,5-TP (Silvex)	93-72-1	200			
		Dicamba	1918-00-9	200			
		picloram	1918-02-1	200			
		dinoseb	88-85-7	400			
		dichlorprop (2,4-DP)	120-36-5	800			
		2,4-D	94-75-7	800			
		bentazon	25057-89-0	1600			
		2,4-DB	94-82-6	1800			
		dalapon	75-99-0	2000			
LK3-206805	Custom Herbicide Mix, Underivatized Acids, 11-8				Methanol	1 ml, Amber Glass	€ 87,00
		pentachlorophenol	87-86-5	200			
		dinoseb	88-85-7	200			
		2,4,5-TP (Silvex)	93-72-1	200			
		Dicamba	1918-00-9	200			
		2,4-D	94-75-7	200			
		2,4-DB	94-82-6	200			
		2,4,5-T	93-76-5	200			
		dichlorprop (2,4-DP)	120-36-5	200			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		dalapon	75-99-0	2000			
		MCPP acid	7085-19-0	20000			
		MCPA acid	94-74-6	20000			
LK3-206305	Chlorinated Herbicide Methyl Derivatives, 11-59				Methanol	1 ml, Amber Glass	€ 79,50
		2,4,5-T methyl ester	1928-37-6	200			
		dicamba methyl ester	6597-78-0	200			
		2,4-D-methyl ester	1928-38-7	200			
		silvex methyl ester	4841-20-7	200			
		dalapon-methyl ester	17640-02-7	2000			
		2,4-DB methyl ester	18625-12-2	200			
		dichlorprop methyl ester	57153-17-0	200			
		dinoseb-methyl ether	6099-79-2	200			
		MCPA-methyl ester	2436-73-9	20000			
		Mecoprop-methyl ester	23844-56-6	20000			
		pentachloroanisole	1825-21-4	200			
LK3-234806	Method 515 Derivatized Herbicide Mix, 12-5				M-t-BE	1 ml, Amber Glass	€ 103,50
		pentachloroanisole	1825-21-4	15			
		picloram-methyl ester	14143-55-6	50			
		dicamba methyl ester	6597-78-0	50			
		bentazon-methyl	61592-45-8	50			
		silvex methyl ester	4841-20-7	100			
		dinoseb-methyl ether	6099-79-2	100			
		2,4,5-T methyl ester	1928-37-6	100			
		dichlorprop methyl ester	57153-17-0	350			
		acifluorfen methyl ester	50594-67-7	50			
		dalapon methyl ester	17640-02-7	500			
		2,4-D-methyl ester	1928-38-7	500			
		2,4-DB methyl ester	18625-12-2	1000			
LK3-008806	Method 515 Herbicide Mix, 12-pentachlorophenol		87-86-5	15	M-t-BE	1 ml, Amber Glass	€ 103,50

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		acifluorfen	50594-66-6	50			
		picloram	1918-02-1	50			
		Dicamba	1918-00-9	50			
		bentazon	25057-89-0	50			
		2,4,5-TP (Silvex)	93-72-1	100			
		dinoseb	88-85-7	100			
		2,4,5-T	93-76-5	100			
		dichlorprop (2,4-DP)	120-36-5	350			
		dalapon	75-99-0	500			
		2,4-D	94-75-7	500			
		2,4-DB	94-82-6	1000			
LK3-234203	Chlorinated Herbicides Mixture, 13-1, 100 mg/L				M-t-BE	1 ml, Amber Glass	€ 55,50
		acifluorfen	50594-66-6	100			
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	100			
		Dicamba	1918-00-9	100			
		3,5-dichlorobenzoic acid	51-36-5	100			
		pentachlorophenol	87-86-5	100			
		picloram	1918-02-1	100			
		2,4,5-TP (Silvex)	93-72-1	100			
		2,4,5-T	93-76-5	100			
		bentazon	25057-89-0	100			
		2,4-dichlorophenoxyacetic acid	94-75-7	100			
		dinoseb	88-85-7	100			
		dichlorprop (2,4-DP)	120-36-5	100			
		2,4-DB	94-82-6	100			
LK3-234303	Herbicide Mix, Derivitized, EPA Method 515, 13-1				Hexane:Toluene: 1:1	1 ml, Amber Glass	€ 72,00
		Chlorthal (Dacthal)	2905-67-1	100			
		2,4,5-T methyl ester	1928-37-6	100			
		dicamba methyl ester	6597-78-0	100			
		2,4-D-methyl ester	1928-38-7	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		silvex methyl ester	4841-20-7	100			
		2,4-DB methyl ester	18625-12-2	100			
		dichlorprop-methyl ester	17640-02-7	100			
		dinoseb-methyl ether	6099-79-2	100			
		pentachloroaniso	1825-21-4	100			
		picloram-methyl ester	14143-55-6	100			
		methyl 3,5-dichlorobenzoate	2905-67-1	100			
		acifluorfen methyl ester	50594-67-7	100			
		bentazon-methyl	61592-45-8	100			
LK3-021703	Custom Chlorinated Herbicide Mix, 14-21				Acetonitrile	1 ml, Amber Glass	€ 103,50
LK3-021501	*Custom Chlorinated Phenoxyacide Herbicides, 14-7					1 ml, Amber Glass	€ 118,50
		2,4-D	94-75-7	100			
		2,4-DB	94-82-6	100			
		2,4,5-T	93-76-5	100			
		2,4,5-TP (Silvex)	93-72-1	100			
		dalapon	75-99-0	100			
		Dicamba	1918-00-9	100			
		dichlorprop (2,4-DP)	120-36-5	100			
		dinoseb	88-85-7	100			
		MCPA acid	94-74-6	100			
		MCPA acid	7085-19-0	100			
		2,4-D butoxyethyl ester	1929-73-3	100			
		2,4-D ethylhexyl ester	1928-43-4	100			
		2,4,5-T butoxyethyl ester	2545-59-7	100			
		2,4,5-T butyl ester	93-79-8	100			
LK3-021703-5pck	Custom Chlorinated Herbicide Mix, 14-21				Acetonitrile	1 ml, Amber Glass	€ 346,50
		2,4-DB	94-82-6	100			
		2,4-dichlorophenoxyacetic acid	94-75-7	100			
		picloram	1918-02-1	100			
		dalapon	75-99-0	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,4,5-trichlorophenoxyacetic acid	93-76-5	100			
		2,4,5-TP (Silvex)	93-72-1	100			
		Dicamba	1918-00-9	100			
		dinoseb	88-85-7	100			
		MCPA acid	94-74-6	100			
		MCPP (Mecoprop)	93-65-2	100			
		dichlorprop (2,4-DP)	120-36-5	100			
		acifluorfen	50594-66-6	100			
		bentazon	25057-89-0	100			
		imazapyr	81334-34-1	100			
LK3-021403	Herbicide Mix, 14-27	2,4-D	94-75-7	100	Acetonitrile	1 ml, Amber Glass	€ 118,50
		2,4-DB	94-82-6	100			
		2,4,5-T	93-76-5	100			
		2,4,5-TP (Silvex)	93-72-1	100			
		Dicamba	1918-00-9	100			
		dichlorprop (2,4-DP)	120-36-5	100			
		dinoseb	88-85-7	100			
		MCPA acid	94-74-6	100			
		MCPP acid	7085-19-0	100			
		picloram	1918-02-1	100			
		trifluralin	1582-09-8	100			
		diclofop-methyl	51338-27-3	100			
		brominal	1689-84-5	100			
		bendiocarb	22781-23-3	100			
LK3-021602	Custom Herbicide Mix, 15-14	2,4-D	94-75-7	2000	M-t-BE	1 ml, Amber Glass	€ 157,50
		2,4,5-T	93-76-5	2000			
		2,4,5-TP (Silvex)	93-72-1	2000			
		Dicamba	1918-00-9	2000			
		dichlorprop (2,4-DP)	120-36-5	2000			
		dinoseb	88-85-7	2000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		MCPA acid	94-74-6	2000			
		MCPP acid	7085-19-0	2000			
		picloram	1918-02-1	2000			
		trifluralin	1582-09-8	2000			
		diclofop-methyl	51338-27-3	2000			
		brominal	1689-84-5	2000			
		bendiocarb	22781-23-3	2000			
		dalapon	75-99-0	2000			
		2,4-D butoxyethyl ester	1929-73-3	2000			
LK3-207504	Chlorinated Herbicides Mixture, 100mg/L	16-1			M-t-BE	1 ml, Amber Glass	€ 55,50
LK3-207503	Chlorinated Herbicides Mixture, 16-1, Various Concentrations					1 ml, Amber Glass	€ 118,50
		acifluorfen	50594-66-6				
		chlormamben	133-90-4				
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0				
		Dicamba	1918-00-9				
		3,5-dichlorobenzoic acid	51-36-5				
		4-nitrophenol	100-02-7				
		pentachlorophenol	87-86-5				
		picloram	1918-02-1				
		2,4,5-TP (Silvex)	93-72-1				
		2,4,5-T	93-76-5				
		bentazon	25057-89-0				
		2,4-dichlorophenoxyacetic acid	94-75-7				
		dinoseb	88-85-7				
		dichlorprop (2,4-DP)	120-36-5				
		2,4-DB	94-82-6				
		dalapon	75-99-0				
LK3-234403	Herbicide Mix, Derivitized, EPA Method 515,	16-1			Hexane:Toluene: 1:1	1 ml, Amber Glass	€ 79,50
		chlorthal-dimethyl (dacthal)	1861-32-1	100			
		2,4,5-T methyl ester	1928-37-6	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		dicamba methyl ester	6597-78-0	100			
		2,4-D-methyl ester	1928-38-7	100			
		silvex methyl ester	4841-20-7	100			
		dalapon methyl ester	17640-02-7	100			
		2,4-DB methyl ester	18625-12-2	100			
		dichlorprop methyl ester	57153-17-0	100			
		dinoseb-methyl ether	6099-79-2	100			
		pentachloroanisole	1825-21-4	100			
		picloram-methyl ester	14143-55-6	100			
		4-nitroanisole	100-17-4	100			
		methyl 3,5-dichlorobenzoate	2905-67-1	100			
		acifluorfen methyl ester	50594-67-7	100			
		chloramben methyl ester	7286-84-2	100			
		bentazon-methyl	61592-45-8	100			
LK3-207401	Herbicide Mix, Underivatized Acids, EPA Method 515, 18-1				Acetone	1 ml, Amber Glass	€ 126,00
		MCPP acid	7085-19-0	100			
		MCPA acid	94-74-6	100			
		3,5-dichlorobenzoic acid	51-36-5	100			
		Dicamba	1918-00-9	100			
		pentachlorophenol	87-86-5	100			
		2,4,5-TP (Silvex)	93-72-1	100			
		2,4,5-T	93-76-5	100			
		picloram	1918-02-1	100			
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	100			
		4-nitrophenol	100-02-7	100			
		dichlorprop (2,4-DP)	120-36-5	100			
		2,4-dichlorophenoxyacetic acid	94-75-7	100			
		chloramben	133-90-4	100			
		dinoseb	88-85-7	100			
		bentazone	25057-89-0	100			
		acifluorfen	50594-66-6	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		dalapon	75-99-0	100			
		2,4-DB	94-82-6	100			
LK3-234605	Herbicide Mix, Underivatized Acids, EPA Method 515, 18-26				Methanol	1 ml, Amber Glass	€ 150,00
		MCPP (Mecoprop)	93-65-2	5000			
		MCPA acid	94-74-6	5000			
		3,5-dichlorobenzoic acid	51-36-5	100			
		Dicamba	1918-00-9	100			
		pentachlorophenol	87-86-5	100			
		2,4,5-TP (Silvex)	93-72-1	100			
		2,4,5-T	93-76-5	100			
		picloram	1918-02-1	100			
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	100			
		dalapon	75-99-0	100			
		2,4-DB	94-82-6	100			
		acifluorfen	50594-66-6	100			
		bentazon	25057-89-0	100			
		chloramben	133-90-4	100			
		2,4-dichlorophenoxyacetic acid	94-75-7	100			
		dichlorprop (2,4-DP)	120-36-5	100			
		dinoseb	88-85-7	100			
		4-nitrophenol	100-02-7	100			
LK3-207001	Custom Herbicide Mix, Underivatized Acids, EPA Method 515, 20-1				Tol.:Hex.:MtBE 12	1 ml, Amber Glass	€ 118,50
		MCPA acid	94-74-6	6400			
		MCPP acid	7085-19-0	6400			
		3,5-dichlorobenzoic acid	51-36-5	6,4			
		pentachlorophenol	87-86-5	6,4			
		2-methyl-4,6-dinitrophenol	534-52-1	12,8			
		2,4,5-T	93-76-5	6,4			
		2,4,5-TP (Silvex)	93-72-1	6,4			
		Dicamba	1918-00-9	6,4			

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
picloram	1918-02-1	6,4			
tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	6,4			
4-nitrophenol	100-02-7	16			
2,4-dichlorophenylacetic acid	19719-28-9	32			
dichlorprop (2,4-DP)	120-36-5	32			
2,4-dichlorophenoxyacetic acid	94-75-7	32			
chloramben	133-90-4	32			
dinoseb	88-85-7	16			
bentazon	25057-89-0	32			
acifluorfen	50594-66-6	16			
dalapon	75-99-0	64			
2,4-DB	94-82-6	64			

LK3-208201	Custom Herbicide Mix, Underivatized Acids, EPA Method 515, 20-9		Methanol	1 ml, Amber Glass	€ 111,00
------------	---	--	----------	-------------------	----------

dalapon	75-99-0	200
2,4-DB	94-82-6	200
3,5-dichlorobenzoic acid	51-36-5	200
4-nitrophenol	100-02-7	200
pentachlorophenol	87-86-5	200
picloram	1918-02-1	200
2,4-dichlorophenylacetic acid	19719-28-9	200
Dicamba	1918-00-9	200
dichlorprop (2,4-DP)	120-36-5	200
2,4-dichlorophenoxyacetic acid	94-75-7	200
2,4,5-TP (Silvex)	93-72-1	200
2-methyl-4,6-dinitrophenol	534-52-1	200
2,4,5-T	93-76-5	200
chloramben	133-90-4	200
dinoseb	88-85-7	200
bentazon	25057-89-0	200
tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	200
acifluorfen	50594-66-6	200

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Triclopyr	55335-06-3	200			
		Lontrel	1702-17-6	200			
LK3-206906	Custom Herbicide Mix, Underivatized Acids, EPA Method 515, 22-9				Methanol	1 ml, Amber Glass	€ 126,00
		MCPA acid	94-74-6	2000			
		MCPP acid	7085-19-0	2000			
		dalapon	75-99-0	200			
		2,4-DB	94-82-6	200			
		Dicamba	1918-00-9	200			
		dichlorprop (2,4-DP)	120-36-5	200			
		2,4,5-TP (Silvex)	93-72-1	200			
		2-methyl-4,6-dinitrophenol	534-52-1	200			
		2,4,5-T	93-76-5	200			
		dinoseb	88-85-7	200			
		2, 4 Dichlorophenylacetic acid	94-82-6	200			
		acifluorfen	50594-66-6	200			
		bentazon	25057-89-0	200			
		chloramben	133-90-4	200			
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	200			
		3,5-dichlorobenzoic acid	51-36-5	200			
		4-nitrophenol	100-02-7	200			
		pentachlorophenol	87-86-5	200			
		picloram	1918-02-1	200			
		2,4-dichlorophenoxyacetic acid	94-75-7	200			
		Lontrel	1702-17-6	200			
		Triclopyr	55335-06-3	200			
LK3-006703-4pck	Custom Herbicide Mix, 23-29	p,p'-DDT	50-29-3	100	Hexane	1 ml, Amber Glass	€ 552,00
		o,p'-DDD	53-19-0	100			
		o,p'-DDT	789-02-6	100			
		a-BHC	319-84-6	100			
		b-BHC	319-85-7	100			

	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
	p,p'-DDE	72-55-9	100			
	p,p'-DDD	72-54-8	100			
	o,p'-DDE	3424-82-6	100			
	chlorpyrifos	2921-88-2	100			
	chlorpyrifos-methyl	5598-13-0	100			
	methyl parathion	298-00-0	100			
	malathion	121-75-5	100			
	chlordane (mixture of isomers)	309-00-2	100			
	dieldrin	60-57-1	100			
	aldrin	309-00-2	100			
	endrin	72-20-8	100			
	isodrin	465-73-6	100			
	coumaphos	56-72-4	100			
	hexachlorobenzene	118-74-1	100			
	disulfoton	298-04-4	100			
	sulfotep	3689-24-5	100			
	pirimiphos methyl	29232-93-7	100			
	2,4,5,6-tetrachloro-m-xylene	877-09-8	100			
LK3-421503-5pck Custom Herbicide Mix, 24-29	p,p'-DDT	50-29-3	100	Methylene Chloric 1 ml, Amber Glass		€ 622,50
	o,p'-DDD	53-19-0	100			
	o,p'-DDT	789-02-6	100			
	a-BHC	319-84-6	100			
	b-BHC	319-85-7	100			
	p,p'-DDE	72-55-9	100			
	p,p'-DDD	72-54-8	100			
	o,p'-DDE	3424-82-6	100			
	chlorpyrifos	2921-88-2	100			
	chlorpyrifos-methyl	5598-13-0	100			
	methyl parathion	298-00-0	100			
	malathion	121-75-5	100			
	chlordane (mixture of isomers)	309-00-2	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		dieldrin	60-57-1	100			
		aldrin	309-00-2	100			
		endrin	72-20-8	100			
		isodrin	465-73-6	100			
		coumaphos	56-72-4	100			
		hexachlorobenzene	118-74-1	100			
		disulfoton	298-04-4	100			
		sulfotep	3689-24-5	100			
		pirimiphos methyl	29232-93-7	100			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	100			
		g-BHC	58-89-9	100			
LK3-030502	GC/ECD Performance Check Solution				M-t-BE	1 ml, Amber Glass	€ 28,50
LK3-030502-5pck	GC/ECD Performance Check Solution					5 x 1 ml, Amber	€ 114,00
		d-BHC	319-86-8	0,4			
		chlorothalonil	1897-45-6	0,5			
		chlorpyrifos	2921-88-2	0,02			
		Chlorthal (Dacthal)	1861-32-1	0,5			
LK3-030501	GC/ECD Performance Check Solution, 2/40/50 mg/L				M-t-BE	1 ml, Amber Glass	€ 24,00
		d-BHC	319-86-8	40			
		chlorothalonil	1897-45-6	50			
		chlorpyrifos	2921-88-2	2			
		Chlorthal (Dacthal)	1861-32-1	50			
LK3-209003	GC-MS N/P Solution, 19-1	p-terphenyl-d14	1718-51-0	50	Methylene Chloric	5 ml, Amber Glass	€ 198,00
		2-fluorobiphenyl	321-60-8	50			
		nitrobenzene-d5	4165-60-0	50			
		1,3-dimethyl-2-nitrobenzene	81-20-9	50			
		methyl parathion	298-00-0	50			
		dichlorvos	62-73-7	50			
		diazinon	333-41-5	50			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		disulfoton	298-04-4	50			
		malathion	121-75-5	50			
		Carbophenothion	786-19-6	50			
		Guthion®	86-50-0	50			
		vernolate	1929-77-7	50			
		atrazine	1912-24-9	50			
		chlorpyrifos	2921-88-2	50			
		dimethoate	60-51-5	50			
		Alachlor	15972-60-8	50			
		phorate	298-02-2	50			
		fenchlorphos	299-84-3	50			
		tetrachlorvinphos (Rabon)	961-11-5	50			
LK3-210102	Custom GC-MS Pesticide Solution, 23-1, Various Concentrations				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 87,00
		a-BHC	319-84-6	100			
		b-BHC	319-85-7	100			
		d-BHC	319-86-8	100			
		g-BHC	58-89-9	100			
		aldrin	309-00-2	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		a-chlordane	5103-71-9	100			
		g-chlordane	5103-74-2	100			
		p,p'-DDD	72-54-8	100			
		p,p'-DDE	72-55-9	100			
		p,p'-DDT	50-29-3	100			
		dieldrin	60-57-1	100			
		endrin	72-20-8	100			
		endrin aldehyde	7421-93-4	100			
		endrin ketone	53494-70-5	100			
		endosulfan I	959-98-8	100			
		endosulfan II	33213-65-9	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endosulfan sulfate	1031-07-8	100			
		methoxychlor	72-43-5	100			
		p-terphenyl-d14	1718-51-0	200			
		2-fluorobiphenyl	321-60-8	200			
		nitrobenzene-d5	4165-60-0	200			
LK3-210002	Custom GC-MS Pesticide Solution, 1,2,8 & 24 mg/L , 25-1				Methylene Chloric	1 ml, Amber Glass	€ 55,50
LK3-210004	Custom GC-MS Pesticide Solution, 1,2,8 & 24 mg/L , 25-1				M-t-BE	1 ml, Amber Glass	€ 55,50
LK3-210001	Custom GC-MS Pesticide Solution, 5,10,40 & 120 mg/L , 25-1				Methylene Chloric	5 ml, Amber Glass	€ 198,00
LK3-210003	Custom GC-MS Pesticide Solution, 5,10,40 & 120 mg/L , 25-1				M-t-BE	5 ml, Amber Glass	€ 198,00
		a-BHC	319-84-6	8			
		b-BHC	319-85-7	8			
		d-BHC	319-86-8	8			
		g-BHC	58-89-9	8			
		aldrin	309-00-2	8			
		heptachlor	76-44-8	8			
		heptachlor epoxide	1024-57-3	8			
		a-chlordane	5103-71-9	8			
		g-chlordane	5103-74-2	8			
		p,p'-DDD	72-54-8	8			
		p,p'-DDE	72-55-9	8			
		p,p'-DDT	50-29-3	8			
		dieldrin	60-57-1	8			
		endrin	72-20-8	8			
		endrin aldehyde	7421-93-4	8			
		endrin ketone	53494-70-5	8			
		endosulfan I	959-98-8	8			
		endosulfan II	33213-65-9	8			
		endosulfan sulfate	1031-07-8	8			
		methoxychlor	72-43-5	8			
		p-terphenyl-d14	1718-51-0	24			
		2-fluorobiphenyl	321-60-8	24			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		nitrobenzene-d5	4165-60-0	24			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	1			
		decachlorobiphenyl	2051-24-3	2			
LK3-016701	Acephate and Methamidophos	Acephate	30560-19-1	1000	Acetone	1 ml, Amber Glass	€ 48,00
		methamidophos	10265-92-6	1000			
LK3-20130405	Alachlor & Endosulfan Solution				Hexane	5 ml, BR, Amber	€ 72,00
LK3-20130625	Alachlor & Endosulfan Solution				M-t-BE	25 ml ,BR, Amber	€ 198,00
		alachlor	15972-60-8	400			
		endosulfan (alpha + beta)	115-29-7	400			
LK3-026102	alpha- & gamma-Chlordane Sol	a-chlordane	5103-71-9	1000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 55,50
		g-chlordane	5103-74-2	1000			
LK3-220008	Pesticide High Level Calibration Check Standard, 50 & 500 ug/L,	4-1			Methanol:Aceton	1 ml, Amber Glass	€ 40,50
LK3-220001	Pesticide High Level Calibration Check Standard, 100 & 500 mg/L,	4-1			Isooctane	1 ml, Amber Glass	€ 19,50
		Aroclor 1016	12674-11-2	500			
		Aroclor 1260	11096-82-5	500			
		decachlorobiphenyl	2051-24-3	50			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	50			
LK3-107703	Aroclor 1016, 1248, 1232, & 12	Aroclor 1016	12674-11-2	200	Hexane	1 ml, Amber Glass	€ 55,50
		Aroclor 1248	12672-29-6	200			
		Aroclor 1232	11141-16-5	200			
		Aroclor 1260	11096-82-5	200			
LK3-220216	Aroclor 1221 + 1254 with Surr	Aroclor 1221	11104-28-2	200	Isooctane	1 ml, Amber Glass	€ 19,50
		Aroclor 1254	11097-69-1	100			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	8			
		decachlorobiphenyl (BZ 209)	2051-24-3	16			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-229101	Aroclor 1221 Intermediate Soln	Aroclor 1221	11104-28-2	10	Hexane	25 ml, BR, Amber	€ 48,00
		decachlorobiphenyl (BZ 209)	2051-24-3	1			
LK3-220506	Aroclor 1232 & 1262 Solution	Aroclor 1232	11141-16-5	200	Hexane	1 ml, Amber Glass	€ 24,00
		Aroclor 1262	37324-23-5	200			
LK3-220416	Aroclor 1232 + 1262 with Surrogate	Aroclor 1232	11141-16-5	100	Isooctane	1 ml, Amber Glass	€ 19,50
		Aroclor 1262	37324-23-5	100			
		decachlorobiphenyl (BZ 209)	2051-24-3	16			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	8			
LK3-226001	Aroclor 1232 Intermediate Soln	Aroclor 1232	11141-16-5	10	Hexane	25 ml, BR, Amber	€ 48,00
		decachlorobiphenyl (BZ 209)	2051-24-3	1			
LK3-106510	Aroclor 1232 Soil Standard, 0.1 mg/Kg, 100 gr				Loam	100 ml Amber Glass	€ 772,50
		Aroclor 1232	11141-16-5	0,1			
LK3-001305	Aroclor 1242 + 1260 Solution	Aroclor 1260	11096-82-5	50	Transformer Oil	1 ml, Amber Glass	€ 51,00
		Aroclor 1242	53469-21-9	50			
LK3-220616	Aroclor 1242 + 1268 with TCMX	Aroclor 1242	53469-21-9	100	Isooctane	1 ml, Amber Glass	€ 19,50
		Aroclor 1268	11100-14-4	100			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	8			
LK3-007505-5pck	Chordane, DCB, & TCMX Solution, Various Concentrations				Isooctane	1 ml, Amber Glass	€ 72,00
		chlordane (mixture of isomers)	57-74-9	2			
		decachlorobiphenyl (BZ 209)	2051-24-3	0,64			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	0,32			
LK3-029902	DBC and TCMX Solution	dibutyl chlorendate	1770-80-5	10	Hexane	25 ml, BR, Amber	€ 55,50
		2,4,5,6-tetrachloro-m-xylene	877-09-8	10			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-227001	Aroclor 1242 Intermediate Soln	Aroclor 1242	53469-21-9	10	Hexane	25 ml, BR, Amber	€ 48,00
		decachlorobiphenyl (BZ 209)	2051-24-3	1			
LK3-228001	Aroclor 1248 Intermediate Soln	Aroclor 1248	12672-29-6	10	Hexane	25 ml, BR, Amber	€ 48,00
		decachlorobiphenyl (BZ 209)	2051-24-3	1			
LK3-220816	Aroclor 1248 with Surrogates, 1:1	Aroclor 1248	12672-29-6	100	Isooctane	1 ml, Amber Glass	€ 19,50
		decachlorobiphenyl (BZ 209)	2051-24-3	16			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	8			
LK3-229001	Aroclor 1254 Intermediate Soln	Aroclor 1254	11097-69-1	10	Hexane	25 ml, BR, Amber	€ 48,00
		decachlorobiphenyl (BZ 209)	2051-24-3	1			
LK3-229201	Aroclor 1262 Intermediate Soln	decachlorobiphenyl (BZ 209)	2051-24-3	1	Hexane	25 ml, BR, Amber	€ 48,00
		Aroclor 1262	37324-23-5	10			
LK3-229301	Aroclor 1268 Intermediate Soln	decachlorobiphenyl (BZ 209)	2051-24-3	1	Hexane	25 ml, BR, Amber	€ 48,00
		Aroclor 1268	11100-14-4	10			
LK3-230001	Aroclor Mix, 4-1, 500 mg/L				Hexane	25 ml, BR, Amber	€ 157,50
LK3-230002	Aroclor Mix, 4-1, 1,000 mg/L					25 ml, BR, Amber	€ 237,00
		Aroclor 1242	53469-21-9				
		Aroclor 1254	11097-69-1				
		Aroclor 1248	12672-29-6				
		Aroclor 1262	37324-23-5				
LK3-224103	Carbamate Pesticides Mix, 1:1	aldicarb	116-06-3	100	Methanol	1 ml, Amber Glass	€ 103,50
		aldicarb sulfone	1646-88-4	100			
		aldicarb sulfoxide	1646-87-3	100			
		Carbaryl	63-25-2	100			
		carbofuran	1563-66-2	100			
		3-hydroxycarbofuran	16655-82-6	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Methiocarb	2032-65-7	100			
		methomyl	16752-77-5	100			
		1-naphthol	90-15-3	100			
		oxamyl	23135-22-0	100			
		propoxur	114-26-1	100			
LK3-252503-5pck	*Carbonyl DNPH Derivatives M	acrolein-DNPH	888-54-0	100	Acetonitrile	1 ml, Amber Glass	€ 748,50
		acetone-DNPH	1567-89-1	100			
		benzaldehyde-DNPH	1157-84-2	100			
		butanal-DNPH	1527-98-6	100			
		propionaldehyde-DNPH	725-00-8	100			
		acetaldehyde-DNPH	1019-57-4	100			
		formaldehyde 2,4-dinitro-phenylhydrazone	1081-15-8	100			
		2,5-dimethylbenzaldehyde-DNPH	152477-96-8	100			
		crotonaldehyde-DNPH	1527-96-4	100			
		hexanal-DNPH	1527-97-5	100			
		isovaleraldehyde-DNPH	590-86-3	100			
		M-TOLUALDEHYDE-DNPH	2880-05-9	100			
		o-tolualdehyde-DNPH	1773-44-0	100			
		p-tolualdehyde-DNPH	2571-00-8	100			
		valeraldehyde-DNPH	110-62-3	100			
LK3-252003	Carbonyl DNPH Derivatives Mix, 100 mg/L				Acetonitrile	1 ml, Amber Glass	€ 31,50
LK3-252002	*Carbonyl DNPH Derivatives Mix, Various Concentrations					1 ml, Amber Glass	€ 31,50
		acrolein-DNPH	888-54-0	500			
		acetone-DNPH as acetone	1567-89-1	100			
		benzaldehyde-DNPH	1157-84-2	100			
		butanal-DNPH as butanal	1527-98-6	100			
		propionaldehyde-DNPH	725-00-8	100			
		acetaldehyde-DNPH	1019-57-4	100			
		formaldehyde 2,4-DNPH	1081-15-8	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-684001-4pck	Chloral Hydrate & Chloropicrin	trichloroacetaldehyde, hydrate	302-17-0	1000	Methanol	1 ml, Amber Glass	€ 198,00
		chloropicrin	76-06-2	1000			
LK3-119403	Chlorothalonil & Permethrin Mixture				Methanol	1 ml, Amber Glass	€ 40,50
LK3-119403-4pck	Chlorothalonil & Permethrin Mixture					1 ml, Amber Glass	€ 135,00
		chlorothalonil	1897-45-6	100			
		Permethrin (mixture of isomers)	52645-53-1	100			
LK3-212301	Chlorothalonil, Benzo(g,h,i)perylene, and Fenamiphos Solution,				Acetone	1 ml, Amber Glass	€ 48,00
		chlorothalonil	1897-45-6	500			
		benzo[ghi]perylene	191-24-2	500			
		fenamiphos	22224-92-6	500			
LK3-211201	Chlorothalonil, Metribuzin, and Benzo(g,h,i)perylene Solution				Acetone/MeCl ₂ 1: 1	1 ml, Amber Glass	€ 40,50
		chlorothalonil	1897-45-6	1000			
		metribuzin	21087-64-9	1000			
		benzo[ghi]perylene	191-24-2	1000			
LK3-211903	Chlorothalonil, Benzo(g,h,i)perylene, Fenamiphos and Methyl Paraoxon				Acetone	1 ml, Amber Glass	€ 55,50
		chlorothalonil	1897-45-6	500			
		benzo[ghi]perylene	191-24-2	500			
		fenamiphos	22224-92-6	500			
		dimethyl-p-nitrophenylphosphate	950-35-6	500			
LK3-211901	Chlorothalonil, Metribuzin, Benzo(g,h,i)perylene, Fenamiphos and Methyl Paraoxon				Acetone	1 ml, Amber Glass	€ 55,50
		chlorothalonil	1897-45-6	500			
		metribuzin	21087-64-9	500			
		benzo[ghi]perylene	191-24-2	500			
		fenamiphos	22224-92-6	500			
		dimethyl-p-nitrophenylphosphate	950-35-6	500			
LK3-212903	Chlorothalonil, Benzo(g,h,i)perylene, Fenamiphos, Mevinophos, and Methyl Paraoxon				Acetone	1 ml, Amber Glass	€ 55,50

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		chlorothalonil	1897-45-6	500			
		benzo[ghi]perylene	191-24-2	500			
		fenamiphos	22224-92-6	500			
		methyl paraoxon	950-35-6	500			
		mevinphos (Phosdrin (TM))	7786-34-7	500			
LK3-246401	Chlorpyrifos and Chlorfenvinphos	chlorpyrifos	2921-88-2	100	Acetone	1 ml, Amber Glass	€ 48,00
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	100			
LK3-050205	cis- & trans-Permethrin Calibration Standard, 50 µg/L				Hexane	25 ml, BR, Amber	€ 48,00
LK3-050204	cis- & trans-Permethrin Calibration Standard, 100 µg/L					25 ml, BR, Amber	€ 48,00
LK3-050203	cis- & trans-Permethrin Calibration Standard, 200 µg/L					25 ml, BR, Amber	€ 48,00
LK3-050202	cis- & trans-Permethrin Calibration Standard, 300 µg/L					25 ml, BR, Amber	€ 48,00
LK3-050201	cis- & trans-Permethrin Calibration Standard, 500 µg/L					25 ml, BR, Amber	€ 48,00
LK3-050298	cis- & trans-Permethrin Calibration Set, 50, 100, 200, 300, 500 µg/L ,					25 ml, each	€ 237,00
		cis-permethrin	54774-45-7				
		trans-permethrin	551877-74-8				
		decachlorobiphenyl	2051-24-3				
		2,4,5,6-tetrachloro-m-xylene	877-09-8				
LK3-050302	cis- & trans-Permethrin Solution, 1 mg/L				Hexane	25 ml, BR, Amber	€ 40,50
LK3-050301	cis- & trans-Permethrin Solution, 10 mg/L					10 ml, BR, Amber	€ 79,50
		cis-permethrin	54774-45-7				
		trans-permethrin	551877-74-8				
LK3-200601	Custom Chlorobenzene Mix, 3- 1,2,4,5-tetrachlorobenzene		95-94-3	100	Methylene Chloric	1 ml, Amber Glass	€ 40,50
		pentachlorobenzene	608-93-5	100			
		hexachlorobenzene	118-74-1	100			
LK3-252103	Custom Derivatized Aldehydes (DNPH) Mix, 13-19				Acetonitrile	1 ml, Amber Glass	€ 103,50
		acetaldehyde-DNPH	1019-57-4	100			
		furfural-DNPH	98-01-1	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		acrolein-DNPH	888-54-0	100			
		propionaldehyde-DNPH	725-00-8	100			
		crotonaldehyde-DNPH	1527-96-4	100			
		benzaldehyde-DNPH	1157-84-2	100			
		heptanal-DNPH	111-71-7	100			
		hexanal-DNPH	1527-97-5	100			
		isovaleraldehyde-DNPH	590-86-3	100			
		valeraldehyde-DNPH	110-62-3	100			
		isobutyraldehyde-DNPH	78-84-2	100			
		butanal-DNPH	1527-98-6	100			
		gluteraldehyde-DNPH		100			
LK3-264301	Custom Glycol Mix, 5-50	2-methoxyethanol	109-86-4	10000	Methanol	1 ml, Amber Glass	€ 118,50
		2-ethoxyethanol	110-80-5	10000			
		ethylene glycol	107-21-1	10000			
		di(ethylene glycol)	111-46-6	10000			
		tri(ethylene glycol)	112-27-6	10000			
LK3-035401	Custom Internal Standard Solution, 1 & 0.1 mg/L				Isooctane	1 ml, Amber Glass	€ 2.284,50
LK3-035410	Custom Internal Standard Solution, 3 & 0.3 mg/L					1 ml, Amber Glass	€ 3.150,00
		2,2',4,5',6-pentachlorobiphenyl (BZ# 103)	60145-21-3	1			
		2,2',3,3',4,5,5',6-octachlorobiphenyl (BZ# 198)	68194-17-2	1			
		13C12-2,2',4,4',6-pentabromodiphenyl ether (BDE-100)		1			
		4,4'-DDE (ring-13C12)		1			
		hexachlorobenzene-13C6		1			
		g-chlordane-13C10		1			
		13C6-4-chlorobiphenyl (BZ# 3)	2051-62-9	0,1			
		13C12-4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	0,1			
		13C12-2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-1	0,1			
		13C12-2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	0,1			
		13C12-2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	0,1			
		13C12-2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	0,1			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		13C12-2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	0,1			
		13C12-2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	0,1			
		13C12-2,2',3,3',4,5,5',6,6'-nonachlorobiphenyl (BZ# 2)	52663-77-1	0,1			
		13C12-decachlorobiphenyl	2051-24-3	0,1			
LK3-010601	Custom Mix: Chlorobenzene, 2,4' - DDE, 4,4' - DDE, 2,4' - DDD, 4,4' - DDD, 2,4' - DDT, 4,4' - DDT				Methylene Chloric	1 ml, Amber Glass	€ 78,00
		chlorobenzene	108-90-7	2000			
		o,p'-DDD	53-19-0	2000			
		p,p'-DDD	72-54-8	2000			
		o,p'-DDE	3424-82-6	2000			
		p,p'-DDE	72-55-9	2000			
		o,p'-DDT	789-02-6	2000			
		p,p'-DDT	50-29-3	2000			
LK3-010701	Custom Mix: Chlorobenzene, 4,4' - DDE, 2,4' - DDD, 4,4' - DDD, 4,4' - DDT				Methylene Chloric	1 ml, Amber Glass	€ 118,50
		chlorobenzene	108-90-7	10000			
		o,p'-DDD	53-19-0	10000			
		p,p'-DDD	72-54-8	10000			
		p,p'-DDE	72-55-9	10000			
		p,p'-DDT	50-29-3	10000			
LK3-211610	Custom NP Pesticide Mix 1- 90 isophorone		78-59-1	100	Acetone	1 ml, Amber Glass	€ 315,00
		chlorobenzilate	510-15-6	100			
		propyzamide (pronamide)	23950-58-5	100			
		prometryn	7287-19-6	100			
		propazine	139-40-2	100			
		trifluralin	1582-09-8	100			
		atrazine	1912-24-9	100			
		ethoprophos (prophos)	13194-48-4	100			
		chlorpyrifos	2921-88-2	100			
		alachlor	15972-60-8	100			
		chlorthal-dimethyl (dacthal)	1861-32-1	100			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
propachlor	1918-16-7	100			
dichlorvos	62-73-7	100			
tetrachlorvinphos (Rabon)	961-11-5	100			
acetochlor	34256-82-1	100			
Tebuthiuron	34014-18-1	100			
simetryn	1014-70-6	100			
metolachlor	51218-45-2	100			
EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	100			
ametryne	834-12-8	100			
terbutryne	886-50-0	100			
cis-permethrin	54774-45-7	100			
trans-permethrin	51877-74-8	100			
hexachlorobenzene	118-74-1	100			
hexachlorocyclopentadiene	77-47-4	100			
butachlor	23184-66-9	100			
molinate	2212-67-1	100			
terbacil	5902-51-2	100			
Cycloate	1134-23-2	100			
Fenarimol	60168-88-9	100			
butylate	2008-41-5	100			
vernolate	1929-77-7	100			
pebulate	1114-71-2	100			
etridiazole	2593-15-9	100			
chloroneb	2675-77-6	100			
chlorpropham	101-21-3	100			
bromacil	314-40-9	100			
triadimefon	43121-43-3	100			
MGK-264 - isomer a	113-48-4	100			
MGK-264 - isomer b	113-48-4	100			
diphenamid	957-51-7	100			
napropamide	15299-99-7	100			
beam™	41814-78-2	100			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
norflurazon	27314-13-2	100			
hexazinon	51235-04-2	100			
fluridone	59756-60-4	100			
mevinphos (Phosdrin (TM))	7786-34-7	100			
2,4-dinitrotoluene	121-14-2	100			
2,6-dinitrotoluene	606-20-2	100			
acenaphthylene	208-96-8	100			
benzo[b]fluoranthene	205-99-2	100			
dibenz[a,h]anthracene	53-70-3	100			
p,p'-DDT	50-29-3	100			
p,p'-DDD	72-54-8	100			
p,p'-DDE	72-55-9	100			
a-BHC	319-84-6	100			
g-BHC	58-89-9	100			
dieldrin	60-57-1	100			
methoxychlor	72-43-5	100			
endosulfan II	33213-65-9	100			
endosulfan sulfate	1031-07-8	100			
endrin	72-20-8	100			
endrin aldehyde	7421-93-4	100			
aldrin	309-00-2	100			
b-BHC	319-85-7	100			
d-BHC	319-86-8	100			
heptachlor	76-44-8	100			
heptachlor epoxide isomer B	1024-57-3	100			
a-chlordane	5103-71-9	100			
g-chlordane	5103-74-2	100			
endosulfan I	959-98-8	100			
trans-nonachlor	39765-80-5	100			
anthracene	120-12-7	100			
benzo[a]anthracene	56-55-3	100			
benzo[a]pyrene	50-32-8	100			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		benzo[k]fluoranthene	207-08-9	100			
		chrysene	218-01-9	100			
		indeno[1,2,3-cd]pyrene	193-39-5	100			
		phenanthrene	85-01-8	100			
		pyrene	129-00-0	100			
		2-chlorobiphenyl (BZ# 1)	2051-60-7	100			
		2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	100			
		2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	100			
		2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	100			
		2,2',3',4,6-pentachlorobiphenyl (BZ# 98)	60233-25-2	100			
		2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	100			
		2,2',3,3',4,4',6-heptachlorobiphenyl (BZ# 171)	52663-71-5	100			
		2,2',3,3',4,5',6'-octachlorobiphenyl (BZ# 201)	40186-71-8	100			
		simazine	122-34-9	100			
		fluorene	86-73-7	100			
LK3-211611	Custom NP Pesticide Mix 2- 27	benzo[ghi]perylene	191-24-2	100	Acetone	1 ml, Amber Glass	€ 79,50
		Benfluralin (Benefin)	1861-40-1	100			
		desethyl atrazine	6190-65-4	100			
		Atraton	1610-17-9	100			
		pentachlorophenol	87-86-5	100			
		profluralin	26399-36-0	100			
		terbufos	13071-79-9	100			
		fonofos	944-22-9	100			
		diazinon	333-41-5	100			
		chlorothalonil	1897-45-6	100			
		disulfoton	298-04-4	100			
		chlordene	3734-48-3	100			
		metribuzin	21087-64-9	100			
		methyl parathion	298-00-0	100			
		malathion	121-75-5	100			
		bladex	21725-46-2	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Isopropalin	33820-53-0	100			
		pendimethalin	40487-42-1	100			
		fenamiphos	22224-92-6	100			
		oxadiazon	19666-30-9	100			
		Oxyfluorfen	42874-03-3	100			
		carboxine	5234-68-4	100			
		endrin ketone	53494-70-5	100			
		mirex	2385-85-5	100			
		cis-nonachlor	5103-73-1	100			
		prometon	1610-18-0	100			
		methyl paraoxon	950-35-6	100			
LK3-114503	Custom NP Pesticide Mix, 10-1	butylate	2008-41-5	500	Acetone	25 ml, BR, Amber	€ 198,00
		simazine	122-34-9	500			
		atrazine	1912-24-9	500			
		diazinon	333-41-5	500			
		acetochlor	34256-82-1	500			
		prometryn	7287-19-6	500			
		metolachlor	51218-45-2	500			
		diphenamid	957-51-7	500			
		butachlor	23184-66-9	500			
		chlorobenzilate	510-15-6	500			
LK3-214903	Custom OPP Mix, 5-40	chloropyrifos-methyl	5598-13-0	100	Methanol	1 ml, Amber Glass	€ 87,00
		dimethoate	60-51-5	100			
		malathion	121-75-5	100			
		parathion	56-38-2	100			
		thionazine (zinophos)	297-97-2	100			
LK3-670002	Custom Phenol Mix, 14-162	4-chloro-3-methylphenol	59-50-7	2000	Acetone	1 ml, Amber Glass	€ 63,00
		2-chlorophenol	95-57-8	2000			
		4-chlorophenol	106-48-9	2000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,4-dichlorophenol	120-83-2	2000			
		2,4-dimethylphenol	105-67-9	2000			
		2,5-dimethylphenol	95-87-4	2000			
		3,5-dimethylphenol	108-68-9	2000			
		2-methylphenol	95-48-7	2000			
		3-methylphenol	108-39-4	2000			
		2,4-dinitrophenol	51-28-5	2000			
		pentachlorophenol	87-86-5	2000			
		phenol	108-95-2	2000			
		2,4,5-trichlorophenol	95-95-4	2000			
		2,4,6-trichlorophenol	88-06-2	2000			
LK3-002806	Decachlorobiphenyl, 2,4,5,6-Tetrachloro-m-xylene, & Sulfur Solution				Hexane	25 ml, BR, Amber	€ 103,50
		decachlorobiphenyl (BZ 209)	2051-24-3	0,1			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	0,1			
		sulfur	7704-34-10	100			
LK3-002807	Decachlorobiphenyl, 2,4,5,6-Tetrachloro-m-xylene, & Sulfur Solution				Hexane	25 ml, BR, Amber	€ 118,50
		decachlorobiphenyl (BZ 209)	2051-24-3	0,1			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	0,1			
		sulfur	7704-34-10	2000			
LK3-013106	Decomposition Solution, 10 mg/L				Hexane	1 ml, Amber Glass	€ 31,50
LK3-013105	Decomposition Solution, 500 mg/L				Methanol	1 ml, Amber Glass	€ 24,00
LK3-013103	Decomposition Solution, 200, 100 mg/L				Isooctane	1 ml, Amber Glass	€ 28,50
		p,p'-DDT	50-29-3				
		endrin	72-20-8				
LK3-252203	Derivatized Aldehyde (DNPH as Aldehyde) Mix, 13-9				Acetonitrile	1 ml, Amber Glass	€ 75,00
		acetaldehyde-DNPH	1019-57-4	50			
		furfural-DNPH as furfural	98-01-1	50			
		acrolein-DNPH	888-54-0	50			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		propionaldehyde-DNPH	725-00-8	50			
		crotonaldehyde-DNPH	1527-96-4	50			
		benzaldehyde-DNPH as benzaldehyde	1157-84-2	50			
		heptanal-DNPH as heptanal	111-71-7	50			
		hexanal-DNPH as hexanal	1527-97-5	50			
		isovaleraldehyde-DNPH as isovaleraldehyde	2256-01-1	50			
		valeraldehyde-DNPH as valeraldehyde	2057-84-3	50			
		isobutyraldehyde-DNPH	2057-82-1	50			
		butanal-DNPH as butanal	1527-98-6	50			
		gluteraldehyde-DNPH as gluteraldehyde		50			
LK3-246803	Diazinon, Methyl Parathion, Fenitrothion, & Malathion, 100 mg/L				Acetone	1 ml, Amber Glass	€ 48,00
LK3-246803-5pck	Diazinon, Methyl Parathion, Fenitrothion, & Malathion, 100 mg/L					5 x 1 ml, Amber	€ 189,00
		diazinon	333-41-5	100			
		methyl parathion	298-00-0	100			
		fenitrothion	122-14-5	100			
		malathion	121-75-5	100			
LK3-246905	Diazinon, Parathion, & Fenitrot diazinon		333-41-5	200	Hexane	1 ml, Amber Glass	€ 55,50
		fenitrothion	122-14-5	200			
		parathion	56-38-2	200			
LK3-205203	Dichlofluanid, Heptachlor Epoxide Isomer A, & Tau-fluvalinate Solution				Acetone	1 ml, Amber Glass	€ 48,00
		Dichlofluanid	1085-98-9	100			
		heptachlor epoxide isomer A	28044-83-9	100			
		tau-fluvalinate	102851-06-9	100			
LK3-221704	Dicofol & Mirex Solution, 10 mg/ml, Kelthane® (Dicofol)		115-32-2	10	Hexane	1 ml, Amber Glass	€ 40,50
		mirex	2385-85-5	10			
LK3-257201	Dieldrin & 4,4-DDD Solution	dieldrin	60-57-1	1000	Methylene Chloric	1 ml, Amber Glass	€ 40,50
		p,p'-DDD	72-54-8	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-257401	Dieldrin & 4,4-DDE Soultion	dieldrin	60-57-1	1000	Methylene Chloric	1 ml, Amber Glass	€ 40,50
		p,p'-DDE	72-55-9	1000			
LK3-257601	Dieldrin & 4,4-DDT Soultion	dieldrin	60-57-1	1000	Methylene Chloric	1 ml, Amber Glass	€ 40,50
		p,p'-DDT	50-29-3	1000			
LK3-251504	Dimethoate and Terbufos-sulfc	terbufos-sulfone	56070-16-7	40	EthylAcetate	1 ml, Amber Glass	€ 31,50
		dimethoate	60-51-5	70			
		terbufos-sulfone	56070-16-7	40			
		dimethoate	60-51-5	70			
LK3-200302-3pck	ECD Check Solution, 0.033% (w/v)				Isooctane	3 x 0.5 ml, Amber	€ 72,00
LK3-200303-3pck	ECD Check Solution, 33 pg/uL					1 ml, Amber Glass	€ 58,50
		g-BHC	58-89-9	33			
		aldrin	309-00-2	33			
LK3-001517	Pesticide MDL Solution, 10 ug/L				Methanol	10 ml, BR, Amber	€ 87,00
LK3-001510	CD Pesticide MDL Intermediate Solution, 2, 4, 20 mg/L					25 ml, BR, Amber	€ 157,50
LK3-001511	CD Pesticide MDL Intermediate Solution, 20, 40, 200 mg/L				Tol.:Hex. 1:1	10 ml, BR, Amber	€ 55,50
LK3-001513	CD Pesticide MDL Intermediate Solution, 10, 20, 100 mg/L					25 ml, BR, Amber	€ 72,00
		b-BHC	319-85-7	2			
		d-BHC	319-86-8	2			
		g-BHC	58-89-9	2			
		a-BHC	319-84-6	2			
		a-chlordane	5103-71-9	2			
		g-chlordane	5103-74-2	2			
		heptachlor	76-44-8	2			
		heptachlor epoxide	1024-57-3	2			
		aldrin	309-00-2	4			
		p,p'-DDE	72-55-9	4			
		p,p'-DDT	50-29-3	4			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		p,p'-DDD	72-54-8	4			
		endrin	72-20-8	4			
		endrin aldehyde	7421-93-4	4			
		endrin ketone	53494-70-5	4			
		dieldrin	60-57-1	4			
		endosulfan I	959-98-8	2			
		endosulfan sulfate	1031-07-8	4			
		endosulfan II	33213-65-9	4			
		methoxychlor	72-43-5	20			
LK3-100101	Endosulfan I, II, and Sulfate Sol	endosulfan I	959-98-8	2000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 48,00
		endosulfan II	33213-65-9	2000			
		endosulfan sulfate	1031-07-8	2000			
LK3-200001	Endosulfan Solution	endosulfan I	959-98-8	5	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 63,00
		endosulfan II	33213-65-9	5			
		endosulfan sulfate	1031-07-8	5			
		endosulfan ether	3369-52-6	5			
		endosulfan lacton	3868-61-9	5			
LK3-004502	Floricil Pesticide Evaluation Sta	aldrin	309-00-2	8	Hexane	1 ml, Amber Glass	€ 48,00
		b-BHC	319-85-7	8			
		d-BHC	319-86-8	8			
		a-BHC	319-84-6	8			
		g-BHC	58-89-9	8			
		heptachlor	76-44-8	8			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	8			
		p,p'-DDE	72-55-9	16			
		decachlorobiphenyl	2051-24-3	16			
		endosulfan II	33213-65-9	16			
		endosulfan sulfate	1031-07-8	16			
		endrin aldehyde	7421-93-4	16			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endrin ketone	53494-70-5	16			
LK3-004401	Florilc Pesticide Evaluation Sta	endosulfan I	959-98-8	5	Hexane	25 ml, BR, Amber	€ 79,50
		g-BHC	58-89-9	5			
		a-BHC	319-84-6	5			
		heptachlor	76-44-8	5			
		dieldrin	60-57-1	10			
		endrin	72-20-8	10			
		p,p'-DDT	50-29-3	10			
		p,p'-DDD	72-54-8	10			
		methoxychlor	72-43-5	50			
		decachlorobiphenyl	2051-24-3	20			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	20			
LK3-005501	GAPA MDL Solution, 4 ug/L				M-t-BE	25 ml, BR, Amber	€ 87,00
LK3-005504	GAPA MDL Solution, 100 ug/L					25 ml, BR, Amber	€ 355,50
LK3-005502	GAPA MDL Intermediate Solution, 20 mg/L					1 ml, Amber Glass	€ 135,00
LK3-005503	GAPA MDL Intermediate Solution, 20 mg/L					25 ml, BR, Amber	€ 135,00
		hexachlorobenzene	118-74-1	100			
		o,p'-DDD	53-19-0	100			
		o,p'-DDE	3424-82-6	100			
		o,p'-DDT	789-02-6	100			
		mirex	2385-85-5	100			
		cis-nonachlor	5103-73-1	100			
		trans-nonachlor	39765-80-5	100			
		oxychlorane	27304-13-8	100			
LK3-004602	GAPA Solution				Hexane	1 ml, Amber Glass	€ 355,50
LK3-004603	GAPA Solution				Methanol	10 ml, BR, Amber	€ 591,00
		hexachlorobenzene	118-74-1	200			
		o,p'-DDD	53-19-0	200			
		o,p'-DDE	3424-82-6	200			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		o,p'-DDT	789-02-6	200			
		mirex	2385-85-5	200			
		cis-nonachlor	5103-73-1	200			
		trans-nonachlor	39765-80-5	200			
		oxychlordane	27304-13-8	200			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	200			
		decachlorobiphenyl (BZ 209)	2051-24-3	200			
LK3-004101	GAPA Solution, 11-01	hexachlorobenzene	118-74-1	200	Hexane	1 ml, Amber Glass	€ 315,00
		o,p'-DDD	53-19-0	200			
		o,p'-DDE	3424-82-6	200			
		o,p'-DDT	789-02-6	200			
		mirex	2385-85-5	200			
		cis-nonachlor	5103-73-1	200			
		trans-nonachlor	39765-80-5	200			
		oxychlordane	27304-13-8	200			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	200			
		decachlorobiphenyl	2051-24-3	400			
		tris(2,3-dibromopropyl)phosphate	126-72-7	400			
LK3-201109	Hexachlorobenzene, Chlorpyrifos, Dicofof, & Mirex Calibration Check				Isooctane	1 ml, Amber Glass	€ 24,00
		hexachlorobenzene	118-74-1	80			
		chlorpyrifos	2921-88-2	80			
		Kelthane ® (Dicofof)	115-32-2	80			
		mirex	2385-85-5	80			
LK3-102904	Hexachlorocyclohexanes Mix	a-BHC	319-84-6	10	Acetonitrile	1 ml, Amber Glass	€ 79,50
		b-BHC	319-85-7	10			
		g-BHC	58-89-9	10			
		d-BHC	319-86-8	10			
LK3-020801	Lontrel and Triclopyr	Lontrel	1702-17-6	50	Methanol	1 ml, Amber Glass	€ 48,00

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Tryclopyr	55335-06-3	50			
LK3-207904	MCPA & MCPP Derivatized Sol	MCPA-methyl ester		20000	Hexane	1 ml, Amber Glass	€ 181,50
		Mecoprop-methyl ester		20000			
LK3-253701	Triazine Herbicides Solution,	1alachlor	15972-60-8	1000	Acetone	1 ml, Amber Glass	€ 103,50
		ametryne	834-12-8	1000			
		atrazine	1912-24-9	1000			
		bladex	21725-46-2	1000			
		diuron	330-54-1	1000			
		metolachlor	51218-45-2	1000			
		metribuzin	21087-64-9	1000			
		prometon	1610-18-0	1000			
		prometryn	7287-19-6	1000			
		propazine	139-40-2	1000			
		simazine	122-34-9	1000			
		terbutryne	886-50-0	1000			
LK3-221301	Mirex & Kepone Solution	mirex	2385-85-5	1000	Methanol	1 ml, Amber Glass	€ 48,00
		kepone	143-50-0	1000			
LK3-004801	Mirex/2,4'-DDT Solution	mirex	2385-85-5	1000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 24,00
		o,p'-DDT	789-02-6	1000			
LK3-211501	Nitrogen/ Phosphorus Pesticide	alachlor	15972-60-8	100	Acetone	1 ml, Amber Glass	€ 252,00
		Ametryne	834-12-8	100			
		Atraton	1610-17-9	100			
		atrazine	1912-24-9	100			
		bromacil	314-40-9	100			
		butachlor	23184-66-9	100			
		butylate	2008-41-5	100			
		chlorpropham	101-21-3	100			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
chlorpyrifos	2921-88-2	100			
Cycloate	1134-23-2	100			
bladex	21725-46-2	100			
dichlorvos	62-73-7	100			
diphenamid	957-51-7	100			
EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	100			
ethoprophos (prophos)	13194-48-4	100			
Fenarimol	60168-88-9	100			
fluridone	59756-60-4	100			
hexazinon	51235-04-2	100			
methyl paraoxon	950-35-6	100			
metolachlor	51218-45-2	100			
simazine	122-34-9	100			
mevinphos (Phosdrin (TM))	7786-34-7	100			
MGK-264 Mix of Isomers	113-48-4	100			
molinate	2212-67-1	100			
napropamide	15299-99-7	100			
norflurazon	27314-13-2	100			
pebulate	1114-71-2	100			
prometon	1610-18-0	100			
prometryn	7287-19-6	100			
propyzamide (pronamide)	23950-58-5	100			
propachlor	1918-16-7	100			
propazine	139-40-2	100			
Simetryn	1014-70-6	100			
tetrachlorvinphos (Rabon)	961-11-5	100			
Tebuthiuron	34014-18-1	100			
Terbacil	5902-51-2	100			
Terbutryne	886-50-0	100			
triadimefon	43121-43-3	100			
tricyclazole	41814-78-2	100			
trifluralin	1582-09-8	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		vernolate	1929-77-7	100			
LK3-010902	OC Pesticide Degradation Check Solution, 50 ug/L				Hexane	10 ml, BR, Amber	€ 16,50
LK3-010901	OC Pesticide Degradation Check Solution, 100 mg/L				M-t-BE	1 ml, Amber Glass	€ 24,00
		endrin	72-20-8				
		p,p'-DDT	50-29-3				
LK3-020001	OCC Pesticide Standard, 100 mg/L				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 94,50
LK3-020002	OCC Pesticide Standard, 100 mg/L (bis(2-ethylhexyl)phthalate @ 2,000 mg/L)					1 ml, Amber Glass	€ 94,50
		alachlor	15972-60-8				
		Benfluralin (Benefin)	1861-40-1				
		captan	133-06-2				
		carbophenothion	786-19-6				
		chlorothalonil	1897-45-6				
		chlorthal-dimethyl (dacthal)	1861-32-1				
		2,6-Dichlorobenzonitrile (Dichlobenil)	1194-65-6				
		Kelthane ® (Dicofol)	115-32-2				
		Nitrofen	1836-75-5				
		oxadiazon	19666-30-9				
		Oxyfluorfen	42874-03-3				
		propachlor	1918-16-7				
		o,p'-DDD	53-19-0				
		o,p'-DDE	3424-82-6				
		o,p'-DDT	789-02-6				
		hexachlorobenzene	118-74-1				
		dioctyl phthalate	117-81-7				
		pentachloronitrobenzene	82-68-8				
LK3-025201	OCC2 Pesticide Standard, 500 & 100 mg/L				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 48,00
LK3-025202	OCC2 Pesticide Standard, 5,000 & 1,000 mg/L					1 ml, Amber Glass	€ 48,00
		Perthane (TM)	72-56-0	500			
		Folpet	133-07-3	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Tetradifon	116-29-0	100			
LK3-222001	OEXT Pesticide Spiking Solution, 21-1				Hexane	5 ml, Amber Glass	€ 24,00
LK3-222002	OEXT Pesticide Spiking Solution, 21-1					2 x 25 ml, BR, Amber	€ 135,00
		b-BHC	319-85-7	0,16			
		d-BHC	319-86-8	0,16			
		g-BHC	58-89-9	0,16			
		a-BHC	319-84-6	0,16			
		a-chlordane	5103-71-9	0,16			
		g-chlordane	5103-74-2	0,16			
		heptachlor	76-44-8	0,16			
		heptachlor epoxide isomer B	1024-57-3	0,16			
		aldrin	309-00-2	0,16			
		p,p'-DDE	72-55-9	0,16			
		p,p'-DDT	50-29-3	0,16			
		p,p'-DDD	72-54-8	0,16			
		endrin	72-20-8	0,16			
		endrin aldehyde	7421-93-4	0,16			
		endrin ketone	53494-70-5	0,16			
		dieldrin	60-57-1	0,16			
		endosulfan I	959-98-8	0,16			
		endosulfan sulfate	1031-07-8	0,16			
		endosulfan II	33213-65-9	0,16			
		methoxychlor	72-43-5	0,4			
		mirex	2385-85-5	0,4			
LK3-016001	Organonitrogen Pesticide Solut EPTC (s-ethyl dipropylthiocarbamate)		759-94-4	1000	Hexane	1 ml, Amber Glass	€ 72,00
		metolachlor	51218-45-2	1000			
		atrazine	1912-24-9	1000			
LK3-021001	Organonitrogen Pesticide Solut propazine		139-40-2	200	Hexane	1 ml, Amber Glass	€ 118,50
		prometon	1610-18-0	200			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		simetryn	1014-70-6	200			
		terbuthylazine	5915-41-3	200			
		prometryn	7287-19-6	200			
LK3-015901	Custom Organonitrogen Pesticide Mix, 6 Compounds				Acetone	1 ml, Amber Glass	€ 51,00
		ametryne	834-12-8	200			
		Atraton	1610-17-9	200			
		atrazine	1912-24-9	200			
		bladex	21725-46-2	200			
		terbutryne	886-50-0	200			
		simazine	122-34-9	200			
LK3-101101	Organophosphorus Pesticide Sol	disulfoton	298-04-4	2500	Acetone:Hexane 1	0.5 oz, BR, Amber	€ 55,50
		ethoprophos (prophos)	13194-48-4	2500			
		methyl parathion	298-00-0	2500			
		mevinphos (Phosdrin (TM))	7786-34-7	2500			
		phorate	298-02-2	2500			
LK3-101603	Organophosphorus Pesticide Sol	dimethoate	60-51-5	400	Hexane	1 ml, Amber Glass	€ 48,00
		famphur	52-85-7	400			
		anilazine	101-05-3	400			
LK3-100901	Organophosphorus Pesticide Sol	fenthion	55-38-9	5000	Tol.:Hex. 1:1	0.5 oz, BR, Amber	€ 87,00
		fenchlorphos	299-84-3	5000			
		Tokuthion®	34643-46-4	5000			
		tributylphosphoro-trithioite (Merphos)	150-50-5	5000			
LK3-100801	Organophosphorus Pesticide Sol	trichloronate	327-98-0	5000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 87,00
		sulprofos (Bolstar®)	35400-43-2	5000			
		coumaphos	56-72-4	5000			
		Demeton O&S (24.4% O & 67.2% S)	8065-48-3	5000			
		fensulfothion	115-90-2	5000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-101001	Organophosphorus Pesticide S	tetrachlorvinphos (Rabon)	961-11-5	5000	Tol.:Hex. 1:1	0.5 oz, BR, Amber	€ 87,00
		dibrom	300-76-5	5000			
		diazinon	333-41-5	5000			
		Guthion®	86-50-0	5000			
		dichlorvos	62-73-7	5000			
LK3-240001	Organophosphorus Pesticide S	dimethoate	60-51-5	200	Acetone:Hexane 1:1	1 ml, Amber Glass	€ 118,50
		EPN	2104-64-5	200			
		malathion	121-75-5	200			
		parathion	56-38-2	200			
		sulfotep	3689-24-5	200			
		famphur	52-85-7	200			
		thionazine (zinophos)	297-97-2	200			
LK3-100203	Organophosphorous Pesticide	O,O,O-triethylphosphorothioate	126-68-1	200	Acetone	1 ml, Amber Glass	€ 55,50
		thionazine (zinophos)	297-97-2	200			
		phorate	298-02-2	200			
		dimethoate	60-51-5	200			
		disulfoton	298-04-4	200			
		methyl parathion	298-00-0	200			
		parathion	56-38-2	200			
		famphur	52-85-7	200			
LK3-004701	Organophosphorus Pesticide N	monocrotophos	6923-22-4	200	Methanol	1 ml, Amber Glass	€ 237,00
		mevinphos (Phosdrin (TM))	7786-34-7	200			
		malathion	121-75-5	200			
		methyl parathion	298-00-0	200			
		EPN	2104-64-5	200			
		methidathion	950-37-8	200			
		ethoprophos (prophos)	13194-48-4	200			
		methamidophos	10265-92-6	200			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		dimethoate	60-51-5	200			
LK3-042803	Organophosphorus Solution	chlorpyrifos	2921-88-2	100	Acetone	1 ml, Amber Glass	€ 72,00
		Guthion®	86-50-0	100			
		malathion	121-75-5	100			
		Demeton O&S (24.4% O & 67.2% S)	8065-48-3	100			
		diazinon	333-41-5	100			
		disulfoton	298-04-4	100			
		parathion	56-38-2	100			
		phorate	298-02-2	100			
		tetraethyl dithiopyrophosphate	3689-24-5	100			
LK3-242001	Organophosphorous Pesticide 'Guthion®		86-50-0	2500	Methanol	1 ml, Amber Glass	€ 103,50
		methamidophos	10265-92-6	100			
		diazinon	333-41-5	100			
		malathion	121-75-5	100			
		chlorpyrifos	2921-88-2	100			
		parathion	56-38-2	100			
		ethion	563-12-2	100			
		Carbophenothion	786-19-6	100			
LK3-242201	Organophosphorous Mix, 200 mg/L , 8-13				Acetone	1 ml, Amber Glass	€ 118,50
LK3-242204	Organophosphorous Mix, 1,000 mg/L , 8-13				Acetone:Hexane 1	1 ml, Amber Glass	€ 48,00
		Guthion®	86-50-0	1000			
		Demeton O&S (31.3% O & 57.6% S)	8065-48-3	1000			
		diazinon	333-41-5	1000			
		disulfoton	298-04-4	1000			
		ethion	563-12-2	1000			
		malathion	121-75-5	1000			
		parathion	56-38-2	1000			
		methyl parathion	298-00-0	1000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-100201	Organophosphorous Pesticide Solution 1, 9-1				Hexane	1 ml, Amber Glass	€ 135,00
LK3-100202	Organophosphorous Pesticide Solution 1, 9-1				Methylene Chloric	1 ml, Amber Glass	€ 135,00
LK3-100205	Organophosphorous Pesticide Solution 1, 9-1				Hexane:Acetone: 1	1 ml, Amber Glass	€ 103,50
		O,O,O-triethylphosphorothioate	126-68-1	2000			
		thionazine (zinophos)	297-97-2	2000			
		phorate	298-02-2	2000			
		dimethoate	60-51-5	2000			
		disulfoton	298-04-4	2000			
		methyl parathion	298-00-0	2000			
		parathion	56-38-2	2000			
		famfur	52-85-7	2000			
		sulfotep	3689-24-5	2000			
LK3-035101	Organophosphate Pesticied Mi monocrotophos		6923-22-4	200	Methanol	1 ml, Amber Glass	€ 142,50
		mevinphos (Phosdrin (TM))	7786-34-7	200			
		ethoprophos (prophos)	13194-48-4	200			
		methyl parathion	298-00-0	200			
		EPN	2104-64-5	200			
		malathion	121-75-5	200			
		methidathion	950-37-8	200			
		methamidophos	10265-92-6	200			
		dimethoate	60-51-5	200			
LK3-035901	Organophosphate Pesticied Mi monocrotophos		6923-22-4	200	Methanol	1 ml, Amber Glass	€ 142,50
		mevinphos (Phosdrin (TM))	7786-34-7	200			
		malathion	121-75-5	200			
		methyl parathion	298-00-0	200			
		EPN	2104-64-5	200			
		methidathion	950-37-8	200			
		ethoprophos (prophos)	13194-48-4	200			
		methamidophos	10265-92-6	200			
		dimethoate	60-51-5	200			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		phosfolan	947-02-4	200			
		profenofos	41198-08-7	200			
		triazophos	24017-47-8	200			
		chlorpyrifos	2921-88-2	200			
LK3-035801	Organophosphate Pesticide Mix (15 Compounds)				Methanol	1 ml, Amber Glass	€ 135,00
		monocrotophos	6923-22-4	1000			
		mevinphos (Phosdrin (TM))	7786-34-7	1000			
		malathion	121-75-5	1000			
		methyl parathion	298-00-0	1000			
		EPN	2104-64-5	1000			
		methidathion	950-37-8	1000			
		ethoprophos (prophos)	13194-48-4	1000			
		methamidophos	10265-92-6	1000			
		dimethoate	60-51-5	1000			
		chlorpyrifos	2921-88-2	1000			
		parathion	56-38-2	1000			
		Demeton O&S (24.4% O & 67.2% S)	8065-48-3	1000			
		profenofos	41198-08-7	1000			
		triazophos	24017-47-8	1000			
		phosalone	2310-17-0	1000			
LK3-213001	Organophosphorous Pesticide Solution, 31-1				Hexane	1 ml, Amber Glass	€ 315,00
		chlorpyrifos	2921-88-2	100			
		trichloronate	327-98-0	100			
		sulprofos (Bolstar®)	35400-43-2	100			
		coumaphos	56-72-4	100			
		Dementon O&S (31.3% O & 57.6% S)	8065-48-3	100			
		fensulfothion	115-90-2	100			
		fenthion (Abate®)	3383-96-8	100			
		fenchlorphos	299-84-3	100			
		Tokuthion®	34643-46-4	100			

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
tributylphosphoro-trithioite (Merphos)	150-50-5	100			
tetrachlorvinphos	961-11-5	100			
dibrom	300-76-5	100			
diazinon	333-41-5	100			
Guthion®	86-50-0	100			
dichlorvos	62-73-7	100			
disulfoton	298-04-4	100			
ethoprophos (prophos)	13194-48-4	100			
methyl parathion	298-00-0	100			
mevinphos (Phosdrin (TM))	7786-34-7	100			
phorate	298-02-2	100			
tetraethyldithio pyrophosphate	3689-24-5	100			
dichlofenthion	97-17-6	100			
EPN	2104-64-5	100			
dimethoate	60-51-5	100			
parathion	56-38-2	100			
famphur	52-85-7	200			
malathion	121-75-5	100			
monocrotophos	6923-22-4	100			
thionazine (zinophos)	297-97-2	200			
ethion	563-12-2	100			
TEPP	107-49-3	100			

LK3-101503	Organophosphorus Pesticide Soltuion, 31-03		Hexane	1 ml, Amber Glass	€ 103,50
	trichloronate	327-98-0	400		
	sulprofos (Bolstar®)	35400-43-2	400		
	coumaphos	56-72-4	400		
	Dementon O&S (31.3% O & 57.6% S)	8065-48-3	400		
	fensulfothion	115-90-2	400		
	fenthion	55-38-9	400		
	fenchlorphos	299-84-3	400		
	Tokuthion®	34643-46-4	400		

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
tributylphosphoro-trithioite (Merphos)	150-50-5	400			
tetrachlorvinphos (Rabon)	961-11-5	400			
dibrom	300-76-5	400			
diazinon	333-41-5	400			
Guthion®	86-50-0	400			
dichlorvos	62-73-7	400			
O,O,O-triethylphosphorothioate	126-68-1	400			
thionazine (zinophos)	297-97-2	400			
phorate	298-02-2	400			
disulfoton	298-04-4	400			
methyl parathion	298-00-0	400			
parathion	56-38-2	400			
sulfotep	3689-24-5	400			
atrazine	1912-24-9	400			
propazine	139-40-2	400			
simazine	122-34-9	400			
EPN	2104-64-5	400			
malathion	121-75-5	400			
monocrotophos	6923-22-4	400			
TEPP	107-49-3	400			
chlorpyrifos	2921-88-2	400			
ethoprophos (prophos)	13194-48-4	400			
mevinphos (Phosdrin (TM))	7786-34-7	400			

LK3-101303	Organophosphorus Pesticide Soltuion, 34-4		Acetone:Hexane 1 1 ml, Amber Glass	€ 118,50
	trichloronate	327-98-0	400	
	sulprofos (Bolstar®)	35400-43-2	400	
	coumaphos	56-72-4	400	
	Demeton O&S (24.4% O & 67.2% S)	8065-48-3	400	
	fensulfothion	115-90-2	400	
	fenthion	55-38-9	400	
	fenchlorphos (Ronnel)	299-84-3	400	

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
Tokuthion®	34643-46-4	400			
tributylphosphoro-trithioite (Merphos)	150-50-5	400			
tetrachlorvinphos (Rabon)	961-11-5	400			
dibrom	300-76-5	400			
diazinon	333-41-5	400			
Guthion®	86-50-0	400			
dichlorvos	62-73-7	400			
O,O,O-triethylphosphorothioate	126-68-1	400			
thionazine (zinophos)	297-97-2	400			
phorate	298-02-2	400			
dimethoate	60-51-5	400			
disulfoton	298-04-4	400			
methyl parathion	298-00-0	400			
parathion	56-38-2	400			
famphur	52-85-7	400			
sulfotep	3689-24-5	400			
anilazine	101-05-3	400			
atrazine	1912-24-9	400			
propazine	139-40-2	400			
simazine	122-34-9	400			
EPN	2104-64-5	400			
malathion	121-75-5	400			
monocrotophos	6923-22-4	400			
TEPP	107-49-3	400			
chlorpyrifos	2921-88-2	400			
ethoprophos (prophos)	13194-48-4	400			
mevinphos (Phosdrin (TM))	7786-34-7	400			

LK3-122501	Oxychlordane, trans-Nonachlor, & cis-Nonachlor Solution		Tol.:Hex. 1:1	1 ml, Amber Glass	€ 433,50
	oxychlordane	27304-13-8	1000		
	trans-nonachlor	39765-80-5	1000		
	cis-nonachlor	5103-73-1	1000		

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-004301	P2, CLP Pesticide Mix B, 11-1	a-BHC	319-84-6	5	Tol.:Hex.:Ethyl Ac	1 ml, Amber Glass	€ 48,00
		aldrin	309-00-2	10			
		b-BHC	319-85-7	10			
		d-BHC	319-86-8	10			
		p,p'-DDE	72-55-9	10			
		endrin	72-20-8	10			
		a-chlordane	5103-71-9	10			
		g-chlordane	5103-74-2	10			
		p,p'-DDD	72-54-8	20			
		endosulfan sulfate	1031-07-8	20			
		endrin ketone	53494-70-5	20			
LK3-010401	PCNB/B2EHPhth Solution	dioctyl phthalate	117-81-7	10000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 31,50
		pentachloronitrobenzene	82-68-8	10000			
LK3-036298-5pck	Pesticide Calibration Set (0.2, 1, 2, 5, 15 mg/L, Pesticide IS Only)				Isooctane	1 ml, Amber Glass	€ 513,00
		methoxychlor	72-43-5	15			
		endosulfan I	959-98-8	15			
		endosulfan sulfate	1031-07-8	15			
		endosulfan II	33213-65-9	15			
		a-chlordane	5103-71-9	15			
		g-chlordane	5103-74-2	15			
		oxychlordane	27304-13-8	15			
		trans-nonachlor	39765-80-5	15			
		cis-nonachlor	5103-73-1	15			
		o,p'-DDE	3424-82-6	15			
		o,p'-DDD	53-19-0	15			
		o,p'-DDT	789-02-6	15			
		atrazine	1912-24-9	15			
		hexachlorobenzene	118-74-1	15			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		heptachlor epoxide	1024-57-3	15			
		mirex	2385-85-5	15			
		p,p'-DDE	72-55-9	15			
		p,p'-DDD	72-54-8	15			
		p,p'-DDT	50-29-3	15			
		a-BHC	319-84-6	15			
		b-BHC	319-85-7	15			
		g-BHC	58-89-9	15			
		aldrin	309-00-2	15			
		dieldrin	60-57-1	15			
		4,4'-DDE (ring-13C12)	60145-21-3	0,75			
		hexachlorobenzene-13C6	68194-17-2	0,75			
		g-chlordane-13C10	135-98-8	0,75			
LK3-006607	Pesticide Calibration Solution, 4 ug/L				Hexane	25 ml BR, Amber	€ 79,50
LK3-006606	Pesticide Calibration Solution, 8 ug/L					25 ml BR, Amber	€ 79,50
LK3-006605	Pesticide Calibration Solution, 12 ug/L					25 ml BR, Amber	€ 79,50
LK3-006604	Pesticide Calibration Solution, 16 ug/L					25 ml BR, Amber	€ 79,50
LK3-006603	Pesticide Calibration Solution, 20 ug/L					25 ml BR, Amber	€ 79,50
LK3-006602	Pesticide Calibration Solution, 40 ug/L					25 ml BR, Amber	€ 79,50
LK3-006698	Pesticide Calibration Set, 4, 8, 12, 16, 20, & 40 µg/L					25 ml, each	€ 472,50
		a-BHC	319-84-6				
		b-BHC	319-85-7				
		d-BHC	319-86-8				
		g-BHC	58-89-9				
		aldrin	309-00-2				
		heptachlor	76-44-8				
		heptachlor epoxide isomer B	1024-57-3				
		p,p'-DDD	72-54-8				
		p,p'-DDE	72-55-9				
		p,p'-DDT	50-29-3				
		dieldrin	60-57-1				

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endrin	72-20-8				
		endrin aldehyde	7421-93-4				
		endrin ketone	53494-70-5				
		a-chlordane	5103-71-9				
		g-chlordane	5103-74-2				
		endosulfan I	959-98-8				
		endosulfan II	33213-65-9				
		endosulfan sulfate	1031-07-8				
		methoxychlor	72-43-5				
		decachlorobiphenyl (BZ 209)	2051-24-3				
		2,4,5,6-tetrachloro-m-xylene	877-09-8				
LK3-006903	Pesticide Calibration Solution, 20 ug/L				Hexane	25 ml BR, Amber	€ 79,50
LK3-006902	Pesticide Calibration Solution, 40 ug/L					25 ml BR, Amber	€ 79,50
LK3-006998	Pesticide Calibration Set, 4, 8, 12, 16, 20, & 40 µg/L		319-84-6			25 ml each	€ 472,50
		a-BHC	319-84-6				
		b-BHC	319-85-7				
		d-BHC	319-86-8				
		g-BHC	58-89-9				
		aldrin	309-00-2				
		heptachlor	76-44-8				
		heptachlor epoxide isomer B	1024-57-3				
		p,p'-DDD	72-54-8				
		p,p'-DDE	72-55-9				
		p,p'-DDT	50-29-3				
		dieldrin	60-57-1				
		endrin	72-20-8				
		endrin aldehyde	7421-93-4				
		endrin ketone	53494-70-5				
		a-chlordane	5103-71-9				
		g-chlordane	5103-74-2				

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endosulfan I	959-98-8				
		endosulfan II	33213-65-9				
		endosulfan sulfate	1031-07-8				
		methoxychlor	72-43-5				
		alachlor	15972-60-8				
		decachlorobiphenyl (BZ 209)	2051-24-3				
		2,4,5,6-tetrachloro-m-xylene	877-09-8				
LK3-233602	Pesticide Calibration Standard, 8-5				EthylAcetate	1 ml, Amber Glass	€ 118,50
LK3-233603	Pesticide Calibration Standard, 8-6				Methanol	1 ml, Amber Glass	€ 118,50
		Atrazine desisopropyl	1007-28-9	1000			
		desethyl atrazine	6190-65-4	1000			
		atrazine	1912-24-9	1000			
		alachlor	15972-60-8	1000			
		acetochlor	34256-82-1	1000			
		metolachlor	51218-45-2	1000			
		Cyanazine	21725-46-2	1000			
		simazine	122-34-9	1000			
LK3-005004	Pesticide Performance Evaluation Mix CLP Standard, 1, 2, 5, 10 & 25mg/L, 8-12				Hexane	1 ml, Amber Glass	€ 48,00
LK3-005005	Pesticide Performance Evaluation Mix CLP Standard, 10, 20, 50, 100 & 250mg/L, 8-12				Isooctane	1 ml, Amber Glass	€ 36,00
		a-BHC	319-84-6	1			
		b-BHC	319-85-7	1			
		g-BHC	58-89-9	1			
		decachlorobiphenyl	2051-24-3	2			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	2			
		endrin	72-20-8	5			
		p,p'-DDT	50-29-3	10			
		methoxychlor	72-43-5	25			
LK3-002022	Pesticide CLP GPC Check Solution, 200 & 500 mg/L				Methanol	1 ml, Amber Glass	€ 40,50
LK3-002020	Pesticide CLP GPC Check Solution, 25 & 50 mg/L				Acetone	1 ml, Amber Glass	€ 31,50

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-002021	Pesticide CLP GPC Check Solution, 0,1 & 0,2 mg/L				Acetone	25 ml, BR, Amber	€ 55,50
		g-BHC	58-89-9	0,1			
		heptachlor	76-44-8	0,1			
		aldrin	309-00-2	0,1			
		p,p'-DDT	50-29-3	0,2			
		endrin	72-20-8	0,2			
		dieldrin	60-57-1	0,2			
LK3-034905	Pesticide Custom Calibration Standard, 0.2 mg/L				Isooctane	1 ml, Amber Glass	€ 103,50
LK3-034904	Pesticide Custom Calibration Standard, 1 mg/L					1 ml, Amber Glass	€ 103,50
LK3-034903	Pesticide Custom Calibration Standard, 2 mg/L					1 ml, Amber Glass	€ 103,50
LK3-034902	Pesticide Custom Calibration Standard, 5 mg/L					1 ml, Amber Glass	€ 103,50
LK3-034901	Pesticide Custom Calibration Standard, 15 mg/L					1 ml, Amber Glass	€ 513,00
LK3-034998-5pck	Pesticide Custom Calibration Curve, 0,2, 1, 2, 5, 15 mg/L					1 ml, each	€ 513,00
		methoxychlor	72-43-5				
		endosulfan I	959-98-8				
		endosulfan II	33213-65-9				
		endosulfan sulfate	1031-07-8				
		a-chlordane	5103-71-9				
		g-chlordane	5103-74-2				
		oxychlordane	27304-13-8				
		trans-nonachlor	39765-80-5				
		cis-nonachlor	5103-73-1				
		o,p'-DDE	3424-82-6				
		o,p'-DDD	53-19-0				
		o,p'-DDT	789-02-6				
		atrazine	1912-24-9				
		hexachlorobenzene	118-74-1				
		heptachlor epoxide	1024-57-3				
		mirex	2385-85-5				
		p,p'-DDE	72-55-9				
		p,p'-DDD	72-54-8				

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		p,p'-DDT	50-29-3				
		a-BHC	319-84-6				
		b-BHC	319-85-7				
		g-BHC	58-89-9				
		aldrin	309-00-2				
		dieldrin	60-57-1				
		2,2',4,5',6-pentachlorobiphenyl (BZ# 103)	60145-21-3	0,75			
		2,2',3,3',4,5,5',6-octachlorobiphenyl (BZ# 198)	68194-17-2	0,75			
		13C12-2,2',4,4',6-pentabromodiphenyl ether (BDE-10)	135-98-8	0,75			
		4,4'-DDE (ring-13C12)	98-06-6	0,75			
		hexachlorobenzene-13C6	98-82-8	0,75			
		g-chlordane-13C10	99-87-6	0,75			
		13C6-4-chlorobiphenyl (BZ# 3)	2051-62-9	0,075			
		13C12-4,4'-dichlorobiphenyl (BZ# 15)	2050-68-2	0,075			
		13C12-2,4,4'-trichlorobiphenyl (BZ# 28)	7012-37-1	0,075			
		13C12-2,2',5,5'-tetrachlorobiphenyl (BZ# 52)	35693-99-3	0,075			
		13C12-2,3',4,4',5-pentachlorobiphenyl (BZ# 118)	31508-00-6	0,075			
		13C12-2,2',4,4',5,5'-hexachlorobiphenyl (BZ# 153)	35065-27-1	0,075			
		13C12-2,2',3,4,4',5,5'-heptachlorobiphenyl (BZ# 180)	35065-29-3	0,075			
		13C12-2,2',3,3',4,4',5,5'-octachlorobiphenyl (BZ# 194)	35694-08-7	0,075			
		13C12-2,2',3,3',4,5,5',6,6'-nonachlorobiphenyl (BZ# 2)	52663-77-1	0,075			
		13C12-decachlorobiphenyl	2051-24-3	0,075			
LK3-221508	Pesticide Intermediate Solution, 1, 2, & 3 mg/L , 23-1				Tol.:Hex. 1:1	10 ml, BR, Amber	€ 94,50
LK3-221506	Pesticide Intermediate Solution, 10, 20, & 30 mg/L , 23-1					10 ml, BR, Amber	€ 118,50
		a-BHC	319-84-6	1			
		b-BHC	319-85-7	1			
		d-BHC	319-86-8	1			
		g-BHC	58-89-9	1			
		aldrin	309-00-2	1			
		heptachlor	76-44-8	1			
		heptachlor epoxide isomer B	1024-57-3	1			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		a-chlordane	5103-71-9	1			
		g-chlordane	5103-74-2	1			
		p,p'-DDD	72-54-8	1			
		p,p'-DDE	72-55-9	1			
		p,p'-DDT	50-29-3	1			
		dieldrin	60-57-1	1			
		endrin	72-20-8	1			
		endrin aldehyde	7421-93-4	1			
		endrin ketone	53494-70-5	1			
		endosulfan I	959-98-8	1			
		endosulfan II	33213-65-9	1			
		endosulfan sulfate	1031-07-8	1			
		methoxychlor	72-43-5	3			
		mirex	2385-85-5	3			
		decachlorobiphenyl (BZ 209)	2051-24-3	2			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	2			
LK3-221002	Pesticide Intermediate Solution, 10, 20, 30 & 40 mg/L , 23-1				Tol.:Hex. 1:1	10 ml, BR, Amber	€ 79,50
LK3-221001	Pesticide Intermediate Solution, 40, 80, 120 & 160 mg/L , 23-1				Hexane	25 ml, BR, Amber	€ 118,50
LK3-217001	Pesticide Intermediate Solution, 40, 80, 120 & 160 mg/L , 23-1				Methanol	25 ml, BR, Amber	€ 315,00
		a-BHC	319-84-6	10			
		b-BHC	319-85-7	10			
		d-BHC	319-86-8	10			
		g-BHC	58-89-9	10			
		aldrin	309-00-2	10			
		heptachlor	76-44-8	10			
		heptachlor epoxide isomer B	1024-57-3	10			
		a-chlordane	5103-71-9	10			
		g-chlordane	5103-74-2	10			
		p,p'-DDE	72-55-9	10			
		p,p'-DDT	50-29-3	10			
		p,p'-DDD	72-54-8	10			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endrin	72-20-8	10			
		endrin aldehyde	7421-93-4	10			
		endrin ketone	53494-70-5	10			
		endosulfan I	959-98-8	10			
		endosulfan II	33213-65-9	10			
		endosulfan sulfate	1031-07-8	10			
		mirex	2385-85-5	10			
		dieldrin	60-57-1	20			
		decachlorobiphenyl (BZ 209)	2051-24-3	20			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	30			
		methoxychlor	72-43-5	40			
LK3-001502	Pesticide LCS Solution in Acetonitrile				Acetonitrile	25 ml, BR, Amber	€ 198,00
LK3-001504	Pesticide LCS Solution in Methanol				Methanol	25 ml, BR, Amber	€ 198,00
		a-BHC	319-84-6	0,5			
		b-BHC	319-85-7	0,5			
		d-BHC	319-86-8	0,5			
		g-BHC	58-89-9	0,5			
		aldrin	309-00-2	0,5			
		heptachlor	76-44-8	0,5			
		heptachlor epoxide isomer B	1024-57-3	0,5			
		a-chlordane	5103-71-9	0,5			
		g-chlordane	5103-74-2	0,5			
		p,p'-DDD	72-54-8	1,25			
		p,p'-DDE	72-55-9	1,25			
		p,p'-DDT	50-29-3	1,25			
		dieldrin	60-57-1	1,25			
		endrin	72-20-8	1,25			
		endrin aldehyde	7421-93-4	1,25			
		endrin ketone	53494-70-5	1,25			
		endosulfan I	959-98-8	0,5			
		endosulfan sulfate	1031-07-8	1,25			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endosulfan II	33213-65-9	1,25			
		methoxychlor	72-43-5	5			
LK3-100301	Pesticide LCS Intermediate Sol	a-BHC	319-84-6	100	Tol.:Hex. 1:1	40 ml VOA Vial, Amber	€ 79,50
		b-BHC	319-85-7	100			
		d-BHC	319-86-8	100			
		g-BHC	58-89-9	100			
		aldrin	309-00-2	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
		a-chlordane	5103-71-9	100			
		g-chlordane	5103-74-2	100			
		p,p'-DDD	72-54-8	250			
		p,p'-DDE	72-55-9	250			
		p,p'-DDT	50-29-3	250			
		dieldrin	60-57-1	250			
		endrin	72-20-8	250			
		endrin aldehyde	7421-93-4	250			
		endrin ketone	53494-70-5	250			
LK3-001601	Pesticide LCS Solution, 400 ug/L				Hexane	25 ml, BR, Amber	€ 103,50
LK3-001610	Pesticide LCS Solution, 10 mg/L					25 ml, BR, Amber	€ 157,50
		a-BHC	319-84-6				
		b-BHC	319-85-7				
		d-BHC	319-86-8				
		g-BHC	58-89-9				
		aldrin	309-00-2				
		heptachlor	76-44-8				
		heptachlor epoxide	1024-57-3				
		p,p'-DDD	72-54-8				
		p,p'-DDE	72-55-9				

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		p,p'-DDT	50-29-3				
		dieldrin	60-57-1				
		endrin	72-20-8				
		endrin aldehyde	7421-93-4				
		endrin ketone	53494-70-5				
		a-chlordane	5103-71-9				
		g-chlordane	5103-74-2				
		endosulfan I	959-98-8				
		endosulfan II	33213-65-9				
		endosulfan sulfate	1031-07-8				
		methoxychlor	72-43-5				
		decachlorobiphenyl	2051-24-3				
LK3-001921	Pesticide LCS Solution, 400 ug/	a-BHC	319-84-6	400	M-t-BE	25 ml, BR, Amber	€ 87,00
		b-BHC	319-85-7	400			
		d-BHC	319-86-8	400			
		g-BHC	58-89-9	400			
		aldrin	309-00-2	400			
		heptachlor	76-44-8	400			
		heptachlor epoxide isomer B	1024-57-3	400			
		p,p'-DDD	72-54-8	400			
		p,p'-DDE	72-55-9	400			
		p,p'-DDT	50-29-3	400			
		dieldrin	60-57-1	400			
		endrin	72-20-8	400			
		endrin aldehyde	7421-93-4	400			
		endrin ketone	53494-70-5	400			
		a-chlordane	5103-71-9	400			
		g-chlordane	5103-74-2	400			
		endosulfan I	959-98-8	400			
		endosulfan II	33213-65-9	400			
		endosulfan sulfate	1031-07-8	400			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		methoxychlor	72-43-5	400			
		alachlor	15972-60-8	400			
LK3-004001	Pesticide Performance Evaluati	aldrin	309-00-2	100	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 40,50
		p,p'-DDT	50-29-3	100			
		endrin	72-20-8	100			
LK3-016801	Pesticide Primary Stock Mix	Guthion®	86-50-0	2000	Acetone	1 ml, Amber Glass	€ 198,00
		chlorpyrifos	2921-88-2	2000			
		diazinon	333-41-5	2000			
		dimethoate	60-51-5	2000			
		disulfoton	298-04-4	2000			
		malathion	121-75-5	2000			
		methidathion	950-37-8	2000			
		molinate	2212-67-1	2000			
		methyl parathion	298-00-0	2000			
		phorate	298-02-2	2000			
		benthocarb	28249-77-6	2000			
		atrazine	1912-24-9	2000			
		bladex	21725-46-2	2000			
		tributyl phosphate	126-73-8	2000			
		triphenyl phosphate	115-86-6	2000			
		Demeton O&S (29.0% O & 62.6% S)	8065-48-3	2000			
		dichlorvos	62-73-7	2000			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	2000			
		parathion	56-38-2	2000			
		pendimethalin	40487-42-1	2000			
		trifluralin	1582-09-8	2000			
LK3-051201	Pesticide Resolution Check Mix	g-chlordane	5103-74-2	10	Hexane	1 ml, Amber Glass	€ 45,00
		endosulfan I	959-98-8	10			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		p,p'-DDE	72-55-9	20			
		dieldrin	60-57-1	20			
		endosulfan sulfate	1031-07-8	20			
		endrin ketone	53494-70-5	20			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	20			
		decachlorobiphenyl (BZ 209)	2051-24-3	20			
		methoxychlor	72-43-5	100			
LK3-010101	Pesticide Soltuion	alachlor	15972-60-8	2000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 55,50
		Benfluralin (Benefin)	1861-40-1	2000			
		captan	133-06-2	2000			
		carbophenothion	786-19-6	2000			
LK3-010501	Pesticide Soltuion	chlorothalonil	1897-45-6	2000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 55,50
		chlorthal-dimethyl (dacthal)	1861-32-1	2000			
		2,6-Dichlorobenzonitrile (Dichlobenil)	1194-65-6	2000			
		Kelthane ® (Dicofol)	115-32-2	2000			
LK3-010301	Pesticide Soltuion	Nitrofen	1836-75-5	2000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 55,50
		oxadiazon	19666-30-9	2000			
		Oxyfluorfen	42874-03-3	2000			
		propachlor	1918-16-7	2000			
LK3-028101	Pesticide Solution 10-1	Captafol	2425-06-1	200	Acetone	1 ml, Amber Glass	€ 94,50
		hexachlorobenzene	118-74-1	20			
		Dichloran	99-30-9	40			
		Ethalfuralin	55283-68-6	20			
		Fenarimol	60168-88-9	40			
		Vinclozolin	50471-44-8	20			
		Propiconazol (mixture of isomers)	60207-90-1	100			
		Systhane (Myclobutanil)	88671-89-0	200			
		Baythroid (mixture of isomers)	68359-37-5	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Bifenthrin	82657-04-3	100			
LK3-070501	Pesticide Solution 13-1	Carbaryl	63-25-2	100	Acetone	1 ml, Amber Glass	€ 472,50
		dimethoate	60-51-5	100			
		Pirimiphos-ethyl	23505-41-1	100			
		methyl paraoxon	950-35-6	100			
		Terbufos	13071-79-9	100			
		trans-permethrin	51877-74-8	100			
		Imazalil	35554-44-0	100			
		furathiocarb	65907-30-4	100			
		mecarbam	2595-54-2	100			
		Acephate	30560-19-1	100			
		Dichloran	99-30-9	100			
		Kelthane® (Dicofol)	115-32-2	100			
		acetochlor	34256-82-1	100			
LK3-010801	Pesticide Solution 15-1	dichlorvos	62-73-7	1000	Methylene Chloric	1 ml, Amber Glass	€ 103,50
		diazinon	333-41-5	1000			
		disulfoton	298-04-4	1000			
		methyl parathion	298-00-0	1000			
		malathion	121-75-5	1000			
		Carbophenothion	786-19-6	1000			
		Guthion®	86-50-0	1000			
		phorate	298-02-2	1000			
		fenchlorphos	299-84-3	1000			
		tetrachlorvinphos	961-11-5	1000			
		Alachlor	15972-60-8	1000			
		dimethoate	60-51-5	1000			
		vernolate	1929-77-7	1000			
		chlorpyrifos	2921-88-2	1000			
		atrazine	1912-24-9	1000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-070001	Pesticide Solution 25-1	dichlorvos	62-73-7	100	Acetone	1 ml, Amber Glass	€ 945,00
		mevinphos (Phosdrin (TM))	7786-34-7	100			
		trichlorfon	52-68-6	100			
		Omethoate	1113-02-6	100			
		diphenylamine	122-39-4	100			
		monocrotophos	6923-22-4	100			
		Demeton S	126-75-0	100			
		Ethoxyquin	91-53-2	100			
		pentachloronitrobenzene	82-68-8	100			
		disulfoton	298-04-4	100			
		Vinclozolin	50471-44-8	100			
		paraoxon ethyl	311-45-5	100			
		Methiocarb	2032-65-7	100			
		parathion	56-38-2	100			
		captan	133-06-2	100			
		Folpet	133-07-3	100			
		methidathion	950-37-8	100			
		Vamidothion	2275-23-2	100			
		fenamiphos	22224-92-6	100			
		bupirimate	41483-43-6	100			
		tebuconazol	107534-96-3	100			
		Guthion®	86-50-0	100			
		Cypermethrin	52315-07-8	100			
		fenvalerate	51630-58-1	100			
		Deltamethrin	52918-63-5	100			
LK3-070101	Pesticide Solution 25-2	chlorpropham	101-21-3	100	Acetone	1 ml, Amber Glass	€ 945,00
		phorate	298-02-2	100			
		methyl paraoxon	950-35-6	100			
		Tolclofos-Methyl	57018-04-9	100			
		fenitrothion	122-14-5	100			
		Dichlofluanid	1085-98-9	100			

			Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
			fenthion	55-38-9	100			
			chlorpyrifos	2921-88-2	100			
			penconazol	66246-88-6	100			
			Isofenphos	25311-71-1	100			
			Phenthoate	2597-03-7	100			
			Procymidone	32809-16-8	100			
			endosulfan I	959-98-8	100			
			dieldrin	60-57-1	100			
			Systhane (TM)	88671-89-0	100			
			endosulfan II	33213-65-9	100			
			p,p'-DDD	72-54-8	100			
			ethion	563-12-2	100			
			endosulfan sulfate	1031-07-8	100			
			p,p'-DDT	50-29-3	100			
			nuarimol	63284-71-9	100			
			Pyridaphenthion	119-12-0	100			
			Tetradifon	116-29-0	100			
			Fenarimol	60168-88-9	100			
			Guthion Ethyl	2642-71-9	100			
LK3-070401	Pesticide Solution	25-4	chlorotoluron	15545-48-9	100	Acetone	1 ml, Amber Glass	€ 945,00
			Atrazine desisopropyl	1007-28-9	100			
			desethyl atrazine	6190-65-4	100			
			trifluralin	1582-09-8	100			
			Benfluralin (Benefin)	1861-40-1	100			
			simazine	122-34-9	100			
			prometon	1610-18-0	100			
			Terbumeton	33693-04-8	100			
			propyzamide (pronamide)	23950-58-5	100			
			secbumethon	26259-45-0	100			
			Dimethachlon	24096-53-5	100			
			alachlor	15972-60-8	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Metalaxyl	57837-19-1	100			
		linuron	330-55-2	100			
		metolachlor	51218-45-2	100			
		fenson	80-38-6	100			
		nitroal isopropyl	10552-74-6	100			
		Isopropalin	33820-53-0	100			
		thiabendazole	148-79-8	100			
		pendimethalin	40487-42-1	100			
		tetrachlorvinphos (Rabon)	961-11-5	100			
		Methoprottryne	841-06-5	100			
		o,p'-DDT	789-02-6	100			
		Oxadixyl	77732-09-3	100			
		Prochloraz	67747-09-5	100			
LK3-028701	Pesticide Solution, 6-7	eserine	57-47-6	1000	Acetone	1 ml, Amber Glass	€ 118,50
		eserine salicylate	57-64-7	1000			
		formetanate HCL	23422-53-9	1000			
		prosulfocarb	52888-80-9	1000			
		metolcarb	1129-41-5	1000			
		thiophanate methyl	23564-05-8	1000			
LK3-028803	Pesticide Solution, 8-8	prosulfocarb	52888-80-9	100	Acetone	1 ml, Amber Glass	€ 72,00
		fenitrothion	122-14-5	100			
		malathion	121-75-5	100			
		metolachlor	51218-45-2	100			
		fenthion	55-38-9	100			
		chlorpyrifos	2921-88-2	100			
		parathion	56-38-2	100			
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	100			
LK3-035005	Pesticide Spike Solution, 0.5 mg/L					10 ml, BR, Amber	€ 315,00
LK3-035002	Pesticide Spike Solution, 1 mg/L					1 ml, Amber Glass	€ 237,00

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-035003	Pesticide Spike Solution, 20 mg/L	methoxychlor	72-43-5		Isooctane	1 ml, Amber Glass	€ 237,00
		endosulfan I	959-98-8				
		endosulfan II	33213-65-9				
		endosulfan sulfate	1031-07-8				
		a-chlordane	5103-71-9				
		g-chlordane	5103-74-2				
		oxychlordane	27304-13-8				
		trans-nonachlor	39765-80-5				
		cis-nonachlor	5103-73-1				
		o,p'-DDE	3424-82-6				
		o,p'-DDD	53-19-0				
		o,p'-DDT	789-02-6				
		atrazine	1912-24-9				
		hexachlorobenzene	118-74-1				
		heptachlor epoxide isomer B	1024-57-3				
		mirex	2385-85-5				
		p,p'-DDE	72-55-9				
		p,p'-DDD	72-54-8				
		p,p'-DDT	50-29-3				
		a-BHC	319-84-6				
		b-BHC	319-85-7				
		g-BHC	58-89-9				
		aldrin	309-00-2				
		dieldrin	60-57-1				
LK3-037104	Pesticide Spike Solution, 26-3	methoxychlor	72-43-5	1	Isooctane	1 ml, Amber Glass	€ 118,50
		endosulfan I	959-98-8	1			
		endosulfan II	33213-65-9	1			
		endosulfan sulfate	1031-07-8	1			
		oxychlordane	27304-13-8	1			
		cis-nonachlor	5103-73-1	1			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		o,p'-DDT	789-02-6	1			
		atrazine	1912-24-9	1			
		hexachlorobenzene	118-74-1	1			
		mirex	2385-85-5	1			
		dieldrin	60-57-1	1			
		chlorpyrifos	2921-88-2	1			
		a-BHC	319-84-6	1			
		b-BHC	319-85-7	1			
		o,p'-DDE	3424-82-6	1			
		p,p'-DDE	72-55-9	1			
		p,p'-DDD	72-54-8	1			
		g-BHC	58-89-9	1			
		heptachlor	76-44-8	1			
		heptachlor epoxide isomer B	1024-57-3	1			
		aldrin	309-00-2	1			
		a-chlordane	5103-71-9	1			
		g-chlordane	5103-74-2	1			
		trans-nonachlor	39765-80-5	1			
		p,p'-DDT	50-29-3	1			
		o,p'-DDD	53-19-0	1			
LK3-104001	Pesticide Standard Mix A, 9-1	a-BHC	319-84-6	5	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 63,00
		dieldrin	60-57-1	10			
		endosulfan I	959-98-8	5			
		endrin	72-20-8	10			
		g-BHC	58-89-9	5			
		heptachlor	76-44-8	5			
		methoxychlor	72-43-5	50			
		p,p'-DDD	72-54-8	10			
		p,p'-DDT	50-29-3	10			
LK3-104101	Pesticide Standard Mix B, 11-1				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 54,00

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-104102	Pesticide Standard Mix B, 11-1	aldrin	309-00-2	5	Hexane:Toluene 1:1	1 ml, Amber Glass	€ 54,00
		a-chlordane	5103-71-9	5			
		b-BHC	319-85-7	5			
		d-BHC	319-86-8	5			
		endosulfan II	33213-65-9	10			
		endosulfan sulfate	1031-07-8	10			
		endrin aldehyde	7421-93-4	10			
		endrin ketone	53494-70-5	10			
		g-chlordane	5103-74-2	5			
		heptachlor epoxide	1024-57-3	5			
		p,p'-DDE	72-55-9	10			
LK3-002306	Pesticide Surrogate Solution, 10 mg/L				Hexane	25 ml, BR, Amber	€ 60,00
LK3-002309	Pesticide Surrogate Solution, 10 mg/L				Methanol	25 ml, BR, Amber	€ 60,00
LK3-002312	Pesticide Surrogate Solution, 1,000 ug/L	decachlorobiphenyl (BZ 209)	2051-24-3		Tol.:Hex. 1:1	1 ml, Amber Glass	€ 40,50
		2,4,5,6-tetrachloro-m-xylene	877-09-8				
LK3-007001	Pesticide Surrogate Solution, 1,000 mg/L					1 ml, Amber Glass	€ 40,50
LK3-007002	Pesticide Surrogate Solution, 5,000 mg/L	dibutyl chlorendate	1770-80-5	1000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 79,50
		decachlorobiphenyl	2051-24-3	1000			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	1000			
LK3-245002	TCL Pesticides Mix, 18-1				Methanol	1 ml, Amber Glass	€ 87,00
LK3-245001	TCL Pesticides Mix, 18-1				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 118,50
LK3-245001-5pck	TCL Pesticides Mix, 18-1				Tol.:Hex. 1:1	5 x 1 ml, Amber	€ 118,50
		a-BHC	319-84-6	1000			
		b-BHC	319-85-7	1000			
		d-BHC	319-86-8	1000			
		g-BHC	58-89-9	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		aldrin	309-00-2	1000			
		heptachlor	76-44-8	1000			
		heptachlor epoxide isomer B	1024-57-3	1000			
		endosulfan I	959-98-8	1000			
		endosulfan II	33213-65-9	1000			
		endosulfan sulfate	1031-07-8	1000			
		p,p'-DDD	72-54-8	1000			
		p,p'-DDE	72-55-9	1000			
		p,p'-DDT	50-29-3	1000			
		dieldrin	60-57-1	1000			
		endrin	72-20-8	1000			
		endrin aldehyde	7421-93-4	1000			
		endrin ketone	53494-70-5	1000			
		methoxychlor	72-43-5	1000			
LK3-054101	TCLP Herbicide Spiking Solution, 2 mg/L				Methanol	25 ml, BR, Amber	€ 24,00
LK3-054102	TCLP Herbicide Spiking Solution, 1000 mg/L					1 ml, Amber Glass	€ 31,50
		2,4-D	94-75-7				
		2,4,5-TP (Silvex)	93-72-1				
LK3-001801	TCLP Pesticide Solution				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 40,50
LK3-001802-4pck	TCLP Pesticide Solution				Methanol	4 x 1 ml, Amber	€ 94,50
		endrin	72-20-8	100			
		methoxychlor	72-43-5	1000			
		g-BHC	58-89-9	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide isomer B	1024-57-3	100			
LK3-001803	TCLP Pesticide Spiking Solution, 2,000 mg/L				Hexane	1 ml, Amber Glass	€ 31,50
		endrin	72-20-8	2000			
		heptachlor	76-44-8	2000			
		heptachlor epoxide	1024-57-3	2000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		g-BHC	58-89-9	2000			
		methoxychlor	72-43-5	2000			
LK3-003901	TCLP Pesticide Spiking Solution	chlordane	57-74-9	1000	Methanol	1 ml, Amber Glass	€ 36,00
		endrin	72-20-8	1000			
		heptachlor	76-44-8	1000			
		heptachlor epoxide	1024-57-3	1000			
		g-BHC	58-89-9	1000			
		methoxychlor	72-43-5	1000			
		toxaphene	8001-35-2	2000			
		chlordane	57-74-9	1000			
		endrin	72-20-8	1000			
		heptachlor	76-44-8	1000			
		heptachlor epoxide	1024-57-3	1000			
		g-BHC	58-89-9	1000			
		methoxychlor	72-43-5	1000			
		toxaphene	8001-35-2	2000			
LK3-023805	Technical Chlordane & Toxaphene	chlordane, technical mixture	12789-03-6	2000	Methanol	1 ml, Amber Glass	€ 31,50
		toxaphene	8001-35-2	4000			
LK3-020102	Technical Chlordane Constituent Solution				Isooctane	1 ml, Amber Glass	€ 16,50
LK3-020101	Technical Chlordane Constituent Solution					5 ml, Amber Glass	€ 31,50
		chlordane	3734-48-3	0,2			
		cis-nonachlor	5103-73-1	0,4			
		trans-nonachlor	39765-80-5	0,4			
		2,4,5,6-tetrachloro-m-xylene	877-09-8	0,4			
		decachlorobiphenyl	2051-24-3	0,8			
LK3-025001	Toxaphene with Surrogates Solution	toxaphene	8001-35-2	1000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 40,50
		2,4,5,6-tetrachloro-m-xylene	877-09-8	500			
		decachlorobiphenyl (BZ 209)		50			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-253503	Triazine Solution, 11-25	terbutryne	886-50-0	100	Acetone	1 ml, Amber Glass	€ 103,50
		prometon	1610-18-0	100			
		simazine	122-34-9	100			
		atrazine	1912-24-9	100			
		propazine	139-40-2	100			
		metribuzin	21087-64-9	100			
		alachlor	15972-60-8	100			
		ametryne	834-12-8	100			
		prometryn	7287-19-6	100			
		metolachlor	51218-45-2	100			
		bladex	21725-46-2	100			
LK3-016101	Tributyl- and Triphenylphosphate Soltuion				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 55,50
LK3-016102	Tributyl- and Triphenylphosphate Soltuion				Acetone	1 ml, Amber Glass	€ 36,00
		tributyl phosphate	126-73-8	1000			
		triphenyl phosphate	115-86-6	1000			
LK3-103801	1-Chloronaphthalene and 1,4-Dichloronaphthalene				Methylene Chloric	1 ml, Amber Glass	€ 40,50
		1-chloronaphthalene	90-13-1	100			
		1,4-dichloronaphthalene	1825-31-6	100			
LK3-255202	2,4-D Methyl Ester & Silvex Methyl Ester Solution				Hexane	1 ml, Amber Glass	€ 31,50
		2,4-D-methyl ester	1928-38-7	2000			
		silvex methyl ester	4841-20-7	2000			
LK3-122001	2,4-DDE, 2,4-DDD, & 2,4-DDT Solution				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 103,50
		o,p'-DDE	3424-82-6	1000			
		o,p'-DDD	53-19-0	1000			
		o,p'-DDT	789-02-6	1000			
LK3-255301	2,4,5-TP (Silvex) & 2,4-Dichlorophenoxyacetic acid (2,4-D) Mix				Acetone	1 ml, Amber Glass	€ 103,50

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,4,5-TP (Silvex)	93-72-1	500			
		2,4-dichlorophenoxyacetic acid	94-75-7	500			
LK3-248801	4,4'-DDE, 2,4'-DDT, & 4,4'-DDD Solution				Toluene	1 ml, Amber Glass	€ 55,50
		p,p'-DDE	72-55-9	1000			
		o,p'-DDT	789-02-6	1000			
		p,p'-DDD	72-54-8	1000			
LK3-002104	4,4'-Dibromooctafluorobiphenyl & Decachlorobiphenyl (BZ 209) Solution				Toluene	1 ml, Amber Glass	€ 63,00
		4,4'-dibromooctafluorobiphenyl	10386-84-2	10000			
		decachlorobiphenyl (BZ 209)	2051-24-3	10000			
LK3-005801	Method 505 Organohalide Pesticide Solution				Methanol	1 ml, Amber Glass	€ 72,00
		Alachlor	15972-60-8	100			
		aldrin	309-00-2	10			
		trans-nonachlor	39765-80-5	10			
		a-chlordane	5103-71-9	10			
		g-chlordane	5103-74-2	10			
		dieldrin	60-57-1	10			
		endrin	72-20-8	10			
		heptachlor	76-44-8	10			
		heptachlor epoxide	1024-57-3	10			
		hexachlorobenzene	118-74-1	10			
		hexachlorocyclopentadiene	77-47-4	10			
		g-BHC	58-89-9	10			
		methoxychlor	72-43-5	50			
		cis-nonachlor	5103-73-1	10			
LK3-005702	Method 505 Organohalide Pesticide Solution, 100 mg/L				Acetone	1 ml, Amber Glass	€ 48,00
LK3-005703	Method 505 Organohalide Pesticide Solution, 200 mg/L				Methanol	1 ml, Amber Glass	€ 55,50
		Alachlor	15972-60-8	100			
		aldrin	309-00-2	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		trans-nonachlor	39765-80-5	100			
		a-chlordane	5103-71-9	100			
		g-chlordane	5103-74-2	100			
		dieldrin	60-57-1	100			
		endrin	72-20-8	100			
		heptachlor	76-44-8	100			
		heptachlor epoxide	1024-57-3	100			
		hexachlorobenzene	118-74-1	100			
		hexachlorocyclopentadiene	77-47-4	100			
		g-BHC	58-89-9	100			
		methoxychlor	72-43-5	100			
		cis-nonachlor	5103-73-1	100			
		atrazine	1912-24-9	100			
		Simazine	122-34-9	100			
LK3-104501	EPA Method 507 Solution, 6-1				M-t-BE	1 ml, Amber Glass	€ 40,50
LK3-104501-5pck	EPA Method 507 Solution, 6-1					5 x 1 ml, Amber	€ 157,50
		Ametryne	834-12-8	1000			
		Cycloate	1134-23-2	1000			
		disulfoton	298-04-4	1000			
		fenamiphos	22224-92-6	1000			
		tributylphosphoro-trithioite (Merphos)	150-50-5	1000			
		Prometon	1610-18-0	1000			
LK3-208501	Custom 507 Standard A, 9-10	atrazine	1912-24-9	1000	Acetone	1 ml, Amber Glass	€ 72,00
		dimethoate	60-51-5	1000			
		disulfoton	298-04-4	1000			
		prometon	1610-18-0	1000			
		benthiocarb	28249-77-6	1000			
		prometryn	7287-19-6	1000			
		butachlor	23184-66-9	1000			
		metribuzin	21087-64-9	1000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		propyzamide (pronamide)	23950-58-5	1000			
LK3-208701	Custom 507 Standard B	mevinphos (Phosdrin (TM))	7786-34-7	1000	Acetone	1 ml, Amber Glass	€ 87,00
		molinate	2212-67-1	1000			
		metolachlor	51218-45-2	1000			
		Demeton O&S (27.7% O & 63% S)	8065-48-3	1000			
		diazinon	333-41-5	1000			
		simetryn	1014-70-6	1000			
		terbutryne	886-50-0	1000			
		bromacil	314-40-9	1000			
		simazine	122-34-9	1000			
		alachlor	15972-60-8	1000			
		tributyl phosphate	126-73-8	1000			
		triphenyl phosphate	115-86-6	1000			
LK3-211706	Method 507 Pesticides Mix, 32-3				M-t-BE	1 ml, Amber Glass	€ 394,50
		ametryne	834-12-8	2500			
		disulfoton sulfone	2497-06-5	2500			
		disulfoton sulfoxide	2497-07-6	2500			
		methyl paraoxon	950-35-6	2500			
		Atraton	1610-17-9	1000			
		bromacil	314-40-9	1000			
		carboxine	5234-68-4	1000			
		chlorpropham	101-21-3	1000			
		Cycloate	1134-23-2	1000			
		diazinon	333-41-5	1000			
		dichlorvos	62-73-7	1000			
		diphenamid	957-51-7	1000			
		disulfoton	298-04-4	1000			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000			
		ethoprophos (prophos)	13194-48-4	1000			
		fenamiphos	22224-92-6	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Fenarimol	60168-88-9	1000			
		fluridone	59756-60-4	1000			
		hexazinon	51235-04-2	1000			
		tributylphosphoro-trithioite (Merphos)	150-50-5	1000			
		mevinphos (Phosdrin (TM))	7786-34-7	1000			
		MGK-264 Mix of Isomers (65%a,35%b)	113-48-4	1000			
		napropamide	15299-99-7	1000			
		norflurazon	27314-13-2	1000			
		pebulate	1114-71-2	1000			
		prometon	1610-18-0	1000			
		prometryn	7287-19-6	1000			
		propyzamide (pronamide)	23950-58-5	1000			
		simetryn	1014-70-6	1000			
		tetrachlorvinphos (Stirophos)	961-11-5	1000			
		Tebuthiuron	34014-18-1	1000			
		terbacil	5902-51-2	1000			
LK3-053401	Method 508.1 Solution 2-3	atrazine	1912-24-9	1000	Ethyl Acetate	1 ml, Amber Glass	€ 31,50
		simazine	122-34-9	1000			
LK3-053101	Method 508.1 Solution 5-1	chlorothalonil	1897-45-6	500	Ethyl Acetate	1 ml, Amber Glass	€ 48,00
		etridiazole	2593-15-9	500			
		hexachlorobenzene	118-74-1	500			
		hexachlorocyclopentadiene	77-47-4	500			
		trifluralin	1582-09-8	500			
LK3-053301	Method 508.1 Solution 7-1	alachlor	15972-60-8	5000	Ethyl Acetate	1 ml, Amber Glass	€ 94,50
		chloroneb	2675-77-6	5000			
		Cyanazine	21725-46-2	5000			
		metolachlor	51218-45-2	5000			
		propachlor	1918-16-7	5000			
		butachlor	23184-66-9	5000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Permethrin (mixture of isomers)	52645-53-1	5000			
LK3-233501	508.1 Chlorinated Pesticides Calibration Mix				EthylAcetate	1 ml, Amber Glass	€ 63,00
		alachlor	15972-60-8	500			
		atrazine	1912-24-9	500			
		chlorothalonil	1897-45-6	500			
		Cyanazine	21725-46-2	500			
		hexachlorocyclopentadiene	77-47-4	500			
		metolachlor	51218-45-2	500			
		metribuzin	21087-64-9	500			
		simazine	122-34-9	500			
LK3-013005	Method 508-Chlorinated Pesticide Solution, 2 mg/L , 17-1				EthylAcetate	1 ml, Amber Glass	€ 40,50
LK3-013004	Method 508-Chlorinated Pesticide Solution, 500 mg/L , 17-1				Ethyl Acetate	1 ml, Amber Glass	€ 94,50
LK3-013002	Method 508-Chlorinated Pesticide Solution, 1,000 mg/L , 17-1				M-t-BE	1 ml, Amber Glass	€ 126,00
LK3-013006	Method 508-Chlorinated Pesticide Solution, 1,000 mg/L , 17-1				Hexane	1 ml, Amber Glass	€ 150,00
		aldrin	309-00-2	1000			
		a-BHC	319-84-6	1000			
		b-BHC	319-85-7	1000			
		d-BHC	319-86-8	1000			
		g-BHC	58-89-9	1000			
		p,p'-DDD	72-54-8	1000			
		p,p'-DDE	72-55-9	1000			
		p,p'-DDT	50-29-3	1000			
		dieldrin	60-57-1	1000			
		endrin	72-20-8	1000			
		heptachlor	76-44-8	1000			
		heptachlor epoxide isomer B	1024-57-3	1000			
		endrin aldehyde	7421-93-4	1000			
		endosulfan I	959-98-8	1000			
		endosulfan II	33213-65-9	1000			
		endosulfan sulfate	1031-07-8	1000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		methoxychlor	72-43-5	1000			
LK3-214406	Custom Method 508 Mix, 17-8	endrin	72-20-8	10	Dimethyl Formam	1 ml, Amber Glass	€ 79,50
		a-chlordane	5103-71-9	20			
		g-chlordane	5103-74-2	20			
		g-BHC	58-89-9	20			
		heptachlor epoxide isomer B	1024-57-3	20			
		heptachlor	76-44-8	40			
		simazine	122-34-9	70			
		atrazine	1912-24-9	100			
		alachlor	15972-60-8	200			
		butachlor	23184-66-9	100			
		metolachlor	51218-45-2	100			
		metribuzin	21087-64-9	100			
		propachlor	1918-16-7	100			
		dieldrin	60-57-1	100			
		hexachlorobenzene	118-74-1	100			
		hexachlorocyclopentadiene	77-47-4	100			
		methoxychlor	72-43-5	100			
LK3-023101	Custom Method 508 Solution, 19-12				Acetone	1 ml, Amber Glass	€ 229,50
		alachlor	15972-60-8	1000			
		atrazine	1912-24-9	1000			
		hexachlorobenzene	118-74-1	1000			
		hexachlorocyclopentadiene	77-47-4	1000			
		butachlor	23184-66-9	1000			
		methoxychlor	72-43-5	1000			
		metolachlor	51218-45-2	1000			
		propachlor	1918-16-7	1000			
		bis(2-ethylhexyl) adipate	103-23-1	1000			
		bis(2-ethylhexyl)phthalate	117-81-7	1000			
		chlordane, technical mixture	12789-03-6	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		endrin	72-20-8	1000			
		aldrin	309-00-2	1000			
		dieldrin	60-57-1	1000			
		heptachlor	76-44-8	1000			
		heptachlor epoxide isomer B	1024-57-3	1000			
		g-BHC	58-89-9	1000			
		benzo[a]pyrene	50-32-8	1000			
		simazine	122-34-9	1000			
LK3-023301	Custom Method 508 Solution, 19-12				Acetone	1 ml, Amber Glass	€ 229,50
		alachlor	15972-60-8	1000			
		atrazine	1912-24-9	1000			
		hexachlorobenzene	118-74-1	1000			
		hexachlorocyclopentadiene	77-47-4	1000			
		butachlor	23184-66-9	1000			
		methoxychlor	72-43-5	1000			
		metolachlor	51218-45-2	1000			
		propachlor	1918-16-7	1000			
		bis(2-ethylhexyl) adipate	103-23-1	1000			
		bis(2-ethylhexyl)phthalate	117-81-7	1000			
		endrin	72-20-8	1000			
		aldrin	309-00-2	1000			
		dieldrin	60-57-1	1000			
		heptachlor	76-44-8	1000			
		heptachlor epoxide isomer B	1024-57-3	1000			
		g-BHC	58-89-9	1000			
		benzo[a]pyrene	50-32-8	1000			
LK3-023201	Method 508-Chlorinated Pesticide Solution, 20-5				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 174,00
		aldrin	309-00-2	1000			
		a-BHC	319-84-6	1000			
		b-BHC	319-85-7	1000			

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
d-BHC	319-86-8	1000			
g-BHC	58-89-9	1000			
p,p'-DDD	72-54-8	1000			
p,p'-DDE	72-55-9	1000			
p,p'-DDT	50-29-3	1000			
dieldrin	60-57-1	1000			
endrin	72-20-8	1000			
heptachlor	76-44-8	1000			
heptachlor epoxide isomer B	1024-57-3	1000			
endrin aldehyde	7421-93-4	1000			
endosulfan I	959-98-8	1000			
endosulfan II	33213-65-9	1000			
endosulfan sulfate	1031-07-8	1000			
methoxychlor	72-43-5	1000			
alachlor	15972-60-8	1000			
atrazine	1912-24-9	1000			
chlordane (mixture of isomers)	57-74-9	1000			

LK3-023006	Custom Method 508 Solution, 20-5			Dimethyl Formam 1 ml, Amber Glass	€ 103,50
	alachlor	15972-60-8	200		
	aldrin	309-00-2	100		
	atrazine	1912-24-9	100		
	butachlor	23184-66-9	100		
	Dichloran	99-30-9	100		
	dieldrin	60-57-1	100		
	hexachlorobenzene	118-74-1	100		
	hexachlorocyclopentadiene	77-47-4	100		
	methoxychlor	72-43-5	100		
	metolachlor	51218-45-2	100		
	metribuzin	21087-64-9	100		
	propachlor	1918-16-7	100		
	simazine	122-34-9	70		

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		heptachlor	76-44-8	40			
		trifluralin	1582-09-8	20			
		heptachlor epoxide isomer B	1024-57-3	20			
		g-BHC	58-89-9	20			
		a-chlordane	5103-71-9	20			
		g-chlordane	5103-74-2	20			
		endrin	72-20-8	10			
LK3-013201-5pck	515 Performance Check Solution 3-1, Various Concentration (ppb)				M-t-BE	1 ml, Amber Glass	€ 126,00
		dinoseb-methyl ether	6099-79-2	4			
		methyl 3,5-dichlorobenzoate	2905-67-1	600			
		4-nitroanisole	100-17-4	1600			
LK3-200701	Method 525 Perchlorinate Mix, 3-24				Acetone	1 ml, Amber Glass	€ 36,00
		hexachlorobenzene	118-74-1	100			
		hexachlorocyclopentadiene	77-47-4	100			
		pentachlorophenol	87-86-5	400			
LK3-104305	Method 525 Surrogate Solution, 10 mg/L, 4-1				Cyclohexane	1 ml, Amber Glass	€ 31,50
LK3-104302	Method 525 Surrogate Solution, 500 mg/L, 4-1				EthylAcetate	1 ml, Amber Glass	€ 31,50
LK3-104301	Method 525 Surrogate Solution, 2,000 mg/L, 4-1				EthylAcetate	1 ml, Amber Glass	€ 72,00
		acenaphthene-d10	15067-26-2				
		phenanthrene-d10	1517-22-2				
		chrysene-d12	1719-03-5				
		perylene-d12	1520-96-3				
LK3-023506	Custom 525 Mix, 24-1	2-chlorobiphenyl (BZ# 1)	2051-60-7	20	Acetone	1 ml, Amber Glass	€ 72,00
		hexachlorobutadiene	87-68-3	20			
		aldrin	309-00-2	20			
		2,2',3',4,6-pentachlorobiphenyl (BZ# 98)	60233-25-2	20			
		butachlor	23184-66-9	20			
		2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 200)	40186-71-8	20			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		benzo[a]pyrene	50-32-8	20			
		dieldrin	60-57-1	26			
		fluorene	86-73-7	40			
		propachlor	1918-16-7	40			
		phenanthrene	85-01-8	40			
		pyrene	129-00-0	40			
		trans-nonachlor	39765-80-5	40			
		indeno[1,2,3-cd]pyrene	193-39-5	40			
		metribuzin	21087-64-9	60			
		dioctyl phthalate	117-81-7	100			
		metolachlor	51218-45-2	200			
		dimethyl phthalate	131-11-3	320			
		bis(2-ethylhexyl) adipate	103-23-1	320			
		diethyl phthalate	84-66-2	400			
		di-n-butyl phthalate	84-74-2	400			
		butyl benzyl phthalate	85-68-7	400			
		molinate	2212-67-1	400			
		benthiocarb	28249-77-6	400			
LK3-023606	Custom 525 Mix, 24-4	2-chlorobiphenyl (BZ# 1)	2051-60-7	20	Acetone	1 ml, Amber Glass	€ 72,00
		2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	20			
		2,2',3',4,6-pentachlorobiphenyl (BZ# 98)	60233-25-2	20			
		2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 201)	40186-71-8	20			
		dieldrin	60-57-1	26			
		metolachlor	51218-45-2	200			
		benzo[a]pyrene	50-32-8	20			
		aldrin	309-00-2	20			
		butachlor	23184-66-9	20			
		fluorene	86-73-7	40			
		indeno[1,2,3-cd]pyrene	193-39-5	40			
		phenanthrene	85-01-8	40			
		pyrene	129-00-0	40			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		trans-nonachlor	39765-80-5	40			
		propachlor	1918-16-7	40			
		metribuzin	21087-64-9	60			
		dimethyl phthalate	131-11-3	320			
		bis(2-ethylhexyl) adipate	103-23-1	320			
		bis(2-ethylhexyl)phthalate	117-81-7	400			
		butyl benzyl phthalate	85-68-7	400			
		diethyl phthalate	84-66-2	400			
		di-n-butyl phthalate	84-74-2	400			
		benthocarb	28249-77-6	400			
		molinate	2212-67-1	400			
LK3-212402	525 Custom Pesticide Solution, 44-3				Acetone	1 ml, Amber Glass	€ 315,00
		chlorobenzilate	510-15-6	1000			
		propyzamide (pronamide)	23950-58-5	1000			
		prometryn	7287-19-6	1000			
		propazine	139-40-2	1000			
		trifluralin	1582-09-8	1000			
		atrazine	1912-24-9	1000			
		ethoprophos (prophos)	13194-48-4	1000			
		chlorpyrifos	2921-88-2	1000			
		alachlor	15972-60-8	1000			
		chlorthal-dimethyl (dacthal)	1861-32-1	1000			
		propachlor	1918-16-7	1000			
		dichlorvos	62-73-7	1000			
		tetrachlorvinphos (Rabon)	961-11-5	1000			
		acetochlor	34256-82-1	1000			
		Tebuthiuron	34014-18-1	1000			
		simetryn	1014-70-6	1000			
		metolachlor	51218-45-2	1000			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000			
		ametryne	834-12-8	1000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		terbutryne	886-50-0	1000			
		cis-permethrin	54774-45-7	500			
		trans-permethrin	51877-74-8	500			
		butachlor	23184-66-9	1000			
		molinate	2212-67-1	1000			
		terbacil	5902-51-2	1000			
		Cycloate	1134-23-2	1000			
		Fenarimol	60168-88-9	1000			
		butylate	2008-41-5	1000			
		vernolate	1929-77-7	1000			
		pebulate	1114-71-2	1000			
		etridiazole	2593-15-9	1000			
		chloroneb	2675-77-6	1000			
		chlorpropham	101-21-3	1000			
		bromacil	314-40-9	1000			
		triadimefon	43121-43-3	1000			
		MGK-264 - isomer a	113-48-4	500			
		MGK-264 - isomer b	113-48-4	500			
		diphenamid	957-51-7	1000			
		napropamide	15299-99-7	1000			
		beam™	41814-78-2	1000			
		norflurazon	27314-13-2	1000			
		hexazinon	51235-04-2	1000			
		fluridone	59756-60-4	1000			
		simazine	122-34-9	1000			
LK3-212402	525 Custom Pesticide Solution, 44-3				Acetone	1 ml, Amber Glass	€ 315,00
		chlorobenzilate	510-15-6	1000			
		propyzamide (pronamide)	23950-58-5	1000			
		prometryn	7287-19-6	1000			
		propazine	139-40-2	1000			
		trifluralin	1582-09-8	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
atrazine	1912-24-9	1000			
ethoprophos (prophos)	13194-48-4	1000			
chlorpyrifos	2921-88-2	1000			
alachlor	15972-60-8	1000			
chlorthal-dimethyl (dacthal)	1861-32-1	1000			
propachlor	1918-16-7	1000			
dichlorvos	62-73-7	1000			
tetrachlorvinphos (Rabon)	961-11-5	1000			
acetochlor	34256-82-1	1000			
Tebuthiuron	34014-18-1	1000			
simetryn	1014-70-6	1000			
metolachlor	51218-45-2	1000			
EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000			
ametryne	834-12-8	1000			
terbutryne	886-50-0	1000			
cis-permethrin	54774-45-7	500			
trans-permethrin	51877-74-8	500			
butachlor	23184-66-9	1000			
molinate	2212-67-1	1000			
terbacil	5902-51-2	1000			
Cycloate	1134-23-2	1000			
Fenarimol	60168-88-9	1000			
butylate	2008-41-5	1000			
vernolate	1929-77-7	1000			
pebulate	1114-71-2	1000			
etridiazole	2593-15-9	1000			
chloroneb	2675-77-6	1000			
chlorpropham	101-21-3	1000			
bromacil	314-40-9	1000			
triadimefon	43121-43-3	1000			
MGK-264 - isomer a	113-48-4	500			
MGK-264 - isomer b	113-48-4	500			

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
diphenamid	957-51-7	1000			
napropamide	15299-99-7	1000			
beam™	41814-78-2	1000			
norflurazon	27314-13-2	1000			
hexazinon	51235-04-2	1000			
fluridone	59756-60-4	1000			
simazine	122-34-9	1000			

LK3-212201	525 Custom Pesticide Solution, 45-3		Acetone	1 ml, Amber Glass	€ 315,00
	chlorobenzilate	510-15-6	1000		
	propyzamide (pronamide)	23950-58-5	1000		
	prometryn	7287-19-6	1000		
	propazine	139-40-2	1000		
	trifluralin	1582-09-8	1000		
	atrazine	1912-24-9	1000		
	mevinphos (Phosdrin (TM))	7786-34-7	1000		
	ethoprophos (prophos)	13194-48-4	1000		
	chlorpyrifos	2921-88-2	1000		
	alachlor	15972-60-8	1000		
	chlorthal-dimethyl (dacthal)	1861-32-1	1000		
	propachlor	1918-16-7	1000		
	dichlorvos	62-73-7	1000		
	tetrachlorvinphos (Rabon)	961-11-5	1000		
	acetochlor	34256-82-1	1000		
	Tebuthiuron	34014-18-1	1000		
	simetryn	1014-70-6	1000		
	metolachlor	51218-45-2	1000		
	EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000		
	ametryne	834-12-8	1000		
	terbutryne	886-50-0	1000		
	cis-permethrin	54774-45-7	1000		
	trans-permethrin	51877-74-8	1000		

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		butachlor	23184-66-9	1000			
		molinate	2212-67-1	1000			
		terbacil	5902-51-2	1000			
		Cycloate	1134-23-2	1000			
		Fenarimol	60168-88-9	1000			
		butylate	2008-41-5	1000			
		vernolate	1929-77-7	1000			
		pebulate	1114-71-2	1000			
		etridiazole	2593-15-9	1000			
		chloroneb	2675-77-6	1000			
		chlorpropham	101-21-3	1000			
		bromacil	314-40-9	1000			
		triadimefon	43121-43-3	1000			
		MGK-264 - isomer a	113-48-4	1000			
		MGK-264 - isomer b	113-48-4	1000			
		diphenamid	957-51-7	1000			
		napropamide	15299-99-7	1000			
		beam™	41814-78-2	1000			
		norflurazon	27314-13-2	1000			
		hexazinon	51235-04-2	1000			
		fluridone	59756-60-4	1000			
		simazine	122-34-9	1000			
LK3-212101	525 Custom Solution, 56-1	acenaphthene	83-32-9	1000	Acetone	1 ml, Amber Glass	€ 433,50
		acenaphthylene	208-96-8	1000			
		anthracene	120-12-7	1000			
		benzo[a]anthracene	56-55-3	1000			
		benzo[b]fluoranthene	205-99-2	1000			
		benzo[k]fluoranthene	207-08-9	1000			
		benzo[a]pyrene	50-32-8	1000			
		chrysene	218-01-9	1000			
		dibenz[a,h]anthracene	53-70-3	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
fluoranthene	206-44-0	1000			
fluorene	86-73-7	1000			
indeno[1,2,3-cd]pyrene	193-39-5	1000			
naphthalene	91-20-3	1000			
phenanthrene	85-01-8	1000			
pyrene	129-00-0	1000			
bis(2-ethylhexyl)phthalate	117-81-7	1000			
butyl benzyl phthalate	85-68-7	1000			
diethyl phthalate	84-66-2	1000			
dimethyl phthalate	131-11-3	1000			
di-n-butyl phthalate	84-74-2	1000			
hexachlorobenzene	118-74-1	1000			
hexachlorocyclopentadiene	77-47-4	1000			
2-methylnaphthalene	91-57-6	1000			
2,4-dinitrotoluene	121-14-2	1000			
2,6-dinitrotoluene	606-20-2	1000			
isophorone	78-59-1	1000			
1-methylnaphthalene	90-12-0	1000			
bis(2-ethylhexyl) adipate	103-23-1	1000			
p,p'-DDT	50-29-3	1000			
p,p'-DDD	72-54-8	1000			
p,p'-DDE	72-55-9	1000			
a-BHC	319-84-6	1000			
g-BHC	58-89-9	1000			
dieldrin	60-57-1	1000			
methoxychlor	72-43-5	1000			
endosulfan II	33213-65-9	1000			
endosulfan sulfate	1031-07-8	1000			
endrin	72-20-8	1000			
endrin aldehyde	7421-93-4	1000			
aldrin	309-00-2	1000			
b-BHC	319-85-7	1000			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		d-BHC	319-86-8	1000			
		heptachlor	76-44-8	1000			
		heptachlor epoxide isomer B	1024-57-3	1000			
		a-chlordane	5103-71-9	1000			
		g-chlordane	5103-74-2	1000			
		endosulfan I	959-98-8	1000			
		trans-nonachlor	39765-80-5	1000			
		2-chlorobiphenyl (BZ# 1)	2051-60-7	1000			
		2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	1000			
		2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	1000			
		2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	1000			
		2,2',3',4,6-pentachlorobiphenyl (BZ# 98)	60233-25-2	1000			
		2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	1000			
		2,2',3,3',4,4',6-heptachlorobiphenyl (BZ# 171)	52663-71-5	1000			
		2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 201)	40186-71-8	1000			
LK3-211101	Custom 525 Spike Mix 94-1 (No PCB Congeners)				Acetone	2 ml, Amber Glass	€ 709,50
		isophorone	78-59-1	1000			
		naphthalene	91-20-3	1000			
		dichlorvos	62-73-7	1000			
		1-methylnaphthalene	90-12-0	1000			
		2-methylnaphthalene	91-57-6	1000			
		hexachlorocyclopentadiene	77-47-4	1000			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000			
		mevinphos (Phosdrin (TM))	7786-34-7	1000			
		butylate	2008-41-5	1000			
		vernolate	1929-77-7	1000			
		dimethyl phthalate	131-11-3	1000			
		pebulate	1114-71-2	1000			
		etridiazole	2593-15-9	1000			
		2,6-dinitrotoluene	606-20-2	1000			
		acenaphthylene	208-96-8	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
acenaphthene	83-32-9	1000			
chloroneb	2675-77-6	1000			
Tebuthiuron	34014-18-1	1000			
2,4-dinitrotoluene	121-14-2	1000			
molinate	2212-67-1	1000			
diethyl phthalate	84-66-2	1000			
fluorene	86-73-7	1000			
propachlor	1918-16-7	1000			
ethoprophos (prophos)	13194-48-4	1000			
Cycloate	1134-23-2	1000			
chlorpropham	101-21-3	1000			
trifluralin	1582-09-8	1000			
a-BHC	319-84-6	1000			
hexachlorobenzene	118-74-1	1000			
simazine	122-34-9	1000			
atrazine	1912-24-9	1000			
propazine	139-40-2	1000			
b-BHC	319-85-7	1000			
g-BHC	58-89-9	1000			
propyzamide (pronamide)	23950-58-5	1000			
phenanthrene	85-01-8	1000			
anthracene	120-12-7	1000			
methyl paraoxon	950-35-6	1000			
d-BHC	319-86-8	1000			
terbacil	5902-51-2	1000			
acetochlor	34256-82-1	1000			
simetryn	1014-70-6	1000			
alachlor	15972-60-8	1000			
ametryne	834-12-8	1000			
heptachlor	76-44-8	1000			
prometryn	7287-19-6	1000			
terbutryne	886-50-0	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
di-n-butyl phthalate	84-74-2	1000			
bromacil	314-40-9	1000			
aldrin	309-00-2	1000			
metolachlor	51218-45-2	1000			
chlorpyrifos	2921-88-2	1000			
triadimefon	43121-43-3	1000			
chlorthal (dacthal)	1861-32-1	1000			
MGK-264 - isomer a	113-48-4	1000			
diphenamid	957-51-7	1000			
MGK-264 - isomer b	113-48-4	1000			
heptachlor epoxide isomer B	1024-57-3	1000			
fluoranthene	206-44-0	1000			
g-chlordane	5103-74-2	1000			
pyrene	129-00-0	1000			
tetrachlorvinphos (Rabon)	961-11-5	1000			
endosulfan I	959-98-8	1000			
butachlor	23184-66-9	1000			
a-chlordane	5103-71-9	1000			
trans-nonachlor	39765-80-5	1000			
napropamide	15299-99-7	1000			
Beam	41814-78-2	1000			
p,p'-DDE	72-55-9	1000			
dieldrin	60-57-1	1000			
endrin	72-20-8	1000			
chlorobenzilate	510-15-6	1000			
p,p'-DDD	72-54-8	1000			
endrin aldehyde	7421-93-4	1000			
endosulfan sulfate	1031-07-8	1000			
endosulfan II	33213-65-9	1000			
butyl benzyl phthalate	85-68-7	1000			
norflurazon	27314-13-2	1000			
p,p'-DDT	50-29-3	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		hexazinon	51235-04-2	1000			
		bis(2-ethylhexyl) adipate	103-23-1	1000			
		benzo[a]anthracene	56-55-3	1000			
		chrysene	218-01-9	1000			
		methoxychlor	72-43-5	1000			
		dioctyl phthalate	117-81-7	1000			
		Fenarimol	60168-88-9	1000			
		cis-permethrin	54774-45-7	1000			
		trans-permethrin	51877-74-8	1000			
		benzo[b]fluoranthene	205-99-2	1000			
		benzo[k]fluoranthene	207-08-9	1000			
		benzo[a]pyrene	50-32-8	1000			
		fluridone	59756-60-4	1000			
		indeno[1,2,3-cd]pyrene	193-39-5	1000			
		dibenz[a,h]anthracene	53-70-3	1000			
LK3-104401	Custom 525 Spike Mix 97-1 (No PCB Congeners)				Acetone	2 ml, Amber Glass	€ 709,50
		isophorone	78-59-1	1000			
		naphthalene	91-20-3	1000			
		dichlorvos	62-73-7	1000			
		1-methylnaphthalene	90-12-0	1000			
		2-methylnaphthalene	91-57-6	1000			
		hexachlorocyclopentadiene	77-47-4	1000			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000			
		mevinphos (Phosdrin (TM))	7786-34-7	1000			
		butylate	2008-41-5	1000			
		vernolate	1929-77-7	1000			
		dimethyl phthalate	131-11-3	1000			
		pebulate	1114-71-2	1000			
		etridiazole	2593-15-9	1000			
		2,6-dinitrotoluene	606-20-2	1000			
		acenaphthylene	208-96-8	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
acenaphthene	83-32-9	1000			
chlorneb	2675-77-6	1000			
Tebuthiuron	34014-18-1	1000			
2,4-dinitrotoluene	121-14-2	1000			
Molinate	2212-67-1	1000			
diethyl phthalate	84-66-2	1000			
fluorene	86-73-7	1000			
Propachlor	1918-16-7	1000			
ethoprophos (prophos)	13194-48-4	1000			
Cycloate	1134-23-2	1000			
chlorpropham	101-21-3	1000			
trifluralin	1582-09-8	1000			
a-BHC	319-84-6	1000			
hexachlorobenzene	118-74-1	1000			
Simazine	122-34-9	1000			
atrazine	1912-24-9	1000			
propazine	139-40-2	1000			
b-BHC	319-85-7	1000			
g-BHC	58-89-9	1000			
propyzamide (pronamide)	23950-58-5	1000			
phenanthrene	85-01-8	1000			
anthracene	120-12-7	1000			
methyl paraoxon	950-35-6	1000			
d-BHC	319-86-8	1000			
Terbacil	5902-51-2	1000			
Chlorothalonil	1897-45-6	1000			
metribuzin	21087-64-9	1000			
acetochlor	34256-82-1	1000			
Simetryn	1014-70-6	1000			
Alachlor	15972-60-8	1000			
Ametryne	834-12-8	1000			
heptachlor	76-44-8	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
prometryne	7287-19-6	1000			
Terbutryne	886-50-0	1000			
di-n-butyl phthalate	84-74-2	1000			
bromacil	314-40-9	1000			
aldrin	309-00-2	1000			
metolachlor	51218-45-2	1000			
chlorpyrifos	2921-88-2	1000			
triademefon	43121-43-3	1000			
Chlorthal (Dacthal)	1861-32-1	1000			
MGK-264 - isomer a	113-48-4	1000			
diphenamid	957-51-7	1000			
MGK-264 - isomer b	113-48-4	1000			
heptachlor epoxide	1024-57-3	1000			
fluoranthene	206-44-0	1000			
g-chlordane	5103-74-2	1000			
pyrene	129-00-0	1000			
tetrachlorvinphos	961-11-5	1000			
endosulfan I	959-98-8	1000			
Butachlor	23184-66-9	1000			
a-chlordane	5103-71-9	1000			
trans-nonachlor	39765-80-5	1000			
napropamide	15299-99-7	1000			
tricyclazole	41814-78-2	1000			
p,p'-DDE	72-55-9	1000			
dieldrin	60-57-1	1000			
endrin	72-20-8	1000			
chlorobenzilate	510-15-6	1000			
p,p'-DDD	72-54-8	1000			
endrin aldehyde	7421-93-4	1000			
endosulfan sulfate	1031-07-8	1000			
endosulfan II	33213-65-9	1000			
butyl benzyl phthalate	85-68-7	1000			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		norflurazon	27314-13-2	1000			
		p,p'-DDT	50-29-3	1000			
		hexazinon	51235-04-2	1000			
		bis(2-ethylhexyl) adipate	103-23-1	1000			
		benzo[a]anthracene	56-55-3	1000			
		chrysene	218-01-9	1000			
		methoxychlor	72-43-5	1000			
		bis(2-ethylhexyl)phthalate	117-81-7	1000			
		Fenarimol	60168-88-9	1000			
		cis-permethrin	54774-45-7	1000			
		trans-permethrin	51877-74-8	1000			
		benzo[b]fluoranthene	205-99-2	1000			
		benzo[k]fluoranthene	207-08-9	1000			
		benzo[a]pyrene	50-32-8	1000			
		benzo[ghi]perylene	191-24-2	1000			
		fluridone	59756-60-4	1000			
		indeno[1,2,3-cd]pyrene	193-39-5	1000			
		dibenz[a,h]anthracene	53-70-3	1000			
LK3-211405	Custom 525 Spike Mix 100-1, Various Concentrations				Ethyl Acetate	10 ml BR, Amber	€ 355,50
LK3-211405-5pck	Custom 525 Spike Mix 100-1, Various Concentrations					5 x 2 ml, Amber	€ 355,50
		acenaphthene	83-32-9	2			
		acenaphthylene	208-96-8	2			
		anthracene	120-12-7	2			
		benzo[a]anthracene	56-55-3	2			
		benzo[b]fluoranthene	205-99-2	2			
		benzo[k]fluoranthene	207-08-9	2			
		benzo[ghi]perylene	191-24-2	2			
		benzo[a]pyrene	50-32-8	2			
		chrysene	218-01-9	2			
		fluoranthene	206-44-0	2			
		fluorene	86-73-7	2			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
indeno[1,2,3-cd]pyrene	193-39-5	2			
naphthalene	91-20-3	2			
phenanthrene	85-01-8	2			
pyrene	129-00-0	2			
dibenz[a,h]anthracene	53-70-3	2			
isophorone	78-59-1	2			
dichlorvos	62-73-7	2			
2-methylnaphthalene	91-57-6	2			
butylate	2008-41-5	2			
vernolate	1929-77-7	2			
pebulate	1114-71-2	2			
etridiazole	2593-15-9	2			
propachlor	1918-16-7	2			
Cycloate	1134-23-2	2			
chlorpropham	101-21-3	2			
EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	2			
molinate	2212-67-1	2			
trifluralin	1582-09-8	2			
hexachlorobenzene	118-74-1	2			
propyzamide (pronamide)	23950-58-5	2			
a-BHC	319-84-6	2			
b-BHC	319-85-7	2			
d-BHC	319-86-8	2			
g-BHC	58-89-9	2			
aldrin	309-00-2	2			
heptachlor	76-44-8	2			
heptachlor epoxide isomer B	1024-57-3	2			
atrazine	1912-24-9	2			
propazine	139-40-2	2			
metribuzin	21087-64-9	2			
alachlor	15972-60-8	2			
chlorothalonil	1897-45-6	2			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
fenamiphos	22224-92-6	2			
metolachlor	51218-45-2	2			
chlorthal-dimethyl (dacthal)	1861-32-1	2			
diphenamid	957-51-7	2			
g-chlordane	5103-74-2	2			
MGK-264 - isomer a	113-48-4	2			
MGK-264 - isomer b	113-48-4	2			
p,p'-DDE	72-55-9	2			
2-chlorobiphenyl (BZ# 1)	2051-60-7	5			
2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	2			
2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	2			
2,2',3',4,6-pentachlorobiphenyl (BZ# 98)	60233-25-2	2			
2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	2			
2,2',3,3',4,4',6-heptachlorobiphenyl (BZ# 171)	52663-71-5	2			
2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 201)	40186-71-8	2			
2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	5			
mevinphos (Phosdrin (TM))	7786-34-7	5			
ethoprophos (prophos)	13194-48-4	5			
dimethyl-p-nitrophenylphosphate	23184-66-9	5			
terbacil	5902-51-2	5			
napropamide	15299-99-7	5			
norflurazon	27314-13-2	5			
hexazinon	51235-04-2	5			
acetochlor	34256-82-1	5			
bromacil	314-40-9	5			
chlorpyrifos	2921-88-2	5			
tetrachlorvinphos (Rabon)	961-11-5	5			
endosulfan I	959-98-8	5			
endosulfan II	33213-65-9	5			
endosulfan sulfate	1031-07-8	5			
chloroneb	2675-77-6	5			
simazine	122-34-9	5			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		butachlor	23184-66-9	5			
		a-chlordane	5103-71-9	5			
		trans-nonachlor	39765-80-5	5			
		dieldrin	60-57-1	5			
		endrin	72-20-8	5			
		p,p'-DDD	72-54-8	5			
		p,p'-DDT	50-29-3	5			
		methoxychlor	72-43-5	5			
		cis-permethrin	54774-45-7	5			
		trans-permethrin	51877-74-8	5			
		2,6-dinitrotoluene	606-20-2	7,5			
		2,4-dinitrotoluene	121-14-2	7,5			
		butyl benzyl phthalate	85-68-7	7,5			
		fluridone	59756-60-4	7,5			
		triadimefon	43121-43-3	10			
		endrin aldehyde	7421-93-4	10			
		dimethyl phthalate	131-11-3	15			
		diethyl phthalate	84-66-2	15			
		di-n-butyl phthalate	84-74-2	15			
		beam™	41814-78-2	15			
		bis(2-ethylhexyl) adipate	103-23-1	15			
		hexachlorocyclopentadiene	77-47-4	20			
		bis(2-ethylhexyl)phthalate	117-81-7	20			
		Fenarimol	60168-88-9	50			
		chlorobenzilate	510-15-6	67			
LK3-104201	Custom 525 Spike Mix 105-1	isophorone	78-59-1	1000	Isooctane	2 ml, Amber Glass	€ 787,50
		naphthalene	91-20-3	1000			
		dichlorvos	62-73-7	1000			
		1-methylnaphthalene	90-12-0	1000			
		2-methylnaphthalene	91-57-6	1000			
		hexachlorocyclopentadiene	77-47-4	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000			
mevinphos (Phosdrin (TM))	7786-34-7	1000			
butylate	2008-41-5	1000			
vernolate	1929-77-7	1000			
dimethyl phthalate	131-11-3	1000			
pebulate	1114-71-2	1000			
etridiazole	2593-15-9	1000			
2,6-dinitrotoluene	606-20-2	1000			
acenaphthylene	208-96-8	1000			
acenaphthene	83-32-9	1000			
chlorneb	2675-77-6	1000			
Tebuthiuron	34014-18-1	1000			
2,4-dinitrotoluene	121-14-2	1000			
Molinate	2212-67-1	1000			
diethyl phthalate	84-66-2	1000			
fluorene	86-73-7	1000			
Propachlor	1918-16-7	1000			
ethoprophos (prophos)	13194-48-4	1000			
Cycloate	1134-23-2	1000			
chlorpropham	101-21-3	1000			
trifluralin	1582-09-8	1000			
a-BHC	319-84-6	1000			
hexachlorobenzene	118-74-1	1000			
Simazine	122-34-9	1000			
atrazine	1912-24-9	1000			
propazine	139-40-2	1000			
b-BHC	319-85-7	1000			
g-BHC	58-89-9	1000			
propyzamide (pronamide)	23950-58-5	1000			
phenanthrene	85-01-8	1000			
anthracene	120-12-7	1000			
methyl paraoxon	950-35-6	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
d-BHC	319-86-8	1000			
Terbacil	5902-51-2	1000			
Chlorothalonil	1897-45-6	1000			
metribuzin	21087-64-9	1000			
acetochlor	34256-82-1	1000			
Simetryn	1014-70-6	1000			
Alachlor	15972-60-8	1000			
Ametryne	834-12-8	1000			
heptachlor	76-44-8	1000			
prometryne	15972-60-8	1000			
Terbutryne	886-50-0	1000			
di-n-butyl phthalate	84-74-2	1000			
bromacil	314-40-9	1000			
aldrin	309-00-2	1000			
metolachlor	51218-45-2	1000			
chlorpyrifos	2921-88-2	1000			
triademefon	43121-43-3	1000			
Chlorthal (Dacthal)	1861-32-1	1000			
MGK-264 - isomer a	113-48-4	1000			
diphenamid	957-51-7	1000			
MGK-264 - isomer b	113-48-4	1000			
heptachlor epoxide	1024-57-3	1000			
fluoranthene	206-44-0	1000			
g-chlordane	5103-74-2	1000			
pyrene	129-00-0	1000			
tetrachlorvinphos	961-11-5	1000			
endosulfan I	959-98-8	1000			
Butachlor	23184-66-9	1000			
a-chlordane	5103-71-9	1000			
trans-nonachlor	39765-80-5	1000			
napropamide	15299-99-7	1000			
tricyclazole	41814-78-2	1000			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
p,p'-DDE	72-55-9	1000			
dieldrin	60-57-1	1000			
endrin	72-20-8	1000			
chlorobenzilate	510-15-6	1000			
p,p'-DDD	72-54-8	1000			
endrin aldehyde	7421-93-4	1000			
endosulfan sulfate	1031-07-8	1000			
endosulfan II	33213-65-9	1000			
butyl benzyl phthalate	85-68-7	1000			
norflurazon	27314-13-2	1000			
p,p'-DDT	50-29-3	1000			
hexazinon	51235-04-2	1000			
bis(2-ethylhexyl) adipate	103-23-1	1000			
benzo[a]anthracene	56-55-3	1000			
chrysene	218-01-9	1000			
methoxychlor	72-43-5	1000			
bis(2-ethylhexyl)phthalate	117-81-7	1000			
Fenarimol	60168-88-9	1000			
cis-permethrin	54774-45-7	1000			
trans-permethrin	51877-74-8	1000			
benzo[b]fluoranthene	205-99-2	1000			
benzo[k]fluoranthene	207-08-9	1000			
benzo[a]pyrene	50-32-8	1000			
benzo[ghi]perylene	191-24-2	1000			
fluridone	59756-60-4	1000			
indeno[1,2,3-cd]pyrene	193-39-5	1000			
dibenz[a,h]anthracene	53-70-3	1000			
2-chlorobiphenyl (BZ# 1)	2051-60-7	1000			
2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	1000			
2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	1000			
2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	1000			
2,2',3',4,6-pentachlorobiphenyl (BZ# 98)	60233-25-2	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	1000			
		2,2',3,3',4,4',6-heptachlorobiphenyl (BZ# 171)	52663-71-5	1000			
		2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 200)	40186-71-8	1000			
LK3-105401	EPA Method 525.1-Capillary Solution, 4-1				Acetone	1 ml, Amber Glass	€ 72,00
		2-chlorobiphenyl (BZ# 1)	2051-60-7	1000			
		2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	1000			
		2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	1000			
		2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	1000			
LK3-105501	EPA Method 525.1-Capillary Solution, 4-2				Acetone	0.5 oz BR, Amber	€ 72,00
		2,2',3,3',4,4',6-heptachlorobiphenyl (BZ# 171)	52663-71-5	1000			
		2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	1000			
		2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 200)	40186-71-8	1000			
		2,2',3,4,6-pentachlorobiphenyl (BZ# 88)	55215-17-3	1000			
LK3-114301	Method 525.2 Nitrogen/Phosphorus Pesticide Mx, 6-13				Acetone	1 ml, Amber Glass	€ 94,50
		carboxine	5234-68-4	1000			
		diazinon	333-41-5	1000			
		disulfoton	298-04-4	1000			
		fenamiphos	22224-92-6	1000			
		tributylphosphoro-trithioite (Merphos)	150-50-5	1000			
		Terbufos	13071-79-9	1000			
LK3-104601	EPA Method 525.1-Capillary Solution, 8-1				Acetone	1 ml, Amber Glass	€ 72,00
		2-chlorobiphenyl (BZ# 1)	2051-60-7	500			
		2,3-dichlorobiphenyl (BZ# 5)	16605-91-7	500			
		2,2',4,4'-tetrachlorobiphenyl (BZ# 47)	2437-79-8	500			
		2,4,5-trichlorobiphenyl (BZ# 29)	15862-07-4	500			
		2,2',3,3',4,4',6-heptachlorobiphenyl (BZ# 171)	52663-71-5	500			
		2,2',4,4',5,6'-hexachlorobiphenyl (BZ# 154)	60145-22-4	500			
		2,2',3,3',4,5',6,6'-octachlorobiphenyl (BZ# 200)	40186-71-8	500			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		2,2',3,4,6-pentachlorobiphenyl (BZ# 88)	55215-17-3	500			
LK3-114601	EPA Method 525.2-Capillary Solution, 8-9				M-t-BE	1 ml, Amber Glass	€ 48,00
LK3-114601-5pck	EPA Method 525.2-Capillary Solution, 8-9					5 x 1 ml, Amber	€ 189,00
		napropamide	15299-99-7	1000			
		dichlorvos	62-73-7	1000			
		Fenarimol	60168-88-9	1000			
		pebulate	1114-71-2	1000			
		Simetryn	1014-70-6	1000			
		fluridone	59756-60-4	1000			
		Tebuthiuron	34014-18-1	1000			
		Terbacil	5902-51-2	1000			
LK3-205705	Method 525.2 Organochlorine Pesticide Mix, 8-21				Acetone	1 ml, Amber Glass	€ 55,50
		chlorobenzilate	510-15-6	500			
		chloroneb	2675-77-6	500			
		chlorothalonil	1897-45-6	500			
		chlorthal-dimethyl (dacthal)	1861-32-1	500			
		heptachlor epoxide isomer A	28044-83-9	500			
		trans-nonachlor	39765-80-5	500			
		cis-permethrin	54774-45-7	500			
		trans-permethrin	51877-74-8	500			
LK3-104801	EPA Method 525.2-Capillary Solution, 9-1				M-t-BE	1 ml, Amber Glass	€ 63,00
		atrazine	1912-24-9	1000			
		4-nitroaniline	100-01-6	1000			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000			
		ethoprophos (prophos)	13194-48-4	1000			
		mevinphos (Phosdrin (TM))	7786-34-7	1000			
		prometryn	7287-19-6	1000			
		propazine	139-40-2	1000			
		Terbutryne	886-50-0	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		triadimefon	43121-43-3	1000			
LK3-114801	EPA Method 525.2-Capillary So 1ml				M-t-BE	1 ml, Amber Glass	€ 55,50
LK3-114801-5pck	EPA Method 525.2-Capillary So 5 X 1ml					5 x1 ml, Amber	€ 213,00
		atrazine	1912-24-9	1000			
		diphenamid	957-51-7	1000			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	1000			
		ethoprophos (prophos)	13194-48-4	1000			
		mevinphos (Phosdrin (TM))	7786-34-7	1000			
		prometryn	7287-19-6	1000			
		propazine	139-40-2	1000			
		Terbutryne	886-50-0	1000			
		triadimefon	43121-43-3	1000			
LK3-114401	EPA Method 525.2 Capillary So 1ml				M-t-BE	1 ml, Amber Glass	€ 103,50
LK3-114401-5pck	EPA Method 525.2 Capillary So 5 X 1ml					5 x1 ml, Amber	€ 220,50
		butachlor	23184-66-9	1000			
		carboxine	5234-68-4	1000			
		diazinon	333-41-5	1000			
		metolachlor	51218-45-2	1000			
		metribuzin	21087-64-9	1000			
		MGK-264 Mix of Isomers	113-48-4	1000			
		norflurazon	27314-13-2	1000			
		Terbufos	13071-79-9	1000			
		vernolate	1929-77-7	1000			
LK3-104701	EPA Method 525.2-Capillary Solution, 10-1					1 ml, Amber Glass	€ 55,50
LK3-104701-5pck	EPA Method 525.2-Capillary Solution, 10-1					5 x 1 ml, Amber	€ 244,50
		alachlor	15972-60-8	1000	M-t-BE		
		Atraton	1610-17-9	1000			
		bromacil	314-40-9	1000			
		butylate	2008-41-5	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		chlorpropham	101-21-3	1000			
		hexazinon	51235-04-2	1000			
		molinate	2212-67-1	1000			
		propyzamide (pronamide)	23950-58-5	1000			
		tetrachlorvinphos (Rabon)	961-11-5	1000			
		tricyclazole	41814-78-2	1000			
LK3-105101	EPA Method 525.2-Capillary Solution, 11-1				M-t-BE	1 ml, Amber Glass	€ 40,50
		bis(2-ethylhexyl) adipate	103-23-1	500			
		bis(2-ethylhexyl)phthalate	117-81-7	500			
		butyl benzyl phthalate	85-68-7	500			
		di-n-butyl phthalate	84-74-2	500			
		diethyl phthalate	84-66-2	500			
		dimethyl phthalate	131-11-3	500			
		2,4-dinitrotoluene	121-14-2	500			
		2,6-dinitrotoluene	606-20-2	500			
		hexachlorocyclopentadiene	77-47-4	500			
		isophorone	78-59-1	500			
		pentachlorophenol	87-86-5	2000			
LK3-105002	EPA Method 525.2-Capillary Solution, 100 mg/L , 13-1				Acetone	1 ml, Amber Glass	€ 48,00
LK3-105003	EPA Method 525.2-Capillary Solution, 500 mg/L , 13-1				Acetone	1 ml, Amber Glass	€ 55,50
LK3-105001	EPA Method 525.2-Capillary Solution, 500 mg/L , 13-1				M-t-BE	1 ml, Amber Glass	€ 55,50
		acenaphthylene	208-96-8				
		anthracene	120-12-7				
		benzo[a]anthracene	56-55-3				
		benzo[b]fluoranthene	205-99-2				
		benzo[k]fluoranthene	207-08-9				
		benzo[a]pyrene	50-32-8				
		benzo[ghi]perylene	191-24-2				
		chrysene	218-01-9				
		dibenz[a,h]anthracene	53-70-3				

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		fluorene	86-73-7				
		indeno[1,2,3-cd]pyrene	193-39-5				
		phenanthrene	85-01-8				
		pyrene	129-00-0				
LK3-104901	EPA Method 525.2-Capillary Solution, 15-1				M-t-BE	1 ml, Amber Glass	€ 63,00
		a-chlordane	5103-71-9	1000			
		g-chlordane	5103-74-2	1000			
		chlorobenzilate	510-15-6	1000			
		chloroneb	2675-77-6	1000			
		Chlorothalonil	1897-45-6	1000			
		chlorpyrifos	2921-88-2	1000			
		bladex	21725-46-2	1000			
		tetrachloroterephthalic acid (DCPA Diacid)	2136-79-0	1000			
		etridiazole	2593-15-9	1000			
		hexachlorobenzene	118-74-1	1000			
		trans-nonachlor	39765-80-5	1000			
		cis-permethrin	54774-45-7	1000			
		trans-permethrin	51877-74-8	1000			
		Propachlor	1918-16-7	1000			
		trifluralin	1582-09-8	1000			
LK3-119505	Method 525.2 Organonitrogen Pesticide Mix, 37-1				Acetone	1 ml, Amber Glass	€ 292,50
		alachlor	15972-60-8	500			
		ametryne	834-12-8	500			
		Atraton	1610-17-9	500			
		atrazine	1912-24-9	500			
		bromacil	314-40-9	500			
		butachlor	23184-66-9	500			
		butylate	2008-41-5	500			
		chlorpropham	101-21-3	500			
		bladex	21725-46-2	500			

		Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		Cycloate	1134-23-2	500			
		diphenamid	957-51-7	500			
		EPTC (s-ethyl dipropylthiocarbamate)	759-94-4	500			
		etridiazole	2593-15-9	500			
		Fenarimol	60168-88-9	500			
		fluridone	59756-60-4	500			
		hexazinon	51235-04-2	500			
		metolachlor	51218-45-2	500			
		metribuzin	21087-64-9	500			
		MGK-264 Mix of Isomers	113-48-4	500			
		molinate	2212-67-1	500			
		napropamide	15299-99-7	500			
		norflurazon	27314-13-2	500			
		pebulate	1114-71-2	500			
		prometon	1610-18-0	500			
		prometryn	7287-19-6	500			
		propyzamide (pronamide)	23950-58-5	500			
		propachlor	1918-16-7	500			
		propazine	139-40-2	500			
		simazine	122-34-9	500			
		simetryn	1014-70-6	500			
		Tebuthiuron	34014-18-1	500			
		terbacil	5902-51-2	500			
		terbutryne	886-50-0	500			
		triadimefon	43121-43-3	500			
		beam™	41814-78-2	500			
		trifluralin	1582-09-8	500			
		vernolate	1929-77-7	500			
LK3-040101	Method 532 Mix 1	Diflubenzuron	35367-38-5	2000	Methanol	1 ml, Amber Glass	€ 66,00
		diuron	330-54-1	2000			
		fluometuron	2164-17-2	2000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		linuron	330-55-2	2000			
		Propanil	709-98-8	2000			
		Siduron	1982-49-6	2000			
LK3-040201	Method 532 Mix 2	Tebuthiuron	34014-18-1	2000	Acetone	1 ml, Amber Glass	€ 66,00
		Thidiazuron	51707-55-2	2000			
LK3-040001	Method 532 Surrogate Solution	carbazole	86-74-8	5000	Methanol	1 ml, Amber Glass	€ 51,00
		Monuron	150-68-5	5000			
LK3-013704	Method 552.2 Lab Performance Check Solution, 1mg/L				M-t-BE	1 ml, Amber Glass	€ 63,00
		methyl-2,3-dibromopropionate	1729-67-5	1			
		methyl bromochloroacetate	20428-74-4	1			
		methyl chloroacetate	96-34-4	1			
		methyl chlorodibromoacetate	20428-75-5	1			
LK3-252701	Method 614.1 Mix, 4-1	dioxathion (mixture of isomers)	78-34-2	1000	Acetone:Hexane 1:1	1 ml, Amber Glass	€ 37,50
		EPN	2104-64-5	1000			
		ethion	563-12-2	1000			
		Terbufos	13071-79-9	1000			
LK3-252801	*Method 619 - Triazine Pesticides Solution, 500 mg/L				Acetone	1 ml, Amber Glass	€ 135,00
LK3-252805	Method 619 - Triazine Pesticides Solution, 1000 mg/L					1 ml, Amber Glass	€ 213,00
		ametryne	834-12-8	1000			
		Atraton	1610-17-9	1000			
		atrazine	1912-24-9	1000			
		prometon	1610-18-0	1000			
		prometryn	7287-19-6	1000			
		propazine	139-40-2	1000			
		secbumethon	26259-45-0	1000			
		Simetryn	1014-70-6	1000			
		simazine	122-34-9	1000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		terbuthylazine	5915-41-3	1000			
		terbutryne	886-50-0	1000			
LK3-252904	*Method 619 - Triazine Pesticides Solution, 10-16				Acetone	1 ml, Amber Glass	€ 51,00
		ametryne	834-12-8	500			
		Atraton	1610-17-9	500			
		atrazine	1912-24-9	500			
		propazine	139-40-2	500			
		prometryn	7287-19-6	500			
		simetryn	1014-70-6	500			
		terbuthylazine	5915-41-3	500			
		terbutryne	886-50-0	500			
		prometon	1610-18-0	500			
		simazine	122-34-9	500			
LK3-100001	608 Pesticide Mix A	aldrin	309-00-2	2000	Tol.:Hex. 1:1	1 ml, Amber Glass	€ 79,50
		a-BHC	319-84-6	2000			
		b-BHC	319-85-7	2000			
		d-BHC	319-86-8	2000			
		g-BHC	58-89-9	2000			
		p,p'-DDD	72-54-8	2000			
		p,p'-DDE	72-55-9	2000			
		p,p'-DDT	50-29-3	2000			
		dieldrin	60-57-1	2000			
		endrin	72-20-8	2000			
		heptachlor	76-44-8	2000			
		heptachlor epoxide isomer B	1024-57-3	2000			
LK3-101202	608 Pesticide Mix, 13-2	p,p'-DDT	50-29-3	100	Isooctane:Toluene	1 ml, Amber Glass	€ 103,50
		o,p'-DDT	789-02-6	100			
		p,p'-DDE	72-55-9	100			
		o,p'-DDE	3424-82-6	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		p,p'-DDD	72-54-8	100			
		o,p'-DDD	53-19-0	100			
		a-BHC	319-84-6	100			
		b-BHC	319-85-7	100			
		g-BHC	58-89-9	100			
		d-BHC	319-86-8	100			
		aldrin	309-00-2	100			
		dieldrin	60-57-1	100			
		endrin	72-20-8	100			
LK3-001412	608 Pesticide Solution, 20 mg/L				Methanol	1 ml, Amber Glass	€ 63,00
LK3-001413	608 Pesticide Solution, 2000 mg/L				Methanol	1 ml, Amber Glass	€ 111,00
LK3-001407	608 Pesticide Solution, 10 mg/L				Methylene Chloric	1 ml, Amber Glass	€ 63,00
LK3-001402	608 Pesticide Solution, 100 mg/L				Methylene Chloric	1 ml, Amber Glass	€ 103,50
LK3-001408	608 Pesticide Check Solution, 40 mg/L				Methylene Chloric	1 ml, Amber Glass	€ 79,50
LK3-001410	608 Pesticide Solution, 200 mg/L				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 78,00
		aldrin	309-00-2				
		a-BHC	319-84-6				
		b-BHC	319-85-7				
		d-BHC	319-86-8				
		g-BHC	58-89-9				
		p,p'-DDD	72-54-8				
		p,p'-DDE	72-55-9				
		p,p'-DDT	50-29-3				
		dieldrin	60-57-1				
		endrin	72-20-8				
		heptachlor	76-44-8				
		heptachlor epoxide	1024-57-3				
		endrin aldehyde	7421-93-4				
		endosulfan I	959-98-8				
		endosulfan II	33213-65-9				
		endosulfan sulfate	1031-07-8				

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-200101	608 Pesticide Performance Check Solution				Methanol	1 ml, Amber Glass	€ 72,00
		aldrin	309-00-2	20			
		a-BHC	319-84-6	20			
		b-BHC	319-85-7	20			
		d-BHC	319-86-8	20			
		g-BHC	58-89-9	20			
		p,p'-DDD	72-54-8	100			
		p,p'-DDE	72-55-9	20			
		p,p'-DDT	50-29-3	100			
		dieldrin	60-57-1	20			
		endrin	72-20-8	20			
		heptachlor	76-44-8	20			
		heptachlor epoxide	1024-57-3	20			
		endrin aldehyde	7421-93-4	100			
		endosulfan I	959-98-8	20			
		endosulfan II	33213-65-9	100			
		endosulfan sulfate	1031-07-8	100			
		methoxychlor	72-43-5	20			
LK3-012201	Method 680 Retention Time Calibration Solution				Hexane	1 ml, Amber Glass	€ 24,00
		3,3',4,4'-tetrachlorobiphenyl (BZ# 77)	32598-13-3	2,5			
		2,2',4,6,6'-pentachlorobiphenyl (BZ# 104)	56558-16-8	2,5			
		2,3,3',4,4',5,5'-heptachlorobiphenyl (BZ# 189)	39635-31-9	2,5			
		2,2',3,3',5,5',6,6'-octachlorobiphenyl (BZ# 202)	2136-99-4	2,5			
LK3-110401	Method 8080 Organochlorine Pesticide Mix, 6-5				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 103,50
		chlorobenzilate	510-15-6	1000			
		1,2-dibromo-3-chloropropane	96-12-8	1000			
		hexachlorobenzene	118-74-1	1000			
		hexachlorocyclopentadiene	77-47-4	1000			
		isodrin	465-73-6	1000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		di-allyl	2303-16-4	1000			
LK3-110402	*8080 Organochlorine Pesticide Mix, 7-5				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 63,00
		chlorobenzilate	510-15-6	200			
		cis-di-allyl	17708-57-5	200			
		trans-di-allyl	17708-58-6	200			
		1,2-dibromo-3-chloropropane	96-12-8	200			
		hexachlorobenzene	118-74-1	200			
		hexachlorocyclopentadiene	77-47-4	200			
		isodrin	465-73-6	200			
LK3-306605	Custom Organochlorine Solution, Method 8080A/8081				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 118,50
		hexachlorobenzene	118-74-1	200			
		a-BHC	319-84-6	200			
		b-BHC	319-85-7	200			
		g-BHC	58-89-9	200			
		d-BHC	319-86-8	200			
		p,p'-DDE	72-55-9	200			
		p,p'-DDD	72-54-8	200			
		o,p'-DDT	789-02-6	200			
		p,p'-DDT	50-29-3	200			
		heptachlor	76-44-8	200			
		aldrin	309-00-2	200			
		heptachlor epoxide isomer B	1024-57-3	200			
		dieldrin	60-57-1	200			
		endrin	72-20-8	200			
LK3-110601	Method 8081 Organochlorine Pesticide Mix, 7-57				Tol.:Hex. 1:1	1 ml, Amber Glass	€ 63,00
		chlorobenzilate	510-15-6	1000			
		1,2-dibromo-3-chloropropane	96-12-8	1000			
		hexachlorobenzene	118-74-1	1000			
		hexachlorocyclopentadiene	77-47-4	1000			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		isodrin	465-73-6	1000			
		di-allate (mixture of isomers)	2303-16-4	1000			
		kepone	143-50-0	1000			
LK3-005603	Method 8081 Pesticide Mix, 100 mg/L , 8-42				Tol.:Isooct. 1:1	1 ml, Amber Glass	€ 79,50
LK3-005601	Method 8081 Pesticide Mix, 1,000 mg/L , 8-43					1 ml, Amber Glass	€ 213,00
		hexachlorobenzene	118-74-1	100			
		oxychlorane	27304-13-8	100			
		o,p'-DDE	3424-82-6	100			
		cis-nonachlor	5103-73-1	100			
		o,p'-DDD	53-19-0	100			
		trans-nonachlor	39765-80-5	100			
		o,p'-DDT	789-02-6	100			
		hexachlorobutadiene	87-68-3	100			
LK3-255401	Method 8081 Pesticide Mix, 9-42				Tol.:Isooct. 1:1	1 ml, Amber Glass	€ 213,00
		hexachlorobenzene	118-74-1	1000			
		hexachlorobutadiene	87-68-3	1000			
		o,p'-DDD	53-19-0	1000			
		cis-nonachlor	5103-73-1	1000			
		trans-nonachlor	39765-80-5	1000			
		mirex	2385-85-5	1000			
		oxychlorane	27304-13-8	1000			
		o,p'-DDT	789-02-6	1000			
		o,p'-DDE	3424-82-6	1000			
LK3-251203	Method 8141 Pesticides Mix, 4-14				Acetone:Hexane 1:1	1 ml, Amber Glass	€ 55,50
		dimethoate	60-51-5	100			
		malathion	121-75-5	100			
		parathion	56-38-2	100			
		thionazine (zinophos)	297-97-2	100			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
LK3-020508	8140/8141 Organophosphorous Pesticide Solution A, 40 mg/L				Hexane:Acetone: 1 ml, Amber Glass		€ 118,50
LK3-020501	8140/8141 Organophosphorous Pesticide Solution A, 200 mg/L				Hexane:Acetone: 1 ml, Amber Glass		€ 181,50
LK3-02050305	8140/8141 Organophosphorous Pesticide Solution A, 100 mg/L				Methylene Chloric 5 ml, Amber Glass		€ 748,50
		chlorpyrifos	2921-88-2				
		trichloronate	327-98-0				
		sulprofos (Bolstar®)	35400-43-2				
		coumaphos	56-72-4				
		Dementon O&S (31.3% O & 57.6% S)	8065-48-3				
		fensulfothion	115-90-2				
		fenthion	55-38-9				
		fenchlorphos	299-84-3				
		Tokuthion®	34643-46-4				
		tributylphosphoro-trithioite (Merphos)	150-50-5				
		tetrachlorvinphos (Rabon)	961-11-5				
		dibrom	300-76-5				
		diazinon	333-41-5				
		Guthion®	86-50-0				
		dichlorvos	62-73-7				
		disulfoton	298-04-4				
		ethoprophos (prophos)	13194-48-4				
		methyl parathion	298-00-0				
		mevinphos (Phosdrin (TM))	7786-34-7				
		phorate	298-02-2				
LK3-238001	8141 Organophosphorous Pesticide Solution B, 200 mg/L				Hexane:Acetone: 1 ml, Amber Glass		€ 118,50
LK3-23800305	8141 Organophosphorous Pesticide Solution B, 100 mg/L				Hexane	5 ml, Amber Glass	€ 189,00
		dimethoate	60-51-5				
		EPN	2104-64-5				
		malathion	121-75-5				
		monocrotophos	6923-22-4				
		parathion	56-38-2				
		tetraethyl dithiopyrophosphate	3689-24-5				

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		TEPP	107-49-3				
LK3-248001	Method 8141 Solution A	atrazine	1912-24-9	100	Acetone	1 ml, Amber Glass	€ 276,00
		Guthion Ethyl	2642-71-9	100			
		Guthion®	86-50-0	100			
		Carbophenothion	786-19-6	100			
		chlorfenvinphos (mixture of z and E isomers)	470-90-6	100			
		chlorpyrifos	2921-88-2	100			
		coumaphos	56-72-4	100			
		Demeton O&S (31.3% O & 57.6% S)	8065-48-3	100			
		diazinon	333-41-5	100			
		ethion	563-12-2	100			
		fenchlorphos	299-84-3	100			
		tributylphosphoro-trithioite (Merphos)	150-50-5	100			
		mevinphos (Phosdrin (TM))	7786-34-7	100			
		dibrom (Naled)	300-76-5	100			
		parathion	56-38-2	100			
		phosmet	732-11-6	100			
		phosphamidon	13171-21-6	100			
		simazine	122-34-9	100			
		sulfotep	3689-24-5	100			
		TEPP	107-49-3	100			
LK3-249001	Method 8141 Solution B	sulprofos (Bolstar®)	35400-43-2	100	Acetone	1 ml, Amber Glass	€ 315,00
		dichlorvos	62-73-7	100			
		dimethoate	60-51-5	100			
		dioxathion	78-34-2	100			
		disulfoton	298-04-4	100			
		EPN	2104-64-5	100			
		ethoprophos (prophos)	13194-48-4	100			
		famphur	52-85-7	100			
		fensulfothion	115-90-2	100			

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		fenthion	55-38-9	100			
		leptophos	21609-90-5	100			
		malathion	121-75-5	100			
		methyl parathion	298-00-0	100			
		monocrotophos	6923-22-4	100			
		phorate	298-02-2	100			
		tetrachlorvinphos (Rabon)	961-11-5	100			
		Terbufos	13071-79-9	100			
		Tokuthion®	34643-46-4	100			
		Tris(o-cresyl)phosphate	78-30-8	100			
		trichloronate	327-98-0	100			
		thionazine (zinophos)	297-97-2	100			
LK3-016605	Method 8141 Pesticides Mix	Guthion®	86-50-0	1000	M-t-BE:Acetone 8: 1 ml,	Amber Glass	€ 292,50
		carbophenothion	786-19-6	1000			
		chlorpyrifos	2921-88-2	1000			
		coumaphos	56-72-4	1000			
		diazinon	333-41-5	1000			
		dichlorvos	62-73-7	1000			
		dimethoate	60-51-5	1000			
		disulfoton	298-04-4	1000			
		ethion	563-12-2	1000			
		famphur	52-85-7	1000			
		malathion	121-75-5	1000			
		mevinphos (Phosdrin (TM))	7786-34-7	1000			
		methyl parathion	298-00-0	1000			
		phosmet	732-11-6	1000			
		Demeton O&S (24.4% O & 67.2% S)	8065-48-3	1000			
		ethoprophos (prophos)	13194-48-4	1000			
		fensulfothion	115-90-2	1000			
		monocrotophos	6923-22-4	2000			
		parathion	56-38-2	2000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		phorate	298-02-2	2000			
		atrazine	1912-24-9	3000			
LK3-110502	8270C Pesticides Mix, 6-1	aramite	140-57-8	2000	Methylene Chloric	1 ml, Amber Glass	€ 103,50
		chlorobenzilate	510-15-6	2000			
		di-allate	2303-16-4	2000			
		isodrin	465-73-6	2000			
		kepone	143-50-0	2000			
		propyzamide (pronamide)	23950-58-5	2000			
LK3-040301	8270 Organophosphorous Pesticide Mix 11-1				Methylene Chloric	1 ml, Amber Glass	€ 111,00
		alachlor	15972-60-8	1000			
		atrazine	1912-24-9	1000			
		Guthion®	86-50-0	1000			
		carbophenothion	786-19-6	1000			
		chlorpyrifos	2921-88-2	1000			
		diazinon	333-41-5	1000			
		dichlorvos	62-73-7	1000			
		dimethoate	60-51-5	1000			
		malathion	121-75-5	1000			
		methyl parathion	298-00-0	1000			
		vernolate	1929-77-7	1000			
LK3-110302	8270C Calibration Mix, 11-2	di-allate (mixture of isomers)	2303-16-4	2000	Methylene Chloric	1 ml, Amber Glass	€ 181,50
		dimethoate	60-51-5	2000			
		disulfoton	298-04-4	2000			
		famphur	52-85-7	2000			
		methyl parathion	298-00-0	2000			
		parathion	56-38-2	2000			
		phorate	298-02-2	2000			
		propyzamide (pronamide)	23950-58-5	2000			
		thionazine (zinophos)	297-97-2	2000			

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		O,O,O-triethylphosphorothioate	126-68-1	2000			
LK3-040901	8270 Organophosphorous Pesticide Mix, 17-16					Methylene Chloric 1 ml, Amber Glass	€ 174,00
		Guthion®	86-50-0	1000			
		atrazine	1912-24-9	1000			
		vernolate	1929-77-7	1000			
		tetrachlorvinphos (Rabon)	961-11-5	1000			
		alachlor	15972-60-8	1000			
		carbophenothion	786-19-6	1000			
		dichlorvos	62-73-7	1000			
		phorate	298-02-2	1000			
		dimethoate	60-51-5	1000			
		diazinon	333-41-5	1000			
		methyl parathion	298-00-0	1000			
		fenchlorphos	299-84-3	1000			
		malathion	121-75-5	1000			
		chlorpyrifos	2921-88-2	1000			
		p-terphenyl-d14	1718-51-0	1000			
		2-fluorobiphenyl	321-60-8	1000			
		nitrobenzene-d5	4165-60-0	1000			
LK3-001498-7pck	8270 Pesticide Calibration Set, 10, 20, 40, 50, 80, 100, 120 mg/L					Methylene Chloric 1 ml, Amber Glass	€ 237,00
		aldrin	309-00-2				
		a-BHC	319-84-6				
		b-BHC	319-85-7				
		d-BHC	319-86-8				
		g-BHC	58-89-9				
		p,p'-DDD	72-54-8				
		p,p'-DDE	72-55-9				
		p,p'-DDT	50-29-3				
		dieldrin	60-57-1				
		endrin	72-20-8				

BR=Boston Round

Catalog Number	Description	Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
		heptachlor	76-44-8				
		heptachlor epoxide	1024-57-3				
		endrin aldehyde	7421-93-4				
		endosulfan I	959-98-8				
		endosulfan II	33213-65-9				
		endosulfan sulfate	1031-07-8				
LK3-028502	Custom 8321 Internal Standard Pesticide Mix, 100 mg/L , 4-8					Methylene Chloric 1 ml, Amber Glass	€ 40,50
LK3-028502-10	Custom 8321 Internal Standard Pesticide Mix, 100 mg/L , 4-8					Methylene Chloric 10 x 1 ml, Amber	€ 237,00
LK3-028501	Custom 8321 Internal Standard Pesticide Mix, 1,000 mg/L , 4-8					Methylene Chloric 1 ml, Amber Glass	€ 55,50
LK3-028501-5pck	Custom 8321 Internal Standard Pesticide Mix, 1,000 mg/L , 4-8				M-t-BE	5 x 1 ml, Amber	€ 198,00
		Metalaxyl	57837-19-1				
		monocrotophos	6923-22-4				
		benthiocarb	28249-77-6				
		dioxacarb	6988-21-2				
LK3-028402	Custom 8321 Pesticide Mix, 100, 200 & 500 mg/L , 29-8,					Methylene Chloric 1 ml, Amber Glass	€ 103,50
LK3-028401	Custom 8321 Pesticide Mix, 1,000, 2,000 & 5,000 mg/L , 29-8,					1 ml, Amber Glass	€ 198,00
		aldicarb	116-06-3	100			
		aldicarb sulfone	1646-88-4	100			
		aminocarb	2032-59-9	100			
		benomyl	17804-35-2	100			
		bendiocarb	22781-23-3	100			
		bromacil	314-40-9	100			
		Carbaryl	63-25-2	100			
		carbofuran	1563-66-2	100			
		chloroxuron	1982-47-4	100			
		diuron	330-54-1	100			
		Fenarimol	60168-88-9	100			
		fenuron	101-42-8	100			
		fluometuron	2164-17-2	100			
		linuron	330-55-2	100			



Catalog Number Description

BR=Boston Round

Compound	CAS#	Concentr. mg/L	Matrix	Container	List Price
Methiocarb	2032-65-7	100			
methomyl	16752-77-5	100			
Monuron	150-68-5	100			
neburon	555-37-3	100			
oryzalin	19044-88-3	100			
oxamyl	23135-22-0	100			
propachlor	1918-16-7	100			
propoxur	114-26-1	100			
Siduron	1982-49-6	100			
Tebuthiuron	34014-18-1	100			
aldicarb sulfoxide	1646-87-3	100			
chlorpropham	101-21-3	200			
mexacarbate	315-18-4	200			
propham	122-42-9	200			
barban	101-27-9	500			