

## Safety data sheet according to 1907/2006/EC, Article 31

Printing date 13.01.2023

Version number 1

Revision: 13.01.2023

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

· **1.1 Product identifier**

· **Trade name:** Rubidium Rb - 1000 mg/l in diluted HNO<sub>3</sub> for ICP (LK1-00370205)

· **Registration number**

A registration number is not available for this substance as the substance or its uses are exempted for registration, the annual tonnage does not require a registration or the registration is envisaged for a later registration deadline.

· **1.2 Relevant identified uses of the substance or mixture and uses advised against**

No further relevant information available.

· **Application of the substance / the mixture** Laboratory Reagent

· **1.3 Details of the supplier of the safety data sheet**

· **Manufacturer/Supplier:**

CPAchem Ltd.  
2 Ivanka Terzieva Str.  
Bogomilovo 6065  
Stara Zagora, BULGARIA  
info@cpachem.com

LabKings  
Utrechtseweg 5,  
1213 TK Hilversum  
The Netherlands  
info@labkings.com

· **Further information obtainable from:** Product safety department

· **1.4 Emergency telephone number:**

**EMERGENCY HEALTH INFORMATION:**

Austria +43 1 31304 5620, Belgium +32022649636, Bulgaria +359 2 9154 409, Croatia +38514686910, Cyprus +3572240561, Czech Republic +420267082257, Denmark +45 72 54 40 00, Estonia +3726943384, Finland +358 5052 000, France +33 3 85 21 92, Germany +49-30-18412-0, Greece +302106479250, Hungary +34 (1) 476 1136, Ireland +35318092566, Italy +390649906140, Latvia +371 67032600, Lithuania +370 70662008, Luxembourg +352 24785551, Netherland +31 88 75 585 61, Norway +47 21 07 70 00, Poland +48 42 2530 400, Portugal +351213303271, Romania +40213183606, Slovakia +421 2 5465 2307, Slovenia +38614006039, Spain +34 917689800, Sweden +46104566750, United Kingdom (England or Wales) 0845 46 47 or Scotland 08454 24 24 24 (UK only).

### SECTION 2: Hazards identification

· **2.1 Classification of the substance or mixture**

· **Classification according to Regulation (EC) No 1272/2008**



GHS07

Skin Irrit. 2 H315 Causes skin irritation.

Eye Irrit. 2 H319 Causes serious eye irritation.

· **2.2 Label elements**

· **Labelling according to Regulation (EC) No 1272/2008**

The product is classified and labelled according to the CLP regulation.

· **Hazard pictograms**



GHS07

· **Signal word** Warning

· **Hazard statements**

H315 Causes skin irritation.

(Contd. on page 2)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 13.01.2023




Version number 1

Revision: 13.01.2023

**Trade name: Rubidium Rb - 1000 mg/l in diluted HNO<sub>3</sub> for ICP (LK1-00370205)**

(Contd. of page 1)

*H319 Causes serious eye irritation.***· Precautionary statements***P264 Wash thoroughly after handling.**P280 Wear eye protection / face protection.**P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.**P332+P313 If skin irritation occurs: Get medical advice/attention.**P362+P364 Take off contaminated clothing and wash it before reuse.**P337+P313 If eye irritation persists: Get medical advice/attention.***· 2.3 Other hazards****· Results of PBT and vPvB assessment****· PBT:** Not applicable.**· vPvB:** Not applicable.**SECTION 3: Composition/information on ingredients****· 3.2 Mixtures****· Description:** Mixture: consisting of the following components.**· Dangerous components:**

|                            |  |      |
|----------------------------|--|------|
| CAS: 7697-37-2             | nitric acid  | 2.0% |
| EINECS: 231-714-2          |  Ox. Liq. 3, H272;  Acute Tox. 3, H331;  Skin Corr. 1A, H314; |      |
| Index number: 007-030-00-3 | Eye Dam. 1, H318, EUH071   |      |
|                            | ATE: LC50/4 h inhalative: > 20 mg/l  |      |
|                            | Specific concentration limits: Ox. Liq. 3; H272: C ≥ 65 %  |      |
|                            | Skin Corr. 1A; H314: C ≥ 20 %  |      |
|                            | Skin Corr. 1B; H314: 5 % ≤ C < 20 %  |      |

**· Additional information:** For the wording of the listed hazard phrases refer to section 16.**SECTION 4: First aid measures****· 4.1 Description of first aid measures****· After inhalation:** Supply fresh air; consult doctor in case of complaints.**· After skin contact:** Generally the product does not irritate the skin.**· After eye contact:** Rinse opened eye for several minutes under running water.**· After swallowing:** If symptoms persist consult doctor.**· 4.2 Most important symptoms and effects, both acute and delayed** No further relevant information available.**· 4.3 Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

**SECTION 5: Firefighting measures****· 5.1 Extinguishing media****· Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.**· 5.2 Special hazards arising from the substance or mixture** No further relevant information available.**· 5.3 Advice for firefighters****· Protective equipment:** No special measures required.**SECTION 6: Accidental release measures****· 6.1 Personal precautions, protective equipment and emergency procedures** Not required.**· 6.2 Environmental precautions:** Do not allow to enter sewers/ surface or ground water.**· 6.3 Methods and material for containment and cleaning up:**

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

(Contd. on page 3)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 13.01.2023

Version number 1

Revision: 13.01.2023

**Trade name: Rubidium Rb - 1000 mg/l in diluted HNO<sub>3</sub> for ICP (LK1-00370205)**

(Contd. of page 2)

### · 6.4 Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

## SECTION 7: Handling and storage

· **7.1 Precautions for safe handling** No special precautions are necessary if used correctly.

· **Information about fire - and explosion protection:** No special measures required.

· **7.2 Conditions for safe storage, including any incompatibilities**

· **Storage:**

· **Requirements to be met by storerooms and receptacles:** No special requirements.

· **Information about storage in one common storage facility:** Not required.

· **Further information about storage conditions:** None.

· **7.3 Specific end use(s)** No further relevant information available.

## SECTION 8: Exposure controls/personal protection

### · 8.1 Control parameters

· **Ingredients with limit values that require monitoring at the workplace:**

7697-37-2 nitric acid

IOELV Short-term value: 2.6 mg/m<sup>3</sup>, 1 ppm

· **Additional information:** The lists valid during the making were used as basis.

### · 8.2 Exposure controls

· **Appropriate engineering controls** No further data; see item 7.

· **Individual protection measures, such as personal protective equipment**

· **General protective and hygienic measures:** Wash hands before breaks and at the end of work.

· **Respiratory protection:** Not required.

· **Hand protection**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye/face protection** Goggles recommended during refilling

## SECTION 9: Physical and chemical properties

### · 9.1 Information on basic physical and chemical properties

· **General Information**

· **Physical state**

Fluid

· **Colour:**

According to product specification

· **Odour:**

Characteristic

· **Odour threshold:**

Not determined.

(Contd. on page 4)

**Safety data sheet**  
according to 1907/2006/EC, Article 31

Printing date 13.01.2023

Version number 1

Revision: 13.01.2023

Trade name: Rubidium Rb - 1000 mg/l in diluted HNO<sub>3</sub> for ICP (LK1-00370205)

(Contd. of page 3)

|  |   |
|--|---|
| · <b>Melting point/freezing point:</b>   | Undetermined.                                 |
| · <b>Boiling point or initial boiling point and boiling range</b>                      | Undetermined.                                 |
| · <b>Flammability</b>  | Not applicable.                               |
| · <b>Lower and upper explosion limit</b>   |   |
| · <b>Lower:</b>  | Not determined.                               |
| · <b>Upper:</b>  | Not determined.                               |
| · <b>Flash point:</b>  | Not applicable.                               |
| · <b>Decomposition temperature:</b>  | Not determined.                               |
| · <b>pH</b>  | Not determined.                               |
| · <b>Viscosity:</b>  |   |
| · <b>Kinematic viscosity</b>   | Not determined.                               |
| · <b>Dynamic:</b>  | Not determined.                               |
| · <b>Solubility</b>  |   |
| · <b>water:</b>  | Not miscible or difficult to mix.             |
| · <b>Partition coefficient n-octanol/water (log value)</b>                             | Not determined.                               |
| · <b>Vapour pressure:</b>  | Not determined.                               |
| · <b>Density and/or relative density</b>   |   |
| · <b>Density:</b>  | Not determined.                               |
| · <b>Relative density</b>  | Not determined.                               |
| · <b>Vapour density</b>  | Not determined.                               |
| · <b>9.2 Other information</b>   |   |
| · <b>Appearance:</b>   |   |
| · <b>Form:</b>   | Fluid   |
| · <b>Important information on protection of health and environment, and on safety.</b> |   |
| · <b>Auto-ignition temperature:</b>  | Product is not selfigniting.                  |
| · <b>Explosive properties:</b>   | Product does not present an explosion hazard. |
| · <b>Solvent content:</b>  |   |
| · <b>Organic solvents:</b>   | 0.0 %   |
| · <b>Solids content:</b>   | 0.2 %   |
| · <b>Change in condition</b>   |   |
| · <b>Evaporation rate</b>  | Not determined.                               |
| · <b>Information with regard to physical hazard classes</b>                            |   |
| · <b>Explosives</b>  | Void  |
| · <b>Flammable gases</b>   | Void  |
| · <b>Aerosols</b>  | Void  |
| · <b>Oxidising gases</b>   | Void  |
| · <b>Gases under pressure</b>  | Void  |
| · <b>Flammable liquids</b>   | Void  |
| · <b>Flammable solids</b>  | Void  |
| · <b>Self-reactive substances and mixtures</b>   | Void  |
| · <b>Pyrophoric liquids</b>  | Void  |
| · <b>Pyrophoric solids</b>   | Void  |
| · <b>Self-heating substances and mixtures</b>  | Void  |
| · <b>Substances and mixtures, which emit flammable gases in contact with water</b>     | Void  |
| · <b>Oxidising liquids</b>   | Void  |
| · <b>Oxidising solids</b>  | Void  |
| · <b>Organic peroxides</b>   | Void  |
| · <b>Corrosive to metals</b>   | Void  |
| · <b>Desensitised explosives</b>   | Void  |

## SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.

(Contd. on page 5)

# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 13.01.2023

Version number 1

Revision: 13.01.2023

**Trade name: Rubidium Rb - 1000 mg/l in diluted HNO<sub>3</sub> for ICP (LK1-00370205)**

(Contd. of page 4)

- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known.
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

### SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **Skin corrosion/irritation** Causes skin irritation.
- **Serious eye damage/irritation** Causes serious eye irritation.
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **11.2 Information on other hazards**

- **Endocrine disrupting properties**

None of the ingredients is listed.

### SECTION 12: Ecological information

- **12.1 Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability** No further relevant information available.
- **12.3 Bioaccumulative potential** No further relevant information available.
- **12.4 Mobility in soil** No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **12.6 Endocrine disrupting properties**
- **The product does not contain substances with endocrine disrupting properties.**
- **12.7 Other adverse effects**
- **Additional ecological information:**
- **General notes:**
- **Water hazard class 1 (German Regulation) (Self-assessment):** slightly hazardous for water
- **Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.**

### SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation**
- **Must not be disposed together with household garbage. Do not allow product to reach sewage system.**
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

EU

(Contd. on page 6)

**Safety data sheet**  
according to 1907/2006/EC, Article 31

Printing date 13.01.2023


Version number 1

Revision: 13.01.2023

Trade name: Rubidium Rb - 1000 mg/l in diluted HNO<sub>3</sub> for ICP (LK1-00370205)

(Contd. of page 5)

**SECTION 14: Transport information**

|   |   |
|---|---|
| · 14.1 UN number or ID number<br>· ADR, IMDG, IATA  | UN1760  |
| · 14.2 UN proper shipping name<br>· ADR<br>· IMDG, IATA   | 1760 CORROSIVE LIQUID, N.O.S. (NITRIC ACID)<br>CORROSIVE LIQUID, N.O.S. (NITRIC ACID)   |
| · 14.3 Transport hazard class(es)<br>· ADR, IMDG, IATA  | <br>· Class<br>· Label<br>8 Corrosive substances.<br>8 |
| · 14.4 Packing group<br>· ADR, IMDG, IATA   | III   |
| · 14.5 Environmental hazards:<br>· Marine pollutant:  | No  |
| · 14.6 Special precautions for user<br>· Hazard identification number (Kemler code):<br>· EMS Number:<br>· Segregation groups<br>· Stowage Category<br>· Stowage Code | Warning: Corrosive substances.<br>80<br>F-A,S-B<br>(SGG1) Acids<br>A<br>SW2 Clear of living quarters.                                   |
| · 14.7 Maritime transport in bulk according to IMO instruments  | Not applicable.   |
| · Transport/Additional information:   |   |
| · ADR<br>· Limited quantities (LQ)<br>· Excepted quantities (EQ)  | 5L<br>Code: E1<br>Maximum net quantity per inner packaging: 30 ml<br>Maximum net quantity per outer packaging: 1000 ml                  |
| · Transport category<br>· Tunnel restriction code   | 3<br>E  |
| · IMDG<br>· Limited quantities (LQ)<br>· Excepted quantities (EQ)   | 5L<br>Code: E1<br>Maximum net quantity per inner packaging: 30 ml<br>Maximum net quantity per outer packaging: 1000 ml                  |
| · UN "Model Regulation":  | UN 1760 CORROSIVE LIQUID, N.O.S. (NITRIC ACID), 8, III  |

**SECTION 15: Regulatory information**

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture  
 · Labelling according to Regulation (EC) No 1272/2008  
 The product is classified and labelled according to the CLP regulation.

(Contd. on page 7)

EU



# Safety data sheet

## according to 1907/2006/EC, Article 31

Printing date 13.01.2023

Version number 1

Revision: 13.01.2023

Trade name: Rubidium Rb - 1000 mg/l in diluted HNO<sub>3</sub> for ICP (LK1-00370205)

(Contd. of page 6)

· **Hazard pictograms**

GHS07

· **Signal word** Warning· **Hazard statements**

H315 Causes skin irritation.

H319 Causes serious eye irritation.

· **Precautionary statements**

P264 Wash thoroughly after handling.

P280 Wear eye protection / face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P332+P313 If skin irritation occurs: Get medical advice/attention.

P362+P364 Take off contaminated clothing and wash it before reuse.

P337+P313 If eye irritation persists: Get medical advice/attention.

· **Directive 2012/18/EU**· **Named dangerous substances - ANNEX I** None of the ingredients is listed.· **REGULATION (EC) No 1907/2006 ANNEX XVII** Conditions of restriction: 3· **DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II**

None of the ingredients is listed.

· **REGULATION (EU) 2019/1148**· **Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3))**

None of the ingredients is listed.

· **Annex II - REPORTABLE EXPLOSIVES PRECURSORS**

None of the ingredients is listed.

· **Regulation (EC) No 273/2004 on drug precursors**

None of the ingredients is listed.

· **Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors**

None of the ingredients is listed.

· **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Relevant phrases**

H272 May intensify fire; oxidiser.

H314 Causes severe skin burns and eye damage.

H318 Causes serious eye damage.

H331 Toxic if inhaled.

EUH071 Corrosive to the respiratory tract.

· **Department issuing SDS:** Product safety department· **Contact:** Mrs. Taralova· **Date of previous version:** 24.02.2022· **Abbreviations and acronyms:**

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

(Contd. on page 8)

**Safety data sheet**  
**according to 1907/2006/EC, Article 31**

Printing date 13.01.2023

Version number 1

Revision: 13.01.2023

**Trade name: Rubidium Rb - 1000 mg/l in diluted HNO<sub>3</sub> for ICP (LK1-00370205)**

(Contd. of page 7)

*ICAO: International Civil Aviation Organisation**ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)**IMDG: International Maritime Code for Dangerous Goods**IATA: International Air Transport Association**GHS: Globally Harmonised System of Classification and Labelling of Chemicals**EINECS: European Inventory of Existing Commercial Chemical Substances**ELINCS: European List of Notified Chemical Substances**CAS: Chemical Abstracts Service (division of the American Chemical Society)**PBT: Persistent, Bioaccumulative and Toxic**vPvB: very Persistent and very Bioaccumulative**Ox. Liq. 3: Oxidizing liquids – Category 3**Acute Tox. 3: Acute toxicity – Category 3**Skin Corr. 1A: Skin corrosion/irritation – Category 1A**Skin Irrit. 2: Skin corrosion/irritation – Category 2**Eye Dam. 1: Serious eye damage/eye irritation – Category 1**Eye Irrit. 2: Serious eye damage/eye irritation – Category 2*

—EU—