

Ice Pack Required?

Quote Specification Sheet

Page 1.000

Catalog No: LK2- CB101842	Storage: -18°C +/- 4°C	Solvent: Acetone	Container: 1 mL Ampule, Amber Glass	Description: Nicosulfuron Solution, 50 ug/mL, 1 ml
--	----------------------------------	----------------------------	--	--

*Estimated Days to Expiration: 730

Compound

CAS No

Concentration

nicosulfuron

111991-09-4

50 +/- 2.5 ug/mL

*We will attempt to prepare solutions as requested, however, on new solutions please note concentrations and/or solvent may change due to chemistry or solubility issues. The expiration period on custom solutions is estimated based on available chemical knowledge and/or data. The actual expiration period may change once a stability study is completed on this particular product, which may take up to three years. Please review the specifications carefully as the product will be made as listed above. By ordering this product, you are approving these specifications. Client signature may still be requested.

Client Specifications Approval: _____

Date: _____